

09 How do I review readings?


Ref: CORP-059-01 Version 1 7-Dec-20

Task: View readings in a graph view

Previous actions taken: Activate DBm-Health app / Open the app / Submit readings

Step 1

From the hub screen, tap on the **Complete readings** circle and you will get to the graphical view.



Step 2


On this graph you can view your readings. The dots are as follows:

- Red:** High reading
- Green:** Reading within normal limits
- Blue:** Low reading

Step 3

The graph can be filtered by meal type.

To do this, click on the **Filter** icon in the top right-hand side of the screen.



Step 4

Select the meal type which you would like to view the details.

You will then see the graph filtered by requested meal type.

Task: View readings in a diary format

Previous actions taken: Activate DBm-Health app / Open the app / Submit readings

Step 1

From the hub screen, tap on the **Complete readings** circle and you will get to the graphical view.

NOTE: Before proceeding to **Step 2**, ensure your **Screen Orientation** setting is unlocked on your phone.

Step 2

You will then see a message which states **Rotate phone for diary view**.


If you do this you will see your readings in a diary format.



	BREAKFAST		LUNCH		EVENING MEAL	
	BEFORE	AFTER	BEFORE	AFTER	BEFORE	AFTER
16/12/18	4.8	---	4.7	8.1	5.0	8.1
	08:30	12:30	12:30	18:30	18:30	18:54
UNTAGGED 8.5 at 12:38 Wed 8 Nov, 2 units Nov						
16/12/18	5.0	8.0	5.7	7.6	5.4	8.6
	08:30	09:00	12:30	13:30	18:30	19:00
18/12/18	3.7	6.4	5.6	---	---	7.6
	08:30	09:00	12:30	13:30	19:30	19:00
18/12/18	4.2	7.9	5.7	7.7	5.5	7.4
	08:30	09:00	12:30	13:30	19:30	19:00

Step 3

Once in this view, you can look at each reading in more detail. Tap on the reading you are interested in and more detail will appear.



	BEFORE	AFTER
16/12/18	4.8	8.1
	08:30	18:54
16/12/18	5.0	8.6
	08:30	19:00
18/12/18	3.7	7.6
	08:30	19:00
18/12/18	4.2	7.4
	08:30	19:00

Comment

LAST UPDATED WED 15 NOV 2017 @ 08:30 PM

Didn't eat very much this morning as wasn't feeling great.

MEDICATION

Metformin 2 units

Novorapid 5 units

Step 4

To exit this screen press on the **X** in the top right corner of the screen.

To leave the diary view, rotate your screen back to portrait