ASI Evolution

Automated Nontreponemal RPR Syphilis Analyzer

FDA CLEARED



For Diagnostic, Blood Donor and Cadaveric Tissue Screening


- Reduce Labor
- Improve Efficiency
- Better Quality
- Lower Costs

ARLINGTON SCIENTIFIC

Your Syphilis Authority®

Arlington Scientific, Inc.® (ASI®) is an industry leader in the development, manufacturing and commercialization of *in vitro* diagnostic products with an emphasis in syphilis. For over 36 years ASI has been a pioneer in delivering new products and technologies to meet the growing demands of diagnostic, blood donor and cadaveric screening laboratories worldwide.


Syphilis Screening Algorithms



Nontreponemal vs. Treponemal

A 2015 CDC study determined that nontreponemal tests can detect infection up to 14 days earlier than treponemal tests.

Courtesy of Centers for Disease Control and Prevention



The Automated RPR Advantage

Clinical Value

- Standardized methodology
- Fewer false-positives
- Detect, diagnose and monitor therapy
- Uses the Nontreponemal (RPR) screening method
- Accurate and reproducible results

Operational Value

- · Walk-a-way system
- Run 190 tests/hour/analyzer
- Titers up to 1:2048
- Single operator can run multiple ASI Evolutions simultaneously
- Data Management (Analyze, Archive, Retrieve)

Financial Value

- · Lowest cost syphilis automated test
- · Minimal consumables required
- · Lower total healthcare cost
- Reduce labor expense

The ASI Evolution Automated Syphilis Analyzer

The ASI Evolution® is the world's first and only fully-automated nontreponemal RPR syphilis analyzer, FDA cleared for diagnostic and blood donor screening and cadaveric (Non-heart beating) tissue screening. Simply load your reagents and samples, then walk-away. The Evolution automates the processing, analysis, reporting and archiving of results for RPR screens

End Point Titers 1:2048

The ASI Evolution® uses a proprietary software algorithm that analyzes the agglutination of each microtiter well and reports a qualitative test result (reactive or nonreactive). In addition the ASI Evolution can serially dilute samples up to 1:2048 for each reactive microtiter well and reports semiquantitative endpoint titers for those wells.

Automation Time Saving Advantage

The ASI Evolution reduces the amount of hands on time required to run an RPR test; freeing up valuable lab resources.



In a study comparing the hands-on time required to run 900 tests on the ASI Evolution vs manually, it was found that the ASI Evolution required 92% less hands-on time vs the manual process. A lab tech using 30 well cards was able to perform 900 samples manually in 3 hours, 2 minutes and 3 seconds. For the Evolution®, samples were placed into 5 racks (4x192 + 1x132) prior to running the samples. Total hands-on time to run 5 racks of samples was 15 minutes and 28 seconds, 1/12 the hands-on time required to perform the test manually.

As a fully automated system, the ASI Evolution can process and analyze up to 190 samples per hour. One technician can easily run three Evolutions simultaneously,



Digital Results

The ASI Evolutions HD camera, allows the proprietary algorithm to digitally analyze particles to determine reactivity.





The Nontreponemal Advantage

Competitive Matrix				
	Trep	Nontrep		
Can monitor response to therapy		\checkmark		
Lowest cost per result		\checkmark		
Can provide titers of sample		\checkmark		
Can be automated for high thruput	✓	\checkmark		
Earliest detection of infection method		✓		
Recommended for prenatal screening*		√		
Determines active infection		√		

^{*}American Academy of Pediatrics, American College of Obstetricians and Gynecologists, and March of Dimes Birth Defects Foundation. Guidelines for Perinatal Care. 6th ed. Elk Grove Village, III.: American Academy of Pediatrics; 2008.

https://www.aafp.org/afp/2010/0115/od1.html#afp20100115od1-b3

ASI Evolution® System (For Technical Specs go to: http://www.arlingtonscientific.com/evolutionspecs.html)

ASI Evolution® Automated Syphilis Analyzer	1	Operations Manual
Laptop	_ 1	Power Cord
ASI Evolution Sample Racks	_ 3	Data Cord
Bar Code Reader	_ 1	

Ordering Information

ASI Evolution® Automated Syphilis Analyzer	#2800-000	ASI Automated RPR Control Set 2.5ml	_#905002.5A
ASI Automated 480 Test Kit	_ #900480A	ASI Automated RPR Control Set 5ml	_#905005A
ASI Automated 768 Test Kit	_ #900768A	ASI Evolution® Waste System (Stand and waste container)	#2800-019
ASI Automated 4800 Test Kit	_ #9004800A	ASI Evolution Sample Rack —	_ #2800-021



(800) 654-0146 (801) 489-8911
Fax (801) 489-5552
Springville, UT, USA
info@arlingtonscientific.com
www.arlingtonscientific.com, www.automatedrpr.com, www.asievolution.com