

FLOUR DISCRIMINATION FOR WAFER PRODUCTION

APPLICATION FAQs – 2019-037

Industrial processes usually do not appreciate changes in raw material quality. There is a strong need for controlling the properties of wheat flour to ensure it will perform well during processing and allow for a consistent final product which meets customer specifications.

The study presented below was done in partnership with an industrial wafer producer. The objective is to discriminate 2 “good” flours from a “medium” flour and a “bad” flour. The study demonstrates that the Alveolab complies with this objective.



Q1. ARE QUANTITATIVE ANALYSES ENOUGH TO DISCRIMINATE BETWEEN THE SAMPLES?

NO. Methods such as protein content or Falling Number are not sufficient to segregate good and bad flour (Table 1).

Q2. ARE THE STANDARD ALVEOGRAPHIC RESULTS DISCRIMINATING THE SAMPLES?

YES. The standard Alveographic curves allow discrimination between the “bad” sample and the others. Nevertheless, the “medium” sample has a very similar behavior to one of the good flours (Figure 1). The Alveolab integrates a new feature allowing one to measure dough consistency during mixing. Using these tools it was possible to discriminate the “medium” sample as the one having the longest tolerance to mixing (Figures 2 and 3)

Table 1: Quantitative results and process behavior for the samples to discriminate

Sample	Quality	Moisture content	Protein content	Falling Number	Comments
1	Good	14.5%	11.35%	342 s	Standard production
2	Good	14.1%	10.70%	411 s	Standard production
3	Medium	14.3%	10.83%	336 s	Gluten network too tenacious during mixing
4	Bad	13.2%	10.30%	395 s	Liquid dough, wafers are cracked or even broken




Figure 1: Standard Alveographic curves




Figure 2: Mixing behavior




Figure 3: Tolerance during mixing

This document is a part of a more comprehensive study which is accessible upon request at: labo.application@chopin.fr

SERVICES FOR CEREAL INDUSTRY PROFESSIONALS



AT YOUR SERVICE :

- **Technical support** for development of protocols and calibrations **adapted to your needs**
- a fully equipped laboratory to conduct your own studies and researches.
- customized training sessions **to fine tune** your **maintenance skills, analysis and data interpretation.**

We provide you with these services.

Contact us at labo.application@chopin.fr