



Graver Technologies

Citadel™ Series
All-fluoropolymer
Filter Cartridges

Performance
Qualification Guide

ISO 9001:2015

Preface

Each section of this Validation Guide represents only the summary portion of the actual test. If your company has a need for expanded detail on any particular test method or the actual data, please contact Graver Technologies Liquid Filter Group for assistance at 800-249-1990.



Table of Contents

| | |
|--------------------------------------|----|
| <i>Preface</i> | 2 |
| <i>Table of Contents</i> | 3 |
| Introduction | 4 |
| Nomenclature & Construction | 5 |
| Product Traceability | 6 |
| Cartridge Integrity Test | 7 |
| Flow Rate Test | 8 |
| Cartridge Extraction | 10 |
| Total Organic Carbon | |
| Non-volatile Residue | |
| Particle Shedding | |
| Resistivity Recovery | |
| Particle Retention | 11 |
| Leachables (Metals) | 12 |
| 5% Hydrochloric Acid | 13 |
| Ultra Pure Water (UPW) | 15 |
| Endotoxin Test | 17 |
| European Regulation No 1935/2004 and | 18 |
| European Regulation 10/2011 | |
| Notes | 19 |

Introduction

The Graver Technologies CITADEL™ (0.05 µm, 0.1 µm, 0.2 µm, 0.45 µm, 1.0 µm and 10 µm) pleated filter cartridges are designed as an exceptionally clean, non-leaching, non-shedding barrier for membrane filtration. These filters offer reliable performance in removing minute contaminants larger than their rated pore sizes. All CITADEL filter cartridges incorporate PTFE (polytetrafluoroethylene) membrane with fluoropolymer components (cage, core, end caps and support layers), either PFA (perfluoroalkoxy), PTFE or FEP (Fluorinated ethylene propylene). CITADEL filters are manufactured in an ISO Class 7 cleanroom environment. Each filter is rinsed with 18 megohm-cm resistivity de-ionized water to single digit TOC readings and integrity tested before release from manufacture. The Citadel filter products are fabricated in an ISO 9001, Rev. 2015 Registered manufacturing facility.

This report contains results of laboratory tests performed on Graver Technologies' CITADEL membrane cartridge filters.

Nomenclature & Construction

| CITADEL NOMENCLATURE INFORMATION | | | | | | | |
|----------------------------------|----------------------------|------|-------------------------|-----------------------------|------------------|--------------------------------|------------|
| Filter Type | Retention Rating (microns) | | Nominal Length (inches) | End Configuration | Gasket or O-Ring | | Option |
| Citadel Series Filters | 0.05 | 0.45 | -10 | P2 226/Flat Single Open End | C | Chemraz | -W Pre-Wet |
| | 0.1 | 1 | -20 | | K | Kalrez | |
| | 0.2 | 10 | -30 | P3 222/Flat Single Open End | T | Teflon encap. Viton (Standard) | |
| | | | -40 | | | | |
| Example: CTL 0.05-20P3T-W | | | | | | | |
| CTL | 0.05 | | -20 | P3 | T | | -W |

Materials of Construction

Membrane: Single layer expanded polytetrafluoroethylene (e-PTFE) – 8.6 ft² (0.8 m²)
 Drainage Layer: FEP
 Core: PFA
 Cage/Outer Sleeve: PFA
 End Caps: PFA
 O-Rings: Teflon Encapsulated Viton (standard)
 Chemraz®
 Kalrez®

Product Traceability

CITADEL Filter Elements are manufactured in conformance with established current Good Manufacturing Practice (cGMP) standards. The filter elements are produced and distributed according to a Quality Management System that is registered for compliance to EN ISO 9001:2015. All filters are non-destructively integrity tested and flushed with Purified Water with a maximum conductivity of 1.1 $\mu\text{S}/\text{cm}$ @ 20°C (68°F) and a maximum TOC (Total Organic Carbon) content of 0.5 mg (500 ppb) of carbon per liter. They are then dried using HEPA filtered air and sealed in a protective polyethylene bag within the cleanroom. Filters that are provided as pre-wet are not dried, but rather are subjected to a further rinse of ozonated water at 0.2 – 0.3 ppm and then vacuum sealed in a cleanroom grade polyethylene bag. Each filter module is marked with an individual serial number, a lot number, product code which is also shown on both the bag label and on the outer product box, therefore all data concerning materials used and production data are documented, accessible and fully traceable.

Cartridge Integrity Test

Graver Technologies, as part of its quality process, integrity tests all CITADEL filter cartridges before release from manufacturing. The specific test used is a Diffusion Test. If further details for the test method are required, please contact Graver Technologies Technical Service. For an integral cartridge, the air diffusion rate which is a measure of the rate at which air diffuses through the water-filled pores of the membrane, must be below a specified value at the Integrity Test pressure. A cartridge with even a minor defect will exhibit much higher airflow rates when measured by this test.

Test Procedure

- 1) The filter cartridge is first wetted by completely immersing in 60/40 IPA/de-ionized water mixture for 5 minutes.
- 2) The water flow is shut off and a pressure of 5 psid (0.34 bard) of compressed air is applied upstream of the filter. Any excess water in the housing passes through the filter and is drained from the downstream side of the housing.
- 3) The air pressure is increased to the value shown in Table below, “Diffusion Pressure” and the system is allowed to stabilize for 2 minutes
- 4) The diffusive air flow through the filter system is measured and the filter passes the integrity test only if the diffusion flow value is less than the “Maximum Diffusion” shown in the table below.

| Pore Size | Diffusion Test Pressure (psig) | Maximum Diffusion (cc/min) per 10-Inch Cartridge Length | Bubble Point (PSIG) |
|--------------------|--------------------------------|---|---------------------|
| 0.05 μm | 22 (1.5 bar) | ≤ 90 | ≥ 28 (1.9 bar) |
| 0.1 μm | 18 (1.2 bar) | ≤ 50 | ≥ 22 (1.5 bar) |
| 0.2 μm | 12 (0.8 bar) | ≤ 50 | ≥ 15 (1.0 bar) |
| 0.45 μm | 5 (0.34 bar) | ≤ 50 | ≥ 6 (0.41 bar) |
| 1.0 μm | 3 (0.2 bar) | ≤ 50 | ≥ 4 (0.27 bar) |
| 10 μm | N/A* | N/A | N/A |

** The 10 μm media is quite open and thus no practical bubble point or diffusion value is available.*

Flow Rate Testing

To contribute to the overall operating economics of an existing filter system, it is important that process filter cartridges offer high flow rates at low-pressure drops. For new systems, this can also allow a smaller filter housing to be used with a resultant savings in capital cost.

Test Procedure

- 1) A filter cartridge is wetted with 60/40 IPA/water mixture. An integrity test is performed and the results are recorded. (See Page 4 for Integrity Test Procedure).
- 2) The cartridge is then flushed with clean water for 15 minutes to remove any residual IPA.
- 3) The filter system is connected to a source of clean water. The pressure of water can be regulated and was adjusted to 18 psi (1.2 bar).
- 4) The flow through the filter is adjusted 3 GPM and the pressure recorded. This is repeated for other values up to 10 GPM.
- 5) The values are averaged over the range and reported out as GPM/PSID.

Results

The filter cartridges at each pore size tested showed flow rates as summarized below, meeting the minimum specifications for that pore size.

| Average Flux Data for Citadel 0.05 micron | | | |
|---|--------------------------------------|----------------------------|-----------------|
| Cartridge ID: C22, 23, 3 & 27 | | | |
| Flowrate (GPM) | Temperature Adjusted flow rate (GPM) | DP (Total - system) (PSID) | Flux (PSID/GPM) |
| 3 | 2.458 | 2.200 | 0.90 |
| 5 | 4.097 | 3.900 | 0.95 |
| 7 | 5.736 | 5.700 | 0.99 |

| Average Flux Data for Citadel 0.1 micron | | | |
|--|--------------------------------------|----------------------------|-----------------|
| Cartridge ID # 75, 76, 77 & 78 | | | |
| Flowrate (GPM) | Temperature Adjusted flow rate (GPM) | DP (Total - system) (PSID) | Flux (PSID/GPM) |
| 3 | 2.55 | 0.875 | 0.34 |
| 5 | 4.25 | 1.525 | 0.35 |
| 7 | 5.95 | 2.222 | 0.37 |
| 10 | 8.5 | 3.625 | 0.42 |

| Average Flux Data for Citadel 0.2 micron | | | |
|--|--------------------------------------|----------------------------|-----------------|
| Cartridge ID: C1, 11, 13 & 15 | | | |
| Flowrate (GPM) | Temperature Adjusted flow rate (GPM) | DP (Total - system) (PSID) | Flux (PSID/GPM) |
| 3 | 2.480 | 0.575 | 0.23 |
| 5 | 4.133 | 1.150 | 0.28 |
| 7 | 5.786 | 1.675 | 0.29 |
| 9 | 7.464 | 2.300 | 0.31 |
| 10 | 8.266 | 2.500 | 0.30 |

| Average Flux Data for Citadel 0.45 micron | | | |
|---|--------------------------------------|----------------------------|-----------------|
| Cartridge ID: CTL 0.45 5546 - 002 & 003 | | | |
| Flowrate (GPM) | Temperature Adjusted flow rate (GPM) | DP (Total - system) (PSID) | Flux (PSID/GPM) |
| 3 | 2.444 | 0.10 | 0.04 |
| 5 | 4.074 | 0.40 | 0.10 |
| 7 | 5.704 | 0.65 | 0.11 |
| 9 | 7.333 | 0.80 | 0.11 |

| Average Flux Data for Citadel 1.0 micron | | | |
|--|--------------------------------------|----------------------------|-----------------|
| Cartridge ID: CTL 1.0 5492 - 002 & 003 | | | |
| Flowrate (GPM) | Temperature Adjusted flow rate (GPM) | DP (Total - system) (PSID) | Flux (PSID/GPM) |
| 3 | 2.444 | 0.10 | 0.04 |
| 5 | 4.074 | 0.10 | 0.02 |
| 7 | 5.704 | 0.20 | 0.04 |
| 9 | 7.333 | 0.20 | 0.03 |

Conclusions

Based on this testing, the typical flow rate/pressure drop characteristics of CITADEL cartridges per 10-inch cartridge length are:

0.05 µm: 1.0 gpm/psid
 0.1 µm: 2.7 gpm/psid
 0.2 µm: 3.5 gpm/psid
 0.45 µm: 9.1 gpm/psid
 1.0 µm: >15 gpm/psid*
 10 µm: >15 gpm/psid*

** The media is quite open and thus flow resistance is low making. Value far exceeds practical recommendations.*

Cartridge Extraction

Testing was conducted by CT Associates of Eden Prairie MN on a CITADEL cartridge (pre-wet) by flushing with 18 megohm-cm DI water. The cartridge was tested at a flow rate of 2.25 GPM (8.5 LPM).

TOC (Total Organic Carbon)

TOC levels were measured throughout the duration of the test using a GE 500 RL TOC analyzer. The system background showed a TOC level of 0.93 ppb. The test system with the installed cartridge had an initial TOC level of 3.95 ppb and by 2 hours and 30 minutes, the added TOC level was below 1 ppb. Total system TOC levels were below 1 ppb within 8 hours.

NVR (Non-volatile Residue)

Non-volatile residue (NVR) was measured using a Kanomax FMT NRM. The system background showed a NVR level of 0.11 ppb. The test system with the installed cartridge had an initial NVR level of 0.26 ppb and the added NVR was below 0.1 ppb within 78 minutes.

Particle Shedding

Particle shedding from the cartridge was measured over time using Lighthouse NC30+ optical particle counter which is capable of measuring particles as small as 30 nm (0.03 microns).

Particle Cleanliness Summary

| <i>Particle Size (nm)</i> | <i>Rinse Time to Conc. Added (h)</i> | | |
|-------------------------------|--------------------------------------|---------------|-----------------|
| | <i>10 p/mL</i> | <i>1 p/mL</i> | <i>0.1 p/mL</i> |
| 30 | 0.14 | 0.55 | 2.21 |
| 50 | 0.05 | 0.29 | 1.73 |
| 80 | 0.00 | 0.00 | 0.03 |
| 100 | 0.00 | 0.00 | 0.00 |

Resistivity Recovery

Resistivity of the system was measured using a Thornton 200CR conductivity meter. Upon installation, the resistivity measured 18 MΩ-cm. The system returned to baseline levels of 18.2 MΩ-cm within 40 gallons.

During manufacture, each CITADEL filter is wet-out with 60/40 IPA/water and then flushed with semiconductor grade ultrapure water for 15 minutes at a flow rate of 10 GPM, assuring the high purity of each filter shipped. The filters that are provide pre-wet have an additional flush step using ozonated water for 10 minutes. Upon start-up, customers can expect minimal extractable levels (TOC, NVR, resistivity decrease). There are no observable particles at 0.1 micron (100 nm) or above at start-up.


Filter Retention

Testing was conducted by CT Associates of Eden Prairie MN on a 0.05 micron CITADEL cartridge (pre-wet) using a modified version of the test procedure described in the *SEMATECH Provisional Test Method for Determining Particle Contributions and Retention by UPW Distribution System Components* [1]. The cartridges are first flushed with 18 megohm-cm DI water for 72 hours to eliminate/reduce any particles that could shed from the cartridge. The filter is then exposed to 0.1% Triton X-100 solution in DI water to eliminate any potential for charge interactions and further flushed for 2 – 4 hours. The test is performed at a flow rate of 2.6 GPM (10 LPM).

Efficiency Average over 16 hours

| Channel Size | Inlet Particles (#/ml) | Outlet Particles (#/ml) | Retention Efficiency |
|---------------------|-------------------------------|--------------------------------|-----------------------------|
| 30 nm | 165.032 | 116 | 99.93% |
| 50 nm | 67,039 | 6.93 | 99.99% |
| 80 nm | 17,982 | 0.34 | 100.00% |
| 100 nm | 8,464 | 0.19 | 100.00% |

The SEMATECH method measures particle retention in DI water by simulating the quantity of particles that could be seen during the equivalent of 1-year downstream of a mixed bed ion exchange resin tank. The 1-year challenge is simulated using a 16-hour filter test.



The reduction values (LRV – log reduction value) remain relatively stable over the time frame. It should be noted that even though the filter is rated at 0.05 micron, capture efficiency of 0.03 micron particle is quite high.

Leachables

Two factory flushed CITADEL 0.05 µm filter cartridge samples were submitted to Balazs Analytical Services for metals extractables. An extraction was done with 1.5 liters of 5% HCL or 1.5% ultrapure water (UPW) and the residue analyzed.

- Leachable anions using ion chromatography
- Leachable Low Level 30 Elements in UPW

The test samples were prepped and leached in accordance with SEMI F57-0301 and F40-0699 protocols. The F57-0301 standard requires manufacturers to produce semiconductor fluid handling products that follow best practices in material science and do not contaminate processes with extractable ions, metallics or TOCs- areas which can negatively impact yield and reduce the life of semiconductor equipment. In addition, the standard provides test protocols and product performance requirements. SEMI F40-0699 covers the Practice for Preparing Liquid Chemical Distribution Components for Chemical Testing. Results are reported as ppb (part per billion) of each contaminant.

.



Air Liquide Electronics, U.S. LP - Balazs™ Analytical Services
 13468 North Central Expressway, M8 301 Dallas TX 76243
 Telephone (972) 995-2228 FAX (972) 995-3173 Duty Phone (214) 878-4673

ANALYTICAL TEST RESULTS

Customer Information

Name: Graver technologies, Praveen Jana
 Address:
 Contact: Praveen Jana
 Phone: 302-731-3540 Fax:

| | | | |
|---------------------|-------------------------------------|-------------------|---------------------|
| Product | HYDROCHLORIC-EXTERNAL | Sample Number | 1018051 |
| Chemical: | HCL | Status | Authorized |
| Sample Type | | Logged By | AGARCIA |
| Vessel: | | Received Date | 6/8/2015 11:00:01AM |
| Vessel Serial #: | | Completed Date: | 6/12/2015 6:39:58PM |
| Lot Number: | C-0.1-7815-004 | PO Number | 140861 |
| Sample In-Spec: | Yes | Quote # : | None |
| Date Sample Pulled: | None | | |
| IC Bottle ID: | None | Metals Bottle ID: | None |
| Description: | HCL extraction of filter leacheable | | |

DRC_HCL Authorized Approval By: Jerry W. Smith, Production Laboratory Manager

Test Comments: There are no comments associated with this test

| Analysis / Rep | Result Value | Min. | Max. | Units | Entered By | Date Test Completed |
|----------------|--------------|------|------|-------|------------|---------------------|
| Aluminum | 1.788 | | | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Antimony | <0.020 | | | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Arsenic | <0.025 | | | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Barium | 0.030 | | | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Beryllium | <0.005 | | | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Bismuth | 0.005 | | | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Boron | 0.034 | | | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Cadmium | <0.003 | | | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Calcium | 9.326 | | | ppb | BWIDJAJA | 6/12/2015 6:39:57PM |
| Chromium | 0.342 | | | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Cobalt | 0.005 | | | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Copper | 0.112 | | | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Gallium | <0.005 | | | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Germanium | <0.030 | | | ppb | BWIDJAJA | 6/12/2015 6:39:57PM |
| Gold | <0.080 | | | ppb | BWIDJAJA | 6/12/2015 6:39:57PM |
| Indium | <0.004 | | | ppb | LABSTATION | 6/12/2015 6:39:57PM |

Printed On 6/12/2015 10:14:46PM

Sample Revision: 1.00

Sample Summary Report for Sample #: 1018051

Page 1 of 2



| | | | | |
|------------|--------|-----|------------|---------------------|
| Iron | 4.650 | ppb | BWIDJAJA | 6/12/2015 6:39:57PM |
| Lanthanum | <0.004 | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Lead | 0.019 | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Lithium | <0.003 | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Magnesium | 0.437 | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Manganese | 0.058 | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Molybdenum | 0.032 | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Nickel | 0.473 | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Niobium | <0.003 | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Palladium | <0.060 | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Platinum | <0.020 | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Potassium | 1.179 | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Silver | <0.050 | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Sodium | 1.867 | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Strontium | 0.009 | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Tantalum | <0.015 | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Thallium | <0.003 | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Tin | <0.013 | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Titanium | 0.066 | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Tungsten | 0.012 | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Vanadium | 0.127 | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Zinc | 0.830 | ppb | LABSTATION | 6/12/2015 6:39:57PM |
| Zirconium | 0.019 | ppb | LABSTATION | 6/12/2015 6:39:57PM |

Review Notes: There are no review notes associated with this sample

These results were obtained by following standard laboratory procedures and are only representative of the samples as received by the laboratory. The liability of AIR LIQUIDE - BALAZS ANALYTICAL SERVICES ("Balazs") shall not exceed the amount paid for this report. In no event shall Balazs be liable for special or consequential damages. Client agrees not to use Balazs' name in reporting results obtained from tests performed by Balazs without first obtaining Balazs written consent as to such use. Report shall not be reproduced except in full, without the written approval of Balazs.

The information contained in this electronic or original document is confidential information for the intended recipient. If you are not the intended recipient, or the employee or agent responsible for delivering this document to the intended recipient, you are hereby notified that any disclosure, dissemination, distribution, or copying of this communication is strictly prohibited. If you have received this document in error, please immediately notify the sender and destroy this document.

Printed On 6/12/2015 10:14:46PM

Sample Revision: 1.00

Sample Summary Report for Sample #: 1018051

Page 2 of 2





Air Liquide Electronics, U.S. LP - Balazs™ Analytical Services
 13468 North Central Expressway, MS 901 Dallas TX 75243
 Telephone (972) 995-2228 FAX (972) 995-3173 Duty Phone (214) 875-4673

ANALYTICAL TEST RESULTS

Customer Information

Name: Graver technologies, Cris Lemay
 Address:
 Contact:
 Phone: 302-824-3828 Fax:

| | | | |
|---------------------|--|-------------------|---------------------|
| Product | DIW_EXTERNAL-External | Sample Number | 1011857 |
| Chemical: | DIH2O | Status | Authorized |
| Sample Type | | Logged By | AGARCIA |
| Vessel: | | Received Date | 5/4/2015 12:36:33PM |
| Vessel Serial #: | | Completed Date: | 5/15/2015 4:22:27PM |
| Lot Number: | C-0.1-7815-002 | PO Number | N/A |
| Sample In-Spec: | Yes | Quote # : | None |
| Date Sample Pulled: | None | | |
| IC Bottle ID: | None | Metals Bottle ID: | None |
| Description: | UPW extraction of Filter. ICP-MS 30 elements | | |

ICPMS_DIW Authorized Approval By: Jerry W. Smith, Production Laboratory Manager

Test Comments: There are no comments associated with this test

| Analysis / Rep | Result Value | Min. | Max. | Units | Entered By | Date Test Completed |
|----------------|--------------|------|------|-------|------------|---------------------|
| Aluminum | 0.470 | | | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Antimony | <0.004 | | | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Arsenic | <0.007 | | | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Barium | 0.025 | | | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Beryllium | <0.005 | | | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Bismuth | <0.005 | | | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Boron | 0.083 | | | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Cadmium | <0.004 | | | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Calcium | 8.211 | | | ppb | NPHUNG | 5/15/2015 4:22:27PM |
| Calcium | 8.950 | | | ppb | NPHUNG | 5/15/2015 4:22:27PM |
| Chromium | 0.013 | | | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Cobalt | 0.004 | | | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Copper | 0.045 | | | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Gallium | <0.004 | | | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Germanium | 0.010 | | | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Gold | <0.005 | | | ppb | LABSTATION | 5/15/2015 4:22:27PM |

Printed On 5/15/2015 5:40:30PM

Sample Revision: 1.00

Sample Summary Report for Sample #: 1011857

Page 1 of 2



| | | | | |
|------------|--------|-----|------------|---------------------|
| Indium | <0.004 | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Iron | 0.424 | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Lanthanum | <0.003 | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Lead | 0.011 | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Lithium | <0.004 | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Magnesium | 0.307 | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Manganese | 0.020 | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Molybdenum | 0.005 | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Nickel | 0.076 | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Niobium | <0.005 | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Palladium | <0.003 | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Platinum | <0.009 | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Potassium | 0.674 | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Silver | <0.006 | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Sodium | 1.128 | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Sodium | 1.185 | ppb | NPHUNG | 5/15/2015 4:22:27PM |
| Strontium | 0.005 | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Tantalum | <0.008 | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Thallium | <0.002 | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Tin | <0.004 | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Titanium | 0.105 | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Tungsten | 0.006 | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Vanadium | 0.043 | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Zinc | 0.281 | ppb | LABSTATION | 5/15/2015 4:22:27PM |
| Zirconium | <0.012 | ppb | LABSTATION | 5/15/2015 4:22:27PM |

Review Notes: There are no review notes associated with this sample

These results were obtained by following standard laboratory procedures and are only representative of the samples as received by the laboratory. The liability of AIR LIQUIDE - BALAZS ANALYTICAL SERVICES ("Balazs") shall not exceed the amount paid for this report. In no event shall Balazs be liable for special or consequential damages. Client agrees not to use Balazs' name in reporting results obtained from tests performed by Balazs without first obtaining Balazs written consent as to such use. Report shall not be reproduced except in full, without the written approval of Balazs.

The information contained in this electronic or original document is confidential information for the intended recipient. If you are not the intended recipient, or the employee or agent responsible for delivering this document to the intended recipient, you are hereby notified that any disclosure, dissemination, distribution, or copying of this communication is strictly prohibited. If you have received this document in error, please immediately notify the sender and destroy this document.

Printed On 5/15/2015 5:40:30PM

Sample Revision: 1.00

Sample Summary Report for Sample #: 1011857

Page 2 of 2



Endotoxin Test

Endotoxins are complex polysaccharide molecules (LPS) composed of lipid (lipid A) and polysaccharide sides chains and are integral components of the outer membrane of gram negative bacteria. These molecules are not secreted but are released only when the cells are disrupted or destroyed. Above certain levels, endotoxins elicit an antigenic response, resulting in fever and altered resistance to bacterial infections. Because of this sensitivity, it is important to monitor products which may contact fluids that could be administered to humans or animals.

The detection of endotoxins is accomplished using Limulus Amebocyte Lysate (LAL) Kinetic Chromogenic Assay. In this test, a filter element is extracted with non-pyrogenic Water for Injection (WFI). Endotoxin levels in the extracted fluid are then measured spectrophotometrically and compared to standard concentrations. These values are reported as EU/ml (Endotoxin Units/ml).

Results

NAMSA of Northwood Ohio conducted the testing on a sample of Citadel. Levels were reported at 0.005 EU/ml and 5.05 EU/device. This level is well below the criteria established by the US FDA of 0.5 EU/ml and the USP requirement of 20 EU/device.

RESULTS

| | |
|--|---|
| Test Article Extract Dilution: | 1 |
| Positive Product Control Percent Recovery: | 117% (between 50% and 200% is acceptable) |
| Test Article Extract: | < 0.00500 EU/mL (Total Concentration) < 5.05 EU/device |

TEST ACCEPTANCE CRITERIA

| Type of Product | Current FDA Requirements | Current USP Requirements |
|---|-------------------------------------|--------------------------------------|
| Anterior Segment Solid Intraocular Device | Less than or equal to 0.2 EU/device | Not Applicable |
| Medical Device | Less than or equal to 0.5 EU/mL* | Less than or equal to 20.0 EU/device |
| Medical Device Contacting Cerebrospinal Fluid | Less than or equal to 0.06 EU/mL* | Less than or equal to 2.15 EU/device |
| Water for Injection | Not Applicable | Less than or equal to 0.25 EU/mL |

*Based on an extraction volume of 40 mL/device

European Regulation No 1935/2004 and European Regulation 10/2011

The underlying principle of these regulations is to ensure that any material or article intended to come into contact directly or indirectly with food must be sufficiently inert to preclude substances from being transferred to food in quantities large enough to endanger human health or to bring about an unacceptable change or deterioration in the composition or properties of the food. Tests for migration behavior in direct food contact were conducted by Belgium Packaging Institute in a variety of liquids to simulate aqueous, acidic, alcoholic, and fatty foodstuffs.

Results

The test results indicate that the overall migration of all of the individual parts will not exceed the overall migration limit of 10 mg/dm² or 60 mg/kg foodstuffs for simulant A (10% ethanol, representing all aqueous foodstuffs), simulant B (3% acetic acid, representing all foodstuffs with a pH below 4.5) and simulant D2 (95% ethanol and isooctane instead of olive oil representing all fatty foodstuffs) using the given conditions. In accordance with the European Regulation No 10/2011 and amendments, conformity with the overall migration limit for simulants A, B and D2 demonstrates suitability for contact with all kinds of foodstuffs.



NOTES