

A panel of methylation markers for multi-cancer detection from plasma



J.Kinross¹, K. Kruusmaa², M. Bitenc²⁻³, M. Chersicola³, W. Pulverer⁴, A. Weinhaeuser⁴



¹Surgery and Cancer, Imperial College London, London, United Kingdom, ²Research & Development, Universal Diagnostics S.L., Seville, Spain, ³Research & Development, Geneplanet d.o.o., Ljubljana, Slovenia, ⁴Molecular Diagnostics, Austrian Institute of Technology GmbH, Vienna, Austria

Introduction

- Cancer is the leading cause of premature mortality in the world and in dire need of early detection tools (1)
- Measuring cancer-related alterations in the tumour-derived portion of the cell-free DNA (cfDNA) could offer an accurate, non-invasive approach, leading to decreased cancer mortality (2-3)
- We report here a plasma targeted methylation marker panel performance for detection of 4 cancer types – colorectal (CRC), lung (LC), pancreatic (PaC) and breast (BC) cancer

- Random forest feature selection algorithm utilizing Monte-Carlo cross-validation over 50 sub-setting was utilized for building and testing marker panels:
 - Panel 1 was used for overall cancer detection
 - Panel 2 was used for localization of tissue of origin.
- Accuracy was defined as the fraction of correct calls.

Table 1 Sample set

Characteristics	Control (n=71)	Colorectal Cancer (n=20)	Pancreatic Cancer (n=15)	Lung Cancer (n=37)	Breast Cancer (n=28)
Age (years, average (IQR))	60 (32-78)	59 (28-77)	60 (50-71)	59 (42-75)	57 (41-77)
Gender (n (%))					
Female	36 (51%)	10 (50%)	7 (47%)	12 (32%)	28 (100%)
Male	35 (49%)	10 (50%)	8 (53%)	25 (68%)	
Stage					
Stage I		6		15	7
Stage II		9		4	8
Stage III		5	15	17	13
Unknown				1	1

Figure 1 Performance of methylation markers. A. dCt values plotted for 4 markers for control samples (red) and cancer samples (blue). B. ROC curve and AUC for 10-marker panel performance (Panel 1)





Figure 2 Performance of 10 methylation marker panel (Panel 1). A. Accuracy values, where red bar represents the overall sensitivity for cancer, red-outlined bars represent sensitivity by cancer type, and blue bar represents specificity. B. Accuracy values, red-outlined bars represent sensitivity by stage and blue bar represents specificity



- Cancer cases that were correctly separated from control group were further evaluated for their tissue of origin.
- Individual markers showed differential methylation levels between different cancer types (Figure 3)

- 16 methylation marker panel (Panel 2) allowed correct tissue of origin assignment of 80% of CRC (16/20), 78% of LC (21/27), 75% of PaC (9/12) and 62% of BC (13/21) cases (Table 2).

Figure 3 Individual marker performance. A. dCt values plotted for 4 markers for BC (red), CRC (green), LC (dark blue), PaC (light blue)




Table 2 Confusion matrix according to 16-marker panel (Panel 2). Diagonal line indicates the correct classification

	BC	CRC	LC	PaC	class. error	Correct
BC	13	0	4	4	0.38	62%
CRC	0	16	2	2	0.2	80%
LC	1	3	21	2	0.22	78%
PaC	1	1	1	9	0.25	75%

Conclusions

- We show that targeted methylation marker panels have potential for early blood-based detection of multiple cancers with high sensitivity and specificity and further tissue of origin localization.
- This method could serve as the basis for further development of a simple, highly accurate and minimally invasive blood-based multi-cancer screening test.

Methods

- Human Methylation 450K data available from The Cancer Genome Atlas (TCGA) consortium was used for initial methylation marker selection (4)
- Candidate methylation marker regions were evaluated in the plasma samples of 101 patients with cancer and 71 age/gender- matching non-cancer controls (Table 1)
- Methylation-sensitive restriction enzyme-qPCR approach was used to target regions of interest in plasma cfDNA
- dCt (delta cycle threshold) values were used for further analysis

Results

- Differential methylation could be observed on single marker level (Figure 1 A)
- 10 methylation marker panel (Panel 1) showed area under curve (AUC) of 89% (Figure 1B), with 79% (80/101) sensitivity for cancer at 90% (64/71) specificity
- 100% of the CRC (20/20), 80% of the PaC (12/15), 75% of the BC (21/28) and 73% of LC (27/37) were correctly identified as cancer patients (Figure 2A)
- Notably, the sensitivity for stage I cancers was 75% (21/28) (Figure 2B).

Corresponding author contact

<Kristi Kruusmaa>
 <Universal Diagnostics S.L.>
 Email: kristi.kruusmaa@universaldx.com
 Website: https://www.universaldx.com/

Bibliography

1. <https://www.who.int/news-room/fact-sheets/detail/cancer>
2. Fiala C. *et al.* BMC Medicine volume 17, Article number: 159 (2019)
3. Barnaby G, *et al.* J Intern Med 2019 Aug;286(2):118-136.
4. <https://www.cancer.gov/about-nci/organization/ccg/research/structural-genomics/tcga>

Disclosure statement

Dr.Kinross:
 . Funding from NIHR: II-OL-1116-10027, NIH: R01-CA204403-01A1, Horizon H2020: ITN GROWTH.
 . Imperial Biomedical Research Centre, SAGES research grant. Bowel and Cancer Research, CRUK
 . Scientific advisory boards: Verb robotics, Safeheal, LNC therapeutics, Medical iSight.
 . Equity:1 Welbeck Day Surgery (day surgery unit)
 . Directorships: Mangetoo.com (teledietetics), 1Worldmedical, Cerulean health.
 . Consulting: Ethicon (J+J), Medtronic
 . Honoraria: Yakult
 . Study is sponsored by Universal Diagnostics S.L.