



Made for DSI Ltd.  
59 Havoda str.  
Ashdod, 7706300 Israel  
+972-8931-7235  
[www.dsisrael.com](http://www.dsisrael.com)

---

Remineralizing Gel fortified with Bioglass

# DSI Sense-X Bio



## INSTRUCTIONS FOR USE

### DESCRIPTION

**DSI Sense-X Bio** is a re-mineralizing gel that reduces post-whitening sensitivity and effectively enhances the re-mineralization of the enamel.

### COMPOSITION

Glycerine, Potassium Nitrate, Bioglass, Carbopol.

### INDICATIONS

- Remineralization of de-mineralized tooth structure.
- Post whitening sensitivity.
- Sensitivity is caused due to non-carious lesions and sensitivity after oral prophylaxis.

### DIRECTION FOR USE

- Sense-X Bio is presented in a convenient dispenser pen with paint on the brush for easy application. The gel will be dispensed on the brush tip by gently rotating the bottom of the pen clockwise. First-time use requires additional rotations for the gel to appear on the brush tip.
  - Brush your teeth each time before applying Sense-X Bio. Sense-X Bio is more effective on clean teeth.
  - Apply a small amount of Sense-X Bio in a thin layer and cover the entire area under treatment.
-

- Avoid lips from contacting teeth for 30-60 sec.
- Leave the Sense-X Bio on the tooth surface and avoid rinsing, brushing your teeth for 4 to 6 hours.
- Apply Sense-X Bio as often as needed.

## **CAUTION**

Always close the pen with the lid after use to avoid drying of the gel on the brush bristles.

## **CONTRAINDICATIONS**

This product is contra-indicated for use in or by individuals who have known hypersensitivity to any of the constituents of the products.

## **STORAGE**

Store at a temperature between 10 °C to 24 °C. Close the cap carefully after use. Keep away from moisture. Shelf life is not affected after opening if stored in a cool dry place. Shelf life: 3 Years from the date of manufacturing.

## **PACKAGING**

Gel in Dispensing Pen, 2ml

**Ref: SENSE-XBIO**

## **SYMBOLS**



Consult instructions for use



Caution, consult accompanying documents



Temperature limit



Keep away from sunlight



Keep dry



Manufacturer



Catalogue number



Batch code



Use by



Date of Manufacture

*Failure to comply with the conditions of storage leads to a change of the working characteristics of the material and decrease the shelf life of the material.*

*The manufacturer is not responsible for any loss of quality caused by the failure to comply with terms of transportation, storage and use established by the manufacturer for this product. Responsibility for the use of the material for purposes other than those specified by the manufacturer falls on the user.*