



Made for DSI Ltd.
59 haAvoda str.
Ashdod, 7706300 Israel
+972-8931-7235
www.dsisrael.com

Cleaning and degreasing agent

DSI U-Dry



INSTRUCTIONS FOR USE

DESCRIPTION

DSI U-Dry is a fast drying, cleaning and degreasing agent for dental surfaces.

COMPOSITION

Ethyl alcohol and an organic solvent.

INDICATIONS

To remove moisture from dental surfaces.

To dry tooth surface prior to applying pit and fissure sealants and other materials which are adversely affected by the presence of moisture.

DIRECTION FOR USE

- Isolate the operative area in the mouth to maintain a clean dry field (for intra oral applications only).
- Rinse the operative area and air dry with clean, uncontaminated, oil free and moisture free compressed air.
- Shake agent bottle. Dispense a few drops of U-Dry on the desired surface. For maximum control of flow, do not push on the bottle.
- To promote rapid drying, use a gentle flow of clean, uncontaminated air if desired.



Warning:

- U-Dry composition is highly flammable. Handle with care and keep away from open flames.
- Do not place U-Dry directly on dentin as it can desiccate the dentinal tubules.

CONTRAINDICATIONS

This product is contra-indicated for use in or by individuals who have shown hypersensitivity to any of its components.

STORAGE

Store at temperature between 10° C to 24° C. Close the cap carefully after use. Keep away from moisture. Shelf life not affected after opening if stored in a cool dry place. Shelf Life: 2-years from the date of manufacturing.

PACKAGING

Plastic Bottle, liquid 10ml

Ref: UDRY

SYMBOLS



Consult instructions for use



Caution, consult accompanying documents



Temperature limit



Keep away from sunlight



Keep dry



Flammable



Keep away from children



Manufacturer



Catalogue number



Batch code



Use by



Date of Manufacture

Failure to comply with the conditions of storage leads to a change of the working characteristics of the material and decrease the shelf life of the material.

The manufacturer is not responsible for any loss of quality caused by the failure to comply with terms of transportation, storage and use established by the manufacturer for this product. Responsibility for the use of the material for purposes other than those specified by the manufacturer falls on the user.

