


# iST Technology

# Sample Preparation You Can Trust

Setting the Standard for Protein Analysis



## Universal Solution

- From biological fluids, cells and tissues to plants
- Versatile format for protein input from 1 µg to 100 µg
- Diverse workflows: label-free & chemical labeling



## Reliable Results

- Automatable all-in-one solution
- Standardized workflows for reproducibility
- Validated by scientific community





## Focus on Biology

- Minimal hands-on time
- Concentrate on data analysis
- Purified peptides enhance your LC-MS performance

# iST Technology Applications

## Reproducible in-depth coverage of a variety of tissue proteomes



Superb coverage of various mouse tissue proteomes was achieved by combining BeatBox tissue homogenization with the iST-fractionation (iST-frx) workflow.

Outstanding reproducibility (CV ~9%) of tissue homogenization is achieved using the BeatBox – iST workflow.

3-step peptide fractionation increases the proteome depth by ~40% compared to unfractionated samples.

### Workflow:


Mouse tissues (1-2 mg wet weight) were homogenized with the BeatBox (10 min, STANDARD mode).

Homogenates were processed with an iST and iST-fractionation (iST-frx) add-on kit (n=3).

Peptides were measured in a 45 min LC-gradient on an EASY nLC 1200 coupled to a timsTOF Pro.




## Reliable analyses of serum and yeast samples by a sensitive and fully automated proteomic sample processing



### Workflow:

Human serum (1-100  $\mu\text{g}$  of total protein) and yeast pellet (100  $\mu\text{g}$  of total protein) were processed fully automated with the iST kit (n=3) on the PreON platform.

Peptides were measured in a 90 min LC-gradient on an EASY nLC 1000 coupled to a Q Exactive Plus.




Easy-to-use, automatable, and standardized kit-based solution for a wide range of sample types.

Highly confident protein identification rate achieved from 1 to 100  $\mu\text{g}$  of serum protein input.

Excellent proteomics workflow reproducibility using a fully automated iST sample preparation on the PreON instrument.


## Optimized sample preparation workflow for biological fluids



### Workflow:

Human plasma from 3 different donors (100 µg of total protein) was processed with the iST-BCT kit (n=4).

Peptides were measured in a 60 min LC-gradient on an EASY nLC 1200 coupled to a timsTOF Pro.




Easy-to-use 3-step sample preparation to analyze a broad range of biological fluids, from complex plasma to challenging Cerebrospinal fluid.

Nearly 60% of the human plasma proteome is shared between three healthy patient donors.

iST-BCT workflow reduces artificial modifications and improves alkylation rate to enhance proteomic coverage for biological fluids.

## Efficient, reproducible, high-throughput cell lysis on the BeatBox platform




### Workflow:

*E. coli* (2.1E7 and 8.4E7 cells), *S.cerevisiae* (1.1E6 and 4.4E6 cells) and HEK293 (1E5 and 4E5) were lysed with the BeatBox (10 min, STANDARD mode).

Cell lysates were processed with the iST kit (n=3).

Peptides were measured in a 45 min LC-gradient on an EASY nLC 1200 coupled to a timsTOF Pro.




BeatBox enables high-throughput, efficient, and reproducible lysis across many cell types.

Effective cell lysis followed by iST sample preparation for different cell counts resulted in a similar protein group identification rate.

Outstanding reproducibility (CV ~3%) is achieved using the BeatBox – iST workflow.

## Fully automated workflow to enable a highly efficient isobaric labeling




### Workflow:

Cells (50 µg of total protein) were processed with the iST-NHS PreON kit (n=4).

Peptides were labeled with TMT zero, TMT 11-plex and TMTPro 16-plex.

Labeled peptides were measured using a 100 min LC-gradient on an EASY nLC 1200 coupled to a Q Exactive HF-X or Orbitrap Fusion Lumos Tribrid.

Collaboration with Prof. Dr. Bernd Wollscheid's lab (ETH Zurich, CH)



Automation of the iST technology on the PreON enables isobaric labeling of up to 18 samples\* in parallel, from cells to labeled peptides in less than 4 hours.



The fully automated workflow increases the throughput and reproducibility of quantitative proteomics experiments using isobaric tags.

Excellent TMT tag labeling efficiency (>98%) even with a label to peptide ratio of 4:1.

\*With the TMTpro™ 18-plex upgrade set.



# Optimized plasma and serum sample preparation for deeper, high-throughput proteome profiling



## Workflow:

20 µL of plasma and serum (mouse, rat and human, respectively) was processed with the ENRICH-iST (n=3).

Peptides were measured using a 45 min LC-gradient on a nanoElute® 2 coupled to a TimsTOF Pro HT.

A fast, easy-to-use, and standardized protocol for high reproducibility and less hands-on time, automatable high-throughput processing of up to 96 samples in parallel possible in just 5h.

Species-independent and flexible enrichment technology optimized for plasma and serum samples that enhance proteomics depth compared to the iST-BCT workflow.

ENRICH-iST workflow provides superior robustness with overall reproducibility of CV<15%, which is essential for the reliable execution of large-scale studies.



# Selection of publications and products listing

## Selection of publications

- Proteomic analysis distinguishes extracellular vesicles produced by cancerous versus healthy pancreatic organoids, Buenafe et al., 2022, Scientific Reports, doi: <https://doi.org/10.1038/s41598-022-07451-6>
- Supramolecular arrangement of protein in nanoparticle structures predicts nanoparticle tropism for neutrophils in acute lung inflammation, Myerson et al., 2022, Nat. Nanotechnol. 17, 86–97, doi: <https://doi.org/10.1038/s41565-021-00997-y>
- A clinically compatible drug-screening platform based on organotypic cultures identifies vulnerabilities to prevent and treat brain metastasis, Zhu et al., 2022, EMBO Mol Med 14:e14552, doi: <https://doi.org/10.15252/emmm.202114552>
- A lysosomal biogenesis map reveals the cargo spectrum of yeast vacuolar protein targeting pathwa, Eising et al., 2022, J Cell Biol 221 (4): e202107148, doi: <https://doi.org/10.1083/jcb.202107148>

## Products listing

Product	Manufacturer	Product Code
iST 8x	PreOmics GmbH	P.O.00001
iST 96x	PreOmics GmbH	P.O.00027
iST HT 192x	PreOmics GmbH	P.O.00067
iST-BCT 8x	PreOmics GmbH	P.O.00084
iST-BCT 96x	PreOmics GmbH	P.O.00099
iST-BCT-HT 192x	PreOmics GmbH	P.O.00124

Product	Manufacturer	Product Code
iST-NHS 12x	PreOmics GmbH	P.O.00026
iST-NHS 96x	PreOmics GmbH	P.O.00030
iST-NHS-HT 192x	PreOmics GmbH	P.O.00083
ENRICH-iST 8x	PreOmics GmbH	P.O.00163
ENRICH-iST 96x	PreOmics GmbH	P.O.00164
ENRICH-iST 96x HT	PreOmics GmbH	P.O.00165

**Want to know more?**

**Visit [www.preomics.com](http://www.preomics.com)  
to learn more about the iST technology**

**PreOmics GmbH**  
Martinsried | Germany

**PreOmics Inc**  
Billerica, MA | USA