

LAB PROTOCOL: pGREEN BACTERIAL TRANSFORMATION

PREP

WATER BATH at 42°C

Foam tube holder for water bath

CONTAINER OF CRUSHED ICE

pGREEN plasmid

CaCl

1 Eppendorf tube: Label "pGREEN"

ROOM TEMPERATURE

SOC (Luria Broth may be substituted)

E. coli MM294 Culture on LB plate (cultured ~8-16 hours ago)

INCUBATOR AT 37°C

1 LB+AMP Agar Plate (Label "LB + Amp + pGREEN")

SUPPLIES

Trash container for pipettor tips and loops

Pipettors for 100µL, 250µL, 10µL

Sterile Transfer Pipets (100µL, 250µL, 10µL)

5 Glass Beads for spreading on plate

1 Sterile Transfer Loop for E. coli colony

Sharpie

parafilm for sealing plates

RESOURCES

www.laurasplan.com/coalesce

LAB PROTOCOL: pGREEN BACTERIAL TRANSFORMATION

1. Pipette 250µL of CaCl into pGREEN tube

2. Submerge pGREEN tube in ice

3. Use sterile loop to transfer an isolated colony of E. coli from LB plate to pGREEN tube. Try not to transfer agar.

- Immerse loop and vigorously spin in tube to dislodge E. coli cells.
- Hold tube to light to observe that cells have fallen off loop.
- Suspend cells. (Flick tube with middle finger to unclump cells.)
- Avoid creating bubbles or splashing sides.
- Examine against light. There should not be visible clumps of cells.

Solutions should appear milky white.

4. Return pGREEN tube to ice

5. Pipette 10µL of plasmid to pGREEN tube

- Flick to mix in and suspend cells

6. Return pGREEN tube to ice.

- Insert tube into foam tube holder with solution hanging out bottom entirely and re-submerge on ice while in foam holder.

7. Incubate pGREEN tube while in foam holder on ice for 15min.

- Carefully carry ice container over to water bath so you are ready for heat shock step.

LAB PROTOCOL: pGREEN BACTERIAL TRANSFORMATION

8. Heat shock cells in 42°C water bath.

- Place foam holder with pGREEN tube into water bath for 45 seconds.

9. Return pGREEN tube to ice for 1 or more minutes.

10. Pipette 250µL SOC to pGREEN tube. (1:1 ratio for CaCl:SOC).

- Flick tube gently to mix the SOC with the cell suspension solution.
- Place pGREEN tube on rocker in incubator for 30 minute recovery.

11. Pipette 100µL from pGREEN tube to LB+AMP plate. Divvy it out as you pipette by releasing drops around the plate.

12. Spread solution on plate with glass beads

- Slightly open plate lid and shake out 5 glass bead onto plate
- Use back & forth shaking motion to move beads across surface to evenly spread. (not swirling round & round) for 30-60 seconds.
- Remove beads: To Remove glass beads, hold plate vertically over trash container. Clam shell lower part of plate and tap out beads into container.

13. Seal plates with parafilm.

14. Place upside down in incubator at 37°C (98.6°F) for 24-36 hours.

15. Observe florescence with UV Flashlight.