



Data Sheet - UA:18C_{SVT}

Panels:	18 panel
Alcohol:	Yes
SVT Strip:	Yes
Drug names:	Direct Read
Opiates:	300 ng/mL
Product Use:	Forensic Use Only
Package Insert	B22624-03

SECTION 01

Product Overview

The UA:18C rapid drug screen by Labb is a urine point-of-care drug screen that is designed to meet the needs of criminal justice drug testing programs. The UA:18C detects substances like Marijuana, Cocaine, Amphetamine, Benzodiazepines, Barbiturates, Buprenorphine, Methamphetamine, Methadone, Methaqualone, MDMA/MDA, Morphine, Codeine, Oxycodone, Oxymorphone, 6-MAM, PCP Spice (K3), and ETG Alcohol.

This Data Sheet is not a package insert. It is important to note that only the package insert can be used to determine precise product

performance. The package insert is located in the box along with the 25 UA:16C urine drug screen cups. The information contained in this Data Sheet may be outdated.

IMPORTANT: Always refer to the product's official package insert located in box along with 25 rapid drug test cups.

SECTION 03

Drugs & Detection Levels

UA:16C drug test is designed to detect the following drugs and/or metabolites at the following concentration levels.

6-MONOACETYLMORPHINE (6-MAM) - Drug 1

6-Monoacetylmorphine (6-MAM) is one of three active metabolites of heroin (diacetylmorphine), the others being morphine and the much less active 3-acetylmorphine (3-ACM). 6-MAM is rapidly created from heroin in the body, and then is either metabolized into morphine or excreted in the urine. Since 6-ACM is a unique metabolite to heroin, its presence in the urine confirms that heroin was the opioid used. This is significant because on a urine immunoassay drug screen, the test typically tests for morphine, which is a metabolite of a number of legal and illegal opiates/opioids such as codeine, morphine sulphate, and heroin. 6-MAM remains in the urine for no more than 24 hours so a urine specimen must be collected soon after the last heroin use, but the presence of 6-MAM guarantees that heroin was in fact used as recently as within the last day.

6-Monoacetylmorphine	
6-Monoacetylmorphine	10
Heroin HCL	250
Oxymorphone	>100,000
Oxycodone	>100,000
Codeine	>600,000
Norcodine	>100,000
Naltrexone	>100,000
Morphine	>500,000

Hydrocodone	>100,000
Hydromorphone	>100,000

AMPHETAMINE (AMP) - Drug 2

Amphetamine is a Schedule II controlled substance available by prescription (Dexedrine®) and is also available on the illicit market. Amphetamines are a class of potent sympathomimetic agents with therapeutic applications. They are chemically related to the human body's natural catecholamines: epinephrine and norepinephrine. Acute higher doses lead to enhanced stimulation of the central nervous system and induce euphoria, alertness, reduced appetite, and a sense of increased energy and power. Cardiovascular responses to amphetamines include increased blood pressure and cardiac arrhythmias. More acute responses produce anxiety, paranoia, hallucinations, and psychotic behavior. The effects of Amphetamines generally last 2-4 hours following use, and the drug has a half-life of 4-24 hours in the body. About 30% of Amphetamines are excreted in the urine in unchanged form, with the remainder as hydroxylated and deaminated derivatives.

Amphetamine 500	
D-Amphetamine	500
D,L-Amphetamine	500
Methylenedioxy-amphetamine (MDA)	200
L-Amphetamine	25,000
4-Hydroxyamphetamine HCL	300
Diethylstihestrol	2,500
Phentermine	4,000

BARBITURATES (BAR) - Drug 3

Barbiturates are central nervous system depressants. They are used therapeutically as sedatives, hypnotics, and anticonvulsants. Barbiturates are almost always taken orally as capsules or tablets. The effects resemble those of intoxication with alcohol. Chronic use of barbiturates leads to tolerance and physical dependence. Short acting Barbiturates taken at 400 mg/day for 2-3 months can produce a clinically significant degree of

physical dependence. Withdrawal symptoms experienced during periods of drug abstinence can be severe enough to cause death. Only a small amount (less than 5%) of most Barbiturates are excreted unaltered in the urine. The approximate detection time limits for Barbiturates are:
 Short acting (e.g. Secobarbital) 100 mg PO (oral) 4.5 days.
 Long acting (e.g. Phenobarbital) 400 mg PO (oral) 7 days.

Barbiturates 300	
Secobarbital	300
Amobarbital	300
Alphenal	750
Aprobarbital	250
Butabarbital	100
Cyclopentobarbital	500
Phenobarbital	250

BENZODIAZEPINES (BZO) - Drug 4

Benzodiazepines are medications that are frequently prescribed for the symptomatic treatment of anxiety and sleep disorders. They produce their effects via specific receptors involving a neurochemical called gamma aminobutyric acid (GABA). Because they are safer and more effective, Benzodiazepines have replaced barbiturates in the treatment of both anxiety and insomnia. Benzodiazepines are also used as sedatives before some surgical and medical procedures, and for the treatment of seizure disorders and alcohol withdrawal. Risk of physical dependence increases if Benzodiazepines are taken regularly (e.g., daily) for more than a few months, especially at higher than normal doses. Stopping abruptly can bring on such symptoms as trouble sleeping, gastrointestinal upset, feeling unwell, loss of appetite, sweating, trembling, weakness, anxiety and changes in perception. Only trace amounts (less than 1%) of most Benzodiazepines are excreted unaltered in the urine; most of the concentration in urine is conjugated drug. The detection period for the Benzodiazepines in the urine is 3-7 days.

BENZODIAZEPINES	
Alprazolam	200

Bromazepam	200
Chlordiazepoxide HCL	100
Clobazam	100
Clorazepate	2,500
Clorazepate Dipotassium	200
Delorazepam	10
Desalkylflurazepam	400
Diazepam	200
Estazolam	2,500
Flunitrazepam	400
a-Hydroxyalprazolam	1,260
(±) Lorazepam	1,560
Lorazepam D4	500
Nitrazepam	100
Norchlordiazepoxide	200
Nordiazepam	400
Oxazepam	300
Temazepam	100
Triazolam	2500

BUPRENORPHINE (BUP) - Drug 5

Buprenorphine is a semisynthetic opioid analgesic derived from thebain, a component of opium. It has a longer duration of action than morphine when indicated for the treatment of moderate to severe pain, peri-operative analgesia, and opioid dependence. Low doses buprenorphine produces sufficient agonist effect to enable opioid-addicted individuals to discontinue the misuse of opioids without experiencing withdrawal symptoms. Buprenorphine carries a lower risk of abuse, addiction, and side effects compared to full opioid agonists because of the "ceiling effect", which means no longer continue to increase with further increases in dose when reaching a plateau at moderate doses. However, it has also been shown that Buprenorphine has abuse potential and may itself cause dependency. Subutex®, and a Buprenorphine/Naloxone combination product, Suboxone®, are the only

two forms of Buprenorphine that have been approved by FDA in 2002 for use in opioid addiction treatment. Buprenorphine was rescheduled from Schedule V to Schedule III drug just before FDA approval of Suboxone and Subutex.

BUPRENORPHINE	
Buprenorphine	10
Buprenorphine -3-D-Glucuronide 160	160
Cocaine HCl 80,000	80,000
Norbuprenorphine 10	10
Norbuprenorphine-3-D-Glucuronide 200	200

COCAINE (COC) - Drug 6

Cocaine is a potent central nervous system (CNS) stimulant and a local anesthetic. Initially, it brings about extreme energy and restlessness while gradually resulting in tremors, over-sensitivity and spasms. In large amounts, cocaine causes fever, unresponsiveness, difficulty in breathing and unconsciousness. Cocaine is often self-administered by nasal inhalation, intravenous injection and free-base smoking. It is excreted in the urine in a short time primarily as Benzoyllecgonine. 1.2 Benzoyllecgonine, a major metabolite of cocaine, has a longer biological half-life (5-8 hours) than cocaine (0.5-1.5 hours), and can generally be detected for 24-48 hours after cocaine exposure.

Cocaine 150	
Benzoyllecgonine	150
Cocaethylene	150
CocaineHCL	150

FENTANYL (FEN) - Drug 7

Fentanyl is a synthetic opioid. It has the brand names of Sublimaze, Actiq, Durogestic, Fentora and others. The Fentanyl drug is approximately 100 times more potent than morphine, with 100 micrograms of fentanyl approximately equivalent to 10 mg. of morphine or 75 mg. of meperidine in

analgesic activity. The Fentanyl drug is a potent narcotic analgesic with rapid onset and short duration of action. Historically, the fentanyl drug has been used to treat chronic breakthrough pain and is commonly used pre-procedures. Illicit use of pharmaceutical fentanyl drugs first appeared in the mid-1970s. Because the effects of the fentanyl drug last for only a very short time, it is even more addictive than heroin. Regular users may become addicted very quickly. The Fentanyl drug is much more potent than heroin, and tends to produce significantly worse respiratory depression, making it somewhat more dangerous than heroin to users. Overdose of the fentanyl drug has caused death. In the United States, the fentanyl drug is classified as a Schedule II controlled substance.

Fentanyl	
Norfentanyl	10
Fentanyl	50

MARIJUANA (THC) - Drug 8

THC (8.9-tetrahydrocannabinol) is the primary active ingredient in cannabinoids (marijuana). When smoked or orally administered, it produces euphoric effects. Users have impaired short term memory and slowed learning. They may also experience transient episodes of confusion and anxiety. Long term relatively heavy use may be associated with behavioral disorders. The peak effect of smoking marijuana occurs in 20-30 minutes and the duration is 90-120 minutes after one cigarette. Elevated levels of urinary metabolites are found within hours of exposure and remain detectable for 3-10 days after smoking. The main metabolite excreted in the urine is II -nor-A9-tetrahydrocannabinol-9-carboxylic acid (A9-TRC-COOR).

Marijuana 50	
11-nor- Δ^9 -THC-9-COOH	50
11-nor- Δ^8 -THC-9-COOH	50
Δ^8 - Tetrahydrocannabinol	15,000
Δ^9 - Tetrahydrocannabinol	50,000

Δ^8 - THC-COOH	50,000
-----------------------	--------

METHAMPHETAMINE (MET) - Drug 9

Methamphetamine is an addictive stimulant drug that strongly activates certain systems in the brain. Methamphetamine is closely related chemically to amphetamine, but the central nervous system effects of Methamphetamine are greater. Methamphetamine is made in illegal laboratories and has a high potential for abuse and dependence. The drug can be taken orally, injected, or inhaled. Acute higher doses lead to enhanced stimulation of the central nervous system and induce euphoria, alertness, reduced appetite, and a sense of increased energy and power. Cardiovascular responses to Methamphetamine include increased blood pressure and cardiac arrhythmias. More acute responses produce anxiety, paranoia, hallucinations, psychotic behavior, and eventually, depression and exhaustion. The effects of Methamphetamine generally last 2-4 hours and the drug has a half-life of 9-24 hours in the body. Methamphetamine is excreted in the urine as amphetamine and oxidized and delaminated derivatives. However, 10-20% of Methamphetamine is excreted unchanged. Thus, the presence of the parent compound in the urine indicates Methamphetamine use.

Methamphetamine 500	
d-Methamphetamine	500
Methamphetamine	500
3,4-methylenedioxymethamphetamine(MDMA)	1,250
3,4-Methylenedioxymethamphetamine(MDEA)	10,000
Fenfluramine	25,000
p-Hydroxymethamphetamine	5,000

METHADONE (MTD) - Drug 10

Methadone is a narcotic analgesic prescribed for the management of moderate to severe pain and for the treatment of opiate dependence (heroin, Vicodin, Percocet, Morphine). The pharmacology of Oral Methadone is very different from IV Methadone. Oral Methadone is partially stored in the liver for later use. IV Methadone acts more like heroin. In most states you must go to a pain clinic or a Methadone maintenance clinic to be prescribed Methadone. Methadone is a long acting pain reliever producing effects that last from twelve to forty-eight

hours. Ideally, Methadone frees the client from the pressures of obtaining illegal heroin, from the dangers of injection, and from the emotional roller coaster that most opiates produce. Methadone, if taken for long periods and at large doses, can lead to a very long withdrawal period. The withdrawals from Methadone are more prolonged and troublesome than those provoked by heroin cessation, yet the substitution and phased removal of methadone is an acceptable method of detoxification for patients and therapists.

METHADONE (MTD)	
Methadone	300

METHAQUALONE (MQL) - Drug11

Methaqualone (Quaalude, Sopor) is a quinazoline derivative that was first synthesized in 1951 and found clinically effective as a sedative and hypnotic in 1956. It soon gained popularity as a drug of abuse and in 1984 was removed from the US market due to extensive misuse. It is occasionally encountered in illicit form, and is also available in European countries in combination with diphenhydramine (Mandrax.). Methaqualone is extensively metabolized in vivo principally by hydroxylation at every possible position on the molecule. At least 12 metabolites have been identified in the urine.

METHAQUALONE	
Methaqualone	300

ECSTASY (MDMA) - Drug 12

Methylenedioxymethamphetamine (ecstasy) is a designer drug first synthesized in 1914 by a German drug company for the treatment of obesity. Those who take the drug frequently report adverse effects, such as increased muscle tension and sweating. MDMA is not clearly a stimulant, although it has, in common with amphetamine drugs, a capacity to increase blood pressure and heart rate. MDMA does produce some perceptual changes in the form of increased sensitivity to light, difficulty in focusing, and blurred vision in some users. Its mechanism of action is

thought to be via release of the neurotransmitter serotonin. MDMA may also release dopamine, although the general opinion is that this is a secondary effect of the drug (Nichols and Oberlender, 1990). The most pervasive effect of MDMA, occurring in virtually all people who took a reasonable dose of the drug, was to produce a clenching of the jaws.

Ecstasy 500	
3,4-Methylenedioxy-methamphetamine (MDMA)	500
D-Amphetamine	>100,000
3,4-methylenedioxy-amphetamine (MDA)	3,000
3,4-Methylenedioxy-ethylamphetamine (MDEA)	300
l-Methamphetamine	100,000

MORPHINE - Low Level (OP300) - Drug 13

Opiate refers to any drug that is derived from the opium poppy, including the natural products, morphine and codeine, and the semi-synthetic drugs such as heroin. Opioid is more general, referring to any drug that acts on the opioid receptor. Opioid analgesics comprise a large group of substances which control pain by depressing the central nervous system. Large doses of morphine can produce higher tolerance levels, physiological dependence in users, and may lead to substance abuse. Morphine is excreted unmetabolized, and is also the major metabolic product of codeine and heroin. Morphine is detectable in the urine for several days after an opiate dose.⁴

MORPHINE (MOP 300)	
Morphine	300
O6-Acetylmorphine	400
Codeine	300
Dihydrocodeine HCl	1,000
Ethylmorphine	100
Heroin	600

Hydromorphone	10,000
Hydrocodone	10,000
Levorphanol tartrate	10,000
Morphine-3-β-D-glucuronide	1,000
Nalorphine HCl	50,000
Thebaine	6,240

OXYCODONE (OXY) - Drug 14

Oxycodone, [4,5-epoxy-14-hydroxy-3-methoxy-17-methyl-morphinan-6-one, dihydrohydroxycodone] is a semi-synthetic opioid agonist derived from thebaine, a constituent of opium. Oxycodone is a Schedule II narcotic analgesic and is widely used in clinical medicine. The pharmacology of oxycodone is similar to that of morphine, in all respects, including its abuse and dependence liabilities. Pharmacological effects include analgesia, euphoria, feelings of relaxation, respiratory depression, constipation, papillary constriction, and cough suppression. Oxycodone is prescribed for the relief of moderate to high pain under pharmaceutical trade names as OxyContin® (controlled release), OxyIR® , OxyFast® (immediate release formulations), or Percodan® (aspirin) and Percocet® (acetaminophen) that are in combination with other nonnarcotic analgesics. Oxycodone's behavioral effects can last up to 5 hours. The controlled-release product, OxyContin® , has a longer duration of action (8-12 hours).

Oxycodone 100	
Oxycodone	100
Hydrocodone	5,000
Hydromorphone	5,000
Oxymorphone	200

PHENCYCLIDINE (PCP) - Drug 15

Phencyclidine, also known as PCP or Angel Dust, is a hallucinogen that was first marketed as a surgical anesthetic in the 1950s. It was removed from the market because patients receiving it became delirious and

experienced hallucinations. Phencyclidine is used in powder, capsule, and tablet form. The powder is either snorted or smoked after mixing it with marijuana or vegetable matter. Phencyclidine is most commonly administered by inhalation but can be used intravenously, intra-nasally, and orally. After low doses, the user thinks and acts swiftly and experiences mood swings from euphoria to depression. Self-injurious behavior is one of the devastating effects of Phencyclidine. PCP can be found in urine within 4 to 6 hours after use and will remain in urine for 7 to 14 days, depending on factors such as metabolic rate, user's age, weight, activity, and diet. Phencyclidine is excreted in the urine as an unchanged drug (4% to 19%) and conjugated metabolites (25% to 30%).

Phencyclidine 25	
Phencyclidine	25
4-Hydroxy Phencyclidine	90

TRAMADOL (TRA) - Drug 16

Tramadol is a quasi-narcotic analgesic used in the treatment of moderate to severe pain. It is a synthetic analog of codeine, but has a low binding affinity to the mu-opioid receptors. It has been prescribed off-label for the treatment of diabetic neuropathy and restless leg syndrome.² Large doses of Tramadol could develop tolerances and physiological dependency and lead to its abuse. Both 8. (d) and L forms of the isomers are controlled substances. Approximately 30% of the dose is excreted in the urine as unchanged drug, whereas 60% is excreted as metabolites. The major pathways appear to be N- and O- demethylation, glucuronidation or sulfation in the liver.

TRAMADOL (TRA)	
Tramadol	200
N-desmethyl-tramadol	500
O-desmethyl-tramadol	20,000

SYNTHETIC MARIJUANA (K2) - Drug 17

Synthetic Marijuana or K2 is a psychoactive herbal and chemical product that, when consumed, mimics the effects of Marijuana. It is best known by the brand names K2 and Spice, both of which have largely become genericized trademarks used to refer to any synthetic Marijuana product. The studies suggest that synthetic marijuana intoxication is associated with acute psychosis, worsening of previously stable psychotic disorders, and also may have the ability to trigger a chronic (long-term) psychotic disorder among vulnerable individuals such as those with a family history of mental illness. Elevated levels of urinary metabolites are found within hours of exposure and remain detectable for 72 hours after smoking (depending on usage/dosage). As of March 1, 2011, five cannabinoids, JWH-018, JWH-073, CP-47, JWH-200 and cannabicyclo hexanol are now illegal in the US because these substances have the potential to be extremely harmful and, therefore, pose an imminent hazard to the public safety. JWH-018 was developed and evaluated in basic scientific research to study structure activity relationships related to the cannabinoid receptors. JWH-073 has been identified in numerous herbal products, such as "Spice", "K1", "K3" and others. These products may be smoked for their psychoactive effects.

K2 (SYNTHETIC CANNABINOID)	
JWH-018 5-Peotanoic acid metabolite	50
JWH-018 5-Hydroxypentyl metabolite	500
JWH-018 4-Hydroxypentyl metabolite	400
JWH-018 N-(4-hydroxypentyl) metabolite solution	5,000
JWH-019 5-hydroxyhexylmetabolite	10,000
JWH-019 6-Hydroxyhexyl	5,000
JWH-073 4-butanoic acid metabolite	50
JWH-073 4-Hydroxybutyl metabolite	500

ETHYL GLUCURONIDE (ETG) - Drug 18

Ethyl Glucuronide (EtG) is a direct metabolite of ethanol alcohol. The presence of EtG in the urine can be used to detect recent alcohol consumption, even after the ethanol alcohol is no longer measurable. Consequently, the presence of EtG in the urine is a definitive indicator that

alcohol has been ingested. Traditional laboratory practices typically measure the amount of alcohol present in the body. Depending on the amount of alcohol that has been consumed, this method usually reveals alcohol ingestion within the past few hours. The presence of EtG in the urine, on the other hand, demonstrates that ethanol alcohol was ingested within the past three or four days, or roughly 80 hours after the ethanol alcohol has been metabolized by the body. As a result, it can be determined that a urine alcohol test employing EtG is a more accurate indicator of the recent consumption of alcohol as opposed to simply measuring for the existence of ethanol alcohol

ETHYL GLUCURONIDE	
Ethyl-P-D-glucuronide	500
Ethyl- -D-glucuronide-D5	500

SECTION 04

Adulteration SVT

The Adulterant Test Strip contains chemically treated reagent pads. Observation of the color change on the strip compared to the color chart provides a semi-quantitative screen for Oxidants, Specific Gravity, pH, Creatinine, Nitrite and Glutaraldehyde in human urine which can help to assess the integrity of the urine specimen.

Adulteration is the tampering of a urine specimen with the intention of altering the test results. The use of adulterants in the urine specimen can cause false negative results by either interfering with the test and/or destroying the drugs present in the urine. Dilution may also be used to produce false negative drug test results. To determine certain urinary characteristics such as specific gravity and pH, and to detect the presence of oxidants, Nitrite, Glutaraldehyde and Creatinine in urine are considered to be the best ways to test for adulteration or dilution.

Oxidants (OX): Tests for the presence of oxidizing agents such as bleach and peroxide in the urine.

Specific Gravity (S.G): Tests for sample dilution. Normal levels for specific gravity will range from 1.003 to 1.030. Specific gravity levels

of less than 1.003 or higher than 1.030 may be an indication of adulteration or specimen dilution.

pH: tests for the presence of acidic or alkaline adulterants in urine. Normal pH levels should be in the range of 4.0 to 9.0. Values below pH 4.0 or above pH 9.0 may indicate the sample has been altered.

SECTION 05

Collecting a Specimen

Allow the test cup to come to room temperature [15-30° C (59-86° F)] prior to test.

1) Tear the foil bag open, remove the test cup and disposable gloves provided for the donor. Label the device with donor information. (Fig. 1)



2) Open the test cup lid. Urinate directly into the test cup. Be sure to fill up the test cup with the urine specimen between minimum 30ml to maximum 110ml. (Fig. 2)

3) Close the lid securely and place the cup on a flat surface. Start the timer. (Fig. 3)

- 4) Put on the glove provided. Peel off label to reveal test result. (Fig. 4)
- 5) Read the adulteration strip at 2 minutes. Compare the colors on the adulteration strip to the enclosed color chart. If the result indicates adulteration, do not interpret or scan the drug test results. Either retest the urine or collect another specimen from the donor.
- 6) Read the drug strip results after 5 minutes. DO NOT INTERPRET RESULT OR SCAN AFTER 10 MINUTES. (Fig. 5)

NEGATIVE:* Two lines appear. One red line should be in the control region (C), and another apparent red or pink line adjacent should be in the test region (Drug/T). This negative result indicates that the drug concentration is below the detectable level
***NOTE:** The shade of red in the test line region (Drug(f) will vary, but it should be considered negative whenever there is even a faint pink line.

POSITIVE: One red line appears in the control region (C). No line appears in the test region (Drugtr). This positive result indicates that the drug concentration is above the detectable level.
INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test panel. If the problem persists, discontinue using the lot immediately and contact your manufacturer.
 Note: There is no meaning attributed to line color intensity or width.

A preliminary positive test result does not always mean a person took illegal drugs and a negative test result does not always mean a person did not take illegal drugs. There are a number of factors that influence the reliability of drug tests. Certain drugs of abuse tests are more accurate than others. **IMPORTANT:** The result you obtained is called preliminary for a reason. The sample must be tested by laboratory in order to determine if a drug of abuse is actually present. Send any sample which does not give a negative result to a laboratory for further testing.

What Is A False Positive Test?

The definition of a false positive test would be an instance where a substance is identified incorrectly by One Step Multi-Drug Screen Urine Test. The most common causes of a false positive test are cross reactants. Certain foods and medicines, diet plan drugs and nutritional supplements may cause a false positive test result with this product.

What Is A False Negative Test?

The definition of a false negative test is that the initial substance is present but isn't detected by One Step Multi-Drug Screen Urine Test. If the sample is diluted, or the sample is adulterated that may cause false negative results.

SECTION 06

Download Labb Station

The Labb Station platform is a proprietary software technology developed by Labb to make rapid drug testing simple, accurate, & more reliable than traditional point of care tests. To open an account for Labb Station you must be a professional drug testing facility, clinic, occupational health facility, and/or manage an in-house workplace drug testing program with a minimum of 25 tests a month.

If you are already a Labb Station client, you may download and use the Labb Station mobile app found on the Apple iOS app store and the Android Google Play store.

To download on the Apple iOS app store:




[Visit https://apps.apple.com/us/app/id1542054233](https://apps.apple.com/us/app/id1542054233)



Visit https://play.google.com/store/apps/details?id=systems.smartcreens.smartReader.v3.production&hl=en_US&gl=US
 SECTION 07

Labb.com

Labb offers a variety of drug testing products and services. The UA Series rapid drug screen is just on of many drug testing options available by Labb.



The Labb Station mobile app, tablet readers, and UA Series rapid drug screens are manufactured for Labb. Labb is a Virginia corporation located at:



Labb

820 Port Center Parkway
 Portsmouth, VA 23704

(877) 868-9967
 Hello@labb.com

Version3.00