



TESTING KITS FOR DRUGS

Procedure Card for IDT7000T

KETAMINE IDENTIFICATION

Substance Test



Read Before Use

As a precaution, gloves and eye protection should be worn when handling drug substances and using test units. Do not puncture any part of unit. Do not remove probe once a capsule is broken. In the unlikely event of a spill, rinse affected area thoroughly with water. Keep out of reach of children. All results are presumptive. If necessary, further testing by a competent drug testing lab can provide additional identification.

Test Preparation	1	Remove the Smart-Tip™ probe from test unit.	2 Sample suspected substance.	
			Powder	Liquid
		Swipe the paper of the Smart-Tip™ probe across the suspect substance.		Place the absorbent paper of Smart-Tip™ probe into liquid. Make sure paper is wet.

Test Preparation	3	Lock Smart-Tip™ probe into test unit. Press firmly.	STOP Do not continue unless probe is locked below the Flex Pins and cannot be removed from the test unit.	Identification	4	With thumb and forefinger , squeeze chamber 1 until capsule breaks.	5	Gently shake test unit for 2-3 seconds.

Identification	6	With thumb and forefinger , squeeze chamber 2 until capsule breaks.	7	Gently shake test unit for 2-3 seconds.

Results		The appearance of a dark blue-violet <u>particles</u> adhering to the sampler tip, the paper, side of the kit, or in the solution within one minute indicates the presence of Ketamine .
		If no dark blue-violet color develops within one minute, Ketamine is not present.
		* Also Note: If the liquid remains the dark red-purple of the reagent, this does not indicate the presence of Ketamine.



Disposal

The test kit is safe to dispose of in the trash.

0719



Manufactured by IDenta Ltd

www.identa-corp.com

Patented Technology

Made in Israel