







## METHAQUALONE (MANDRAX) IDENTIFICATION



### Substance Test






#### Read Before Use



As a precaution, gloves and eye protection should be worn when handling drug substances and using test units. Do not puncture any part of unit. Do not remove probe once a capsule is broken. In the unlikely event of a spill, rinse affected area thoroughly with water. Keep out of reach of children. All results are presumptive. If necessary, further testing by a competent drug testing lab can provide additional identification.


Test Preparation	1	Remove the Smart-Tip™ probe from test unit.	Sample suspected substance.				
			2	<table border="1"> <thead> <tr> <th>Powder</th> <th>Tablet</th> </tr> </thead> <tbody> <tr> <td>           Touch Smart-Tip™ probe to sample.   </td> <td>           Firmly press Smart-Tip™ probe on the tablet while scrubbing in a circular motion until the tip is completely covered with substance.   </td> </tr> </tbody> </table>	Powder	Tablet	Touch Smart-Tip™ probe to sample. 
Powder	Tablet						
Touch Smart-Tip™ probe to sample. 	Firmly press Smart-Tip™ probe on the tablet while scrubbing in a circular motion until the tip is completely covered with substance. 						

Test Preparation	3	Lock Smart-Tip™ probe into test unit. Press firmly.	Identification	4	5
		<div style="background-color: yellow; padding: 5px; text-align: center;"> <b>STOP</b>            Do not continue unless probe is locked below the Flex Pins and cannot be removed from the test unit.         </div>		With <b>thumb and forefinger</b> , squeeze chamber 1 until capsule breaks. 	Gently shake test unit for 2-3 seconds. 

Preliminary Identification	6	7
	With <b>thumb and forefinger</b> , squeeze chamber 2 until capsule breaks. 	Gently shake test unit for 10 seconds.  <div style="background-color: #e6e6fa; padding: 5px; margin-top: 10px;">             Wait up to 2 minutes for color to fully develop. The color may develop sooner.           </div>

Preliminary Results	<div style="background-color: yellow; width: 30px; height: 20px; display: inline-block;"></div> If a yellow color develops, the presence of Methaqualone (Mandrax) is possibly indicated. <b>Continue to Step 8.</b>
	<div style="border: 1px solid black; padding: 5px; display: inline-block;">            If no yellow color develops, Methaqualone (Mandrax) is <b>not</b> present. Proceed to Disposal.         </div>

Identification	8	9
	With <b>thumb and forefinger</b> , squeeze chamber 3 until capsule breaks. 	Vigorously shake test unit for 20 seconds.  <div style="background-color: #ffe4e1; padding: 5px; margin-top: 10px;">             Hold test unit steady and keep upright, until 2 layers are visibly separated. This may take up to 2 minutes. The process may occur sooner.           </div>

Results	<div style="background: linear-gradient(to bottom, yellow 50%, white 50%); width: 30px; height: 20px; display: inline-block;"></div> The appearance of a yellow color in the <u>lower</u> layer of liquid, while the upper aqueous layers remain colorless, strongly indicates the presence of <b>Methaqualone (Mandrax)</b> .
	<div style="border: 1px solid black; padding: 5px; display: inline-block;">            If the lower layer of liquid does not turn a yellow color, Methaqualone (Mandrax) is <b>not</b> present.         </div>



#### Disposal

The test kit is safe to dispose of in the trash.

01-19