

Caption Health

CAPTION AI™ Echocardiography Guidance Software



Operator's Manual

Part No. 739-00138 rev 15

Caption Health



Caption Health, Inc.
2000 Sierra Point Parkwa, 8th Floor
Brisbane, CA 94005 USA
www.captionhealth.com
415-671-4711



Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands



captionhealth.com/user-documents
Manual is also available on the tablet desktop

Product names are trademarks of their respective owners.

DICOM is a registered trademark of the National Electrical Manufacturers Association.

Patents pending.

This document and the information contained in it is Caption Health's proprietary and confidential information and may not be reproduced, copied in whole or in part, adapted, modified, disclosed to others, or disseminated without the prior written permission of Caption Health. This document is intended to be used either by customers, and is licensed to them as part of their Caption Health product purchase, or to meet regulatory commitments as required by the FDA under 21 CFR 892.2050 (and any amendments to it) and other local regulatory requirements. Use of this document by unauthorized persons is strictly prohibited.

Caption Health provides this document without warranty of any kind, implied or expressed, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose.

Caption Health has taken care to ensure the accuracy of this document. However, Caption Health assumes no liability for errors or omissions and

reserves the right to make changes without further notice to any products herein to improve reliability, function, or design. Caption Health may at any time make improvements or changes in the products or programs described in this document.

Copyright ©2015–2021 by Caption Health, Inc. All rights reserved. No part of this publication may be reproduced, transmitted, transcribed, stored in retrieval systems, or translated into any language or computer language, in any form or by any means, electronic, mechanical, magnetic, optical, chemical, manual, or otherwise, without the prior written permission of Caption Health, Inc. Caption Health, Inc. reserves the right to change its products and services at any time. In addition, this manual is subject to change without notice. Caption Health, Inc. welcomes customer input on corrections or recommendations for improvements to this document. Caption Health, Inc. has attempted to ensure accuracy in this manual. Caption Health, Inc. assumes no liability for any errors or omissions, nor for any damages resulting from the application or use of this information.

Part Number: 739-00138 rev 15

Date of Issue: 09/2021

Language: English

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

WARNING

Installation, use, and operation of this product are subject to the law in the jurisdictions in which the product is used. Install, use, and operate Caption AI™ only in a manner that does not conflict with applicable laws or regulations. Using the product for purposes other than those intended and expressly stated by Caption Health, Inc., as well as incorrect use or operation, may relieve Caption Health or its agents from all or some responsibility for resultant noncompliance, damage, or injury.

License Agreement

All computer programs copyright 2015–2021 by Caption Health, Inc. or its suppliers. Such programs are licensed under the following software agreement:

Caption Health, Inc. or its suppliers retain(s) ownership of and title to any computer program supplied with the Device and to the trade secrets embodied in such computer programs. Subject to the Buyer's acceptance and fulfillment of the obligations in this paragraph, Caption Health, Inc. grants the Buyer a personal, non-transferable, non-exclusive license to use any computer program supplied with the Device that is necessary to operate the Device solely on the medium in which such program is delivered for the purpose of operating the Device in accordance with the instructions set forth in the operator's manuals supplied with the Device and for no other purpose whatsoever. Buyer may not reverse-assemble, reverse-compile or otherwise reverse-engineer such computer programs nor may Buyer make a copy of such program or apply any techniques to derive the trade secrets embodied therein. In the event of a failure by Buyer to comply with the terms of this license, the license granted by this paragraph shall terminate. Further, because unauthorized use of such computer programs will leave Caption Health without an adequate remedy at law, Buyer agrees that injunctive or other equitable relief will be appropriate to restrain such use, threatened or actual. Buyer further agrees that (i) any of Caption Health, Inc. suppliers of software is a direct and intended beneficiary of this end-user sublicense and may enforce it directly against Buyer with respect to software supplied by such supplier, and (ii) NO SUPPLIER OF CAPTION HEALTH SHALL BE LIABLE TO BUYER FOR ANY GENERAL, SPECIAL, DIRECT, INDIRECT, CONSEQUENTIAL, INCIDENTAL OR OTHER DAMAGES ARISING OUT OF THE SUBLICENSE OF THE COMPUTER PROGRAMS SUPPLIED WITH THE DEVICE.

CONTENTS

Chapter 1: Introduction	11
About this manual	11
Conventions used in this manual	11
Cautions	12
Legal considerations	12
Image quality and diagnosis	12
Intended use/Indications for use	12
Other types of studies	13
Electrical shock hazard	13
Acoustic output and ALARA	13
Ultrasound-system hazards	14
Probe compatibility	14
Caption AI accuracy, use, and image saving	14
Changing uSmart Presets	15
Maintain probe contact during AutoCapture	15
Match manually saved clips to the active view	15
Workflow scanning protocol completion	15
Prescriptive Guidance instructions	15
Storing studies	16
PHI protection	16
Cybersecurity precautions and practices	16
Product compatibility	17
Intended audience for Caption AI	17
About Caption AI/UDI information	18
Specifications	19
Frequent tasks using Caption AI	19
Contacting Customer Support	19
Chapter 2: Clinical Studies and Non-Clinical Testing	21
Specialist (Sonographer) Study	21
Pivotal (Nurse) Study	21
Study Design	21
Results	22
Human Factors Validation Testing Study	26
Software Verification and Validation	27
Algorithm Performance Testing	27

Chapter 3: Operator-Safety Considerations	29
Prescription device statement.....	29
Healthy ergonomic practices	29
Rest breaks	30
Body position.....	30
Hand and wrist positions	30
Exercise and stretching.....	30
Chapter 4: Configuring Caption AI	31
System components.....	31
Electrical shock hazard	31
Connecting to the network	31
Starting Caption AI	32
Configuring Caption AI to send studies to the PACS server.....	32
Annotating exported clips with the view name	35
Configuring a Caption AI modality worklist	36
Chapter 5: Setting up for an exam	39
Working with scanning protocols	39
Creating custom scanning protocols	39
Turning off AutoCapture.....	44
Selecting the scanning convention	45
Chapter 6: Scanning with Caption AI	47
Starting a Caption AI study	47
Using modality worklists to populate exam information	50
The Caption AI ultrasound display and controls	51
Setting up your scan	55
Starting and optimizing your scan.....	56
Recording a clip.....	62
AutoCapture.....	62
Using Caption AI's Auto Ejection Fraction (AutoEF) feature	64
Save Best Clip.....	65
Manual recording	66
Non-sequential scanning	67
Recording a lung scan	69
Recording a clip based on Color Doppler.....	70
Recording a free-scan clip.....	71
Changing the scanning preset	73
Scanning with the linear probe	74
Switching the ultrasound probe	74
Reviewing acquired clips during an exam.....	75
Reviewing and sending a study.....	77

Chapter 7: Working with Saved Studies **83**

Displaying the list of saved studies	83
Reviewing and exporting saved studies to the DICOM PACS	85
Updating a patient's information	88
Deleting a study	90

Appendix A: Caption Interpretation AutoEF **A-1**

Introduction	A-1
Symbols and terms	A-1
Customer service.....	A-1
Safety information	A-2
Prescription Device Statement	A-2
Intended use/Indications for use	A-2
Intended audience	A-3
Measurement accuracy.....	A-3
Echocardiography ejection fraction	A-3
Ultrasound systems	A-4
Clip annotation and selection	A-4
Transthoracic and transesophageal echocardiography	A-5
Patient data archival.....	A-6
PHI protection.....	A-6
Cybersecurity precautions and practices	A-6
Product compatibility	A-7
Performance Testing	A-7
Validation testing and performance summary of AutoEF.....	A-7
Clip annotation performance	A-8
Clinical validation	A-9
System Overview and Configuration.....	A-18
System components and workflow	A-18
Hardware and software requirements.....	A-18
Configuration.....	A-18
System Operation.....	A-19
Caption Health AutoEF processing.....	A-19
Single-view EF calculation	A-20
Accompanying information.....	A-21
Reviewing results on the PACS workstation	A-25
Editing results	A-25
System specifications	A-25
Safety and regulatory requirements	A-26

Appendix B: Caption AI Scanning Metadata **B-1**

List of Tables

Table 1.	Product Label Symbols	18
Table 2.	Study results: primary endpoints	22
Table 3.	Study results: secondary endpoints	23
Table 4.	Study results: diagnostic-quality clips.....	24
Table 5.	Prescriptive Guidance Examples	58
Table 6.	Available Caption AI presets	73
Table 7.	Clip Annotator positive predictive value (PPV) and sensitivity	A-9
Table 8.	Caption Interpretation Automated EF performance compared to the reference standard – primary hypothesis testing	A-10
Table 9.	Caption Interpretation Automated EF performance compared to the reference standard – secondary hypothesis testing	A-11
Table 10.	Caption Interpretation Automated EF performance compared to the reference standard – descriptive analysis.....	A-12
Table 11.	EF Qualitative bin assessments	A-12
Table 12.	Paired-difference analysis (example).....	A-13
Table 13.	Caption Interpretation Automated EF qualitative comparison.....	A-14
Table 14.	Caption Interpretation qualitative comparison – hypothesis testing	A-14
Table 15.	Sonographer and cardiologist variability.....	A-15
Table 16.	AutoEF 2.0 PPV (Reduced/ Severely Reduced)	A-16
Table 17.	AutoEF 2.0 sensitivity	A-16
Table 18.	EF assessment outliers – Caption Interpretation AutoEF	A-17
Table 19.	EF assessment outliers – physicians	A-17
Table 20.	Recommendations for cardiac chamber quantification by echocardiography in adults (ASE/AACVI 2015): left ventricular ejection fraction	A-23
Table 21.	AutoEF 2.0 qualitative categorization of EF	A-24
Table 22.	System specifications	A-25

List of Figures

Figure 1.	Caption AI start screen.....	32
Figure 2.	Export Destination dialog: configuring Caption AI	33
Figure 3.	Settings screen: Annotations	35
Figure 4.	Exported clip with annotated view name	35
Figure 5.	Settings screen: Worklist	36
Figure 6.	Modality Worklist Connection dialog	37
Figure 7.	Scanning Protocols screen.....	40
Figure 8.	Scanning Protocols screen: Add custom workflow protocol.....	40
Figure 9.	Create Free Scan dialog.....	41
Figure 10.	Scanning Protocols screen: Indications tab	42
Figure 11.	Scanning Protocols screen: Add custom indications protocol.....	42
Figure 12.	Scanning Protocols screen: Edit protocol	43
Figure 13.	Acquisition screen: AutoCapture	44
Figure 14.	Acquisition screen: Scanning convention	45
Figure 15.	Caption AI start screen: Workflow-based scan	48
Figure 16.	Caption AI start screen: Indication-based scan	48
Figure 17.	Connect Probe message	49
Figure 18.	Modality Worklist screen	50
Figure 19.	Modality Worklist screen: search results	51
Figure 20.	Basic ultrasound display and controls	52
Figure 21.	Details panel	52
Figure 22.	All Views list	53
Figure 23.	End exam dialog	54
Figure 24.	Ultrasound display: features for setting up to scan	55
Figure 25.	Ultrasound display: real-time guidance.....	56
Figure 26.	Quality Meter: basic features	56
Figure 27.	Dual Quality Meter.....	57
Figure 28.	Visualization for Prescriptive Guidance probe-movement icons	58
Figure 29.	Resume scanning dialog.....	62
Figure 30.	Quality Meter: AutoCapture states (home view)	64
Figure 31.	AutoEF estimate: single view	64
Figure 32.	Save Best Clip feature	65
Figure 33.	Record button: manual recording.....	66
Figure 34.	Record button: AutoCapture turned off.....	67
Figure 35.	All Views list	68
Figure 36.	Scanning screen: Lung protocol	69
Figure 37.	Scanning screen: Color Doppler mode	70
Figure 38.	Moving Color Doppler area of focus	71
Figure 39.	Free Scan mode	72
Figure 40.	Free-Scan clip acquired.....	72
Figure 41.	Scanning presets menu	73
Figure 42.	Connect Probe message	74
Figure 43.	Switch Probe message: Linear probe	74
Figure 44.	Switch Probe message: AI-guided cardiac scanning	75

Figure 45. All Views list for active study.....	76
Figure 46. Report screen	77
Figure 47. Lung report screen	79
Figure 48. Lung Zone gallery.....	80
Figure 49. Export Studies dialog: sending a just-completed study	81
Figure 50. Report Not Sent message: sending a just-completed study	81
Figure 51. Saved Studies screen	84
Figure 52. Saved Studies screen: multiple studies selected	85
Figure 53. Export Studies dialog: selected studies	86
Figure 54. Exporting Studies screen	86
Figure 55. Stop Export dialog	87
Figure 56. Patient Details screen	89
Figure 57. Delete study dialog	90
Figure 58. Bland-Altman (left) and Linear Fit (right) plots for Caption Interpretation Automated EF estimation based on best available view(s)	A-11
Figure 59. Detail view with AutoEF	A-20
Figure 60. AutoEF calculation confirmation.....	A-21
Figure 61. Detailed view with single AutoEF calculation	A-21
Figure 62. Report screen with AutoEF estimate: expected error range	A-22
Figure 63. Report screen with AutoEF estimate: global LV function	A-23
Figure 64. Report screen with AutoEF estimate: likelihood of LV function	A-24
Figure 65. Metadata Export screen.....	B-1

CHAPTER 1: INTRODUCTION

About this manual

This manual accompanies **Caption AI™** and provides information on configuring Caption AI and using it to capture diagnostic-quality ultrasound images. This manual is for healthcare professionals who have received appropriate training on ultrasound basics and on using Caption AI, provided either by Caption Health or by a trained clinician using official Caption AI training materials. Please ensure that you do so and that before using Caption AI you read this operator's manual in its entirety.

Note: This manual describes the most extensive configuration of Caption AI, with the maximum number of options and accessories available. Some functions described in this manual may not be included in your product's configuration.

Conventions used in this manual

This manual uses the following symbols and text to indicate specific types of information:



This symbol and term alert you to information regarding patient, operator, or equipment safety. They also indicate information about preventing the loss of patient or product data.

Note: Indicates supplemental information.



Indicates step-by-step instructions.

Cautions

The following is important information regarding patient, operator, or equipment safety, as well as preventing the loss of patient or product data or damage to the Caption AI device or environment.



Keep the Caption AI system outside the MRI Scanner room.



CAUTION

Legal considerations

Do not use the product for purposes other than those intended and expressly stated by Caption Health, Inc. Do not misuse the product, and do not use or operate the product in an incorrect manner.

Installation, use, and operation of this product are subject to the law in the jurisdictions in which the product is used. Install, use, and operate Caption AI only in a manner that does not conflict with applicable laws or regulations, which have the force of law.

Use of the product for purposes other than those intended and expressly stated by Caption Health, Inc., as well as incorrect use or operation, may relieve Caption Health or its agents from all or some responsibility for resultant noncompliance, damage, or injury.



CAUTION

Image quality and diagnosis

Product users are responsible for image quality and diagnosis. The images acquired using Caption AI are to be interpreted only by qualified medical professionals. A qualified medical professional must inspect the data being used for analysis and diagnosis, and ensure that the data is sufficient and appropriate in anatomical correctness and both spatial and temporal resolution for the measurement being employed.



CAUTION

Intended use/Indications for use

The Caption AI software is intended to assist medical professionals in the acquisition of cardiac ultrasound images. Caption AI software is an accessory to compatible general purpose diagnostic ultrasound systems.

The Caption AI software is indicated for use in two-dimensional transthoracic echocardiography (2D-TTE) for adult patients, specifically in the acquisition of the following standard views: Parasternal Long-Axis (PLAX), Parasternal Short-Axis at the Aortic Valve (PSAX-AV), Parasternal Short-Axis at the Mitral Valve (PSAX-MV), Parasternal Short-Axis at the Papillary Muscle (PSAX-PM), Apical 4-Chamber (AP4), Apical 5-Chamber (AP5), Apical 2-Chamber (AP2), Apical 3-Chamber (AP3), Subcostal 4-Chamber (SubC4), and Subcostal Inferior Vena Cava (SubC-IVC).

CAUTION

Other types of studies

Caption AI is not intended for transesophageal echocardiography or any other type of ultrasound study not listed under "Intended Use/Indications for use."

CAUTION

Electrical shock hazard

Shock hazards exist if the AC power adapter for the uSmart 3200T Plus ultrasound system is damaged or is not properly grounded. Use only the supplied medical-grade power adapters and power cords and connect only to receptacles marked "Hospital Only," "Hospital Grade," or equivalent.

CAUTION

Acoustic output and ALARA

The uSmart 3200T Plus ultrasound system, with which Caption AI functions as a software accessory, complies with the Standard for Real-Time Display of Thermal and Mechanical Acoustic Output Indices on Diagnostic Ultrasound Equipment (UD3-2004).

It is important to observe the MI and TI values that are displayed during scanning. When conducting ultrasound studies, follow the ALARA principle: expose the patient to ultrasound energy at a level that is "As Low As Reasonably Achievable." Please read the *Terason uSmart3200T Ultrasound System User Guide Volume 2* for important information about acoustic output and ALARA.

 **CAUTION****Ultrasound-system hazards**

General ultrasound-system hazards not specific to Caption AI are described in the *Terason uSmart3200T Ultrasound System User Guide Volume 2*. You must carefully read the *Terason uSmart3200T Ultrasound System User Guide, Volumes 1 and 2* before using Caption AI.

 **CAUTION****Probe compatibility**

Caption AI operates with the 3200T-compatible Terason 4V2A linear phased array probe (called a “transducer” in the Terason manual). Do not operate Caption AI with any other probe during AI-assisted scanning.

 **CAUTION****Caption AI accuracy, use, and image saving**

Caption AI provides real-time guidance and automated capture during cardiac ultrasound (echocardiographic) examinations for the 10 standard echocardiographic views. The accuracy of Caption AI in classifying correct echocardiographic views and estimating image quality has been verified and validated, but individual patient variations may introduce errors. Also, many echocardiographic views are similar to other views, and this can introduce errors. As a result, the automated saving functions of Caption AI (AutoCapture and Save Best Clip) may occasionally contain errors. It is important to review saved images independently using experienced clinical judgment prior to making a diagnosis. This may be especially important with saved clips noted as being saved through the Save Best Clip function because in these instances Caption AI did not detect a clip of sufficient quality to meet the AutoCapture threshold. It is also possible that some images obtained during scanning may be correct and have sufficient image quality for diagnosis, but Caption AI does not recognize this and does not auto-capture these images, and using the Save Best Clip option does not capture these images either. In these instances, manual saves may be necessary.

Making a diagnosis based solely on Caption AI, without applying clinical judgment regarding view correctness and quality, is not recommended.



Changing uSmart Presets

The uSmart presets should not be changed for guided modes. Only users that have received the appropriate training should make changes to the uSmart presets, and should only do so during free-scan modes.



Maintain probe contact during AutoCapture

It is important to maintain the probe contact and position during AutoCapture to ensure the capture of a good clip. If you move the probe, remove it from the patient, or otherwise interrupt the recording, the clip will not be recorded.



Match manually saved clips to the active view

When you perform a manual clip save, the system will assign the active view label to the saved clip. Confirm that the view you are saving manually is the view listed in order to avoid a mismatch between the label and the actual view captured.



Workflow scanning protocol completion

Caption AI provides workflow scanning protocols that guide you through a standard echocardiographic or lung study, including saving clips for multiple views. If you terminate a study prior to scanning and saving clips for all the views, the diagnostic capability of the study may be compromised.



Prescriptive Guidance instructions

Caption AI provides Prescriptive Guidance that gives you suggestions on how to manipulate the probe to capture the desired images for an echocardiographic study. These suggestions have been verified and validated, but individual patient variations may cause the instructions to perform better in some patients than others. It is recommended that you pay attention to the Quality Meter as you

scan and make adjustments to the probe position to capture diagnostic-quality studies.

 CAUTION

Storing studies

Caption AI is not for permanent storage for patient studies. Caption AI is meant to be a temporary storage for studies you acquire. All patient studies should be archived in your permanent location.

 CAUTION

PHI protection

DICOM studies contain Protected Health Information (PHI). Caption AI can operate within a PACS network and can send patient studies to a DICOM PACS server. To protect PHI, only authorized IT administrators who log into the system using their security credentials should configure Caption AI to communicate with the PACS server.

Caption AI is not intended for long-term storage of studies.

This device adheres to HIPAA security and privacy guidelines. When viewing Caption AI results in your DICOM Viewer, make certain to observe your institution's guidelines and practices for protecting PHI.

Please refer to the *Terason uSmart 3200T Plus User Guide* for important information about PHI protection.

 CAUTION

Cybersecurity precautions and practices

Malware, computer viruses, ransomware, and other cybersecurity threats are an increasing concern in healthcare IT systems.

Please refer to the *Terason uSmart 3200T Plus User Guide* for important information about cybersecurity provisions and management.

For information and guidance on implementing proper cybersecurity in the healthcare IT environment, see "Health IT Privacy and Security Resources for Providers" at <https://www.healthit.gov/topic/privacy-security-and-hipaa/health-it-privacy-and-security-resources-providers>.

 **CAUTION****Product compatibility**

Caption Health, Inc. Caption AI Software Accessory is a software accessory that operates with the Terason uSmart 3200T Plus ultrasound system. It is not designed to operate with any other ultrasound systems. Do not attempt to operate Caption AI with other ultrasound systems.

Do not use your system in combination with other products or components unless Caption Health expressly recognizes those other products or components as compatible. For information about such products and components, see "System components" on page 31 or contact your Caption Health representative.

Changes and additions to the system should be made only by Caption Health or by third parties expressly authorized by Caption Health to do so. Such changes and additions must comply with all applicable laws and regulations that have the force of law within the jurisdictions concerned and with best engineering practices.

Intended audience for Caption AI


Caption AI must be operated by or under the direction and supervision of a licensed physician.


 **CAUTION**

Before using Caption AI, read this manual and strictly observe all warnings and cautions.

About Caption AI/UDI information






► To display the About Us screen and UDI information




1. On the Caption AI **start** screen, tap the **Menu** icon (). The **Main menu** appears.
2. Tap **About**. The **About Us** screen appears.



The symbols on this screen are described in the table below.

Table 1. *Product Label Symbols*

Symbol	Definition
	Consult the instructions for use.
	This is a medical device.
	Operators must read the manual before using the device.
GMDN	The global medical device nomenclature number.
	The product batch code ("lot number").
	The product catalog number ("reference number").
Quantity	A single quantity of this device is used to perform the functions associated with the device.

Symbol	Definition
Rx Only	Federal (United States) law restricts this product to sale by or on the order of a physician.
	The unique device identification number.
	The product manufacturer, including address.
	The date of manufacturer.

Specifications

Please refer to the *Terason uSmart3200T Ultrasound System User Guide Volume 2* for information regarding probe and system specifications. Some of the specifications listed may not be available on your system. The display frame rate for Caption AI is 30 frames per second.

Frequent tasks using Caption AI

Task	Refer to
Set up Caption AI to communicate with the DICOM PACS	Page 32: "Configuring Caption AI to send studies to the PACS server"
Set up Caption AI for scanning	Page 39: <i>Chapter 5: Setting up for an exam</i>
Perform an ultrasound exam with Caption AI	Page 47: <i>Chapter 6: Scanning with Caption AI™</i>
Review and send saved studies	Page 77: "Reviewing and sending a study"
Delete studies from the Caption AI device	Page 90: "Deleting a study"

Contacting Customer Support

Customer service representatives are available to answer questions and provide support to help you use the product. Please contact your local Caption Health representative for assistance or email support@captionhealth.com.

This page intentionally left blank.

CHAPTER 2: CLINICAL STUDIES AND NON-CLINICAL TESTING

Specialist (Sonographer) Study

Use by specialist-users was evaluated in a prospective clinical study, in which 50 patients were scanned by sonographers with the Caption AI™ system (study exam), followed by a reference scan (control exam) on the same patient using the same ultrasound but unassisted by Caption AI. Study and control exams were assessed by three (3) expert cardiologists who graded each clip using the ACEP scale. The data from this study was used to provide descriptive supportive evidence of the use of Caption AI by users with specialized echocardiography training.

The results of the study indicated that sonographers obtained diagnostic-quality images in a high proportion of clips in both study and control exams, demonstrating comparable (non-inferior) image quality in clips acquired using Caption AI compared to unassisted acquisition. Importantly, the high specificity of the AutoCapture feature (97.85% of auto-captured clips were diagnostic) demonstrated that a registered sonographer can confidently rely on the AutoCapture feature when using Caption AI.

Pivotal (Nurse) Study

A prospective non-specialist study was subsequently conducted to evaluate use by medical professionals without specialized echocardiography training.

Study Design

A minimum of eight (8) RNs were to be trained and evaluated on their performance to acquire a 10-view 2D-TTE protocol. Participants were scanned by the RN (study exam) and 10 standard views were obtained using a Terason ultrasound system with Caption AI software: PLAX, PSAX-AV, PSAX-MV, PSAX-PM, AP4, AP5, AP2, AP3, SubC4, and SubC-IVC.

The study continued enrollment until eight RNs had completed scans of 30 patients each. Enrolled patients were stratified into two groups based on cardiac abnormalities to ensure a sufficient number of patients with cardiac abnormalities. In addition, enrolled patients were evenly stratified into three groups based on BMI to ensure a sufficient distribution of patients by BMI. For

the sake of comparison, participants were also scanned by a trained sonographer without Caption AI and the same 10 views were obtained (control exam) using the same Terason ultrasound system.

Following the study and control exams, a panel of five (5) expert cardiologist readers independently provided assessments of whether the patient study, in its totality, provided sufficient information to assess 10 clinical parameters. In addition, a panel of eight (8) expert cardiologists also independently provided assessments of diagnostic image quality per clip using the ACEP scale; each clip was graded by five (5) expert cardiologists. The readers were blinded to assessments from other panel members as well as to which site the images were obtained at and whether the images were obtained by an RN or a sonographer. The results from the expert-panel reads were used for the statistical analysis. To reduce possible sources of bias in the design, the RNs, sonographers, and cardiologists were all blinded to results determined by others.

Four (4) prospectively defined primary endpoints were evaluated sequentially for the study, all of which assessed whether the patient study exam conducted by the RN, taken as a whole, was of sufficient image quality to make these clinical assessments. Specifically, the endpoints assessed whether, in the judgment of expert cardiologists, the studies permitted qualitative visual assessment of left ventricular (LV) size, LV function, right ventricular (RV) size, and non-trivial pericardial effusion.

Results

The four primary endpoints were satisfied and demonstrated the clinical utility of Caption AI for non-specialist users. Specifically, the eight (8) RNs acquired echocardiographic exams of sufficient image quality to make clinical assessments in the proportion of study exams conducted, as shown below.

Table 2. *Study results: primary endpoints*

Endpoint	Percent of diagnostic quality
Qualitative Visual Assessment of LV Size	98.8% (95% CI: 96.7, 100)
Qualitative Visual Assessment of LV Function	98.8% (95% CI: 96.7, 100)
Qualitative Visual Assessment of RV Size	92.5% (95% CI: 88.1, 96.9)
Qualitative Visual Assessment of Non-Trivial Pericardial Effusion	98.8% (95% CI: 96.7, 100)

Note: Secondary endpoints and additional analyses presented below were not evaluated based on a specific hypothesis. Since the evaluation of secondary endpoints and additional analyses did not allow for control of Type I error, the study results are presented as a descriptive demonstration of the use of Caption AI for the specific secondary endpoints and additional analyses.

Additional secondary endpoints were evaluated and demonstrated the robustness of the data, including six (6) additional patient-level clinical parameters were evaluated and each had a high proportion of scans considered to be of sufficient image quality to make each of the six (6) additional patient-level clinical parameter assessments, i.e., qualitative visual assessment of inferior vena cava (IVC) size, RV function, left atrial (LA) size, aortic valve (AV), mitral valve (MV), and tricuspid valve (TV). Specifically, the eight (8) RNs acquired echocardiographic exams of sufficient image quality to make clinical assessments in the proportion of study exams conducted, as shown below.

Table 3. *Study results: secondary endpoints*

Endpoint	Percent of diagnostic quality
Qualitative Visual Assessment of RV Function	91.3% (95% CI: 85.7, 96.8)
Qualitative Visual Assessment of LA Size	94.6% (95% CI: 90.7, 98.5)
Qualitative Visual Assessment of AV	91.7% (95% CI: 88.0, 95.3)
Qualitative Visual Assessment of MV	96.3% (95% CI: 93.9, 98.6)
Qualitative Visual Assessment of TV	83.3% (95% CI: 77.0, 89.7)
Qualitative Visual Assessment of IVC Size	57.5% (95% CI: 41.5, 73.5)

In addition to assessing if image quality was sufficient to make assessments, cardiologists also made specific qualitative visual assessments based on the study and control exams (e.g., presence or absence of non-trivial pericardial effusion). The proportion of scans in which the diagnostic decision was the same between study and control exams was very high, further demonstrating the usability of Caption AI. For primary clinical parameters, the range was 92.5% to 99.6%. Similarly, for secondary clinical parameters, the range was 83.2% to 95.2%.

To provide a robust assessment of the performance of Caption AI, subjects enrolled in the study included a broad range of patient characteristics representative of the intended use population. In particular, effort was made

to include subjects with known cardiac abnormalities at time of enrollment (63.9%), which would be expected to provide a more technically difficult scan. In fact, the standard of care echo revealed a much higher proportion of known cardiac abnormalities (91.4%). In addition, many of the study patients were inpatient or had other challenging characteristics such as high BMI, history of smoking (42.2%), and cardiac implantables (24.6%).

Subgroup analyses were conducted to evaluate the impact of specific baseline and demographic characteristics (i.e., BMI, presence of known cardiac abnormalities, sex, age, scan sequence number within each acquiring nurse, RN user, and study site) on the outcomes of the primary and secondary endpoints. The results demonstrated consistent performance across subgroups.

Furthermore, it was evaluated whether the RN users were able to obtain a high proportion of clips that were considered of diagnostic quality. Specifically, the eight (8) RNs acquired echocardiographic clips of diagnostic-image quality for each of the standard views in the following proportion of study exams conducted.

Table 4. Study results: diagnostic-quality clips¹

View	Percent of diagnostic quality
PLAX	92.1% (95% CI: 87.9, 96.3)
PSAX-AV	66.3% (95% CI: 59.0, 73.5)
PSAX-MV	75.8% (95% CI: 70.7, 80.9)
PSAX-PM	92.9% (95% CI: 89.1, 96.7)
AP4	88.8% (95% CI: 81.5, 96.0)
AP5	78.8% (95% CI: 66.9, 90.6)
AP2	71.3% (95% CI: 61.6, 80.9)
AP3	80.0% (95% CI: 70.4, 89.7)
SC4	76.3% (95% CI: 70.2, 82.3)
SC-IVC	59.2% (95% CI: 43.1, 75.2)

The study also demonstrated the safety profile of Caption AI. No device-related adverse events were reported in the pivotal study.

¹ Median acquisition time per view ranged from 2.04 to 4.03 minutes.

The following two analyses were performed on the Pivotal Study data to assess the performance of Caption AI using objective quantitative metrics:

- **AutoEF:** The exams acquired in the pivotal study were processed by a previous version of the 510(k) cleared Caption Health AutoEF software (K173780). This version of AutoEF provided an automated estimation of left ventricular ejection fraction, and requires an AP4 and AP2 clip of sufficient quality in order to return an estimate. AutoEF returned an EF estimate in 65.5% of study exams and 85.1% of control exams. The nurse-sonographer root-mean-square deviation was 5.19 EF% and mean absolute deviation was 3.96 EF%, indicating a clinically acceptable amount of variability and within what might be expected between experts. Therefore, the results indicate that there is no clinically significant difference in AutoEF estimation for nurse- and sonographer-acquired exams, provided there is an AP2 and AP4 of sufficient quality as assessed by the AutoEF software.

Note: The version of AutoEF included in the Caption AI product is updated from the version used in the pivotal study analysis described above. The updated AutoEF includes new features, such as producing an EF estimate from any combination of the PLAX, AP4, and AP2 views, as well as a more accurate estimation algorithm. Nonetheless, the analysis above provides useful quantitative information supporting the effectiveness of Caption AI when used by non-specialist nurses.

- **PLAX Sonographer Measurements:** Three (3) registered cardiac diagnostic sonographers independently provided measurements for each of the PLAX clips acquired in the pivotal study: septal wall thickness (diastole), posterior wall thickness (diastole), left ventricular internal diameter (diastole), left ventricular internal diameter (systole), and aortic root. Measurability of study exam clips ranged from 89.17% to 92.08%. The study-control clip variability (RMSD) was found to be comparable to the inter-sonographer measurement variability. The results demonstrate that PLAX clips acquired by nurses were highly suitable for linear measurements in clinical use.

Human Factors Validation Testing Study

A Human Factors Validation Study was performed with a total of 28 participants to demonstrate the usability of the device. The Caption AI user interface and training materials were developed through a series of preliminary human factors analyses. The device and training were then implemented and tested during the Human Factors Validation Study.

The study enrolled five (5) user groups:

- Physicians (Hospitalists): 4 users
- Nurse Practitioners and Physician Assistants: 4 users
- Registered Nurses (RNs): 9 users
- Medical Residents: 6 users
- Certified Medical Assistants: 5 users

The human factors usability validation testing on 28 users demonstrated 100% of critical tasks passed across all user types. No use errors were found that could cause harm for the scanner nor the patient. In addition, non-critical tasks that are essential for skills assessment and scanning ability were tested and evaluated, including the AutoCapture rate across 10 views. An average of 76% of views were auto-captured per user across all 5 user groups. The results did not indicate significant differences in performance amongst user types. Given the same level of training, medical professionals across the backgrounds that were tested from certified medical assistants to hospitalists performed consistently during independent scans with Caption AI on a patient model.

Supplemental human-factors validation was performed to validate the usability of the Prescriptive Guidance icons visual system, using participants with varying levels of TTE scanning experience. This study tested the accuracy of movements made by 16 users while following the on-screen guidance. Correct movements were made 89.9% of the time, which validated the usability of the visual system to help guide users to a better-quality image during scanning.

Software Verification and Validation

Caption AI was identified as having a moderate level of concern as defined in the FDA guidance document, "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices." The software documentation included:

1. Software/Firmware Description
2. Device Hazard Analysis
3. Software Requirement Specifications
4. Architecture Design Chart
5. Software Design Specifications
6. Traceability
7. Software Development Environment Description
8. Verification and Validation Documentation
9. Revision Level History
10. Unresolved Anomalies
11. Cybersecurity

A comprehensive risk analysis was provided for the software with detailed description of the hazards, their causes and severity, as well as acceptable methods for control of the identified hazards. Caption Health provided a description, with test protocols including pass/fail criteria and report of results, of acceptable verification and validation activities at the unit, integration, and system level.

Algorithm Performance Testing

Comprehensive non-clinical performance testing of the deep-learning algorithms used in the device was provided to support their clinical performance. Specifically, the performance testing evaluated the performance of the following software functionality:

- Standalone performance
- Feature-level performance:
 - Quality Meter
 - Save Best Clip
 - AutoCapture
 - Prescriptive Guidance

This page intentionally left blank.

CHAPTER 3: OPERATOR-SAFETY CONSIDERATIONS

Prescription device statement

CAUTION

For product usage in the United States of America, the following labeling statement applies:

Federal law restricts this device to use by or on the order of a physician.

Healthy ergonomic practices

The operation of an ultrasound system may be linked to musculoskeletal disorders (MSDs). In ultrasound imaging, ergonomics may be defined as the physical interaction between the operator, the system, and the probe during exams. It is important for the operator of the system to practice good ergonomic techniques to reduce the risk of injury. This section provides guidelines to help you work more comfortably to possibly reduce the risk of musculoskeletal disorders.

When using an ultrasound system, as with many similar physical activities, you may experience occasional discomfort in your hands, fingers, arms, shoulders, eyes, neck, back, or other parts of your body. If you experience symptoms such as constant or recurring discomfort, pain, throbbing, aching, tingling, numbness, burning sensation, or stiffness, you should promptly consult a qualified health professional. These may be symptoms of MSD. MSDs can be painful and may result in potentially disabling injuries to the nerves, muscles, tendons, joints, or other parts of the body. Examples of MSDs include carpal tunnel syndrome and tendinitis.

Described here are steps you can take to guard against discomfort while scanning or the risk of MSDs. It is also recommended that you consult the guidelines of professional medical societies concerned with ultrasound.

Rest breaks

- In between exams, take breaks to rest and give your body a chance to recuperate from the strained positions and repetitive movements of examinations.
- While scanning, avoid maintaining the same body position for extended periods of time by moving and varying your head, neck, body, arm, and leg positions.

Body position

- Avoid bending or stooping.
- Adjust the position of the device so that viewing or reaching the Caption AI controls does not require strained or awkward body positions.
- Whenever possible, use an adjustable chair with good back support, and adjust the chair height to promote good body posture. If possible, adjust the height of the patient bed to optimize body posture.
- Maintain a comfortable and balanced body position with minimal stress on your joints, minimizing bending and twisting.
- Keep elbows close to your side and relax your shoulders in a level position.

Hand and wrist positions

- Do not grasp the probe with excessive force; hold it as lightly as possible.
- Minimize the amount of pressure applied when pressing the probe against the patient.
- Avoid or minimize bending your wrist.

Exercise and stretching

Targeted exercises and stretching may help you avoid the risk of MSDs. Consult with a qualified health professional to define a program suited to your needs.

CHAPTER 4: CONFIGURING CAPTION AI

System components

Caption AI™ is a software-only device for assisting medical professions in performing echocardiographic studies. Caption AI is used in conjunction with the following third-party products:

- Terason uSmart 3200T Plus (3200T) ultrasound system
- Terason 4V2A linear phased array probe
- Terason 15L4A linear array probe
- Terason 3200T NexGen mobile cart
- Bluetooth keyboard for data entry
- Caption AI for performing and recording echocardiographic ultrasound studies

Your customer service representatives can provide support for setting up Caption AI. Please contact your local Caption Health representative for assistance.

CAUTION

Electrical shock hazard

Shock hazards exist if the AC power adapter for the uSmart 3200T Plus ultrasound system is damaged or is not properly grounded. Use only the supplied medical-grade power adapters and power cords and connect only to receptacles marked "Hospital Only," "Hospital Grade," or equivalent.

Connecting to the network

To use DICOM features, setup includes establishing a secure connection between Caption AI and the DICOM PACS network.

CAUTION

Configuration and setup must be done in compliance with your network security policies.

Starting Caption AI

► To start Caption AI

Double-tap the **Caption AI** icon () on the desktop. The Caption AI **start** screen appears.

When the app starts, the Caption AI **start** screen appears.




Figure 1. Caption AI start screen

Configuring Caption AI to send studies to the PACS server

To configure Caption AI to send studies to a PACS server, set up an export destination in the settings. The export destination is the permanent-storage location that the PACS accesses. Please contact Caption Health for the Caption AI DICOM Conformance Statement.

CAUTION


To ensure proper operation of Caption AI for patient scanning, system setup and configuration should be performed by your ultrasound service representative or a qualified biomedical engineer.

CAUTION

Only an authorized IT administrator who logs into the system using their security credentials should configure Caption AI to communicate with the PACS server.

► **To set the export destination**

1. On the Caption AI **start** screen, tap the **Menu** icon (☰), then **Settings**. The **Settings** screen appears.
2. Select the **Export** tab.



3. Next to **Export Destination**, tap **Add**. The **Export Destination** dialog appears. You must fill in all the fields.

 A screenshot of the 'Export Destination' dialog box. The dialog has a title bar with 'Export Destination' and a close button (X). It contains several input fields: 'Tablet AE Title' with the value 'Origin AE Title', 'Remote AE Title' with the value 'Remote AE Title', 'Hostname or IP' with the value 'Hostname', and 'Port' with the value 'Port'. Below these fields is the 'Export Format' set to 'Uncompressed (RGB 24 bpp)'. At the bottom left, there is a checked checkbox labeled 'Make default'. At the bottom right, there is a blue 'Save' button.

Figure 2. *Export Destination dialog: configuring Caption AI*

4. For **Tablet AE** (Application Entity) **Title**, enter the name that Caption AI will identify itself by to the DICOM server. This field can include only alphanumeric characters and accepts a maximum of 32 characters.
5. For **Remote AE Title**, enter the name of the DICOM PACS server where the Caption AI studies will be stored.

6. For **Hostname or IP**, enter the DICOM PACS server hostname or IP address. This will be the name that appears in the **Export Destination** list.
7. For **Port**, enter the port number of the DICOM PACS server.
8. **Export Format** is the image format for the clips sent to the DICOM PACS. The format that Caption AI uses is Uncompressed (RGB 24 bpp).
9. To make the current settings the default whenever a new study is started on the device, check **Make Default**. If only one export destination is defined, it is selected by default.
10. Tap **Save** to save the export destination. It will now be available to operators when they choose to send a study.

Or, if you don't want to save the settings, tap **X** to clear out the fields and return to the **Settings** screen.

11. On the **Settings** screen, to the right of the export destination you just added, tap **Check connection**. The button will change to reflect the DICOM connection's current status. **Connection Error** or **No Response** indicate that the connection details are not accurate or further troubleshooting is required. Contact your IT department for assistance if the error persists.
12. Tap the back arrow to return to the **Main** menu.

► To edit or delete an export destination

1. On the **Export** tab of the **Settings** screen, tap the name of the export destination you want to edit or delete.
2. In the **Export Destination** dialog, make your changes:
 - Change the name, hostname/IP, or port information.
 - Tap **Make Default** to define whether it will be the export destination used to send reports.
 - Delete the export destination.
3. Save or confirm (when deleting) your edits. You will return to the **Settings** screen.

Annotating exported clips with the view name

Exported clips can be annotated with the clip's view name by enabling **Annotations** in the software settings. The setting applies to all export destinations saved in Caption AI.




Figure 3. Settings screen: Annotations

Each exported clip will be labeled with its view name in the scanning protocol used for the study. For more information on scanning protocols, see "Working with scanning protocols" on page 39.




Figure 4. Exported clip with annotated view name

Configuring a Caption AI modality worklist

To configure Caption AI to query a modality worklist for patient information, set up a worklist connection in the settings. You can use the worklist to fill in patient information on the start screen before beginning an exam.


CAUTION

To ensure proper operation of Caption AI for patient scanning, system setup and configuration should be performed by your ultrasound service representative or a qualified biomedical engineer.

CAUTION

Only an authorized IT administrator who logs into the system using their security credentials should configure Caption AI to communicate with the PACS server.

► To add a modality worklist connection

1. On the Caption AI **start** screen, tap the **Menu** icon () , then **Settings**. The **Settings** screen appears.
2. Select the **Worklist** tab.




Figure 5. Settings screen: Worklist

3. Next to **Modality Worklist Connection**, tap **Add**. The **Modality Worklist Connection** dialog appears. You must fill in all the fields.

Figure 6. *Modality Worklist Connection dialog*

4. For **Tablet AE** (Application Entity) **Title**, enter the name that Caption AI will identify itself by to the DICOM server.
5. For **MWL AE Title**, enter the name of the DICOM server where the orders will be queried. This will be the name that appears in the **Modality Worklist Connection** list.
6. For **Hostname or IP**, enter the DICOM server hostname or IP address.
7. For **Port**, enter the port number of the DICOM server.
8. **Date Range** sets the date range for worklist queries. Orders will only be returned if they fall within the specified date range. You can select one the following options:
 - Past 24 hours
 - Past 3 days
 - Past 7 days
 - All orders
9. **Filter to only ultrasound orders** configures Caption AI to search for orders with a "US" modality. This option is deselected by default.
10. To make the connection the default for worklist searches, check **Make Default**. If only one connection is defined, it is selected by default.
11. Tap **Save** to save the worklist connection. It will now be available to operators before starting an ultrasound exam.

Or, if you don't want to save the settings, tap **X** to clear out the fields and return to the **Settings** screen.

12. On the **Settings** screen, tap the back arrow to return to the **Main** menu.

▶ **To edit or delete a modality worklist connection**

1. On the **Worklist** tab of the **Settings** screen, tap the name of the worklist connection destination you want to edit or delete.
2. In the **Modality Worklist Connection** dialog, make your changes:
 - Change the name, hostname/IP, or port information.
 - Change the search criteria.
 - Tap **Make Default** to define whether it will be the connection used to query the worklist.
 - Delete the export destination.
3. Save or confirm (when deleting) your edits. You will return to the **Settings** screen.
4. To the right of the modality worklist you just added, tap **Check Connection**. The button will change to reflect the DICOM connection's current status. **Connection Error** or **No Response** indicate that the connection details are not accurate or further troubleshooting is required. Contact your IT department for assistance if the error persists.

CHAPTER 5: SETTING UP FOR AN EXAM

Working with scanning protocols

When you start an exam, you follow a protocol that specifies the views you'll acquire for the study. Caption AI offers two types of protocols:


- **Workflow-based protocols** are based on the particular views you want to acquire. The default workflow protocol provided with Caption AI is All Views.
- **Indication-based protocols** are based on the patient's symptoms, such as shortness of breath.

When you start a new exam, you first select the scanning protocol you want to use for the exam. You can use a protocol provided with Caption AI or you can create custom protocols specific to your practice.

Creating custom scanning protocols

You can create custom protocol presets on the **Scanning Protocols** screen.

▶ To create a custom workflow protocol

1. On the Caption AI **start** screen, tap the **Menu** icon () , then **Scanning Protocols**. The **Scanning Protocols** screen appears.




Figure 7. Scanning Protocols screen

- In the **Protocols** list on the **Workflow** tab, select **Add Protocol**. An empty field appears.
- Enter the name for your new workflow protocol.





Figure 8. Scanning Protocols screen: Add custom workflow protocol

- In the list to the right, tap a view or echo/lung window to deselect or select it. By default, all views and windows are selected.
- If you want to reorder the views, use the move icon () to drag a view to another location. You can reorder echo views within an echo window, but lung views cannot be reordered within the lung window. You can reorder echo and lung windows between other windows.

- To add a free scan to the custom protocol, tap **Create Free Scan** at the bottom of the views list.






Figure 9. Create Free Scan dialog

- Give the new free-scan view a name, select **Color Doppler** if the view uses color doppler, and choose the scanning preset.



During scanning, you will automatically be put into the settings you indicated for the free-scan view.

Note: Free-scan views and lung views do not have AI Guidance and are not recommended for novice users. The Quality Meter is also not available for these views.

- Tap **Create** to add the free scan and return to the views list.
- To reorder a free-scan view, use the move icon () to drag the view to another location (you cannot move a free scan view into a lung scanning window). You can add up to 50 free scan views to a custom protocol.
- Tap **Save**. The views list closes and the new protocol is selected in the **Workflow Protocols** list.

► **To create a custom Indications protocol**

1. On the Caption AI **start** screen, tap the **Menu** icon (☰), then **Scanning Protocols**.
2. Select the **Indications** tab.




Figure 10. Scanning Protocols screen: Indications tab

3. In the **Protocols** list, tap **Add Protocol** and enter the name for your new indication protocol.




Figure 11. Scanning Protocols screen: Add custom indications protocol

4. In the views list to the right, tap a view to deselect or select it. By default, all views are selected.

5. Tap **Save**. The views list closes and the new protocol is selected in the **Indications Protocols** list.

► To edit a custom protocol

Note: You cannot edit the protocols that come with the Caption AI product.

1. In the **Scanning Protocols** screen, select the type of protocol you want to edit: **Workflow** for a workflow preset or **Indications** for a symptom-based preset.
2. Select the protocol you want to change and tap the edit icon (✎).





Figure 12. Scanning Protocols screen: Edit protocol

3. Make your changes:
 - Rename the protocol
 - Add or remove views (note that you cannot add or remove views within the lung window)
 - Rearrange views within an echo window (note that you cannot rearrange views within the lung window or move a view into a different acoustic-view window)
4. Save your edits.

► To delete a custom protocol

Note: You cannot delete the protocols that come with the Caption AI product.


1. In the **Scanning Protocols** screen, select the type of protocol you want to edit: **Workflow** or **Indications**.
2. Select the protocol you want to delete and tap the delete icon ().
3. Confirm that you want to delete the protocol.

Turning off AutoCapture

Caption AI comes with the AutoCapture feature turned on by default. If, however, you want to acquire clips manually, you can turn off AutoCapture in the settings. For more information on the AutoCapture feature, see “AutoCapture” on page 62.

Note: Manually recording clips should be done only by users who can determine without the assistance of Caption AI that a clip is of sufficient diagnostic quality.

► To turn off AutoCapture

1. On the Caption AI **start** screen, tap the **Menu** icon (), then **Settings**. The **Settings** screen appears.
2. Select the **Acquisition** tab if it is not active.




Figure 13. Acquisition screen: AutoCapture

3. Under **AutoCapture**, tap **Off**. AutoCapture will be turned off for all future studies, unless it is manually turned on again here in the settings.

Selecting the scanning convention

Caption AI provides two scanning conventions: Cardiology (the default) and General.

► To select the scanning convention

1. On the Caption AI **start** screen, tap the **Menu** icon (☰), then **Settings**. The **Settings** screen appears.
2. Select the **Acquisition** tab if it is not active.
3. Under **Scanning Convention**, select the convention you want to use for cardiac scanning.




Figure 14. Acquisition screen: Scanning convention

Caption AI will use this convention for all future studies, unless it is manually changed here in the settings.

This page intentionally left blank.

CHAPTER 6: SCANNING WITH CAPTION AI

CAUTION

If you are not familiar with performing an ultrasound exam using Caption AI™, make sure that you receive the appropriate training before using the system, provided either by Caption Health or by a trained clinician using official Caption AI training materials.


Before you perform a Caption AI exam, it is vital that you are familiar with the Terason uSmart 3200T Plus. See the *Terason uSmart3200T Ultrasound System User Guide, Volumes 1 and 2*.

CAUTION

For important safety information regarding the proper use and application of Caption AI, see *Chapter 1: Introduction* and *Chapter 3: Operator-Safety Considerations* in this manual.

Starting a Caption AI study

► To start Caption AI

1. Double-tap the **Caption AI** icon () on the desktop.
The Caption AI **start** screen appears.




Figure 15. Caption AI start screen: Workflow-based scan




Figure 16. Caption AI start screen: Indication-based scan

2. Pick the exam protocol you want to use:
 - **Workflow-based scan** for a study based on the views to acquire.
 - **Indication-based scan** for a study based on the patient's symptoms. You can select more than one indication; Caption AI will combine all the views for each indication into a single protocol.

3. In the center pane, select the **Workflow** or **Indications** protocol you want to use for the exam. If you do not see an applicable protocol, you may need to create a custom one. See "Creating custom scanning protocols" on page 39.
4. (Optional) Use a modality worklist to fill in the exam information. See "Using modality worklists to populate exam information" on page 50. If you use a worklist to fill in the exam information, you can skip to step 6.
5. (Optional) Enter patient information, provider ID, accession number, and visit ID.

Notes: To enter text, you can use the on-screen keyboard or the portable keyboard that comes with the device. If the on-screen keyboard does not appear, the function may be disabled in the device-operation settings. Refer to the Windows 10 documentation that came with the device or ask your IT support personnel for assistance.

If you do not enter patient information, the system will automatically assign a six-digit alphanumeric ID number for the patient.

6. Tap **Start scan** to go to the first view in the study.

Note: If you do not have a compatible probe attached to the device, the Connect Probe message appears.




Figure 17. Connect Probe message

Connect a compatible probe to the device to begin scanning. For a list of compatible probes, see "System components" on page 31.

If you change your mind about doing the exam, tap **End Exam**. The **start** screen will reappear with the fields cleared out.

Using modality worklists to populate exam information

You can use modality worklists to populate patient and exam information directly from ultrasound orders at your institution. You must select the order from the worklist before starting an exam.

► To choose an order from a worklist


1. On the start screen, tap **Worklist**.

Note: The Worklist button will only be available if a worklist connection has been configured. See “Configuring a Caption AI modality worklist” on page 36.

The screenshot shows a 'Modality Worklist' dialog box. At the top, there is a title 'Modality Worklist' and a close button (X). Below the title, there is a 'SELECT WORKLIST' section with a dropdown menu showing 'Worklist A' and 'DEFAULT', and a 'SHOW ALL' button. Below this, there are six input fields for search criteria: 'PATIENT LAST NAME' (placeholder: Enter last name), 'PATIENT FIRST NAME' (placeholder: Enter first name), 'ACCESSION NUMBER' (placeholder: 0000000000), 'PATIENT ID' (placeholder: Enter ID), 'DATE OF BIRTH' (placeholder: MM/DD/YYYY with a calendar icon), and 'VISIT ID' (placeholder: 0000000000). At the bottom, there are 'CLEAR' and 'SEARCH' buttons.

Figure 18. Modality Worklist screen

2. Select the worklist you want to search. The default connection is selected by default.
3. Tap **Show All** to view all orders in the selected worklist or, to narrow the results or find a specific order, enter the search criteria in the fields. You can enter any combination of criteria.
4. Tap **Search**. When the search has finished, the worklist results appear.



ACCESSION #	PATIENT NAME	PATIENT ID	DOB	RP ID	START DATE
0000000000	Edwards James	0000000000	00/00/0000	00000000	00/00/0000
0000000000	Edwards James	0000000000	00/00/0000	00000000	00/00/0000
0000000000	Edwards James	0000000000	00/00/0000	00000000	00/00/0000
0000000000	Edwards James	0000000000	00/00/0000	00000000	00/00/0000
0000000000	Edwards James	0000000000	00/00/0000	00000000	00/00/0000
0000000000	Edwards James	0000000000	00/00/0000	00000000	00/00/0000
0000000000	Edwards James	0000000000	00/00/0000	00000000	00/00/0000
0000000000	Edwards James	0000000000	00/00/0000	00000000	00/00/0000
0000000000	Edwards James	0000000000	00/00/0000	00000000	00/00/0000
0000000000	Edwards James	0000000000	00/00/0000	00000000	00/00/0000

Figure 19. Modality Worklist screen: search results

You can sort the results using the arrows at the top of each column. Tap **Edit Search** to do another search.

5. Find the relevant order and tap **Select**. You return to the start screen.

You are now ready to start scanning. First, however, let us review the Caption AI ultrasound display and controls. If you are already familiar with the Caption AI real-time display, you can skip to “Starting and optimizing your scan” on page 56.

The Caption AI ultrasound display and controls

The first view you will be prompted to scan is based on the protocol preset you selected. This first view of each echo window is known as the “home” view and provides the foundation for the other clips in that acoustic-window series.





Figure 20. Basic ultrasound display and controls

This screen has the following basic features:

- **Details panel** displays the current view in the sequence you are about to record. The first view in an echo window is called the “home” view. This is the foundational view for the particular acoustic-window series. For example, PLAX is the home view for the Parasternal series.



Current view and location in protocol sequence

Figure 21. Details panel

Beneath the current view is text that tells you which view you are on out of the total number of views in the protocol and what the next view is.

- **Go to All Views** takes you to the **All Views** list, which shows all the views in the protocol, in the order that Caption AI will go through them, unless you choose to go in a different order.




Figure 22. All Views list

Each view shows a reference image, the view name, and, if applicable, the Caption Interpretation AutoEF estimate for that view. Views that Caption Interpretation uses for individual-view AutoEF estimates are PLAX, AP4, and AP2. For more information on Caption Interpretation AutoEF, *Appendix A: Caption Interpretation AutoEF*.

The current view is outlined. Views that have been acquired are marked as follows:

- – Acquired using AutoCapture
- – Recorded using Save Best Clip
- – Manually recorded
- **Free Scan** allows you to go into free-scan mode.
- **CD Color** allows you to go into Color Doppler mode.
- **Patient ID** shows the current patient's identification that you entered on the start screen or a system-generated ID if you started scanning without entering the patient's ID first.

- **MI and TI indices** provide real-time display of the thermal and mechanical acoustic output indices. It is important to observe the MI and TI values during scanning. When conducting ultrasound studies, follow the ALARA principle: expose the patient to ultrasound energy at a level that is “As Low As Reasonably Achievable.” Please read the *Terason uSmart3200T Ultrasound System User Guide Volume 2* for important information about acoustic output and ALARA.
- **Depth** changes the field of view of the image. Changes are in one-centimeter increments. Each view has a default depth setting that provides a good starting point for typical cases. Depth markers are displayed along the right side of the imaging sector (the closer together the dots, the shallower the depth).
- **Gain** controls the overall amplification of the image signal. Setting gain properly, not so high as to introduce extraneous noise or so low as to make targets difficult to distinguish, may optimize the operation of the Quality Meter. Each view has a default gain setting that provides a good starting point for typical cases.
- **AutoEF** displays the overall Caption Interpretation Automated Ejection Fraction estimate.
- **End Exam** closes the current study and saves all the clips that have been recorded for this exam. If you have not recorded all the views for the exam, you will be asked to confirm that you want to stop the exam. You have the option of reviewing the report of the recorded clips or going back to the **start** screen (**Exit to home**).




Figure 23. End exam dialog

Note: When you end an exam, you cannot reopen it later to add new clips or redo existing ones. You can only update patient information, delete the study, or send it to the DICOM PACS. See *Chapter 7: Working with Saved Studies*.

Setting up your scan

On the scanning screen, the following features help you set up your scan.




Figure 24. Ultrasound display: features for setting up to scan

- The **reference ultrasound image** is an example of an optimized, normal image for the view you are acquiring. For each view in the protocol, the system displays the reference image.
- **Patient-positioning guidance** appears for each echo window home view—Parasternal Long-Axis (PLAX), Apical 4-Chamber (AP4), and Subcostal 4-Chamber (SubC4)—and certain lung views. This guidance tells you how to position the patient for those views.
- **Probe-positioning guidance** appears on the reference image for each echo window home view—Parasternal Long-Axis (PLAX), Apical 4-Chamber (AP4), and Subcostal 4-Chamber (SubC4)—and all lung views. This guidance tells you how to position the probe before you start scanning.
- The **probe-orientation diagram** shows where on the patient's body to place the probe and the direction in which to point the indicator for this view. The probe indicator should point in the direction shown by the blue arrow. For echo views, the "clock" numbers indicate the region in which you may need to explore, because the optimal direction will vary based on the patient.

Starting and optimizing your scan

To optimize your scan, Caption AI provides real-time guidance, described next.

Note: Real-time guidance is currently available for echo views only.




Figure 25. Ultrasound display: real-time guidance

The **Quality Meter** provides real-time feedback as you scan. The bars indicate how close you are to the optimal probe position and orientation.

As you get closer, the meter fills up toward the triangular notches, which indicate the diagnostic-quality threshold. When the meter reaches or exceeds this threshold, Caption AI automatically records the clip (provided AutoCapture is turned on; it is on by default).




Figure 26. Quality Meter: basic features

When scanning for subsequent views within an acoustic-window series, two Quality Meters appear. This dual-meter shows the home view to the left and the current view to the right.




Figure 27. Dual Quality Meter

Caption AI displays the dual Quality Meter to help re-orient the operator, where going back to the home view may be helpful in getting to the current view. For example, in getting an Apical 2-Chamber view of the heart, it may be helpful to go back to the Apical 4-Chamber view, establish a good image for that view, and then rotate from there to obtain an Apical 2-Chamber view.

When the dual Quality Meter appears, you may do the following:

- Move the probe to align with the probe position for the home view. The home-view Quality Meter displays the quality bars for the home-view clip.
- Move the probe in very small increments, following any Prescriptive Guidance, to achieve the optimal probe position and orientation for the current view. The current-view Quality Meter will indicate how close you are compared to the reference ideal image.
- **Prescriptive Guidance** may appear while you are scanning. Caption AI presents guidance on how to adjust the probe to optimize the image. This guidance appears as text and, in many cases, as an icon on the active-image area of the screen.

When Prescriptive Guidance appears, follow the suggestions until they disappear. If you make the adjustments correctly, Prescriptive Guidance disappears and the Quality Meter rises as the image quality gets closer to diagnostic quality. Yellow bars in the Quality Meter mean that you should make smaller, fine-tuned movements. Green bars mean the image is of diagnostic quality.

The types of Prescriptive Guidance the system presents are shown in Table 5, below. The icons represent the face of the probe; the white line is the probe indicator.




















Figure 28. Visualization for Prescriptive Guidance probe-movement icons









Note: All probe movements are in relation to the patient. For example, when the guidance says “Tail slowly more lateral,” it means to move the tail away from the patient’s midline. Tailing is always on the longer, non-indicator axis of the probe face, whereas rocking is always on the shorter, indicator axis.

Table 5. Prescriptive Guidance Examples

On-screen	Icon	Movement/Action
Rotate clockwise		

On-screen	Icon	Movement/Action
Rotate counter clockwise		
Tail up		
Tail down		
Tail medial		

On-screen	Icon	Movement/Action
Tail lateral		
Slide medial		
Slide lateral		
Slide up		

On-screen	Icon	Movement/Action
Slide down		
Rock towards indicator		
Rock away from indicator		
Increase depth		Tap highlighted + button
Decrease depth		Tap highlighted - button

If Caption AI does not detect any visual structures, it prompts you with a different message, depending on the view. For the PLAX and AP4 views, it says, "Slow circular sweeps until moving anatomy appears," and for the SubC4 view it says, "Slowly increase pressure until moving anatomy appears." See the troubleshooting tips in your training materials. Some techniques to try include

rolling the patient, using respiratory prompts, or resetting the probe to the correct echo region.

Note: In acquisition mode, if you haven't touched the screen for 30 minutes for a single view, Caption AI will freeze the probe and display a message telling you that the session has been paused. If you exit the dialog or tap **Resume**, you'll go back to the acquisition screen you were on when the probe was frozen. If you tap **Exit to home**, you'll go to the home page.




Figure 29. Resume scanning dialog

Recording a clip

Caption AI provides three ways to capture clips: AutoCapture, Save Best Clip, and manual recording. AutoCapture and Save Best Clip are currently available for guided echo views only. The default method is AutoCapture. Manual recording should be used only by those users who are able to determine, without the assistance of Caption AI, that a clip is of sufficient diagnostic quality.

AutoCapture

The AutoCapture function automatically records a clip if the image quality exceeds the diagnostic-quality Quality Meter threshold and stays above it for at least 2 seconds. If you maintain the image above the diagnostic-quality Quality Meter threshold for more than 2 seconds, Caption AI will capture a longer clip, up to 4 seconds. Recordings will capture 62–124 frames, depending on how long you keep the image quality above the threshold. By capturing 2–4 seconds of frames, Caption AI will capture at least one complete heart cycle, provided the heart rate is not lower than 30 beats per minute.

You can turn AutoCapture on and off in the settings. See *Chapter 5: Setting up for an exam*.

After Caption AI has recorded the clip, it progresses you to the next view to acquire in the protocol.

► To record a clip using AutoCapture

1. Position the patient according to the positioning guidance.
2. Position the probe on the patient according to the probe-position and acoustic-window guidance, with the probe indicator pointing in the direction of the arrow and clock ranges.
3. As you scan, watch the screen and follow any Prescriptive Guidance suggestions that appear, such as rotating the probe or changing the depth (see Table 5 for an explanation of Caption AI's Prescriptive Guidance).

The number of bars in the Quality Meter increases and decreases based on the optimal position and orientation of the probe for the view being scanned. Your goal is to reach the diagnostic-quality threshold, indicated by the two small triangular notches.

Yellow bars mean that you are getting close to the optimal position and should make smaller, more finely tuned movements.

When one or more green bars appear in the Quality Meter, that means the ultrasound image has exceeded the diagnostic-quality threshold and Caption AI immediately starts recording the clip. "Recording..." appears at the top of the screen; the **Record** button turns red, indicating that recording is in progress; and the ring around the button turns blue as the recording progresses.

When the button changes to a circle with two dots, it means that the image quality has been above the threshold long enough for the clip to be saved (at least 2 seconds), however you are encouraged to hold the image above the threshold for an additional 2 seconds, if possible.

CAUTION

To ensure a diagnostic-quality clip, maintain probe contact and position during recording for at least 2 seconds. If you move the probe or otherwise interrupt the recording, the clip will not be recorded.

When the clip has been recorded, you will hear a "beep" and the **Record** button displays a checkmark while Caption AI processes the clip.

When done processing, Caption AI proceeds to the next view in the workflow protocol. "Saving the clip..." appears at the top of the screen

to indicate that Caption AI is processing the saved frames for storage. At this time, you may move the probe.




Figure 30. Quality Meter: AutoCapture states (home view)

Using Caption AI's Auto Ejection Fraction (AutoEF) feature

While scanning, you can use Caption AI's AutoEF feature to get an estimate of clip quality. Caption Health recommends that you acquire at least two relevant views for best results and for this reason, the AutoEF % does not automatically display until you have scanned two relevant views.





Figure 31. AutoEF estimate: single view

After you acquire a single relevant view for EF, you can tap **Calculate** to see the EF estimate. If you choose to calculate a single-view EF, then acquire a second relevant view, the EF estimate will be recalculated automatically, and if you acquire a third relevant view, the EF estimate will automatically be recalculated again.

Save Best Clip

Another way to record a clip is with the **Save Best Clip** feature, described below.




Tap to have Caption AI save the 62 best continuous frames

Figure 32. Save Best Clip feature

The **Save Best Clip** button appears after you have been scanning for some time without getting what Caption AI determines is a diagnostic-quality image. The button displays a quality “score,” which indicates how close your best clip is to the diagnostic-quality threshold. The score is updated at regular intervals as you continue to scan and helps you decide when to save the best clip that the system has seen and move on. If the quality score and image are satisfactory, tap **Save Best Clip**. Otherwise, continue scanning, following the Prescriptive Guidance and general suggestions, to try to get a better quality image.

Manual recording

If you are an experienced user and are able to determine whether a clip is of sufficient diagnostic quality without the assistance of Caption AI, such as a trained sonographer, you might choose to manually record clips. Manual recording is the only way to capture clips for lung views and free scan views.



Tap to manually record a clip

Figure 33. Record button: manual recording

To manually record a clip, tap the record button. You must then hold the probe steady for 2 seconds (for echo views and free scan views) or 3 seconds (for lung views) so Caption AI can capture the requisite number of frames (62 continuous frames for echo views and free scan views, 96 continuous frames for lung views). After the recording is complete, the record button displays a checkmark, and "Saving the clip..." appears at the top of the screen. During this time Caption AI is processing the saved frames for storage. At this time you may move the probe.

After Caption AI has stored the clip, it progresses you to the next view to acquire in the workflow. For echo views, Caption AI progresses you to the next view in the protocol. For lung views, Caption AI allows you to capture multiple clips within a zone: after capturing clip (R) Zone 1 Clip A, you will be able to capture clip (R) Zone 1 Clip B, followed by (R) Zone 1 Clip C, etc. Once you have captured the desired number of clips within a zone, click on the Next arrow below the view name to progress to the next lung zone.

You can manually record a clip at any time while scanning, or you can turn off AutoCapture in the **Settings** screen. A dotted-circle around the record button indicates that AutoCapture is turned off.



Figure 34. Record button: AutoCapture turned off

CAUTION

Caption AI always presents the Quality Meter and Prescriptive Guidance during scanning of guided echo views, but because manual recording is always available during scanning, and manual recording is the only option for lung views and free scan views, Caption AI cannot guarantee the quality of manually captured clips.

CAUTION


When you manually record a clip, Caption AI assigns the active view label to the saved clip. To avoid a mismatch between the label and the actual view captured, confirm that the view you are saving manually is listed at the top of the screen.

Non-sequential scanning

By default, when you record a study, you step through the views according to the protocol you are using for the exam. However, you may decide that you do not want to record a particular view or that you want to record a clip again before moving on to another view. Caption AI allows you to diverge from the protocol to meet the needs of the specific exam and patient. You can move between views, whether or not you have recorded clips for those views, and you can skip views you do not want to scan.

To move between views

- From the **Details panel**, tap the next or previous arrows below the view name.



To go to previous view and next view

- In the **All Views** list (from the **Details panel**, tap **Go to All Views**), tap the view you want to display.




Figure 35. All Views list

If you select a view you have already acquired a clip for, you will have the option to record that clip again. Tap **Re-Record** and acquire a new clip for that view (see “Recording a clip” on page 62).


Note: If there are previous views in the protocol that you skipped over during scanning, Caption AI will not automatically go back to those views, but you can select any view for which you want to record a clip. Before you end the exam, make sure you have acquired all the relevant views for the study.

Recording a lung scan

If you want to record a lung scan, you can use the 8-zone protocol and record clips manually.

► To record a lung scan

1. In protocol setup, select or create a protocol that contains lung views.



2. On the scanning screen, you may use the reference ultrasound image, patient-positioning guidance, probe-positioning guidance, and zone diagram (with probe orientation) to get started. You may also change the scanning preset (see “Changing the scanning preset” on page 73).




Figure 36. Scanning screen: Lung protocol

- When you are satisfied with the view in the image sector, tap Record. The clip is recorded and saved in the format of “(R) Zone 1 Clip A,” where (R) denotes the patient’s side, 1 denotes the zone number, and A denotes the clip within each zone.

Note: Caption AI currently supports an 8-zone protocol. Zones 1-4 cover the patient’s right lung. Zones 5-8 cover the patient’s left lung. Up to 8 clips (A-H) may be captured per zone.

- After a lung clip is captured, you will automatically be taken to the next clip for that zone (up to a maximum of 8 clips per zone). To advance to the next zone, click on the Next arrow below the view name.

Recording a clip based on Color Doppler

If you want to record a clip using Color Doppler, you can turn on **Color** and record the clip manually.

► To record a clip using Color Doppler

- On the scanning screen, tap **CD Color**. The Color Doppler box and options appear.

Note: The Color Doppler clip will be added above your current position in the protocol list.




Figure 37. Scanning screen: Color Doppler mode

2. Select the color velocity scale. The default is Fast Flow for most cardiac imaging. Select Slow Flow for low velocities.
3. To reposition the Doppler box, tap and drag it. While you are repositioning, the box turns into a circle. When you release your finger from the screen, the Doppler box reappears in the new location.




Figure 38. Moving Color Doppler area of focus

4. When you are satisfied with the view in the image sector, tap the **Record** button. The clip is recorded and saved as "Color *n*," where *n* is the number of the Color Doppler clip. For example, if you record three Color Doppler clips, they will be named Color 1, Color 2, and Color 3.

Note: After a Color Doppler clip is recorded, a new Color Doppler clip will be created below it automatically. To exit Color Doppler mode without recording a clip, tap the **X** in the top right corner of the scanning screen.

Recording a free-scan clip

If you want to manually record a clip that is not associated with one of the protocol's views, you can use the Free Scan feature.

► To record a clip using free scan

1. On the scanning screen, tap **Free Scan**. A free scan view will be added to your current position in the protocol.




Figure 39. Free Scan mode

2. When you are satisfied with the view in the image sector, tap **Record**. The clip is recorded and saved as “Free Scan n ,” where n is the number of the free-scan clip.
3. When you are satisfied with the view in the image sector, tap the **Record** button. The clip is recorded and a new free-scan clip will be created below it automatically. To exit Free Scan mode without recording a clip, tap the **X** in the top right corner of the scanning screen.




Figure 40. Free-Scan clip acquired

Changing the scanning preset

Scanning presets optimize imaging parameters for different types of ultrasound scans. You can change the scanning preset at any time while in Lung, Free Scan, or Color Doppler mode using the **Preset** dropdown menu in the top left of the screen.




Figure 41. Scanning presets menu

The list of available scanning presets varies based on scanning mode and ultrasound probe, as shown in the table below.

Table 6. Available Caption AI presets

Lung scans (following Caption AI's 8-zone protocol)	Free Scan or Color Doppler mode
<p>4V2A phased array probe</p> <ul style="list-style-type: none"> • Lung (13cm, optimized for anterior and upper lateral lung zones) • Lung: Deep (16cm, optimized for lower lateral lung zones) • Cardiac <p>15L4A linear array probe</p> <ul style="list-style-type: none"> • Lung 	<p>4V2A phased array probe</p> <ul style="list-style-type: none"> • Cardiac • Lung • Lung: Deep • FAST • Abdominal <p>15L4A linear array probe</p> <ul style="list-style-type: none"> • Vascular: Carotid • Vascular: DVT • Vascular: Access • Lung • Musculoskeletal • Nerve Block

Notes: Caption AI's lung protocol presets are set to Lung for all zones and probes by default, with the exception of the lower lateral views (Zones 4 and 8) for the phased array probe, which are set to Lung: Deep. You can change to another preset using the dropdown menu.

AI guidance is available only in the Cardiac preset. The Preset menu is not available while scanning a view with AI assistance.

Scanning with the linear probe

The linear probe can be used any time while in Lung, Free Scan, or Color Doppler mode. To switch to the linear probe, first unplug the current probe. The Connect Probe message appears.




Figure 42. Connect Probe message

After you plug in the linear probe, the message will disappear automatically within a few seconds and you can begin scanning.

Switching the ultrasound probe

If you come to a free scan or Color Doppler view in the scanning protocol and do not have the appropriate probe plugged in for that view's scanning preset, you will see the Switch Probe message.




Figure 43. Switch Probe message: Linear probe

After you plug in the appropriate probe, the message will disappear automatically within a few seconds. Or, if you want to continue scanning with the current probe, you can dismiss the message by tapping **Continue**. Selecting **Do not show again** will prevent similar messages for the rest of the ultrasound exam.

AI-assisted cardiac scanning is available only for the phased array probe. If you come to an AI-guided view without the phased array probe plugged in, you will be prompted to switch the ultrasound probe.




Figure 44. *Switch Probe message: AI-guided cardiac scanning*

You scan the view once you plug in the phased array probe. Or, if you want to continue scanning with the current probe, tap **Enter Free Scan** to enter Free Scan mode and capture the AI-guided view later.

Reviewing acquired clips during an exam

After you finish recording all clips in a protocol, the **All Views** list appears. On this screen you can review all the clips before sending the exam report. This is where you can choose to record clips for any views that you did not already acquire a clip for or that you want to scan again.




Figure 45. All Views list for active study

- **All Views** shows you all the views in the protocol, including views that you did not acquire clips for. In this list, you can select a view to review or record a clip.
- **Re-Record** puts you in live-scanning mode for the selected view.
- **AutoEF** shows the Caption AI Interpretation AutoEF information.
- **Review Exam** takes you to the Report screen where you can review all the clips in the exam before sending the report to the PACS system.

► To review all the clips in a study

When the last clip in a protocol is saved, the All Views list automatically appears, with the last-recorded clip selected. If you are in Detail view and want to end the exam without acquiring a clip for the last view in the study, tap **Go to All Views** or **End Exam**.

1. In the **All Views** list, select a clip and review it in the scan-image area to determine whether the image is of good-enough quality for the study.
2. If you want to record the clip again, tap **Re-Record**. This takes you back to the live-scanning for that view. (See “Recording a clip” on page 62.)

CAUTION

When you record a clip again, the new clip will replace the previously recorded clip.

3. Repeat steps 1 and 2 until you have reviewed all the clips you plan to send as part of the exam.
4. When you are satisfied with the study, tap **Review Exam**. This takes you to the Report screen. See "Reviewing and sending a study" on page 77.

Reviewing and sending a study

After you are satisfied with the clips you have acquired for a study, you can review and send the study report. You do this on the **Report** screen.






Figure 46. Report screen

- **Re-Record** resumes the study and puts you into live scanning for the current view in the image window. This option is only available immediately after you have captured a study. Once the study is resumed, you can record other clips again by selecting them in the All Views list. To return to the Reports screen after you are satisfied with the images, tap the **Review Exam** button at the bottom of the All views list.

- **Thumbnails** is a visual list of all the clips in the study. A border around a thumbnail indicates the clip that appears in the image window. If the study has more clips than can be displayed at one time in the row of thumbnails, you can swipe left and right to scroll through the set of clips. The thumbnail also indicates how the clip was captured: AutoCapture (), Save Best Clip (the quality score), or manually ().
- The **patient's ID** is shown below the leftmost thumbnail. This is either the ID you entered on the New Patient screen or a system-generated ID that Caption AI created if you did not enter an ID.
- **AutoEF** shows the Caption Interpretation AutoEF results, if the protocol included a PLAX, AP4, or AP2 view. An Ejection Fraction estimate is shown along with a qualitative assessment of global LV function associated with the Ejection Fraction estimate, using a mapping from ASE² and ACEP³ guidelines, shown in the color-coded chart below the overall estimate.

A percentage likelihood is shown, to indicate the probability that the displayed qualitative assessment is in fact the true qualitative assessment of global LV function for the patient. For more information on AutoEF, see *Appendix A: Caption Interpretation AutoEF*.
- **Scanning region menu** allows you to toggle between All views and Lung views only, if you scanned using Caption AI's lung protocol. The Lung report screen is explained in the next section.
- **Close** allows you to close the study. It will be saved to Caption AI, but not sent to the PACS location.
- **Export** transmits the study to the permanent-storage location specified in the Caption AI settings and saves a copy to Caption AI. For more information, see *Chapter 4: Configuring Caption AI*.

2 Lang RM Badano LP Mor-Avi V Afilalo J Armstrong A Ernande L Flachskampf FA Foster E Goldstein SA Kuznetsova T Lancellotti P Muraru D Picard MH Rietzschel ER Rudski L Spencer KT Tsang W Voigt JU. "Recommendations for cardiac chamber quantification by echocardiography in adults: an update from the American society of echocardiography and the European association of cardiovascular imaging." *Eur Heart J Cardiovasc Imaging* 2015;16:233–271.

3 ACEP Emergency Ultrasound Standard Reporting Guidelines. June 2018

⚠ CAUTION

Once you export or close a study, the study is complete. It is “locked” and the clips cannot be modified. Before closing or sending a study, make sure you are satisfied with the quality of the clips.

▶ To view the lung report screen

1. On the **Report** screen, select **Lung** from the **Scanning region** dropdown menu to see the **Lung report** screen.




Figure 47. Lung report screen

2. The **Lung report** screen displays all clips captured using the lung protocol, organized by zone. Each zone contains a thumbnail of the first clip captured in that zone. You can click on the previous and next arrows within the thumbnail to view other clips captured in that zone. Clicking the thumbnail displays the **Lung Zone gallery**.




Figure 48. Lung Zone gallery

The **Lung Zone** gallery displays all clips captured in that zone. From here, you can re-record a clip or go back to the Lung report page to see clips from other zones.

► To export or close out of a study

1. On the **Report** screen, select a clip to view.
2. When you have reviewed the recorded clips and are satisfied with the study, you can do either of the following:
 - Tap **Close**. This closes out of the study but does not send it to the DICOM PACS permanent storage.
 - Tap **Export** to electronically transmit the study to the DICOM PACS storage location.

If an export destination has been selected, the **Export** button changes to reflect that the study is being sent.



If an export destination has not been selected, the **Export Studies** dialog appears, where you can select the PACS system to send the report to.




Figure 49. Export Studies dialog: sending a just-completed study

If no export destination has been set up, you can do so in the **Export Destination** dialog. For more information, see *Chapter 4: Configuring Caption AI* or contact your IT administrator.

If the study is successfully exported, the study is saved and the **Start** screen reappears.

If the study is not sent successfully, you return to the **Report** screen.




Figure 50. Report Not Sent message: sending a just-completed study

Do one of the following:

- Tap **Resend** to have Caption AI make another attempt.
- Tap **Close** to return to the **start** screen. The study is added to the **Saved Studies** list.

Note: If at least one clip has been recorded, the study will be saved.

After a study has been saved to Caption AI, it is available in the **Saved Studies** list. See *Chapter 7: Working with Saved Studies*.

CHAPTER 7: WORKING WITH SAVED STUDIES

The Saved Studies screen is where Caption AI™ lists all the studies currently stored on the device. This list includes all studies in which at least one clip was saved. The studies in this list may or may not have been sent to the DICOM PACS.

From this list you can do the following:

- Edit the patient's name, gender, and date of birth. (The patient ID cannot be changed.)
- Review the clips that were captured for a study.
- Send or resend

Caption AI is meant to be used for *temporary* storage of your studies. The number of studies that can be stored depends on the number and time-length of clips, as well as the storage capacity of the uSmart 3200T Plus.

Displaying the list of saved studies

All studies in which at least one clip has been saved remain in Caption AI until they are deleted.

► To display saved studies

On the **start** screen, tap the **Menu** icon () , then **Saved Studies**. The **Saved Studies** screen appears.





Figure 51. Saved Studies screen

- The **Saved Studies** list displays all the studies currently stored in Caption AI. They are listed by most recent first, based on the date and timestamp of the study. To sort by oldest first, tap on **Created**.
- **Provider or Patient search** allows you to search for a provider by their ID or for a patient by their name or ID.
- **Export** starts the export process.
- **Patient ID** is the patient's ID (entered by the operator or automatically by Caption AI if the operator did not enter a name when starting the study).
- **Patient Name** is the patient's name, if it was entered when setting up the study.
- **Provider ID** is the provider's ID, if it was supplied when setting up the study.
- **Date and Time** is the time the study was saved to the device.
- **Clips** shows how many clips are in the study.
- **Exported** is the status of the study was exported. If the study was successfully exported, this will show the date. If the study was not successfully transmitted, the word "Error" will appear here.

Reviewing and exporting saved studies to the DICOM PACS

You can export one or more saved studies from the **Saved Studies** screen or, if you want to review a study before exporting it, from the **Report** screen.

► To export one or more studies

1. On the **start** screen, tap the **Menu** icon () , then **Saved Studies**.
2. On the **Saved Studies** screen, do one of the following:
 - If you want to export all of today's studies or all the studies that have not yet been transmitted, tap **Export**.
 - If you want to export specific studies, select the studies first, then tap **Export**.




Figure 52. Saved Studies screen: multiple studies selected

3. In the **Export Studies** dialog, confirm the studies to be sent.




Figure 53. *Export Studies dialog: selected studies*

4. If you want to select a different PACS than the one listed in the **Destination** field, tap the down arrow and select the PACS you want to transmit the studies to. (If the PACS system is not listed, see “Configuring Caption AI to send studies to the PACS server” on page 32.)
5. Tap **Export**. The studies are transmitted to the DICOM PACS.




Figure 54. *Exporting Studies screen*

If you do not want to export the studies, tap **Stop Export** and confirm that you want to stop exporting.




Figure 55. Stop Export dialog

If a study was not successfully transmitted, “Error” appears in the **Exported** column and the **Resend** button appears. Tap **Resend** to make another attempt to export all studies that were not sent.

► To review and export a single study

1. On the **start** screen, tap the **Menu** icon (☰), then **Saved Studies**.
2. On the **Saved Studies** screen, tap **Review**. The **Report** screen for that saved study appears (see “Figure 45. Report screen” on page 77).
3. Tap a thumbnail to display the image for that clip. If more clips are available than can be displayed at one time, you can swipe left and right to see other thumbnails.
4. When you are ready to send the study, tap **Send** or **Resend**, whichever appears. Or, if you change your mind about sending the study, tap **Close**.

If the study was not successfully transmitted, “Error” appears in the **Exported** column and the **Resend** button appears (see “Figure 49. Report Not Sent message: sending a just-completed study” on page 81). Tap **Resend** to make another attempt to export the study.


The study remains stored on the device until it is deleted (see “Deleting a study” on page 90).


If a storage location has not been specified, the **Export Destination** dialog appears. For more on what information to enter here, see “Configuring Caption AI to send studies to the PACS server” on page 32.

Updating a patient's information

After a study has been saved, some of the patient's information can be changed. This is especially helpful if the operator who performed the study did not enter this information and only a system-generated ID was assigned to the study.

► To update a patient's information

1. On the **start** screen, tap the **Menu** icon () , then **Saved Studies**.
2. Tap the checkbox next to the patient whose information you want to enter (a checkmark appears to show that it is selected). The **Edit** button appears.



3. Tap **Edit**. The **Patient Details** screen appears.

The screenshot shows a 'Patient Details' form with the following fields and values:

Field	Value
PROVIDER ID	00000000
PATIENT ID	B936F6
FIRST NAME	Kanitha
LAST NAME	Sathien
GENDER	Select
DATE OF BIRTH	July, 21, 1986
ACCESSION NUMBER	00000000
VISIT ID	00000000

Figure 56. Patient Details screen

Note: You cannot change the patient ID once it has been created, therefore, it is grayed out on this screen.

On this screen, you can enter or change the following information:

- **Provider ID:** tap in the field and enter the provider ID.
 - **First and last names:** tap in the field to activate the on-screen keyboard or use the portable keyboard that accompanies the device.
 - **Gender:** tap the down-arrow to display the list, then tap to select the appropriate gender
 - **Date of Birth:** tap the down-arrows to display the year, month, and day, then tap the appropriate value.
 - **Accession number:** tap in the field and enter the accession number.
 - **Visit ID:** tap in the field and enter the visit ID.
4. Tap **Save Changes** to save your entries.

Deleting a study

You can delete studies one at a time or you can delete multiple studies at one time.

CAUTION

Be certain that you want to delete a study because once it is deleted, it cannot be retrieved.

► To delete one or more studies

1. In the **Saved Studies** list, tap the checkbox next to the study or studies you want to delete.

A checkmark appears next to each selected study and the **Delete** button appears.

2. Tap **Delete**. You receive the following confirmation message:




Figure 57. Delete study dialog

3. Tap **Yes, Delete**. The studies are deleted from Caption AI.

APPENDIX A: CAPTION INTERPRETATION AUTOEF

Introduction

This manual is intended to assist you with the safe and effective operation of your Caption Health product. Before attempting to operate the product, read this manual and strictly observe all warnings and cautions. Pay special attention to the information in the section titled “Safety.”

Note: The user information for your Caption Health product may describe the most extensive configuration of the product, with the maximum number of options and accessories. Some functions described may not be included in your product’s configuration.

Symbols and terms

This manual uses various warning, caution, and safety symbols and terms. Before using the product, familiarize yourself with these symbols and carefully read this section to understand their use.



CAUTION

This symbol and term alert you to information regarding patient, operator, or equipment safety. They also indicate information about preventing the loss of patient or product data.

Customer service

Customer service representatives are available to answer questions and provide support for your operation of the product. Please contact your local Caption Health representative for assistance or email support@captionhealth.com.

Safety information

Prescription Device Statement

CAUTION

For product usage in the United States of America, the following labeling statement applies:

Federal law restricts this device to use by or on the order of a physician.

CAUTION

Do not use the product for purposes other than those intended and expressly stated by Caption Health, Inc. Do not misuse the product, and do not use or operate the product in an incorrect manner.

Installation, use, and operation of this product are subject to the law in the jurisdictions in which the product is used. Install, use, and operate Caption Interpretation™ only in a manner that does not conflict with applicable laws or regulations, which have the force of law.

Use of the product for purposes other than those intended and expressly stated by Caption Health, Inc., as well as incorrect use or operation, may relieve Caption Health or its agents from all or some responsibility for resultant noncompliance, damage, or injury.

CAUTION

Product users are responsible for image quality and diagnosis. Inspect the data that is being used for analysis and diagnosis, and ensure that the data is sufficient and appropriate in anatomical correctness and both spatial and temporal resolution for the measurement being employed.

Intended use/Indications for use

This product is intended to be installed, used, and operated only in accordance with the safety procedures and operating instructions given in the product operator information, and only for the purposes for which it was designed. Nothing stated in the operator information reduces the operator's responsibility for informed clinical judgment and best clinical procedure.

The Caption Interpretation™ software is used to process previously acquired transthoracic cardiac ultrasound images, to store images, and to manipulate and

make measurements on images using an ultrasound device, personal computer, or a compatible DICOM-compliant PACS system in order to provide automated estimation of left ventricular ejection fraction. This measurement can be used to assist the clinician in a cardiac evaluation.

The Caption Health AutoEF software application is indicated for use in adult patients.

Intended audience

This manual is intended for healthcare professionals who operate and maintain the Caption Health AutoEF product.

Prior to using this information and the AutoEF product, operators must be familiar with ultrasound techniques. Sonography training and clinical procedures are not described here.

This product must be operated by or under the direction and supervision of a licensed physician.

Measurement accuracy

Measurements can be done on ultrasound images, which may be used with other clinical information to make a diagnosis. The AutoEF product performs measurements on images that have previously been acquired.

Making a diagnosis based solely on ejection-refraction measurements is not recommended. It is also important to understand limitations and factors that affect the accuracy of ultrasound image measurements. These include speed of sound variations in different tissue types, image quality, patient echogenicity, ultrasound scanner design, and operator skill. These factors must be considered when using any quantified data based upon ultrasound images. Consult your ultrasound system operator's manual for further information.

Echocardiography ejection fraction

Left ventricular ejection fraction (EF) by two-dimensional ("B-Mode") imaging is a measurement of the fraction of blood ejected from the left ventricle during systole. It is conventionally measured by dividing the stroke volume (the amount of blood ejected during systole) by the end-diastolic volume. Ejection fraction is stated as a percentage.

With echocardiography, ejection fraction is typically measured using planimetry to measure end-diastole and end-systole endocardial borders from apical views.

The Simpson's method is commonly used to then calculate volumes and ejection fraction.

The Caption Health AutoEF method uses deep learning algorithms that have been trained on very large datasets of echocardiograms with reported Ejection Fraction percentages. The algorithm operates on Parasternal Long Axis (PLAX), Apical Four Chamber (AP4), and Apical Two Chamber (AP2) images. The algorithm calculates an ejection fraction result but does not calculate chamber volumes. In this regard, it may be conceptually likened to the expert echocardiographer using the "eyeball" method to estimate Ejection Fraction.

Ultrasound systems

The development and testing of Caption Health AutoEF employed a large and diverse number of cardiac ultrasound systems. Provided that the other conditions described in the Operator's Manual are met, it is expected that acceptable performance of Caption Health AutoEF may be expected with images acquired by ultrasound systems that meet the following minimum requirements: multi-channel ultrasound systems operating wide bandwidth linear phased array transducers with DICOM storage and sending compatibility. The systems and transducers used to acquire the clips should be cleared by the FDA for cardiac imaging.

When used in Caption AI, Caption Health AutoEF is compatible with the ultrasound systems that are compatible with Caption AI.

Clip annotation and selection

The AutoEF software includes a function that processes video clips in a study in order to automatically classify clips that are PLAX, AP4, and AP2 views. It also includes a function that selects PLAX, AP4, and AP2 clips that are most appropriate for measurements. It is important to understand this function and to pay attention to the results of the Clip Selector when reviewing studies.

Based on verification and validation testing, with a diverse range of patient types, the Clip Selector may reject studies because it fails to find suitable Parasternal Long Axis, Apical 4 Chamber, and Apical 2 Chamber video clips in the study. Of the 166 studies that the Clip Selector was ran on, 162 (97.6%) had a suitable AP4 clip, 155 (93.4%) had a suitable AP2 clip, and 147 (88.6%) had a suitable PLAX clip. The frequency of rejection may be lower with better images and easier patients, and it may be higher with lower quality images, and difficult patients, such as patients with a high body mass index or technically difficult acoustic windows.

The rate is consistent with echocardiography norms where as many as 20% of adult patients are considered to typically present with sufficient technical difficulty that contrast agent enhancement is considered indicated in order to improve diagnostic quality, particularly of left ventricular function. (Lindner, JR, A Practical Approach to Contrast Echocardiography, American College of Cardiology, July 10, 2017 accessed online at <http://www.acc.org/latest-in-cardiology/articles/2017/07/10/09/17/a-practical-approach-to-contrast-echocardiography>.) It is also consistent with echocardiography left ventricular ejection fraction studies, where it is typical that approximately 10-13% of patients are considered too technically difficult to be included for EF analysis. (Malm S, Frigstad S, Sagberg E, Larsson H, Skjaerpe T. Accurate and reproducible measurement of left ventricular volume and ejection fraction by contrast echocardiography. Journal of the American College of Cardiology. 2004 Sep 1;44(5):1030-5.)

For best results, and to minimize the rejection of clips presented to the AutoEF software, submit images that meet the requirements for successful manually traced Biplane Simpson's method of disks measurements: minimal endocardial dropout, and minimal foreshortening of the apical views. (Lang RM, Badano LP, Mor-Avi V, et al. Recommendations for cardiac chamber quantification by echocardiography in adults: an update from the American Society of Echocardiography and the European Association of Cardiovascular Imaging. J Am Soc Echocardiography 2015;28:1-39. e14.)

The clips selected for use are displayed with an indication that they have been used to generate the AutoEF. When reviewing AutoEF results, confirm that PLAX, AP4, and AP2 clips that are appropriate for measuring ejection fraction are being used in the AutoEF estimation.

To produce an AutoEF estimate, the Clip Annotator and Selector requires a minimum number of 16 frames, in order to process enough of the heart cycle to estimate ejection fraction.

CAUTION

Review the video clips used in generating AutoEF measurements to ensure that the clips are appropriate for measuring ejection fraction.

Transthoracic and transesophageal echocardiography

The Caption Health AutoEF software has been trained on transthoracic echocardiograms. It is not intended for use on transesophageal echo images. Because of the similarity of many transthoracic and transesophageal

echocardiographic views, the Clip Annotator and Selector function may occasionally inadvertently select transesophageal images for processing. Ejection fraction results from these images should not be used for diagnosis.

CAUTION

AutoEF is not intended for transesophageal echocardiograms (TEE). Do not operate AutoEF on transesophageal images.

Patient data archival

The Caption Health AutoEF product is not intended to be a permanent archival or storage location of studies. Patient studies should be archived in your permanent location.

CAUTION

The Caption Health product is not a permanent storage device for patient studies.

PHI protection

Patient information contained in DICOM studies contains Protected Health Information (PHI). The Caption Health AutoEF software operates as a microservice within the Caption AI product and users do not have direct access to any data stored in the Caption Health AutoEF software. Instead, users view the results produced by the Caption Health AutoEF software via the Caption AI software and PACS DICOM Viewer used to read studies. Configuration of Caption Health products and PACS DICOM viewer(s) is performed by authorized IT administrators who log into the system using their security credentials.

The Caption Health AutoEF software is not intended for long-term storage of studies.

When viewing AutoEF results in your DICOM Viewer, make certain to observe your institution's guidelines and practices regarding the protection of PHI.

Cybersecurity precautions and practices

Malware, computer viruses, ransomware, and other cybersecurity threats are an increasing concern in healthcare IT systems. The Caption Health AutoEF software operates as a microservice inside Caption AI and behind your

institution's security firewalls. There is no direct connection between the Caption Health AutoEF software and the outside internet.

For information and guidance on implementing proper cybersecurity in the healthcare IT environment, see "Health IT Privacy and Security Resources for Providers" at <https://www.healthit.gov/topic/privacy-security-and-hipaa/health-it-privacy-and-security-resources-providers>.

Product compatibility

Caption Health AutoEF is compatible with Caption Guidance as part of the Caption AI product. Do not use your system in combination with other products or components, unless the Caption Health AutoEF software expressly recognizes those other products or components as compatible. For information about such products and components, see "Specifications" and contact your Caption Health representative.

Changes and additions to the system should be made only by Caption Health or by third parties expressly authorized by Caption Health to do so. Such changes and additions must comply with all applicable laws and regulations that have the force of law within the jurisdictions concerned, and best engineering practices.

Performance Testing

Validation testing and performance summary of AutoEF

Caption Health AutoEF performance, including the integrated operation of the Clip Selector, was validated on a dataset of 166 patient studies. These studies were assembled to represent a broad cross-section of patient conditions, including a diverse range of ejection fraction values, and a significant portion of high body mass index patients as is typically associated with technically difficult studies. The validation dataset also included a diverse and representative mix of ultrasound manufacturers and product models. These included:

Manufacturer/Model

- Acuson/Sequoia
- Siemens/Acuson SC2000
- Philips Medical Systems/CX50
- GEMS US/Vivid i
- GE Vingmed Ultrasound/Vivid7
- Philips Medical Systems/Epiq 7C
- Philips Medical Systems/iE33

For each of these studies a reference ejection fraction by 2D echo using the biplane Modified Simpson's method of disks was available. These had been created by expert registered sonographers and read and reported by expert cardiologists. The test compared the Caption Health AutoEF to the expert produced and reported biplane Modified Simpson's ejection fraction. The expert produced and reported biplane Modified Simpson's ejection fraction was considered the ground truth. The performance target was that the results of the AutoEF software would be within the limits typically represented by inter-observer variability in biplane Modified Simpson's ejection fraction measurements. Root mean squared deviation (RMSD) was used as the endpoint to evaluate this performance. RMSD can be computed from a Bland-Altman analysis. In this test, RMSD is provided in the unit of ejection fraction percentage. The RMSD target was derived from published literature. The performance target was an RMSD of less than or equal to 9.2%. This target was met.

The testing described in this section was performed with Caption Health AutoEF processing operating on a Zotac EN1070 computer, and communication and display of results was with the Osirix MD 9.0 DICOM PACS Server/Workstation software running on an Apple Macintosh computer.

Clip annotation performance

The Clip Annotator was run on 100 studies. For each study, it identified the most probable clip ("best clip") for each of the following four view/mode categories: PLAX, AP4, AP2, and Color Doppler clip. Three (3) expert registered diagnostic cardiac sonographers (RDCS) with more than five (5) years of professional experience were recruited. For a given clip, each expert was independently presented with the best clip and the Clip Annotator's prediction of the view/mode for the clip. The expert was asked whether the Annotator's prediction was correct. This was repeated for all the best clips found for all recruited experts. A consensus paradigm was used to determine ground truth ("reference standard") for each best clip.

The studies for which the Clip Annotator was not able to find any clips for a certain view/mode was presented to the panel of 3 sonographers. Each member of the panel was independently asked whether the study contained at least one clip of the view/mode without use of contrast agent or TEE. This was repeated for all such studies and for all applicable view/modes. The resulting data, combined with the above PPV assessment dataset, was used to assess sensitivity.

The study-level positive predictive value (PPV) and sensitivity were computed, and each was tested for whether it exceeded the pre-specified performance goals of 80% PPV and 80% sensitivity.

The results are shown in Table 6. All endpoints were met, demonstrating acceptable performance of the Caption Interpretation Clip Annotator feature.

Table 7. *Clip Annotator positive predictive value (PPV) and sensitivity*

	View/ Mode	PPV % [95% CI]	p-value	Sensitivity % [95% CI]	p-value
View	PLAX	100% [100, 100]	p < 0.001	97.8% [94.6, 100]	p < 0.001
	AP4	100% [100, 100]	p < 0.001	91.5% [85., 96.8]	p < 0.001
	AP2	98.4% [96, 100]	p < 0.001	94.2% [88.4, 98.8]	p < 0.001
Mode	No color flow present	100% [100, 100]	p < 0.001	100% [100, 100]	p < 0.001

Clinical validation

The Clinical Performance Validation study was a retrospective, multi-center study. The test design is summarized here, and the results are discussed in the Results section.

This study included a selection of 166 patient studies from the following three sites: the Minneapolis Heart Institute (MHI), Duke University (Duke), and Northwestern Medicine (NM). A range of body-mass index (BMI) was included in the dataset, with 38.0% of the patients overweight ($25 \leq \text{BMI} < 30 \text{ kg/m}^2$) and 27.7% obese ($\text{BMI} > 30 \text{ kg/m}^2$). In addition, 43.9% of patients had reduced EF ($30\% \leq \text{EF} < 53\%$), and 16.2% of the patients had severely reduced EF ($\text{EF} < 30\%$).

All studies were traced by three (3) sonographers using the biplane Simpson's method, and overread by three (3) Level III cardiologists to establish the reference standard.

In addition, a panel of ten (10) physicians was recruited to independently provide visual EF assessments based on a single view. These physicians perform visual EF assessments in their clinical practice; the panel included those specializing in cardiology (n = 3), anesthesiology (n = 2), emergency medicine (n = 4), and

hospital medicine (n = 1). For each PLAX, AP4, and AP2 clip in the 166 studies, each physician was asked to provide an assessment of the qualitative categories of left ventricular function based on guidelines from the American Society of Echocardiography and the American College of Emergency Physicians, shown below:

- Hyperdynamic: EF > 73%
- Normal: 53% ≤ EF ≤ 73%
- Reduced: 30% ≤ EF < 53%
- Severely reduced: EF < 30%

There were two sets of primary hypotheses tested in this study:

- **Objective performance goal (best available view(s)):** Caption Interpretation measurements based on the best available view(s) are superior to a literature derived endpoint of 9.2% root-mean-square deviation (RMSD).
- **Comparison to qualitative visual assessment (single views):** Caption Interpretation assessment of qualitative categories of left ventricular function (hyperdynamic, normal, reduced, and severely reduced) based on single views (PLAX, AP4, AP2) are superior to predictions by physicians trained to interpret echocardiographic images.

Results

Caption Interpretation Automated Ejection software produced an EF assessment from at least one view in all 166 (100%) studies. An EF assessment based on at least two views was produced in 163 (98.2%) of the studies.

Objective performance goal—primary hypothesis testing

Results of the hypothesis testing of Caption interpretation Automated EF software evaluating the performance of the algorithm against its objective performance goal of 9.2% RMSD is shown in Table 2 below. The primary endpoint was met.

Table 8. *Caption Interpretation Automated EF performance compared to the reference standard – primary hypothesis testing*

MAD EF% [95% CI]	RMSD EF% [95% CI]	p-value
6.050 [5.595, 6.506]	7.549 [6.986, 8.072]	< 0.001

Bland-Altman and linear fit plots for Automated EF estimates made from the available views (PLAX, AP4, and AP2) in each patient study are shown below. From the plots, it can be observed that Caption Interpretation estimation of EF is highly linear with the reference EF and performs with comparable accuracy across the range of EFs, including in patients with hyperdynamic and severely reduced LV function.




Figure 58. Bland-Altman (left) and Linear Fit (right) plots for Caption Interpretation Automated EF estimation based on best available view(s)

Objective performance goal—secondary hypothesis testing

Results of the hypothesis testing of Caption Interpretation AutoEF software evaluating the performance of the algorithm against its objective performance goal of 9.2% RMSD for combinations of views are shown in Table 3 below. The secondary endpoints were met.

Table 9. Caption Interpretation Automated EF performance compared to the reference standard – secondary hypothesis testing

#	View Combination	MAD EF% [95% CI]	RMSD EF% [95% CI]	p-value
1	AP4 and AP2	6.002 [5.524, 6.480]	7.479 [6.858, 8.052]	< 0.001
2	AP4 and PLAX	6.171 [5.645, 6.696]	7.791 [7.148, 8.384]	< 0.001
3	AP4, AP2, and PLAX	5.685 [5.232, 6.138]	7.094 [6.534, 7.613]	< 0.001
4	AP2 and PLAX	6.193 [5.625, 6.761]	7.812 [7.084, 8.479]	< 0.001

Objective performance goal—descriptive analysis

Results of the hypothesis testing of Caption Interpretation AutoEF software evaluating the performance of the algorithm against its objective performance goal of 9.2% RMSD for single views are shown in Table 9 below.

Table 10. *Caption Interpretation Automated EF performance compared to the reference standard – descriptive analysis*

#	View(s)	MAD EF% [95% CI]	RMSD EF% [95% CI]	p-value
1	AP4 only	6.618 [5.977, 7.258]	8.594 [7.673, 9.426]	0.0794
2	AP2 only	6.906 [6.260, 7.552]	8.790 [7.912, 9.588]	0.1621
3	PLAX only	7.639 [6.817, 8.461]	9.870 [8.562, 11.024]	0.8506

Comparison to qualitative visual assessment (single views)

A paired-difference analysis was performed to demonstrate superiority of Caption Interpretation in assessing EF qualitatively from single views compared to visual assessments by physicians.

For each patient study, three qualitative bin assessments were performed as shown in Table 10 below.

Table 11. *EF Qualitative bin assessments*

Reference Standard	Caption Interpretation	Visual Assessment
Mean biplane-overread assigned to corresponding qualitative bin	Caption Interpretation output assigned to corresponding qualitative bin	Individual physician assessment of qualitative bin

For each patient and each view, the Caption Interpretation bin assessment and each Physician's bin assessment were compared to the reference standard bin assessment. Note that the reference standard is derived from the biplane Simpson's method.

Table 12. Paired-difference analysis (example)

Patient	Reference Standard (biplane)	Caption Interpretation (PLAX)	Physician 1 (PLAX)	Paired Difference
1	Normal	Reduced	Reduced	Caption Equal ("tie")
2	Reduced	Severely Reduced	Normal	Caption Equal ("tie")
3	Reduced	Normal	Reduced	Caption Worse
4	Severely Reduced	Severely Reduced	Reduced	Caption Better

The statistical analysis assessed whether AutoEF assessments were better than visual assessments compared to the reference standard using the sign test for each view. Tested was whether the number of samples with AutoEF being better was significantly more than the samples with AutoEF being worse, excluding the tied samples. Specifically, the empirical probability was computed as follows for each view and tested whether it was greater than 0.5:

$$P(\text{AutoEF Better}) = B / (B + W)$$

where B is the number of samples with AutoEF being better and W is the number of samples with AutoEF being worse.

This qualitative comparison was to proceed considering the variability due to the choice of patient studies first and then due to the choice of studies and individual physician's assessments (fixed sequence procedure).

The raw results of the qualitative comparisons are shown in the table below. For each of the three (3) views, a significant majority (>80%) of AutoEF bin predictions were at least as accurate as the physicians' assessments when an assessment was made (PLAX: 83.7%, AP4: 87.0%, AP2: 86.4%). AutoEF and physician bin assignments could not be compared for a given study/view physician combination when: 1) the physician did not provide an estimate, 2) AutoEF did not provide one, or 3) both. In all three views, AutoEF more accurately assessed a greater percentage of studies than did the physicians' visual assessments.

Table 13. *Caption Interpretation Automated EF qualitative comparison*

	PLAX only n (%)	AP4 only n (%)	AP2 only n (%)
AutoEF Better	242 (14.6%)	340 (20.5%)	323 (19.5%)
AutoEF Tied	933 (56.2%)	1019 (61.4%)	961 (57.9%)
AutoEF Worse	229 (13.8%)	202 (12.2%)	202 (12.2%)
No assessment	256 (15.4%)	99 (6.0%)	174 (10.5%)
Total	1,660 (100%)	1,660 (100%)	1,660 (100%)

The results of the primary hypothesis testing are shown in the table below. When taking into account the variability of studies alone, all three views demonstrate that a statistically significantly greater number of Caption Interpretation bin estimates were closer to the reference bin than those assessed by physicians. When taking into account variability arising from studies and physicians, AP4 and AP2 meet the primary endpoint. “PLAX Only” did not meet its “studies and physicians” primary endpoint, although the point estimate showed that more PLAX bin estimates were closer to the reference bin than those assessed by physicians.

Table 14. *Caption Interpretation qualitative comparison – hypothesis testing*

Source(s) of Variability	#	View(s)	P(AutoEF Better) %, [95% CI]	p-value
Studies	1	AP4 only	0.723 [64.9, 79.7]	< 0.001
	2	PLAX only	0.636 [55.3, 71.9]	< 0.001
	3	AP2 only	0.709 [63.3, 78.5]	< 0.001
Studies and Physicians	4	AP4 only	0.627 [58.7, 66.8]	< 0.001
	5	PLAX only	0.514 [46.9, 55.8]	0.259
	6	AP2 only	0.615 [57.3, 65.7]	N/A

Descriptive analyses

Note: Secondary endpoints and additional analyses presented below were not evaluated based on a specific hypothesis. Since the evaluation of secondary endpoints and additional analyses did not allow for control of Type I error, the study results are presented as a descriptive demonstration using Caption Interpretation for the specific secondary endpoints and additional analyses

Sonographer and cardiologist variability

Human variability was analyzed on all biplane-overread assessments to descriptively compare the experimental value to the literature-derived 9.2% RMSD endpoint from the clinical-validation study. The descriptive results indicate that Caption Interpretation EF estimation based on at least two views performs comparably to a cardiologist biplane-overread. All Caption Interpretation EF estimations were observed to have lower variability than a sonographer biplane trace.

Table 15. *Sonographer and cardiologist variability*

Caption Interpretation Automated EF	RMSD % [95% CI]
PLAX, AP4, and AP2	7.082 [6.561,7.568]
AP4 and AP2	7.460 [6.937,7.949]
PLAX and AP4	7.782 [7.228,8.299]
PLAX and AP2	7.799 [7.226,8.333]
Cardiologist	
PLAX only	9.88 [9.039,10.647]
AP4 only	8.565 [7.906,9.176]
AP2 only	8.759 [8.125,9.350]
Cardiologist Overread	7.907 [7.354, 8.424]
Sonographer Tracing	11.047 [10.057, 11.954]

Sensitivity and PPV of detecting reduced EF

The Sensitivity and PPV of AutoEF in identifying patients with reduced or severely reduced EF (i.e., EF < 53%) was calculated and compared to the sensitivity and PPV of physicians' visual assessments. The results are reported in the tables below.

Table 16. *AutoEF 2.0 PPV (Reduced/ Severely Reduced)*

View(s)	AutoEF PPV [95% CI]	Physician Visual Assessment PPV [95% CI]
PLAX, AP4, and AP2	0.951 [0.902, 0.988]	N/A
PLAX and AP4	0.963 [0.914, 1.000]	N/A
PLAX and AP2	0.918 [0.859, 0.976]	N/A
AP4 and AP2	0.938 [0.885, 0.979]	N/A
PLAX only	0.917 [0.857, 0.976]	0.930 [0.912, 0.947]
AP4 only	0.929 [0.878, 0.980]	0.947 [0.930, 0.961]
AP2 only	0.908 [0.847, 0.959]	0.959 [0.944, 0.972]

Table 17. *AutoEF 2.0 sensitivity*

View(s)	AutoEF Sensitivity (Reduced/ Severely Reduced) [95% CI]	Reader Sensitivity (Reduced/ Severely Reduced) [95% CI]
PLAX, AP4, and AP2	0.780 [0.70, 0.86]	N/A
PLAX and AP4	0.780 [0.700, 0.860]	N/A
PLAX and AP2	0.780 [0.700, 0.860]	N/A
AP4 and AP2	0.900 [0.840, 0.950]	N/A
PLAX only	0.770 [0.690, 0.850]	0.790 [0.760, 0.810]
AP4 only	0.910 [0.850, 0.960]	0.760 [0.730, 0.790]
AP2 only	0.890 [0.830, 0.950]	0.720 [0.690, 0.750]

Outlier assessment

The incidence of outliers in AutoEF assessment relative to the reference standard EF estimate is shown in Table 17 below. The incidence of observed outliers was very low for all EF estimates produced from a combination of views. The observed incidence of outliers was higher for single-view EF estimates, but was lower than the observed outliers of the physicians' estimates.

Table 18. *EF assessment outliers – Caption Interpretation AutoEF*

View(s)	Studies Assessed <i>n</i>	AutoEF error >15% EF <i>n</i> (%)
PLAX, AP4, and AP2	135	0 (0%)
PLAX and AP4	144	2 (1.39%)
PLAX and AP2	137	2 (1.46%)
AP4 and AP2	152	2 (1.32%)
PLAX only	147	11 (7.48%)
AP4 only	162	8 (4.94%)
AP2 only	155	9 (5.81%)

Table 19. *EF assessment outliers – physicians*

View(s)	Studies Assessed <i>n</i>	Physician error >15% EF <i>n</i> (%)
PLAX, AP4, and AP2	–	N/A
PLAX and AP4	–	N/A
PLAX and AP2	–	N/A
AP4 and AP2	–	N/A
PLAX only	1,524	203 (13.32%)
AP4 only	1,562	169 (10.82%)
AP2 only	1,568	224 (14.29%)

System Overview and Configuration

System components and workflow

In operation, the Caption Health AutoEF software provides automated measurements of left ventricular ejection fraction integrated into the Caption AI product. Studies are sent to the Caption Health AutoEF microservice module for processing. When the AutoEF process is completed, the outputs are forwarded to the Caption AI software. Using Caption AI, the reading physician reviews the patient study, which now includes the AutoEF outputs, along with the clips used to generate the Ejection Fraction result. The physician may use the AutoEF number as an input in creating the Patient Report.

Hardware and software requirements

The hardware and software requirements for the operation of the Caption Health AutoEF include the hardware platform for the Caption Health AutoEF software microservice, the operating system software for the Caption Health AutoEF software, and the DICOM-compliant interfaces. These requirements are superseded by the Caption AI hardware requirements. See “System specifications” on page A-25, or contact your Caption Health representative for details.

Configuration

CAUTION

The Caption Health AutoEF product is not a permanent storage device for patient studies.

CAUTION

To ensure proper operation of the Caption Health AutoEF software, configuration and setup of Caption AI and PACS DICOM Viewer(s) should be performed by an IT or PACS Administrator.

System Operation

This chapter describes the operation of Caption Health AutoEF. For important information regarding the proper use and application of AutoEF, see *Chapter 2: Safety*.

Caption Health AutoEF processing

Caption AI automatically forwards the images acquired during a study to the Caption Health AutoEF microservice in parallel to study acquisition. Once the user has acquired two acceptable views, the AutoEF measurement will be displayed, as acquiring at least two views will improve the accuracy of the measurement.

When Caption Health AutoEF receives a study, it automatically processes it to: (1) identify the Parasternal Long Axis, Apical Four Chamber, and Apical Two Chamber view clips; (2) assesses which of the identified clips are suitable for measurement; and (3) of the clips suitable for measurement, selects the best quality clip for each view. At most, three clips will be selected for AutoEF estimation, one for each view. Fewer than three clips may be selected, as few as none, if one or more of the views did not have a clip in the study that was suitable for measurement.

The selected clips are used to produce the AutoEF estimate. An individual Ejection Fraction estimate is produced for each view (PLAX, AP4, AP2) and an overall aggregate Ejection Fraction estimate is produced for the study that averages the Ejection Fraction estimates from each view.

For the aggregate AutoEF estimate produced for the study, the estimate will be accompanied with a list of views used to produce the estimate and an "expected error range" around the estimate to indicate the range in which the true Ejection Fraction is expected to lie. The greater the number of views that the AutoEF estimate is produced from, the more reliable we anticipate the estimate to be. A qualitative assessment of global LV function associated with the Ejection Fraction estimate will also be displayed below using a mapping from ASE guidelines.

- EF greater than 73% shall be mapped to "Hyperdynamic"
- EF greater than or equal to 53% but less than or equal to 73% shall be mapped to "Normal"
- EF greater than or equal to 30% but less than or equal to 52% shall be mapped to "Reduced"
- EF less than 30% shall be mapped to "Severely Reduced"

For the qualitative assessment of LV function displayed, a likelihood will be shown in % to indicate the probability that the qualitative assessment displayed is in fact the true qualitative assessment of global LV function for the patient. This is estimated by modeling the distribution of the true EF as a normal distribution with the mean at the predicted EF and the standard deviation as the Confidence Metric prediction, and then by computing the area under the curve between the bounds of the qualitative bin.

Single-view EF calculation

In certain circumstances, the user may decide that an AutoEF assessment based on a single view is appropriate, such as for patients for whom it is too technically difficult to obtain additional views, or in time-critical situations. To calculate the AutoEF measurement from a single view, the user selects “Calculate” (Figure 44) and confirms their choice (Figure 3). The AutoEF measurement from a single view will be displayed (Figure 4), and will automatically update if the user is able to acquire an acceptable additional view.




Figure 59. Detail view with AutoEF




Figure 60. AutoEF calculation confirmation




Figure 61. Detailed view with single AutoEF calculation

The user should verify that the images include Parasternal Long Axis, Apical Four Chamber, and Apical Two Chamber views that are appropriate for measurements, as applicable, following standard guidelines and their institution's guidelines regarding anatomical presentation, image quality, and the capture of at least one full heart cycle.

Accompanying information

The AutoEF estimate is accompanied by an “expected error range” (Figure 5) around the estimate to indicate the range in which the true Ejection Fraction is generally expected to lie. This error range is computed by modeling the distribution of the observed difference between the AutoEF estimate and the panel mean EF (i.e., the mean of the cardiologist EF estimates, each of which is produced by a single cardiologist in the panel, for a given study), denoted

as $\Delta EF_{\text{AutoEF, Panel}}$ across a large population of echocardiographic studies, as a function of image quality. Specifically, for a given study, this error range represents the expected standard deviation of the difference between the AutoEF and expert cardiologist panel estimates ($\Delta EF_{\text{AutoEF, Panel}}$), based on the image quality of the study. In other words, the error range represents one standard deviation of $\Delta EF_{\text{AutoEF, Panel}}$ that is expected for the image quality. This error range estimation is done with the assumption that the EF difference ($\Delta EF_{\text{AutoEF, Panel}}$) is normally distributed. A descriptive analysis of the Kolmogorov-Smirnov (K-S) Test for Normality was performed, showing a small K-S distance of 0.060 ($n = 166$, for best available views), indicating that the observed difference is normally distributed, supporting the appropriateness of modeling the $\Delta EF_{\text{AutoEF, Panel}}$ distribution.




Figure 62. Report screen with AutoEF estimate: expected error range

By looking at this expected error range, the user can evaluate the confidence in the estimate and any implications that error range has on their clinical assessment. For instance, a large error range may create uncertainty around whether the patient's Ejection Fraction estimate is normal or abnormal. The error range is estimated taking into account the underlying image quality of the clips in the study and is expected to be smaller in studies with higher image quality as well as in studies where the user captured at least two views that can be used for an assessment of Ejection Fraction.

For the aggregate AutoEF estimate produced for the study, the estimate will be accompanied with a list of views used to produce the estimate and an "expected error range" around the estimate to indicate the range in which the true Ejection Fraction is expected to lie. The greater the number of views that the AutoEF estimate is produced from, the more reliable the estimate is anticipated to be.

A qualitative assessment of global LV function associated with the Ejection Fraction estimate will also be displayed (Figure 6) using a mapping from ASE and ACEP guidelines.^{1,2}




Figure 63. Report screen with AutoEF estimate: global LV function

ASE guidelines report the following categories for assessment of left ventricular ejection fraction (EF).

Table 20. Recommendations for cardiac chamber quantification by echocardiography in adults (ASE/AACVI 2015): left ventricular ejection fraction

	Normal (%)	Mildly abnormal (%)	Moderately abnormal (%)	Severely abnormal (%)
Male	52–72	41–51	30–40	< 30
Female	54–74	41–53	30–40	< 30

Caption Health made two adjustments to simplify the guidelines for use at the point of care. First, the minor sex-related differences in the guidelines for the normal-hyperdynamic and normal-mildly abnormal boundaries were averaged so that the guidelines were sex-agnostic. These differences are minor, and the reconciled guidelines differ by only a single EF percentage point at these boundaries. Second, these guidelines were used to determine the categories defined by the ACEP in their reporting guidelines, which recommend the following: hyperdynamic, normal, reduced, and severely reduced. In order to

1 Lang RM, Badano LP, Mor-Avi V, Afilalo J, Armstrong A, Ernande L *et al.* "Recommendations for Cardiac Chamber Quantification by Echocardiography in Adults: An Update from the American Society of Echocardiography and the European Association of Cardiovascular Imaging" *Eur Heart J Cardiovasc Imaging* 2015;16:233–271.

2 ACEP Emergency Ultrasound Standard Reporting Guidelines. June 2018.

reconcile the two guidelines, the ASE “mildly abnormal” and “moderately abnormal” guidelines were collapsed into a single category. Ultimately these changes resulted in the final categorization used by AutoEF, shown in the table below.

Table 21. AutoEF 2.0 qualitative categorization of EF

Hyperdynamic (%)	Normal (%)	Reduced (%)	Severely reduced (%)
> 73	53-73	30-52	< 30

The software maps the point Ejection Fraction estimate to one of the categories of LV function defined above. Taking into consideration the expected error range around the point Ejection Fraction estimate, which may result in an overlap across two different categories of LV function, the software provides a “likelihood” (or probability to be more precise) for the displayed category of LV function (Figure 49). This “likelihood” (in %) indicates to the user how confident the software is that the true Ejection Fraction of the patient lies in the same category of LV function as the point Ejection Fraction estimate predicted by the software.




Figure 64. Report screen with AutoEF estimate: likelihood of LV function

Note: If the study did not contain any clips from any of the three views that were suitable for measurement, an Ejection Fraction estimate will not be produced for the study and the user will be notified that the software did not have sufficient information to produce an Ejection Fraction estimate for the study.

When AutoEF processing is complete, the updated study is forwarded automatically to the destination PACS server defined in the settings.

CAUTION

Review the video clips used in generating AutoEF measurements to ensure that the clips are appropriate for measuring ejection fraction.

Reviewing results on the PACS workstation

The AutoEF results are integrated into the patient study and are viewed via the Caption AI software, and can be forwarded to your DICOM viewer. The AutoEF outputs and the overall study can be viewed via the DICOM viewer. It is important for the physician preparing the Patient Report to carefully inspect the AutoEF result and the images used, as well as other contents of the study, in order to determine a reported ejection fraction estimate.

Editing results

AutoEF results may manually be made a part of a final Patient Report, so the reviewing physician may edit the result as needed. Users may also use the manual measurement tools of the viewing software to create new calculations, such as manual biplane Simpson's ejection fraction on the original versions of the clips.

System specifications

Table 22. *System specifications*

CPU/Motherboard	Core i7 (quad-core, hyperthreaded or higher) or equivalent/better
GPU	GeForce 1070 or 1080 or Titan X or Titan Xp or equivalent/better GPU
Memory	8GB DDR3 or equivalent/higher
Local System Storage	120 GB SSD or equivalent/greater
Power Supply	Thermaltake 500W or higher (depending on the number of GPUs) (or equivalent)
Remote Management	OpenSSH server or equivalent
OS	Ubuntu 16.04 LTS (Linux Long-Term-Support)
Database	Redis (stores async task messages)
Network	10 gigabit network card or equivalent

Safety and regulatory requirements

The system complies with the following voluntary standards:

- ISO 14971, Application of Risk Management to Medical Devices
- IEC 62366, Application of Usability Engineering to Medical Devices
- Digital Imaging and Communications in Medicine (DICOM) (NEMA Standard PS3)

APPENDIX B: CAPTION AI SCANNING METADATA

Caption AI's unique algorithms generate data on the quality of each scan performed with prescriptive guidance. You can use these insights to gauge the performance of scanners and optimize training to ensure diagnostic-quality imaging in the clinical workflow.

► To export scanning metadata

1. On the Caption AI **start** screen, tap the **Menu** icon (☰), then **Settings**. The **Settings** screen appears.
2. Select **Metadata**. The **Metadata Export** screen appears.




Figure 65. Metadata Export screen

Note: Data is anonymized to minimize the risk of compromising patient information. To remove provider IDs, select **Anonymize Provider ID**. Please take steps to secure sensitive information and work with your IT department to configure Caption AI in accordance with your institution's policies.

3. Tap **Export**. The scanning metadata is exported as a JSON file on the device's desktop.

The scanning-metadata file contains information about the use and quality of ultrasound scans on the device. For each study, this data includes:

- The anonymized date on which the study was performed
- The provider who performed the exam (if applicable)
- The scanning protocol used
- Algorithmic interpretations generated by Caption AI

For each view in a study, the following additional information is included:

- The diagnostic quality on a scale from 0 (low quality) to 100 (high quality)
- The amount of time it took to capture the clip
- The recording mode (“AutoCapture”, “Save Best Clip”, or manual)

The metadata file also includes information about the device configuration, such as saved scanning protocols and software settings.

For help analyzing this data or to set up a Quality Review of scanning at your institution, contact your Caption Health representative.