



User manual – Health care platform & Patient app for HCPs

Get insight in the rehabilitation of your patients and provide close follow up

App & Healthcare platform version: 1.34.1 till current

Manual version: 11

Date of Publication: 2023/03/06

Link to e-learning: <https://moveup.talentlms.com/index>

Table of contents

1. Instructions
2. Warnings and precautions
3. Intended use and intended users
4. Indication for use
5. Contact details
6. Patient app
 - a) Patient app functionalities
 - b) How to pair a wearable
7. Healthcare platform
 - a) Healthcare platform functionalities
 - b) Clinical guidelines for physiotherapists
 - c) Registration of a moveUP patient
 - d) Stop the treatment of a patient in moveUP

INSTRUCTIONS



Read carefully this instruction for use. Make sure that you know how to use moveUP, but also understand the limitations of moveUP before you start using our service. If you do have any questions on the use, contact moveUP.



In the final chapter of this document you find our contacts, which you can use during the use of our service, in case you would have technical problems with the application (moveUP dashboard or patient app) or when expertise is required about a care path.



The Healthcare practitioner can't give access to third parties to access the moveUP dashboard. Every other person can, even by accident, send out wrong information or adapt the patient profile, with the result that the advice sent by moveUP no longer corresponds with the needs of the patient.



In case the application is use while moveUP & b.clinic is not asked to do the follow-up of patients, the individual healthcare provider or the health care facility is responsible for the follow-up of the patient and the action that are carried out.



In compliance with the medical device legislation, you are required to report any complaints towards the application, any alteration regarding characteristics / specifications of the application and incidents / complications potentially linked to the use of the application to support@moveup.care. The AFMPS can also be contacted in case of incident via https://www.famhp.be/en/human_use/health_products/medical_devices_accessories/materiovigilance/how_notify



In compliance with the GDPR legislation, please report eventual data leaks linked to the use of moveUP. Please, read the entire privacy policy : <https://www.moveup.care/privacy>

Warnings and precautions



- moveUP is not an emergency tool. In case of emergency please contact your doctor or call 112.
- moveUP is not intended to monitor/treat vital parameters of critical diseases.



- moveUP companion is not a replacement of any treatment you need to follow, but is used as an addition.
- Patients need to be aware if questionnaires are not filled out on frequent basis or not filled out trustworthy, the care team has more difficulty to know their health situation & evolution.



- The moveUP App offers personalized healthcare based on the patient's profile and therefore can not be shared, borrowed or traded between users. The information and advice of the moveUP application do not apply to a patient other than the one whose profile is stored in the application. The patient may not grant access to the moveUP App to third parties. Any other person can, even accidentally, send wrong information or change the profile of the patient, with the result that the moveUP-guided advice no longer corresponds to the needs of the patient.

With software, there is always a residual risk for bugs. If you have the feeling something is wrong or you don't get a response in the app, please reach out to us via support@ or call 0800 88 008.

The moveUP IFU are only supplied in electronic format. If you need a full printed version please contact us at the e-mail address info@moveup.care. This mail is available 24/7



INTENDED USE

For healthcare providers, moveUP provides efficient clinical management of pathologies and treatments. With the use of moveUP, valuable insights are given to the HCP/care team of the status of the patient and its evolution.

For patients, moveUP provides personalised information and instructions to help them managing their symptoms and progress in their rehabilitation. The intensity of follow-up is adapted based on patients' needs and timepoint in the patient pathway:

moveUP companion = monitoring and information, no active follow-up

moveUP companion offers targeted information and evolution reports to patients. Their care team takes the evolution reports into account in their further recovery when appropriate.

moveUP coach = active follow-up by healthcare team

moveUP coach offers targeted information and evolution reports to patients. Their care team is more actively involved and takes the evolution reports into account and can provide advice and exercise suggestions through the digital platform.

moveUP therapy = active follow-up by healthcare team, with data driven validated care protocol

moveUP therapy offers a data-driven validated care protocol with certain category and level of exercises and activities, specifically targeted to the individual patient. Their care team can manually adapt the data driven validated protocol when needed. Patients can fully rehabilitation with moveUP without leaving their home environment.

Optional functionality that can be enabled: interoperability with Class IIa continuous passive motion (CPM) medical device. For knee and hip patients who are using a Class IIa CPM medical device, moveUP can interoperate with the Class IIa CPM software. moveUP acts as a facilitator to easy assign a designated CPM exercise protocol chosen by a physician to a patient and to display the performed CPM exercises in the medical dashboard.



Contra-indications and potential side effects

There are no contra-indications or known side effects.

Human body contact

No human body contact with patient or user, due to nature of the product (software).

Accessories/products used in combination

There are no accessories. If the patient has no compatible tablet/smartphone or activity tracker, moveUP can lease these devices to the patient. The leased devices are CE marked devices that meets the compatibility criteria outlined in the IFU for the app.

Device lifetime

2-years. Depending of the willingness of the user to update the app. We sent out a notification to recommend to update the app, to make sure the app will keep the performance & recent included features.

Claims

moveUP companion/coach/therapy

The intensity of follow-up is adapted based on the needs of the patient, via the moveUP Symptom & QoL monitoring tool.

moveUP enhances the clinical management of the patients, because early detection & management of complications is possible via the symptom & QoL monitoring tool

More efficient clinical management, such as the number of consultations can be reduced

Enforces therapy compliance / adherence

The correct information is provided at the right time

moveUP therapy

With the use of moveUP therapy knee & hip arthroplasty patients can fully rehabilitate via the in-app care team without leaving their home environment.

Intended users

moveUP is intended to be used by patients and healthcare providers. The main user of the mobile app and patient website is the patient.

Inclusion:

- Age: minimum 18 years / maximum no limit
- Health & condition: capable of performing basic activities of daily living
- Language: understanding one of the available languages of the app (Dutch, French, German, English)

Exclusion:

- Patients who are mentally incompetent or having troubles to express what they are feeling (for instance, mentally diseased people, people staying in elderly care centres, ...) are excluded.
- Patients who are not capable of operating a tablet/smartphone and activity tracker.
- Patients who can't understand one of the available languages of the app (Dutch, French, German, English)

The main user of the medical web interface is a healthcare practitioner (group) or clinical researcher (group), named the care team. The care team is able to operate a web interface via web browser on PC/tablet/smartphone. The healthcare practitioner needs to understand one of the available languages of the web interface (today only available in English).


COMPANY DETAILS



Head office Brussel: Parvis sainte-gudule 5,
B-1000 Brussels, Belgium Office Ghent:
Oktrooiplein 1, B-9000 Ghent, Belgium
Tel: 0800 88 008
Email: info@moveUP.care www.moveUP.care

If you need the eIFU in paper format please contact moveUP and we'll supply you a copy in paper in 7 days.

moveUP supplies the IFU in electronic format due to the topology of the product (software) and the availability of the eIFU at any moment that is needed for the user.



+G166G166MOVEUPV1340WS00+\$71341/16D202303100Y

Basic UDI-DI: ++G166MOVEUPV1340WS



In case of questions about the App, the tablet or the smart bracelet, you can contact moveUP via the following channels - in order of priority:

- Via the message system of the App. Please preface the message with **“Technical question:”**
- Via email to support@moveup.care
- Via phone during office hours: **0800 88 008**







moveUP

Patient app


CE



1




2




The patient navigates between 3 different tabs "Overview", "Exercises" or "Questions".

3




4



On top, you will find

- A warm welcome message
- Today's date
- Number of days since the surgery

5



Messaging function.


In the messaging function, patients can ask questions and voice their concerns.

The care team consisting out of specialized doctors and physical therapists to respond and give advice.

All messages are read by the physical therapists. Messages which are not within their professional competence are referred to the doctors.

6






In the messaging function, there is also the possibility to send pictures. For example, to follow up wound or swelling of the operated joint.




7



By clicking the arrow in the top left corner, you return to the overview screen.

8




At specific time points, the patient is asked to complete questionnaires, information modules and video requests in "Surveys and Information".

- **Medical questionnaires** to complete the patient's profile
- **Information modules:** right information on the right time
- **Video requests:** gait analysis, range of motion

All available requests need to be completed as soon as they become available

When the medical questionnaires or video requests are completed, they disappear from the screen.


9



Additional questionnaires: Assess medical status before and after the operation.

- ✓ PROMs
- ✓ Satisfaction, PREMs
- ✓ Psychosocial risk factors
- ✓ ...


10



Video requests: Follow up on gait and range of motion.



11




Every video request has


- a description
- an example video

12







13




14




15



16




17




18





Every day, a personalized step maximum is displayed by a number and a red line. This personal step maximum is not an absolute maximum but avoids overloading the operated joint and makes the patient aware of its own activity.


19



At the bottom of the overview screen you can find the last section.


"More information about your surgery" opens a page with practical information helpful before and after surgery.

20



The mid tab: "Exercises"

21



Therapeutic exercises


In the middle tab you will find the available exercises.

On average, there are 3 to 6 exercises active daily.

Only the exercises that need to be performed that day are shown in random order.

The exercises are adjusted according to the parameters and feedback of the patient.

22




Title and instructions

The name of the exercise and detailed instructions on how to perform the exercise are available.

Additional instructions can be provided via the messaging function.

23



Video


There is a short video available per exercise, which shows how the exercise should be performed.

Click on the triangle ▶ to start the video.

Click on the rectangle ◻ to see the video on a larger screen.

24






Intensity and frequency

Frequency and intensity are provided per exercise. e.g. 6x a day, 1 minute

We recommend not to exceed the number of prescribed series and repetitions.

25




Reporting

The patient reports the number of series performed via the + and - symbol.


In this way, the healthcare team is informed of the number of exercises performed.

26



The right tab: "Questionnaire"

27




Daily questions

By completing the daily questionnaire, the healthcare team is provided with the necessary information to guide and follow up the patient in the best possible way.

The questionnaire is available between 4pm and midnight.

Completing the daily questionnaire takes about 2-3 minutes.

28




Send

Patients must complete the questionnaire on a daily basis to receive proper follow-up.

Click on "send" at the bottom of the questionnaire.

29




Report event

Events can be reported when the recurring questionnaire is not available yet or already filled out.

30






Report event


When reporting an event, an overview of event types is displayed. One or more event categories can be selected.

31



For each selected category, a set of questions is displayed that need to be filled out. Once all questionnaires have been answered, the given answers are sent to the practitioner.


32



App menu

Click on the  drop-down button to display the app menu.


33



moveUP journeys

On this page, all past journeys are visible. This includes surgery details, questionnaires and exercises.

34



Non active-follow up

When a patient is beyond active follow up, the moveUP app will be in a 'non-active follow-up mode'.

In this mode, the navigation has been minimized to what the patient needs at this point with only 2 tabs: An 'overview' and a new 'info page' containing the frequently asked questions.

Once the active follow up for the patient has been stopped, sending chat messages will be available for one more week. The chat messages sent during this week will automatically be converted into an email to the help desk. Our help desk will answer the message in the message function. After this period, it is no longer possible to send messages, but the message history remains accessible.

35




How to pair your Garmin Vivofit 4 & Withings go activity tracker with the moveUP app?


i The following activity trackers are recommended by moveUP for usability and accuracy reasons. You can use any activity tracker of your choice.

i Pairing the activity tracker with the moveUP app is essential to your care trajectory. The pairing only needs to be done once.

To pair your activity tracker click the green button on your homescreen.



Choose "Garmin Vivofit 4" and click continue




Go through all the steps that will appear on your screen.

! Be aware: when you are pairing the activity tracker a 6-digit code will show on the screen of the **activity tracker**.


On some android devices, a pop-up appears on the screen of the phone with a message the code is often 1111 or 1234. Please ignore that message. It is the **6-digit code from the activity tracker** screen you have to enter.





- If everything went well you will end up on this screen.
- You can click on the arrow at the left top to go back to the overview screen.




 If you have problems pairing your activity tracker please send an email to support@moveup.care

Withings Go:

To pair your activity tracker click the green button on your homescreen.


Choose "Personal Activity tracker"



Select the Withings option



Pair the moveUP app with your Withings account. If you have any doubt you can find more info at: <https://support.withings.com/hc/en-us/articles/214484018-Installing-my-Go->



Your device is paired with the moveUP app.

Activity tracker paired!








How to synchronize my steps daily?

We advise you to wear your activity tracker around your wrist day and night. However, using during the day, is more important for us to register your activity. The activity tracker is always on and captures how active you are, hence it does not need to be activated.


You must need to manually sync your activity tracker data **daily** in the moveUP app. This way your care team can monitor how active you are and decide if any adaptations to your exercise scheme are needed.

 **We advise syncing your steps daily in the evening.**
This way your care team will see the next day how active you were.


 Make sure you have at least made 100 steps with your activity tracker before you try to synchronize it.
Don't sync your steps too often during the day. Once a day after filling out the daily questionnaire is more than enough.

 The activity tracker has many options in the menu, but we will not be using them. The instructions below are the only ones you will use during your recovery. At the end of the treatment, if you wish to keep the activity tracker, we will send you a full manual by email.

- Press the synchronizing button marked in red.




- You will see the following screen. Click on “Sync”



- Follow the instructions on the screen. Basically you have to press the silver button on your activity tracker twice for 2 seconds.
- At the end you will automatically be taken back to the homescreen and your steps will be visible in the graph.




- You will see the synced steps data appear on the graph on your home screen.



i If you have problems synchronizing please send an email to support@moveup.care and try to describe exactly in which step you encounter problems and what messages you see or get on your screen.






moveUP

Healthcare Platform

CE




1



2




3



4




5




6






7



8



9



10

Candidate Therapy Coach Companion PROM Completed All

Candidates
Patients who are awaiting a care plan. They downloaded the moveUP app and registered or are created by a hcp.

Therapy
Patients with an active Therapy care plan before or after surgery.

Coach
Patients with an active Coach care plan before or after surgery.


Companion
Patients with an active Companion care plan before or after surgery.

PROM
Patients who will still receive medical questionnaires for long term follow up on specific timepoints (e.g. 3m - 6m, ...).

Completed
Patients who have completed their long term follow up.


All
All patients.

11




12







13




14



15



16



17



18





First line review

- 1 Patient was checked and validated 1 day ago
- ✓ Patient was checked and validated today by your discipline
- ✓ Patient was checked and validated today by other discipline
- 🚩 Patient was flagged for your discipline
- 🚩 Patient was flagged for another discipline

19

Patient compliance is visualized in the patient overview

- Compliant for the last 2 days
- Non compliant > 2 days
- Non compliant > 7 days

20

Contact - support@moveup.care

Need support? Do not hesitate to contact moveUP! →

21

Activated can mean 2 things: either the patient registered him/herself with a code through the app or he has been marked as immediately activated in onboarding.

- Non activated: patients who haven't registered themselves through the app yet and haven't been marked as immediately activated.

22


Data can be extracted by exporting to a csv file

23


Mail templates can be sent via email

24







25




26




27



28




29




30







31




32




33



34




35




36







37




38




39



40



41




42






Patient details - Status


43




44




45



46



47



48





Patient details - Status

Table overview last 3 days: activities, medication and walking aids

Compliance on the prescribed medication
Questions about medication are escalated to the doctors

Compliance on the use of walking aids

Question	Has patient answered?	Has doctor answered?	Priority
Medication	Yes	Yes	Low
Walking aids	Yes	Yes	Low
Medication	Yes	Yes	Low
Walking aids	Yes	Yes	Low
Medication	Yes	Yes	Low
Walking aids	Yes	Yes	Low

49

Patient details - Status

Click, to see all answers on the daily questionnaire

Question	Has patient answered?	Has doctor answered?	Priority
Medication	Yes	Yes	Low
Walking aids	Yes	Yes	Low
Medication	Yes	Yes	Low
Walking aids	Yes	Yes	Low
Medication	Yes	Yes	Low
Walking aids	Yes	Yes	Low

50

Patient details - Status

Scroll down

Plan 2: Medication

Plan 3: Medication

Plan 4: Medication

Plan 5: Medication

Plan 6: Medication

Plan 7: Medication

Plan 8: Medication

Plan 9: Medication

Plan 10: Medication

Plan 11: Medication

Plan 12: Medication

Plan 13: Medication

Plan 14: Medication

Plan 15: Medication

Plan 16: Medication

Plan 17: Medication

Plan 18: Medication

Plan 19: Medication

Plan 20: Medication

Plan 21: Medication

Plan 22: Medication

Plan 23: Medication

Plan 24: Medication

Plan 25: Medication

Plan 26: Medication

Plan 27: Medication

Plan 28: Medication

Plan 29: Medication

Plan 30: Medication

Plan 31: Medication

Plan 32: Medication

Plan 33: Medication

Plan 34: Medication

Plan 35: Medication

Plan 36: Medication

Plan 37: Medication

Plan 38: Medication

Plan 39: Medication

Plan 40: Medication

Plan 41: Medication

Plan 42: Medication

Plan 43: Medication

Plan 44: Medication

Plan 45: Medication

Plan 46: Medication

Plan 47: Medication

Plan 48: Medication

Plan 49: Medication

Plan 50: Medication

Plan 51: Medication

Plan 52: Medication

Plan 53: Medication

Plan 54: Medication

Plan 55: Medication

Plan 56: Medication

Plan 57: Medication

Plan 58: Medication

Plan 59: Medication

Plan 60: Medication

Plan 61: Medication

Plan 62: Medication

Plan 63: Medication

Plan 64: Medication

Plan 65: Medication

Plan 66: Medication

Plan 67: Medication

Plan 68: Medication

Plan 69: Medication

Plan 70: Medication

Plan 71: Medication

Plan 72: Medication

Plan 73: Medication

Plan 74: Medication

Plan 75: Medication

Plan 76: Medication

Plan 77: Medication

Plan 78: Medication

Plan 79: Medication

Plan 80: Medication

Plan 81: Medication

Plan 82: Medication

Plan 83: Medication

Plan 84: Medication

Plan 85: Medication

Plan 86: Medication

Plan 87: Medication

Plan 88: Medication

Plan 89: Medication

Plan 90: Medication

Plan 91: Medication

Plan 92: Medication

Plan 93: Medication

Plan 94: Medication

Plan 95: Medication

Plan 96: Medication

Plan 97: Medication

Plan 98: Medication

Plan 99: Medication

Plan 100: Medication

51

Patient details - Status

Visualized evolution of

- ✓ Activity
- ✓ % executed exercises
- ✓ Reported pain of the operated joint: night, day, during exercises
- ✓ Reported swelling, warmth, stiffness of the operated joint
- ✓ RPE

52

Patiëntgegevens - Status

Preoperative data

Postoperative data

Day of surgery

Current date

53

Patient details - Status

Maximum steps

Minimum steps

Steps / day

Pre-operative days

Post operative days

Day of surgery


Today

Physical activity


Gradual increase in activity is promoted by a step limit

54






55



56



57

58

59

60





Name	Result
Age	67
Sex	M
Height	175
Weight	75
BMI	24

61

ID	Questions	Answers
1	Do you have any pain?	No
2	Do you have any difficulty in walking or any other type of motor activity?	No
3	Do you have any difficulty in sitting or standing for long periods?	No
4	Do you have any difficulty in climbing stairs?	No

62

Questionnaire not (yet) completed

63

See the answers of the (daily) questionnaires

64


Click to see questions and answers

65


Questionnaire	Date	Score	Max Score	Pass
KOOS	Monday, 11 January 2023 @ 14:00:00	0.00	0.00	No
	Monday, 11 January 2023 @ 14:00:00	0.00	0.00	No
Pain	Monday, 11 January 2023 @ 14:00:00	0.00	0.00	No
	Monday, 11 January 2023 @ 14:00:00	0.00	0.00	No
General health	Monday, 11 January 2023 @ 14:00:00	0.00	0.00	No
	Monday, 11 January 2023 @ 14:00:00	0.00	0.00	No
IC completion	Monday, 11 January 2023 @ 14:00:00	0.00	0.00	No
	Monday, 11 January 2023 @ 14:00:00	0.00	0.00	No

66







67




68




69



70



71



72





Patient details - Log

73

Patient details- Log

- ✓ Multidisciplinary communication
- ✓ Adjusted exercise scheme or step limit
- ✓ Automatic mail reminders
- ✓ Patient called
- ✓ Administrative actions
- ✓ Flags

74

Patient details - Log

When you flag someone, this automatically appears in the log of the patient:

GREEN: Solved flag
 RED: Unsolved flag for your discipline
 YELLOW: Unsolved flag for other discipline

75

Patient details - Report

76

Patient details Report

This is a patient/medical professional report used to inform the main milestones and medical scores of the patient, as well as the evolution of the patient during their treatment with moveUP.

77

The report can be downloaded in pdf form by clicking the "Download report" button. This way you can send the report to the patient.

It's also possible to share the link directly to the report, using the button "Copy report link".

78





Patient details Report

The proof of Engagement report gathers the dedication of the care team to the patient on the one hand and the adherence of the patient to the system on the other hand.
 The purpose of the report is to be able to **prove activity for reimbursement**. Therefore two main variables are taken into account:
 1. Time spent by the HCP by being active on the patient specific pages and by (logged) calling to the patient.
 2. Actions by the patient in the app subdivided into: answers on daily questionnaires, sent messages and syncing the activity tracker

An export button copies the displayed table to an excel file

79

Patient details Proof of engagement report

In this slider you can select the time frame of which you want to see the actions

To find out what each column means, you can click on the question mark next to the name above the table.

80

Patient details Proof of engagement report

Below, you can find a weekly/monthly report that summarizes the actions of a given week, for billing purposes.

81

Patient details
-
Careplan
patient's adapted treatment

82

Patient details - Care plan Exercises

83

Patient details - Care plan Exercises

84





Patient details - Care plan Exercises

See the exercises for today
Adapt description, frequency and intensity

85

Patient details - Care plan Exercises

Adapt the daily message displayed in the exercise scheme

86

Patient details - Care plan Exercises

See and adapt the exercise scheme.

The standardized exercise scheme can be extended with a set of exercises that are not part of the exercises defined in the care plan.

87

Patient details - Care plan

A standardized exercise scheme is personalized based on the patient's need

Blue check marks means that the exercise has been activated

Gray check marks means that the exercise has not been activated

88

Patient details - Care plan

A standardized exercise scheme is personalized based on the patient's need

There is a planning feature to schedule exercises over a longer period of time (click on *Advanced settings*)

- Start day/End day:
 - Start day = the day of the cell you clicked on
 - End day = Fill in the day you want this exercise to end
- Repeat input field:
 - Choose how the exercise will repeat till the end day
- Intensity/ Frequency START:
 - This is the value where the patient will start on
- Intensity/Frequency TARGET:
 - This is the value where the patient will end on. E.g. Start day=10, End day=20, Intensity start= 1, Intensity target=10 => On day 15 the intensity will be automatically calculated and in this case, will be set to 5.

89

Patient details - Care plan

A standardized exercise scheme is personalized based on the patient's need

The exercise scheme grid can be filtered using the label filter dropdown at the top.

90





Patient details – Care plan
A standardized exercise scheme is personalized based on the patient's need

When selecting a label in the dropdown the displayed exercise columns will be filtered according to the selected exercise labels. Following filter options are available

- Activated exercises:** Only show exercises that are currently active in the patients scheme
- Standard exercises:** Only show exercises that are part of the care plans core exercise list
- Extended exercises:** Only show exercises that are part of the care plans extended exercise list
- List of labels enabled for the care plan (e.g. Balance, Endurance, ...):** Only show exercises matching the selected labels.

91

Patient details – Care plan Step goal

92

Patient details – Care plan Step goal

Navigate to previous and coming days

93

Patient details – Care plan Step goal

See and adapt the step maximum

94

Patient details – Care plan Step goal

By clicking on "Save", you can adjust the number of steps for that day (example: day 5)

95

Patient details – Care plan step goal

By clicking on "Save & update previous + future days" you can adjust the number of steps for the previous and future days. This way you can change the number of steps for the next 10 days for example.

For example you go 10 days further (in this case you navigate to day 15) and set f.e. step maximum to 4000. Then the number of steps escalate progressively until day 15. (See next slides)

96





Patient details - Care plan Step goal

Fill in '4000'

97

Patient details - Care plan Step goal

All previous and future days are updated automatically


98

Patient details - Care plan Recurring questionnaire

99

Patient details - Care plan Medical questionnaire

100




101


Patient details - Care plan Video analysis

102






103




104




105



106



107




108





Patient details - Documents


109



110

Patient details - More

111




112




113






1




2




3



4



5



6





Direct gain - e.g. visualization of patient reported parameters: Joint warmth, swelling

- Patient's reported warmth and swelling of the operated joint can provide usefull insights on the healing proces.
- See below a decrease in patient's reported swelling (green).

7

Direct gain - e.g. insight in patient's activity captured by a wrist worn activity tracker

Insight in postoperative activity levels

- Which activity level does the patient reach?
- Does the patient gradually increase activity?
- Does the patient follow the advised step maximum?
- ...

Insight in preoperative activity levels

Physical activity

8

Direct gain - e.g. consider daily reported activities, medication-intake and crutches

Table overview last 3 days: activities, medication and crutches

- Check compliance on the prescribed medication
- Questions about medication are escalated to the doctors
- Compliance on the use of walking aids
- Give advice about transition 2 -> 1 -> 0 crutches/walker
- See which activities the patient has already picked up

Question	Yes/No/Not answered	Answer	Escalated
Medication	Yes	Yes	Escalated
Walking aids	Yes	Yes	Escalated
Crutches	Yes	Yes	Escalated

9

Escalation to 2nd line doctors

Flag is automatically reported in patients log:

- Picture of wound
- Alarming or atypical parameters
- High pain levels and possibility to adapt medication
- Non compliance of medication intake
- Question about medication
- Extra motivation
- ...

10

Direct gain - Track your patient's range of motion over time

11

Direct gain - Evaluate specific PROMs answers

Formula Calculations


Name	Result
...	...

KOOS (NL-FR-EN)

Questions	Answers
...	...

12





13



14

Knee arthroplasty

Total knee replacement
Uni knee replacement

15

Preparation	Hospitalization Protection post-surgical knee	Return to daily activity	Return to active lifestyle
Estimated time	Day 0 - 21	Day 7 - 60	Day 95 - ...
Events	Surgery Hospital stay Discharge First week home Removal of staples	First postoperative consult Week 6 PROMs evaluation	Second postoperative consult Month 3 PROMs evaluation Back to work Driving car
Goal	Protection post-surgical knee Prevent complications Pain & symptoms control Regain proprioception Knee movement within pain limits	Gradually decrease in pain <2/10 Gradual stop of medication Gain in pain free basic strength ROM 0-90° week 4 Normalized gait pattern Basic balance / coordination	ROM 0 - 110° week 8 Rough full ROM 0 crutches/walker device Regain functional strength
Actions (See tables)	Educate & support patient Supervised use of crutches/walker Supervise medication intake Provide pain-free exercises Control activity level Check wound Check mobility	Educate & support patient Support transfer 2 crutches > 1/0 crutches Adapt exercise scheme Ask pictures of swelling/haematoma Add non default exercises	Adapt exercise scheme Resume normal progression Add non default exercises Finalize rehabilitation
Criteria -> next phase	Protection & education Pain & symptom control	A minimum of / no pain A minimum feeling of warmth of the knee Painkilling medication stopped or limited to occasional intake of WHO type 1 medication No more NSAIDs	For some /a/ or exercises surgeon approval is advised e.g. driving, biking, lunges...

16

Preparation	Hospitalization Protection post-surgical knee	Return to daily activity	Return to active lifestyle
Isometric quadriceps setting - 3 x 10	Ice - 8 x 15 min. Circulation - 8 x 2 min. Heel raises - 3 x 10	Hiq/Knee flexion/extension standing - 3 x 15 Half rounds stationary bike - 2 x 5 min. Hiq/knee flexion to hip/knee extension - 3 x 10 Hip flexion and knee extension stability exercises on the chair - 3 x 15 Mini squat - 3 x 10 Side toe taps - 3 x 15 Standing march - 3 x 4 min	Stationary bike - 3 x 10 min Squat - 3 x 15 Hip abduction standing - 3 x 15 Step up - 3 x 10 Single leg stance - 3 x 10 Small lunges - 3 x 10 Forward lunge - 3 x 10 Side lunge - 3 x 10 Kneeling - 3 x 10

17

Preparation	Hospitalization Protection post-surgical knee	Return to daily activity	Return to active lifestyle
Activity	1. Preop avg < 4000 + Max steps 1000/day 2. Preop avg 4000 - 5000 Max steps 1500/day 3. Preop avg > 5000 + Max steps 2000/day From day 5: + 100 steps/day	Gradual increase until 50-100% preop steps Positive evolution +100-200 steps / day Negative evolution a. Max steps on hold b. Decrease in Max steps	Gradual increase until 75% - 100% or > preop steps Positive evolution +100-200 steps / day Negative evolution a. Max steps on hold b. Decrease in Max steps Stop PK No NSAID
Medication	PK: WHO 1 - 2 NSAID 2/day - 2 weeks Blood thinning medication	PK: WHO 1 NSAID 1/day - 2 weeks	No NSAID
Mobility	Pain free mobility exercises without forcing	ROM 0° - 90° week 4 Focus on full extension	ROM 0 - 130° week 8 Full ROM
Crutches	2 crutches	2 or 1 crutch	No crutches
Stationary bike	No bike	Stationary bike no resistance No biking	Stationary bike with resistance Biking after surgeon approval
Car		Driving (when feeling safe to do so) *Check the need for surgeon or car insurance approval	

18





Preparation	Hospitalization Protection post-surgical knee	Return to daily activity	Return to active lifestyle
Activity	High activity >50% prep Disappointed unilateral weightbearing Jumping Running Biking Loading activities or sports	High activity >100% prep Jumping Running Biking "High" loading activities or sports	Jumping Running "High" loading activities or sports → Surgeon dependent
Exercises	Pain provoking exercises Exercises against heavy resistance	Pain provoking exercises Exercises inducing pain after execution	Pain provoking exercises Exercises inducing pain after execution
Mobility	Forcing towards flexion / extension	Forcing towards flexion / extension	Forcing towards flexion / extension
Crutches	No use of crutches	No use of crutch(es) for long distances	
Pain	Pain provoking movements or exercises	Pain provoking movements or exercises	Pain provoking movements or exercises
Sleep	Rest/sleep with a fixed knee (e.g. pillow under operated joint)	Rest/sleep with a fixed knee (e.g. pillow under operated joint)	Rest/sleep with a fixed knee (e.g. pillow under operated joint)

! DONT'S
Table knee: Don'ts = contra-indications

moveUP

19

Hospitalization & Protection post surgical knee replacement

is a first postoperative phase.

- Focus on protection of the post-surgical knee and education of the patient on how to do so: the patient should take prescribed medication, use crutches, put ice and remove rest and elevation of the operated leg.
- Symptom control: Aim to get pain, warmth and swelling under control as soon as possible. Persisting swelling, pain and warmth will prohibit a smooth rehab. So, control and prevention of these parameters is key.
- Low load mobility: Motivate the patient to execute low load mobility exercises frequently. Changing position regularly will reduce stiffness.

How?

Exercises: See the exercises in the table. The priority are circulatory and light mobility exercises. Circulatory exercises help drain postoperative blood and swelling and prevent thrombosis. Quadriceps setting stimulates proprioception and maintenance of muscle tone. These exercises are automatically active in the exercise scheme. Exercises against resistance or those which provoke pain, swelling, stiffness and warmth have to be avoided, including those during the first 48 hours.

Activity: Guidance in activity is provided by a daily, personalized step maximum. (see table for reference)

The chat: is used to educate, motivate the patient, address concerns, answer questions and clarify therapeutic actions.

The use of crutches/walkers: should be supervised. Walking aids are important to prevent falling and to lower the load on the knee so soft tissue get time to heal.

Gait pattern: Check on heel-to-toe gait pattern and on neutral foot positioning during video analysis.

Knee range of motion: is recorded by weekly video inquiries. The first mobility check is done at day 2, this first check is a reference point to track improvement. During RDM analysis, also look for movement flunty (e.g. signs of fear or discomfort), early soft knee symptoms and compensation. Reaching full extension without overpressure has high priority in the beginning of the rehabilitation.

What if...

- If a patient doesn't manage to get the pain under control? → In first 48 hours can be triggered to give advice on medication. Most TKA patients are advised to start with third step 3 paracetamol. Step 2 paracetamol can temporarily help high pain persisting patients.
- If a patient reports a 'red flag'? Always flag 2nd line doctors in cases of a bleeding wound, picture of wounds, generic problems, reports of fever, postop nausea, visit to general practitioner, medication advice or any other sign of complication.
- If regarding the wound care: a picture of the wound after removal of electroclips should be asked. Most wounds are closed. The clips are removed starting from 18 days to 21 days postoperatively, until clips are removed a patch should cover the wound. In most cases the patch should be changed until this day (unless specified by surgeon).

20

Common actions Hospitalization & Protection post surgical knee

Exercises (see table)	Messages	Check
Default scheme	Questions	Therapy Compliance
High activity → no more need for circulatory exercises	<ul style="list-style-type: none"> How did surgery go? How do you feel at home? How do the exercises going? 	
	Feedback	Picture Video RDM Early RDM (progression RDM not necessary) Swelling Warmth
	Stimulate	Use of crutches Medication intake Daily activity and step max

Actions Hospitalisation & protection post-surgical knee

moveUP

21

Protection post-surgical knee

Protection & education
Pain & symptom control

Return to daily activity

CRITERIA

- The patient is compliant to the postoperative rehabilitation guidelines.
- The patient found a balance between low load, regular, basic exercises and rest to decrease/control post-operative pain and swelling.
- Transition to the next phase is advised when the patient has control over the postoperative complaints.

Transition phase protection post-surgical knee → Return to daily activity

22

Return to daily activity:

is a second postoperative phase.

- Pursue a continuous decrease in symptoms and pain while reducing WHO type 1-2 medication to a minimum next to increasing load by new exercises and allowing higher activity.
- Support the patient finding the right balance between rest and activity.

How?

Exercises: Find back the exercises in the table. Personalize (advance or delay) them in the exercise scheme based on the daily parameters. A maximum of 4 exercises / day is advised. With mobility, we strive for comfortable and pain-free exercises during the whole rehabilitation. Therefore, when returning exercises from previous phase and/or adding new exercises, a positive emotion should still be noticeable. Performing any of the exercises should not provoke pain, induce warmth and/or swelling, ting and/or bruise of the operated leg is still advised to give these parameters fully under control.

Activity: Guidance in activity is still provided by a daily, personalized step maximum.

- During a positive evolution, the patient is allowed to increase activity by 100-200 steps/day.
- During a negative evolution, the patient should remain on its activity level or have more rest (decrease in activity).
- Avoid pain activity during the whole rehabilitation (see table for reference).

The messaging function: is still commonly used to educate the patient, address concerns, answer questions and clarify therapeutic actions.

Gait pattern: Check on heel-to-toe gait pattern and on neutral foot positioning during video analysis. Check on compensation.

The use of crutches: Transition from 2 to 1 and to no crutches should be motivated when patient is safe.

Knee range of motion: is very important to see continuous improvement in the active ROM. A minimum of 90° active flexion should be achieved before the end of week 4. If not, contact the treating surgeon and provide him/her with this information. Full knee extension is stimulated daily by knee extension posture aided with instructions to also focus on full extension while walking and standing. During active or passive flexion mobilization, it is not advised to force and/or to go into the pain zone.

Besides daily step collection and visualization, during this phase, patient's evolution is tracked by postoperative PROM collection. Compare the pre- to postoperative results and give feedback to the patient. In most cases reassurance of time of progression is sufficient. At the end of this phase, the biggest gains have been made but naturally, patients should still experience improvements. Most patients perceive improvement over one year post-op.

23

Common actions Return to daily activity

Exercises (see table)	Messages	Check
Adapt exercise scheme: advance or delay phase 2 exercises	Questions	Reaction on adaptations in the exercise scheme wound
	Feedback	
Basic strength Proprioception Adaptation	<ul style="list-style-type: none"> Gradual increase in exercises and activity Paracetamol medication lowered to occasional use of NSAIDs type 3 medication Reassurance and supporting when there are ups and downs PROMs results Trend 2 → 1 Crutches Heading towards excellent result 	PROMs: prep vs week 6
Delay exercises which increases pain, warmth, swelling	Stimulate	Knee ROM (RT) before week 4 Gait pattern

Actions Return to daily activity

moveUP

24



Return to daily activity

- ✓ Minimum of pain
- ✓ Painkilling stopped or limited use

Return to active lifestyle

CRITERIA

- ✓ A minimum of / no pain
- ✓ A minimum feeling of warmth of the knee
- ✓ Painkilling medication stopped or limited to occasional intake of WHO type 1 medication
- ✓ No more NSAIDs
- ✓ for some /a/ or exercises surgeon approval e.g. driving, biking, lunges, etc.

Transition Return to daily activity -> Return to active lifestyle

25

Return to active lifestyle

in a 2nd postoperative phase.

- Focus on functional strengthening exercises.

Exercises: Find back the exercises in the table. Gradually add the exercises one by one, while closely following the patient's parameters in the graph or table. Do not forget to remove exercises from exercises phase. A maximum of 4 exercises / day is advised.

Before starting the functional strengthening exercises, pain, warmth and swelling are stable and at a low level. Medication is limited or stopped. During and after the functional strengthening exercises, these parameters should stay stable and/or at a minimum. There should be no need to increase medication intake.

During this phase, leg can still be swollen at the end of a busy day.

Activity: Stimulate a gradual increase in activity. Give recommendations about physical activity in general.

The messaging function: During this phase, reassurance of normal progression is often needed. Point out the fact that most patients reach a positive evolution up to 3 year after surgery.

Gait pattern: Check on heel-toe gait pattern and on neutral foot positioning during video analysis. Check on compensation.

Range of motion: While swelling and pain are reduced to a minimum, mobility should have evolved to 0 - 110° at around week 8-10. Consider that the achievable ROM greatly depends on the patient's preop mobility.

Finalize rehabilitation by checking if all expectations are met. You can check patient's satisfaction and outcome at the month 3 PROM collection.

26

Common actions Return to active lifestyle

Exercises	Messages	Check
Adapt exercise scheme: advance or delay certain exercises	<p>Questions</p> <ul style="list-style-type: none"> • How was the 3m course? • How far from your goals/expectations? • Any movements/activities which are still difficult/painful? 	Evaluate response on higher load
Functional strength	<p>Feedback</p> <ul style="list-style-type: none"> • Right timing 1 -> 0 crutches • PROMs 3 month <p>Stimulate</p> <ul style="list-style-type: none"> • Back to work • Full ACL update 	<p>Videos ROM</p> <p>Videos gait</p> <p>Patient satisfaction</p>
		<ul style="list-style-type: none"> • Functional ROM (110° + full extension) • Normal gait pattern

Actions Return to active lifestyle

27

APPROACH knee arthroplasty

Preop day -4 to day -1: Functional strength & mobility: Abduction, quadriceps setting, isometric gluteus, stationary bike, knee extension posture

Surgery - hospital stay: day 0 to day 5-6: Preop hospital (DCL) start
Elevation and circulation

Symptom control, proprioception and mobility: day 0 to day 21:
Ice circulation, Rest and Regular Low Load Movement
Flexion/extension, quadriceps setting

Basic strength and mobility: day 7 to day 60:
Functional ROM (standing and back-up of ACL)
Basic mobility: Heeling, single foot, extension posture
Basic strength: Side toe taps, mini squat...

Functional strength from day 35:
Low load to higher load
Functional strength Single leg, stairs, jumps, lunges, bridging

Timeline: Day 0 (Surgery) - Day 7-10 (Hospital stay) - Day 10-15 (Return to daily activity) - Day 35 (Return to active lifestyle)

28

Hip arthroplasty

Anterior protocol

- DAA
- ALA

Posterior protocol

29

	Preparation	Hospitalization Protection post-surgical hip	Return to daily activity	Return to active lifestyle
Estimated time	Day 0 - 21		Day 10-60	Day 35 and on...
Events	Surgery Hospital stay Discharge First week home Removal of staples		First postoperative consult	Second postoperative consult Back to work Driving car
Goal	Protection post-surgical hip Prevent complications Pain & symptoms control Proprioception		Gradually have no more pain <2/10 Gradual stop medication Basic strength Normalized gait pattern Prevent/minimize groin pain Basis balance / coordination	Full ROM 0 crutches Functional strength
Actions (See tables)	Educate & support patient Supervise use of crutches/walker Supervise medication intake Provide pain free exercises Control activity level Check wound Gait check		Educate & support patient Advise 2 crutches -> 1/0 crutches Adapt exercise scheme	
Criteria next phase	Pain & symptom control Protection & education		A minimum of / no pain Painkilling medication stopped or limited to occasional intake of WHO type 1 medication No more NSAIDs	For some /a/ or exercises surgeon approval is advised e.g. driving, biking, lunges

Phases hip

30



Preop	Hospitalization Protection post-surgical hip	Return to daily activity	Return to active Lifestyle
Ice - 4 x 15 min. Circulation - 4 x 2 min. Heel raises - 3x10 Isometric gluteus settings - 3 x 15 Isometric quadriceps settings - 3 x 10 Side toe taps - 3 x 15 Triple flexion lying - 3 x 10	Half rounds stationary bike - 3 x 5 min External rotation - 3 x 15 Hip abduction standing - 3 x 15 Hip flexion to hip extension - 3 x 10 Knee curls standing - 3 x 15 Hip extension - 3 x 10 Hip flexion and knee extension mobility exercise on the stairs - 3 x 15 Mini squat - 3 x 15 Balancing - Single leg stance - 3 x 10 Standing march - 3 x 2 min.	Stationary bike - 3 x 10 min Squat - 3 x 15 Single leg stance - knee lift - 3 x 15 Small lunge - 3 x 10 Forward lunge - 3 x 10 Side lunge - 3 x 10 Bridging - 3 x 10 Step up - 3 x 10	

Table hip APPROACH: See differences Anterior - Posterior

31

Differences in surgical approach

Anterior protocol	Posterior protocol
More prone for hip flexor irritation <ul style="list-style-type: none"> Avoid extension past 20 degrees Avoid external rotation past 50 degrees 	More prone for abductor weakness and/or irritation <ul style="list-style-type: none"> Avoid hip flexion past 90 degrees Avoid external rotation past neutral Avoid abduction past neutral
Anterior approach in daily exercise scheme: Less focus on hip flexor exercises in acute phase, e.g.: <ul style="list-style-type: none"> Triple flexion (with assistance) 3x8: d7-d11-d15-d19 ... Hip flexion to extension 3x8: d9-d13-d17-d21 ... 	Posterior approach in daily exercise scheme: More focus on hip flexor exercises, e.g.: <ul style="list-style-type: none"> Triple flexion (without assistance) 3x10: d5-d7-d9-d11 ... Hip flexion to extension 3x10: d6-d8-d10-d12 ...
More focus on abductor and glute exercises in acute phase, e.g.: <ul style="list-style-type: none"> Side toe taps (3x15): d5-d7-d9 Abduction (3x15): d11-d13-d15 ... Mini squat (3x10): d10-d12-d14 ... 	Less focus on abductor and glute exercises in acute phase, e.g.: <ul style="list-style-type: none"> Side toe taps (3x10): d7-d11-d13 Abduction (3x10): d15-d17-d19 ... Mini squat (3x10): d14-d16-d18 ...

32

Preparation	Hospitalization Protection post-surgical hip	Return to daily activity	Return to active Lifestyle
Activity 1. Preop avg < 4000 Max steps 1000/day 2. Preop avg 4000-7000 Max steps 1500/day 3. Preop avg > 7000 Max steps 2000/day From day 5: + 100 steps/day	Gradual increase until 50-100% preop steps Positive evolution: +100-200 steps / day Negative evolution: a. Max steps on hold b. Decrease in Max steps	Gradual increase until 75%- 100% or > preop steps Positive evolution: +200-200 steps / day Negative evolution: a. Max steps on hold b. Decrease in Max steps	
Medication PK WHO 1-2 NSAID 2/day - 2 weeks Blood thinning medication	PK WHO 1 NSAID 1/day - 2 weeks	Stop PK No NSAID	
Mobility Pain free mobility exercises without forcing	Pain free mobility exercises without forcing	Pain free mobility exercises without forcing	
Crutches 2 crutches	2 or 1 crutch	No crutches	
Stationary No stationary bike	Stationary bike no resistance	Stationary with resistance	
Biking Car	Driving (when feeling safe to do so) *Check the need for surgeon or car insurance approval	Biking after surgeon approval	

Table hip: Activity - Medication (*see addendum) - Crutches -

33

Preparation	Hospitalization Protection post-surgical hip	Return to daily activity	Return to active Lifestyle
Activity High activity > 50% preop Biking	High activity > 100% preop Biking	Jumping Running "high" loading activities or sports → Surgeon dependent	
Exercises Pain provoking exercises Exercises against heavy resistance Straight leg raises (SLR)	Pain provoking exercises Exercises inducing pain after execution	Pain provoking exercises Exercises inducing pain after execution	
ANTERIOR don't: Repetitive hip flexion	Excessive ranges: a. POSTERIOR don't: Hip flexion > 90° Endorotation > neutral Adduction > midline	a. Painful excessive ranges	
Mobility a. Excessive ranges b. POSTERIOR don't: Hip flexion > 90° Endorotation > neutral Adduction > midline	No use of crutches with dysfunctional (limping) gait pattern or for long distances	No use of crutches with dysfunctional (limping) gait pattern	
Crutches No use of crutches	Pain provoking movements or exercises	Pain provoking movements or exercises	
Pain Pain provoking movements or exercises	Sleeping on non-op and/or on side (surgeon dependent)	Sleeping on non-op and/or on side (surgeon dependent)	
Sleep Sleeping on non-op and/or on side (surgeon dependent)	No use of crutches with dysfunctional (limping) gait pattern	Pain provoking movements or exercises	
DO NOTs			

Table hip: Don'ts - Contra-indications

34

Hospitalization & Protection post surgical hip

In a first postoperative phase,

- Focus on protection of the post surgical hip and education of the patient on how to do so: the patient should take its prescribed medication, use crutches, put ice and elevate rest and elevation of the operated leg
- Symptom control: Aim to get pain, warmth and swelling under control
- Low load mobility: Low load and light mobility exercises are advised, motivating the patient to move his/her leg regularly without loading and regularly changing position
- Attention:
 - Hip anterior approach: See table don'ts
 - Hip posterior approach: See table don'ts

How?

Exercise: See the exercises in the table. Focus is on circulatory and light mobility exercises. Circulatory exercises help drain postoperative blood and swelling and prevent thrombosis. Warming and stretching exercises help maintain range of motion and muscle tone. These exercises are automatically active in the exercise scheme. Exercises against resistance or that provoke pain, swelling, stiffness and warmth are to be avoided, including during a knee-to-face exercise. Guidance in activity is provided by a diary, personalized step maximum (see table for reference).

The message is used to inform, motivate the patient, address concerns and advise questions.

The use of crutches/walkers should be supervised. Walking aids are important to prevent falling and lower the load on the hip as both require get time to heal and substitution is avoided.

6-8h patient: Check on heel-toe gait pattern and on shower floor posturing during video analysis.

What to do...

- A patient doesn't manage to get his pain under control? And the doctor can be flagged to give advice on medication. Most hip patients are advised to start with 1x10 step 3 painkillers. Step 3 painkillers can temporarily help high pain perceived patients.
- A patient reports a "red flag"? Always flag 2nd line doctors in cases of a bleeding wound, picture of wound, gastric problems, reports of fever, postop nausea, visit to general practitioner, medication advice.

1) Regarding the wound care: a picture of the wound after removal of stitches/staples should be asked. Most wounds are stapled. The staples are removed starting from 24 days to 21 days postoperatively. Small staples are removed a patch above the wound in most cases the patch should not be changed until the day (unless specifically noted to change). Some wounds are

35

Common actions Hospitalization & Protection post surgical hip

Exercises	Messages	Check
Default scheme	Questions How did surgery go? Message when home Removal staples How are the exercises going?	Therapy Compliance Picture hematoma Swelling Warmth
High activity > no more need for circulatory exercises	Feedback Pain is under control You have difficulties at night Warmth, stiffness, swelling Transfer from 2 → 1 crutch: safe and pain free?	Stimulate To move Adherence of step limit Adherence of use of crutches Long & elevation
		Use of crutches Medication intake Daily activity and step max

36





Protection post-surgical hip → Return to daily activity

- ✓ Protection & education
- ✓ Pain & symptom control

CRITERIA

- ✓ The patient is compliant to the postoperative rehabilitation guidelines.
- ✓ The patient found a balance between low load, regular, basic exercises and rest to decrease/control post-operative pain and swelling.
- ✓ Transition to the next phase is advised when the patient has control over its postoperative complaints.

Transition protection post-surgical hip -> Return to daily activity

37

Return to daily activity

In a second postoperative phase.

- Pursue a continuous decrease in symptoms and pain while reducing WHO type 1-2 medication to a minimum, need to increasing load by new exercises and allowing higher activity.
- New having a clear overview on the patient's evolution during the first postoperative period, support the patient finding the right balance between rest and activity.

• **Activities:**

- Hip anterior approach: See table dash's
- Hip posterior approach: See table dash's

How?

Exercise: Find back the exercises in the table. Personalize (advance or delay) them in the exercise scheme based on the daily parameters. With in case, we believe in controlled and passive exercises during the whole rehabilitation. Therefore, adding new exercises, a positive evolution should still be noticeable and performing the exercises should not provoke pain, induce warmth and/or swelling. Stop and evaluate the operation log is still advised to get these parameters fully under control.

• **Hip Anterior:** Focus on hip abductors, provoke anterior region less.

- Groin pain is often seen after hip arthroscopy, temporary, avoid exercises which provoke groin pain.

Activity: Guidance in activity is still provided by a daily, personalized step maximum.

- During a positive evolution, the patient is allowed to increase activity by 100-200 steps/day.
- During a negative evolution, the patient should reduce on its activity level or have more rest (decrease in activity).
- Avoid pain activity during the whole rehabilitation (see table for references).

The messaging function: is still commonly used to motivate the patient, solve concerns and answer questions.

The use of crutches: Transition from 2 to 1 and to no crutches should be motivated as pain free and safe.

Gait pattern: Check on heel-toe gait pattern and on neutral foot positioning during video analysis. Check on abductor strength or compensations.

Under daily data collection and visualization, during this phase, patient's evolution is tracked by postoperative PROM collection. Comparison pre to post can be made. Interpret the results and give feedback to the patient. In most cases reassurance of normal progression is sufficient. Of course there is still room for improvement there most patients still perceive evolution until 1 year postop.

38

Common actions Return to daily activity

Exercises	Messages	Check
<p>Adapt exercise scheme: Advance or delay phase 2 exercises</p> <p>Basic strength Posture Mobilization</p> <p>Delay exercises which increases pain, stiffness, swelling</p>	<p>Questions</p> <ul style="list-style-type: none"> • How was first postoperative consult? <p>Feedback</p> <ul style="list-style-type: none"> • Gradual increase in exercises and activity • Paracetamol medication (avoided) to occasional use of WHO type 1 medication before starting next level exercises • Reassurance and supporting when there are ups and downs • PROMs results • Timing 2-3 x 15 crutches • Heading towards excellent result <p>Stimulate</p> <ul style="list-style-type: none"> • Gradual build up basic strength and ADL without pain • Motivate RCM in ADL 	<p>Check</p> <ul style="list-style-type: none"> • Reaction on adaptations in the exercise scheme • Wound • Gradual increase in activity/load with decrease of symptoms and pain • PROMs: prep vs week 6 • Gait pattern

Actions Return to daily activity

39

Return to daily activity → Return to active lifestyle

- ✓ Minimum of pain
- ✓ Painkilling stopped or limited use

CRITERIA

- ✓ A minimum of / no pain
- ✓ A minimum feeling of stiffness & swelling of the hip
- ✓ Painkilling medication stopped or limited to occasional intake of WHO type 1 medication
- ✓ No more NSAIDs
- ✓ for some /a/ or exercises surgeon approval e.g. driving, biking, lunges, etc.

Transition Return to daily activity -> Return to active lifestyle

40

Return to active lifestyle

In a 3rd postoperative phase.

- focus on functional strengthening exercises.

Exercises: Gradually add the exercises one by one, while closely following the patient's parameters in the graph or table.

Before starting the functional strengthening exercises, pain, warmth and swelling are stable and at a minimum. Medication is limited or stopped. During and after the functional strengthening, these parameters should remain stable and/or at a minimum. There should be no more need for (extra) medication also.

During this phase, icing can still be useful at the end of a busy day.

Activity: Stimulate a gradual increase in activity.

The messaging function: During this phase, reassurance of normal progression is often needed. Point out the fact that most patients notice a positive evolution up till 1 year after surgery.

Gait pattern: Check on heel-toe gait pattern and on neutral foot positioning during video analysis. Check on abductor strength or compensations.

Finalize rehabilitation by checking if all expectations are met. You can check patient's satisfaction and outcome at the month 3 PROM collection.

41


Common actions Return to active lifestyle

Exercises	Messages	Check
<p>Functional strengthening</p>	<p>Questions</p> <ul style="list-style-type: none"> • How was first consult? • How far from your goal/expectations? • Any movements/activities which are still difficult/painful? <p>Feedback</p> <ul style="list-style-type: none"> • Right timing 1-2 x 0 crutches • PROMs 3 month <p>Stimulate</p> <ul style="list-style-type: none"> • Back to work • Full ADL uptake 	<p>Check</p> <ul style="list-style-type: none"> • Video gait • Patient satisfaction • Functional RCM • Normal gait pattern


Actions Return to active lifestyle

42







43



44



45



46





How to onboard a new patient on the moveUP platform?

- ✔ Onboarding is the process of registering and starting a new patient in the moveUP platform. After registering a patient in the platform, an invitation mail can be send. In this mail, the patient will receive all the needed information to install and activate the moveUP application.

1. Create a new patient in the onboarding functionality

- Surf to: <https://patient-onboarding-admin.moveup.care>.
- Fill out your moveUP login credentials and click “Login”.
- Click “New patient”.



The screenshot shows the 'Patient Onboarding Dashboard' for user 'Ellen'. At the top, there is a navigation bar with 'Dashboard' and 'Patients' on the left, and 'Support' and 'Account' on the right. Below the navigation bar, the title 'Patient Onboarding Dashboard' and a welcome message 'Welcome, Ellen' are displayed. A prominent orange warning banner at the top states: 'We're working with sensitive medical data. That's why you should enable two factor authentication. Two factor authentication secures your account better. Enable it in your account settings.' Below the banner, there are five data cards: a blue card for '227 moveUP candidates' with a '+ New patient' button (highlighted by a red arrow), a green card for '3706 patients active', an orange card for '13 patients reminded', a pink card for '0 candidates not interested', and a light green card for '0 filled in anamnesis questionnaire'.

Create new patient

2. Fill out the patient's details

- Fill out the patient's details.
- Click 'save'.

Edit patient

Basic Advanced

Managing Organization *
AZ Alma Eeklo x v
Organization that is the custodian of the patient record.

Name *
A name associated with the patient

Given names * Ellen <small>Given names (not always 'first'). Includes middle names.</small>	Family name * Tset <small>Family name (often called 'Surname')</small>
--	---

+ Add Name

Contact Point *
A contact detail for the individual

System * Email x v	Details * ellen+patient@moveup.care <small>The actual contact point details. Email address, phone number, etc.</small>
---	---

+ Add Contact Point

Edit patient details

[Back to list](#) [Save](#)


Save patient details

! Only one account can be coupled to an e-mail address. It is not possible to reuse an e-mail address that is already known within the moveUP database, e.g.: reusing the same e-mail address for the account of a partner or parent.

3. Add a new care plan to the patient

- ✓ A care plan is a crucial part of the digital care journey of your patient. A care plan determines what will be visible for your patient in the moveUP application, e.g.: type of exercises, information modules, message functionality, etc. A care plan usually consists out of an intervention, a moveUP treatment and (a) care team(s).


- The “Edit Care Plan” tab will automatically open after entering patient’s details.
- Click “Start new care plan”.



Create a new care plan

- Drag and drop an intervention to the care plan.
 - Fill out the intervention date (e.g.: operation date) and the name of the treating doctor (e.g.: surgeon).

Edit care plans for Ellen Test



Drag and drop to create a care plan

Edit care plans for Ellen Test

Patient's care plan

▼ Knee Right - total, TKA General Therapy

Intervention: Knee Right - total

Date 1

Doctor 2

Treatment: TKA General Therapy

Start Date

Code

HCP

Activated care plan (move care plan to active list in medical dashboard)

+ Start new care plan

3

AZ Alma Eeklo Interventions

Bilateral total hip arthroplasty Bilateral total knee arthroplasty Hip left - resurfacing

Hip left - revision Hip left - total Hip right - resurfacing Hip right - revision Hip right - total

Knee Right - total Knee left - revision Knee left - total Knee left - uni Knee right - revision

Knee right - uni

moveUP Treatments

THA Anterior Therapy THA Bilateral Therapy THA General Therapy TKA Bilateral Therapy

TKA General Therapy Unknown procedure date

AZ Alma Eeklo studies

Convention

Add type of intervention, date of intervention and doctor


- Drag and drop a moveUP treatment.
 - REMARK: The treatment's start date is set automatically to today's date if you don't fill it out. In most cases, it is not needed to set a date. This is only done in case the intervention or intervention date is not yet known.
- Drag and drop a study if applicable.
- Drag and drop care team(s) or care team participants.
- Click 'save'.

- ✓ Care team participants determine which healthcare providers will have access to the patient's account on the moveUP platform. The type of moveUP healthcare providers accounts determines which data this healthcare provider can access. Someone with a support role (e.g.: moveUP support staff) will not have access to the same data as the treating surgeon for example.

Care team participants are added automatically based on the type of intervention, however it is possible to manually add or remove a healthcare provider from the patient's care team if needed.


4. Send an invitation email to patients

- Once a care plan is saved, click on “Next”.



Click on 'Next'


- A template of the email(s) that will be sent to the patient is shown.
- Click on “Send selected email(s)” to send the emails to the patient.



Variable	Value
{{ firstname }}	Ellen
{{ lastname }}	Test
{{ code }}	b82w84z2
{{ registrationlink }}	moveup://registrationcode/b82w84z2
{{ organizationName }}	AZ Alma Eeklo
{{ surgeonName }}	profile.firstname 310z042tit profile.lastname 9yf9gpf98
{{ hosp_spec_par_start_NL }}	. Op vraag van uw ziekenhuis neemt u deel aan een nationale studie van het Kenniscentrum Gezondheidszorg (KCE). Deze studie is in kader van uw revalidatie die digitaal ondersteund zal worden met de moveUP applicatie
{{ organizationName }}	AZ Alma
{{ bodyStructure }}	heupprothese

Email template proposal

- A message “mails were sent” will appear for a few seconds so you know the mails are successfully sent.



Confirmation mail(s) sent

5. What happens next?

The patient will receive an **invitation email** with

- **Instructions how to install moveUP**– With a short introduction of moveUP, information link on pricing and instructions on how to install moveUP with login details, and a link to the quick start guide.

In case the moveUP treatment of your patient also uses an **activity tracker**, the patient will receive a second email with :

- **Link to order activity tracker** and confirm his/her subscription– With an order form for the activity tracker.





How to stop the treatment of a patient?

To stop the follow-up of your patient, **there are two mains steps** to follow.

Table of Contents




1. Stop active follow-up in care plan tab
2. Fill out the stop form
3. Optional: Adjust exercises for the patient if you want to personalize them
 - Personalize general exercise message for patient
 - Personalize the available exercises to a patient

1. Stop active follow-up in care plan tab

1. Go to the Care plan tab

2. Click on “Stop active follow up” in the top right of the Care plan tab.



3. The following pop up will appear.

Stop active follow up



Stopping the active follow up will move the patient to the PROM stage. Personalisation of exercises after stop active follow-up can be done in the tab 'care plan'. If exercises are not personalised. The patient will receive an overview with all exercises performed during his/her rehab.

Message

Yes

The message below will be send to the patient using the messaging function of the app.

B I U **Normal** **Sans Serif**

Dear Test,

You have already traveled a fine road in the first phase of your rehabilitation. Given the good progress in your recovery, the active follow up by your medical team will be stopped. Keep in mind that recovering from such an operation takes up to one year or more. The next weeks and months you will certainly continue to experience improvements.

Always take in mind the following important guidelines during next phases of your recovery:

Do not be afraid to try anything, but do not force it.
Always build up gradually and listen to your body.
Feel free to take a break or slow down temporarily if you notice that the joint is giving you trouble.

If you keep doing this, you will certainly continue to evolve positively.


Best regards,

HCP test

Email

Yes

The email below will be send to the patient.



Stop

Cancel

4. A default chat message will be set on “Yes” and below it you may read the chat message that will be sent to the patient when clicking “stop” at the bottom. You can also edit this message if you want to personalize it. If you don’t want to send a chat message you can toggle this option to “No”.



The message is designed for a patient that is doing well at the end of his rehabilitation. In case of opt-out, early end of follow-up, etc., this message really needs to be adapted or toggled off.

Stop active follow up

Stopping the active follow up will move the patient to the PROM stage. Personalisation of exercises after stop active follow-up can be done in the tab 'care plan'. If exercises are not personalised. The patient will receive an overview with all exercises performed during his/her rehab.

Message

The message below will be send to the patient using the messaging function of the app.

B I U ↺ ” ☰ ☷ Normal Sans Serif Ix 🔗 📎

Dear Test,

You have already traveled a fine road in the first phase of your rehabilitation. Given the good progress in your recovery, the active follow up by your medical team will be stopped. Keep in mind that recovering from such an operation takes up to one year or more. The next weeks and months you will certainly continue to experience improvements.

Always take in mind the following important guidelines during next phases of your recovery:

Do not be afraid to try anything, but do not force it.
 Always build up gradually and listen to your body.
 Feel free to take a break or slow down temporarily if you notice that the joint is giving you trouble.

If you keep doing this, you will certainly continue to evolve positively. |



Best regards,

HCP test

5. In addition, a default email with administrative info will be sent to the patient (e.g. how reimbursement and invoicing are arranged, how to send back the activity tracker etc.). If you don't want to send this e-mail, you can toggle it off.

Email

The email below will be send to the patient.




6. If you are ready to stop the active follow-up, you can click on "Stop".



✓ The patient will no longer be listed in the “Therapy” tab. You can now find him in the “PROM” tab or also in the “All” tab.

i In case of opt-out, early end of follow-up or if your patient does not wish to receive any more emails from moveUP, you must complete his follow-up. The patient will appear in the “Completed” tab.



2. Fill out the stop form

To fill out the form:

1. Go to the status tab
2. Click “start consultation”
3. Choose “Stop active therapy patient”

The screenshot displays the 'moveUP' patient management interface. The patient's name is 'Test PATIENT' and the status is 'TestPatient'. The interface includes a 'Status' tab (highlighted with a red box and a red circle '1'), a 'Start consultation' button (highlighted with a red box and a red circle '2'), and a 'Phases' section with three phases: Phase 3: Recovery at home: Protection post surgical knee, Phase 4: Return to daily activity, and Phase 5: Return to active lifestyle. A 'Start consultation' dialog box is open, showing a list of medical questionnaires to complete. The 'Stop active therapy patient' option is highlighted in blue (highlighted with a red box and a red circle '3'). The dialog box also includes a 'Cancel' button and a 'Start questionnaire' button.

moveUP Patients Onboarding

Test PATIENT
TestPatient

Gender
Age years

1 Status

2 Start consultation

3 Start consultation

Select the medical questionnaire to complete.

- ✓ Select a medical questionnaire
- AZMM Pixee
- Onboarding
- Phone call
- Physical therapy session KNEE
- Stop active therapy patient

Cancel Start questionnaire

Phases

Phase 3: Recovery at home: Protection post surgical knee
Starts from discharge: Main focus is to support patient at home AND to prevent and detect complications. The patients searches a balance between low regular, basic exercises and rest to decrease/control post-operative pain and swelling.

Phase 4: Return to daily activity
Support the patient finding the right balance between rest and activity. Pursue a continuous decrease in symptoms and pain while reducing WHO type medication to a minimum next to increasing load by new exercises and allowing higher activity.


Phase 5: Return to active lifestyle
Focus on functional strengthening exercises. Pain, warmth and swelling are stable and at a low level. Medication is limited or stopped. During and after functional strengthening exercises, these parameters should stay stable and/or at a minimum. There should be no need to increase medication intake.

3. Optional: Adjust exercises for the patient if you want to personalize them


Once the active follow-up has been stopped, we still want the patient to have access to the exercises executed during their rehabilitation. By default, the patient will see an overview of these exercises.

Personalize general exercise message for patient

If you want to personalize the message the patient sees in his app, you can do this by adding your message in the medical dashboard in the care plan tab.




The patient will see this message like this in the app:



Personalize the available exercises to a patient


To make sure the patient has the best experience possible, try to personalize which exercises the patient will see (by default it will be all exercises that were active during active follow-up). It is also possible to include a message related to the specific exercises.

To personalize the exercises click “personalize” in the care plan tab.



To change the available exercises you can:

- Add exercises: click on the + (1)
- Remove exercise: click on the x (2)
- Add personal message to an exercise: click edit (3)



In the pop-up menu that appears you can type a personal message that goes with that exercise:

Add exercise



Patient app preview

Ice



1. From a comfortable position, raise the knee and ankle above the hip
2. Place ice (not directly on the skin!)
3. Ice is important to reduce swelling, warmth and pain.

Available exercises



1008

Ice

Message

B **I** **U** Normal

Sans Serif

Type a message

Add

Cancel

The patient will see this message in the exercise overview:

Drop and dangle



1. Sit high enough so your feet are not touching the ground. Sit for a minute or so before starting the exercise.
2. Gently dangle your leg back and forth within your comfort zone.
3. This exercise should be relatively painless



Michiel De Wandelaar

This is a message to the patient!

