

Meril

Endo-Surgery

Oxidized Regenerated Cellulose

The **Gold Standard** in controlling small vessel bleeding

Introducing


Absorbable Hemostat

MERIZELLE™

Oxidized Regenerated Cellulose



Mechanism of Action



Features	Benefits
Rapid in Action	Achieves Hemostasis in 2-5 minutes
Anti-Adhesive	Does not interfere with instrument handling and tissue healing
Adaptable	Applicable for various tissue thicknesses and surgery sites
Absorbable	Within 1-2 weeks if left in place; may be removed prior to closure (at a surgeon's discretion)

Absorbable Hemostat

MERIZELLE™ available as

Oxidized Regenerated Cellulose

STANDARD

Fine Weave of ORC



Indications

Cardiovascular and Thoracic Surgery:

Anastomotic sites

General Surgery:

Cholecystectomy, Liver and Spleen Laceration, Port Site Bleeding

ENT Surgery:

Tonsillectomy

Plastic Surgery:

Skin Grafting

Indications

Neurosurgery:

Epidural Oozing and Capillary Oozing during Craniotomy, Laminectomy, Surgery on Tumor Bed, Head Trauma, Spine and Spinal Cord Surgeries

Cardiovascular and Thoracic Surgery:

Carotid Endarterectomy, Abdominal Aortic Aneurysm, Sternum Closure in CABG and Valve Repair/Replacement

Orthopaedic Surgery:

Femur Shaft Reduction and Fixation

FIBER

Tufts of Soft, Light Weight ORC



WOVEN

Denser Weave of ORC for Heavier Bleeding



Indications

Cardiovascular and Thoracic Surgery:

Cannulation sites, Anastomotic sites in CABG and Valve Repair/Replacement

General Surgery:

Liver Transplant, Resection and Repair, Splenic Injury or Splenorrhaphy



Absorbable Hemostat TM
MERIZELLE
 Oxidized Regenerated Cellulose

Product	Product Code	Dimensions	No. of Foils per box
STANDARD	ORC S214	2 in. x 14 in.	12
	ORC S48	4 in. x 8 in.	
	ORC S23	2 in. x 3 in.	
FIBER	ORC F12	1 in. x 2 in.	12
	ORC F24	2 in. x 4 in.	
	ORC F44	4 in. x 4 in.	
WOVEN	ORC W11	1 in. x 1 in.	12
	ORC W135	1 in. x 3.5 in.	
	ORC W34	3 in. x 4 in.	
	ORC W69	6 in. x 9 in.	



MERIL ENDO-SURGERY PVT. LTD.
Corporate Headquarters, Vapi
 E1-E3, Meril Park, Survey No. 135/2/B & 174/2, Muktanand Marg,
 Chala, Vapi- 396 191, Gujarat, India. Tel.: +91 260 3052 100
 Email: endosurgery@merillife.com, Web: www.merillife.com

India Sales & Marketing, Mumbai
 601-612, Midas, Sahar Plaza Complex, M.V. Road J.B Nagar,
 Andheri (East), Mumbai 400 059, India. Tel.: +91 22 307 32 400
 Email: endosurgery@merillife.com, Web: www.merillife.com

International Sales & Marketing
 301, A Wing, Business Square, Chakala,
 Andheri Kurla Road, Andheri East, Mumbai 400 093.
 Tel.: 022-39350700.