


MIRUS™ Ligating Clip is a sterile, medical grade titanium, non-absorbable, implantable device for general surgical procedures that require vessel or tissue ligation. **MIRUS™** Ligating Clip is intended to securely grip the vessel or tissue to provide complete occlusion, atraumatically.


MIRUS™ Ligating Clip are compatible with **MIRUS™** Ligating clip applicators, allowing the end user to ligate wide range of vessels and tissue structures.

Available as:

- Clip Size : Small, Medium, Medium-Large, Large
- Clips/Cartridges : 6
- Cartridges/Box : 20
- Clips/Box : 120
- Product Code : MLT100, MLT200, MLT300, MLT400



Meril Titanium Clips:



- Lateral & Transverse grooves allow secure fixation on the structure and increase resistance to dislodgement.
- Triangulated cross section maximises surface contact between clip & jaw.
- Precise tip to tip closure with firm vessel grip.

Ligation Clips	Titanium
Procedure Application	Blood vessel ligation for hemostasis

Technical Specifications

Product Code	Clip Size	Length (mm)	Height (mm)	Thickness (mm)
MLT-100	Small	6.80	1.00	0.55
MLT-200	Medium	9.50	2.25	0.80
MLT-300	Medium-Large	15.00	5.00	1.05
MLT-400	Large	17.00	8.50	1.25

Appliers for Open Surgery:



Technical Specifications

Product	Product code	Specification	Length	Applicator size	Clips compatibility
Applicator-Open	FTOA-10015	OPEN APPLICATOR 100-15CM	15 cm	Small	MLT100
Applicator-Open	FTOA-10020	OPEN APPLICATOR 100-20CM	20 cm	Small	MLT100
Applicator-Open	FTOA-10028	OPEN APPLICATOR 100-28CM	28 cm	Small	MLT100
Applicator-Open	FTOA-20015	OPEN APPLICATOR 200-15CM	15 cm	Medium	MLT200
Applicator-Open	FTOA-20020	OPEN APPLICATOR 200-20CM	20 cm	Medium	MLT200
Applicator-Open	FTOA-20028	OPEN APPLICATOR 200-28CM	28 cm	Medium	MLT200
Applicator-Open	FTOA-30015	OPEN APPLICATOR 300-15CM	15 cm	Medium-Large	MLT300
Applicator-Open	FTOA-30020	OPEN APPLICATOR 300-20CM	20 cm	Medium-Large	MLT300
Applicator-Open	FTOA-30028	OPEN APPLICATOR 300-28CM	28 cm	Medium-Large	MLT300
Applicator-Open	FTOA-40015	OPEN APPLICATOR 400-15CM	15 cm	Large	MLT400
Applicator-Open	FTOA-40020	OPEN APPLICATOR 400-20CM	20 cm	Large	MLT400
Applicator-Open	FTOA-40028	OPEN APPLICATOR 400-28CM	28 cm	Large	MLT400

Appliers for Laparoscopic Surgery:



Technical Specifications

Product	Product code	Specification	Length	Applicator size	Compatible Clips	Shaft Diameter	Compatible Trocar
Applicator-Endo	FTEA-00300	LAPAROSCOPIC APPLICATOR 300	34 cm	Medium-Large	MLT300	10 mm	10 mm & above
Applicator-Endo	FTEA-00400	LAPAROSCOPIC APPLICATOR 400	34 cm	Large	MLT400	10 mm	10 mm & above