



Catalogue 2020



BORT. Benefit at your side.

The successful
treatment that offers
a variety of solutions,
qualities and skills.



BORT select®

The first choice
for demanding
people.

BREEZE



BORT Sport:

The sporty side
of BORT.

COOLMAX®

COOLMAX® is a trademark of The Lycra Company.



BORT Generation:

For better mobility
for the aged.



BORT Med:

In the midst
of healthy life

FOR KIDS



BORT **select**® reflects the pursuit of excellence in quality at BORT, where special product features with added value meet the highest user requirements. Here, exclusive design and perfect functionality meet with therapeutic effectiveness, providing wearing comfort at its finest.



Stabilisation and relief for moving moments: the Sport product line combines specific functional requirements with modern, sporty design. Attractive equipment features, the use of high-tech fibre and an innovative fastener system make the supports an ideal accompaniment for all kinds of bodily activity.



Especially easy to handle, particularly comfortable to wear: our special offer for people with restricted freedom of movement or gripping ability is called Generation. The products help to retain independence in everyday life and can considerably improve the quality of living.



From young to old, from the head to the foot - the Med product line encompasses our diverse range of supports and braces. The broad delivery programme enables the care of virtually all joint-specific indications in accordance with the therapy.



BORT - At a glance.

OUR NEWS



Year by year we actively develop new product solutions that meet current treatment needs.

BORT Generation: FOR BETTER MOBILITY FOR THE AGED

Range extension with new Knee Support



Page 17, 133



Page 109



Page 121

OUR MODIFICATIONS




The most important adaptations as a brief overview.

You can find further modifications in the programme contents.



Page 87



Page 91



Page 152

All information without guarantee, error and changes reserved.
Status: January 2020.



OUR PROGRAMME

Product ranges	6
KIDS	23
Overview of indications	33
Neck	34
Shoulder	37
Overview of indications	47
Elbow	48
Overview of indications	59
Fingers	60
Thumb	62
Overview of indications	77
Hand	78
Overview of indications	97
Thoracic Spine	98
Back	99
Overview of indications	117
Torso	118
Groin/Hip	124
Overview of indications	129
Knee	130
Thigh	138
Overview of indications	151
Ankle	152
Foot	166
Shoes	168
Insoles	173
Reha Products	177
Overview OTC Products	181
EasyLife® products for rehabilitation	182
BORT Special Detergent	183
ActiveColor®	184
CLIMACare® warmers	194
PediSoft® foot care	196




The first choice for
demanding people.

select⁺

THE SELECT-BENEFIT

- ⊕ SOFTflex Plus knit fabric in feel-good quality with cooling luxury **BREEZE** fibre
- ⊕ State-of-the-art knitting technology for top wearing comfort and a perfect fit
- ⊕ Attractive product design with high-quality packaging
- ⊕ Made in Germany

**BORT select®
EpiPlus®**
No. 122 700



BREEZE

INDICATIONS

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of the elbow joint, **Epicondylitis** (acute and chronic), tendomyopathy

See Page 48

**BORT select®
ManuZip® volar**
No. 112 220



INDICATIONS

Chronic, post-traumatic or post-operative irritations in the wrist area, e. g. sprain, carpal arthrosis, tenosynovitis

See Page 78

**BORT select®
StabiloGen®**
No. 114 420



BREEZE

INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, gonarthrose, arthritis, e. g. ligament instability, meniscopathy

See Page 130

**BORT select®
GenuZip®**
No. 114 390



INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, gonarthrose, arthritis, e. g. ligament instability, meniscopathy

See Page 131

**BORT select®
Back Support with Pad**
No. 104 650



**BORT select® Lady
Back Support with Pad**
No. 104 640



INDICATIONS

Pain in the lumbar spine, recurrent low back pain, osteochondrosis, spondylarthrosis, spondylolysis without spondylolisthesis, mild lumbar deformity (degenerative spondylolisthesis), dorsalgia

See Page 100

**BORT select® Stabilo®
Back Support with Pad**
No. 104 710



**BORT select® Stabilo® Lady
Back Support with Pad**
No. 104 740




INDICATIONS

Treatment of increased pain in the lumbar spine area, severe low back pain, recurrent low back pain, osteochondrosis, spondylarthrosis, lumbosciatica, mild radicular syndrome, spondylolysis without spondylolisthesis, mild to moderate lumbar deformity (degenerative spondylolisthesis), after intervertebral disc surgery

See Page 101

BORT select® Lumbar Spine Brace for Mobilisation

No. 104 700




INDICATIONS

Moderate lumboischialgia with mild deficits with herniated discs/prolapse, moderate radicular, pseudo-radicular lumbar syndrome, grade I spondylolisthesis with low back pain, moderate lumbar deformity (recurrent treatment, complicated history), i.e. after intervertebral disc surgery with mild to moderate prolapse

See Page 99

BORT select® TaloStabil® Plus

No. 054 100



BREEZE

INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, arthrosis, ligament weakness, after removal of plaster cast, ligament injuries, distortion of upper ankle joint, supination prophylaxis

See Page 153

BORT select® TaloStabil®

No. 054 700



BREEZE


INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, arthrosis, ligament weakness, after removal of plaster cast, ligament injuries, distortion of upper ankle joint, tendomyopathies

See Page 152

BORT select® AchilloStabil® Plus

No. 053 100



BREEZE

INDICATIONS

Chronic or post-traumatic irritation or tenderness in the Achilles tendon area, Haglund syndrome heel irritation, achillodynia, irritation and pain of the Achilles tendon and surrounding areas, after surgery, e.g. Achilles tendon rupture, bursitis subachillea

See Page 152



The sporty side of BORT.


SPORT

THE SPORT-BENEFIT

- ⊕ SOFTflex 3D knitted fabric with a velvety soft surface structure
- ⊕ For more elasticity, stability and very comfortable to wear
- ⊕ High-tech fibre **COOLMAX®**
- ⊕ Innovative: Stabilo® Strap with easyClick system
- ⊕ Modern, sporty design with a reflecting label
- ⊕ Made in Germany

BORT EpiBasic Sport

No. 122 600 SP



COOLMAX®


INDICATIONS

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of the elbow joint, **Epicondylitis** (acute and chronic), tendomyopathy

See Page 48

BORT EpiContur® Sport

No. 022 500 SP



INDICATIONS

Epicondylitis, **chronic or recurrent epicondylopathy**, tendopathy, overstraining caused by work and sport, 'tennis/golf elbow'

See Page 53

BORT SellaFlex Sport Thumb Support

No. 112 030 SP



COOLMAX®

INDICATIONS

Chronic, post-traumatic or post-operative irritations in the area of saddle joint of the thumb, in the case of ulnar and radial ligament instability, **early stages of rhizarthrosis**

See Page 63

BORT ManuBasic® Sport

No. 112 920 SP



COOLMAX®

INDICATIONS

Chronic, post-traumatic or post-operative irritations in the wrist area, e. g. sprain, carpal arthrosis, tenosynovitis

See Page 79

**BORT StabiloBasic Sport
Back Support with Pad**
No. 104 690 SP




COOLMAX®

INDICATIONS

Pain in the lumbar spine, recurrent low back pain, osteochondrosis, spondylarthrosis, spondylolysis without spondylolisthesis, mild lumbar deformity (degenerative spondylolisthesis), ischialgia, dorsalgia

See Page 102

**BORT StabiloBasic Lady Sport
Back Support with Pad**
No. 104 680 SP



COOLMAX®

INDICATIONS

Pain in the lumbar spine, recurrent low back pain, osteochondrosis, spondylarthrosis, spondylolysis without spondylolisthesis, mild lumbar deformity (degenerative spondylolisthesis), ischialgia, dorsalgia

See Page 102

BORT StabiloGen® Sport
No. 114 520 SP




COOLMAX®

INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, e. g. gonarthrose, arthritis, ligament instability, meniscopathy

See Page 132

**BORT Patella
Tendon Strap Sport**
No. 114 010 SP



INDICATIONS

Patellar chondropathy, femoropatellar arthrosis, patellar tendonitis

See Page 139

BORT MyoActiv Sport Thigh Support

No. 116 600 SP



Sport

COOLMAX®

INDICATIONS

Preventive: For highly strained muscle groups, **Acute:** Swollen contusions with haematoma, muscle injury, **Chronic:** following muscle injury, protection against new muscle lesions, reduction of the burden on damaged muscle structures (e.g. scarring after old injuries)

See Page 138

BORT AchilloStabil® Plus Sport

No. 053 100 SP



Sport

COOLMAX®


INDICATIONS

Chronic or post-traumatic irritation or tenderness in the Achilles tendon area, in the case of Haglund syndrome heel irritation, achillodynia, after surgery, e. g. Achilles tendon rupture, bursitis subachillea

See Page 154

BORT TaloStabil® Sport

No. 054 600 SP



Sport

COOLMAX®

INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, sprain, ligament weakness, after removal of plaster cast, ligament injuries, distortion of upper ankle joint, arthrosis

See Page 154

COOLMAX® is a trademark of The Lycra Company

BORT TaloStabil® Plus Sport

No. 054 100 SP



Sport

COOLMAX®

INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, arthrosis, ligament weakness, after removal of plaster cast, ligament injuries, distortion of upper ankle joint, supination prophylaxis

See Page 155



For better mobility for the aged.

Generation



THE GENERATION-BENEFIT

- ⊕ Easy handling for more independence
- ⊕ Safe application due to clear marking
- ⊕ Support without limits
- ⊕ Comfortable to wear



BORT Generation Thumb-Wrist Brace

No. 215 200



Generation



INDICATIONS

Pre-operative, post-traumatic, post-operative, i. e. after thumb ligament reconstruction, immobilisation in the case of rheumatic-inflammatory processes, i. e. **rhizarthrosis**, in the case of degenerative diseases such as carpal tunnel syndrome

See Page 73

BORT Generation Thumb Splint

No. 215 300



Generation



INDICATIONS

Post-operative, post-traumatic, **rhizarthrosis**, **degenerative**, **inflammatory**, i. e. arthritis, ulnar/radial ligament injury of the thumb

See Page 68

BORT Generation Wrist Brace

No. 215 100



Generation



INDICATIONS

Pre-operative, post-operative, post-traumatic, **carpal arthrosis**, **inflammatoryrheumatic**, **degenerative diseases**, e. g. carpal tunnel syndrome (CTS), tenosynovitis, peripheral nerve paralysis

See Page 85

BORT Generation Back Support No. 215 850



INDICATIONS

Degenerative lumbar change, mild to moderate lumboischialgia, slipped disc with mild deficits, spondylosis, spondylolisthesis, spinal canal stenoses, radicular irritation syndrome, condition after inter-vertebral disc surgery

See Page 108

BORT Generation Back Brace No. 215 800



INDICATIONS

Facet syndrome, lumbosciatica, intervertebral disc protrusions/prolapse, radicular irritation syndrome, degenerative lumbar change, condition after inter-vertebral disc surgery, stable fractures, spondylosis, spondylolisthesis, spinal canal stenoses

See Page 108

BORT Generation Osteoporosis Brace No. 215 870



INDICATIONS

Active straightening for relief and correction of the lumbar/thoracic spine associated with restricted movement in the sagittal plane, stable, osteoporotic vertebral fractures, Osteoporosis, cautionary treatment of secondary kyphosis (stable fractures, tumors), muscular insufficiency, pre-operative and post-operative stabilisation

See Page 115

BORT Generation Knee Brace No. 215 400

NEW
DESIGN



INDICATIONS

Degenerative diseases, e. g. mild to moderate ligament instability of the knee joint, gonarthrosis, meniscus injuries

See Page 145

BORT Generation Knee Support No. 215 420

NEW



Generation



INDICATIONS

Gonarthrosis, recurrent articular effusion, arthritis, slight ligament instability

See Page 133

BORT Generation Knee Support Plus No. 114 490 G



Generation



INDICATIONS

Gonarthrosis, recurrent articular effusion, arthritis, slight ligament instability

See Page 133

BORT Generation Ankle Support No. 215 600



Generation



INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, sprain, ligament weakness, after removal of plaster cast, ligament injuries, arthrosis

See Page 156

BORT Generation Hallux Valgus Splint No. 215 930



Generation



INDICATIONS

Hallux valgus, after hallux surgery to ensure the success of treatment, post-operative capsule relief

See Page 173



In the midst
of healthy life.

Med

THE MED-BENEFIT

- ⊕ Full cover of indications
- ⊕ Very practical
- ⊕ Comfortable to wear
- ⊕ Easy to use

BORT Cervical Support^{Eco}

No. 127 460

**INDICATIONS**

Pain in the cervical spine area, minor whiplash injuries, sprain, post-operative

See Page 36

BORT StabiloNec^{®Eco}

No. 127 570

**INDICATIONS**

Severe painful conditions, minor whiplash injuries, cervical syndrome, post-operative

See Page 36

BORT Epi^{Basic}

No. 122 600

**INDICATIONS**

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of the elbow joint, **Epicondylitis** (acute and chronic), tendomyopathy

See Page 49

**BORT SOFT
Thumb Splint Long**

No. 112 720

**INDICATIONS**

Post-operative, post-traumatic, degenerative, inflammatory conditions, such as in case of arthritis, ulnar/radial ligament injury to the thumb, late stage after sinew reconstruction, **rhizarthrosis**

See Page 65

BORT
ManuBasic®
 No. 112 920



Med

BORT
ManuStabil® Short
 No. 112 730

NEW
DESIGN

Med

INDICATIONS

Chronic, post-traumatic or post-operative irritations in the wrist area, e. g. sprain, carpal arthrosis, tenosynovitis

See Page 80

INDICATIONS

Pre-operative, post-operative, post-traumatic, e. g. sprain degenerative, e. g. **carpal tunnel syndrome (CTS)**, tenosynovitis, peripheral nerve paralysis

See Page 86

BORT StabiloGen® Eco
 No. 114 520



Med

KNITTED
Jeans
 EDITION

INDICATIONS


Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, e. g. gonarthrose, arthritis, ligament instability, meniscopathy

See Page 134

BORT StabiloBasic
Back Support with Pad
 No. 104 690



BORT StabiloBasic Lady
Back Support with Pad
 No. 104 680



INDICATIONS

Pain in the lumbar spine, recurrent low back pain, osteochondrosis, spondylarthrosis, spondylolysis without spondylolisthesis, mild lumbar deformity (degenerative spondylolisthesis), dorsalgia

See Page 103

BORT TaloStabil® Eco
 No. 054 600



INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, sprain, ligament weakness, after removal of plaster cast, ligament injuries, arthrosis

See Page 157

BORT AchilloStabil® Eco

No. 054 900



Med

INDICATIONS

Chronic or post-traumatic irritation or tenderness in the Achilles tendon area, Haglund syndrome heel irritation, achillodynia, irritation and pain of the Achilles tendon and surrounding areas, after surgery, e. g. Achilles tendon rupture, bursitis subachillea

See Page 157

BORT MalleoStabil® Brace

No. 100 500



Med

INDICATIONS

Post-operative protection after ligament suture, conservative treatment of ankle distortions of any severity, functional treatment of ligament ruptures ankle fractures (Weber A), chronic ligament instability, shortening of time spent in plaster, ligament injuries, distortion of upper ankle joint, fibular ligament rupture, pronation and/or supination traumas

See Page 165



New leeway: BORT orthopaedic aids



THE KIDS-RANGE

- ⊕ **BORT supports and braces for kids**
Best protection, so that small joints can grow healthy.




PRODUCT-BENEFIT

- Anatomical cervical support for stabilisation of cervical spine
- Soft, open-pored foam
- With recess for larynx
- Adjustable Velcro fastening

COLOUR: ■ blue

THE KIDS-BENEFIT

- Can be worn day and night
- Particularly rounded edges
- Skin-friendly 100% cotton cover



PRODUCT-BENEFIT

- Open shoulder/arm adduction brace for immobilisation of the shoulder joint
- Stabilising strapping system ensures secure fit and immobilisation of the arm in the adduction position
- Soft, tension-resistant straps

COLOUR: ■ blue



THE KIDS-BENEFIT

- Flexible hand and arm adjustment
- Easy access to the surgical area to allow for wound care
- Skin-friendly

INDICATIONS

Pain in the cervical spine area, minor whiplash injuries, sprain, cervical syndrome, post-operative

INDICATIONS

Pre-operative, post-operative, post-traumatic, sprain, dislocation, subcapital humerus fractures, scapular fracture, rotator cuff injury

SIZES Circumference of neck



Children (x-small)_up to 33 cm

Adult sizes see page 34

Chin heights 5.0 cm / 6.0 cm
1 cm = 0,3937 inch

SIZE Length Olecranon to metacarpal knuckle point V




Children _____ 22 – 25 cm

Adult sizes see page 39

right and left symmetrical
1 cm = 0,3937 inch

BORT Clavicle Support

No. 121 700 KIDS



PRODUCT-BENEFIT

- Backpack support
- Soft padded straps
- Three Velcro fastenings, individually adjustable, ensure an ideal blockage

COLOUR: ■ blue

THE KIDS-BENEFIT


- Conforms ideally to the body's contours, especially in the critical armpit area
- Comfortable to wear
- Skin-friendly
- Can be worn day and night

INDICATIONS

Treatment of clavicle fractures

BORT StabiloFix®

No. 104 600 KIDS



PRODUCT-BENEFIT

- Corset for correcting the posture of the thoracic spine
- Soft, non-elastic strap material

COLOUR: ■ blue

THE KIDS-BENEFIT

- Fastening in the chest area makes it possible to put on/remove the item independently
- Conforms perfectly to the body contours, especially in the critical armpit area
- Fitted with tucks to enable straps to be shortened easily
- **Latex-free**

INDICATIONS

Straightening and relief of the thoracic spine, e. g. misalignment of thoracic spine

SIZE Underbust circumference



Children _____ up to 65 cm

Adult sizes see page 37

1 cm = 0,3937 inch


SIZE Underbust circumference



Children _____ 50 – 70 cm

Adult sizes see page 98

1 cm = 0,3937 inch




PRODUCT-BENEFIT

- Finger brace for immobilisation of interphalangeal joints
- Adjustable aluminium splint with soft padding

COLOUR:   grey

THE KIDS-BENEFIT

- Immobilisation of either one or two fingers
- Velcro fasteners, on metacarpals and fingers can be positioned as required



PRODUCT-BENEFIT

- Wrist brace with additional immobilisation of the metacarpophalangeal joint and saddle joint of the thumb
- Palm and dorsal splints to stabilise wrist and metacarpus
- Dorsal thumb splint for additional thumb stabilisation
- Length: approx. 17 cm (6.7 inches)

COLOUR:   grey

THE KIDS-BENEFIT

- Slip-on aid
- Comfortable to wear
- Stabilo® Strap

INDICATIONS


Post-operative, post-traumatic, i. e. redressement, positioning, immobilisation of the interphalangeal joint

INDICATIONS

Pre-operative, post-traumatic, post-operative, i. e. after thumb ligament reconstruction, immobilisation in the case of rheumatic-inflammatory processes

SIZE Hand length
 Children ____ 11 – 14 cm
Adult sizes see page 62

right and left symmetrical
1 cm = 0,3937 inch

SIZE Circumference of wrist
 Children ____ 13 – 15 cm
Adult sizes see page 74

please indicate left or right
1 cm = 0,3937 inch

BORT ManuStabil® short

No. 112 730 KIDS

NEW
DESIGN

PRODUCT-BENEFIT

- Rigid wrist brace with palm and dorsal reinforcement splints, individually shaped
- Short design
- Length: approx. 17 cm (6.7 inches)

COLOUR:  grey 

THE KIDS-BENEFIT

- Stabilo® Strap for effective stabilisation of the wrist
- Breathable material provides optimal comfort
- Can also be worn at night

INDICATIONS

Pre-operative, post-operative, post-traumatic, i. e. sprain, peripheral nerve paralysis, tenosynovitis

SIZE Circumference of wrist



Children — 13 – 15 cm

Adult sizes see page 86

please indicate left or right

1 cm = 0,3937 inch

Osgood-Schlatter Knee Support

No. 114 510 KIDS



PRODUCT-BENEFIT

- Knee support for local decompression of tibial tuberosity
- **SOFTflex knitted fabric** for more elasticity, stability and very comfortable to wear
- Anatomically contoured knit
- Multi-zone silicone pad
- One lateral spiral stay on each side
- Compression class II

COLOUR:  blue

THE KIDS-BENEFIT

- Velvety soft surface structure, creating very pleasant wearing comfort
- Highly elastic bending zone, also protects the popliteal fossa when worn for longer periods
- Adjustable setting using V-shaped strap correction

INDICATIONS

Facilitation of the movement of the patella, **Osgood-Schlatter disease**, tendinosis of the tibial tuberosity ligament, articular effusion, irritations

SIZE d-measure f-measure



Children 23 – 26 cm 28 – 32 cm

Adult sizes see page 137

right and left symmetrical

1 cm = 0,3937 inch




PRODUCT-BENEFIT

- Classical knee support for compression of soft tissues with pad for stabilisation and relief of the knee joint
- **SOFTflex knitted fabric** for more elasticity, stability and very comfortable to wear
- One lateral spiral stay on each side
- Silicone pad
- Highly elastic, soft flexion region in the hollow of the knee
- Compression class II

COLOUR: ■ blue

THE KIDS-BENEFIT

- Velvety soft surface structure
- Highly elastic bending zone, also protects the popliteal fossa when worn for longer periods



PRODUCT-BENEFIT

- Knee support with side joint splints for guidance of the knee joint
- Silicone pad
- Circular fastening straps on the thigh and lower leg
- Anatomically contoured knit
- Compression class II

COLOUR: ■ silver

THE KIDS-BENEFIT

- Extension limit at 0°
- Individually mouldable joint splints


INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, e. g. ligament instability, meniscopathy

INDICATIONS


Mild to moderate lateral ligament instability of the knee joint

SIZE d-measure f-measure

 Children 23 – 26 cm 28 – 32 cm
Adult sizes see page 134

right and left symmetrical
1 cm = 0,3937 inch

SIZE d-measure f-measure

 Children 23 – 26 cm 28 – 32 cm
Adult sizes see page 142

right and left symmetrical
1 cm = 0,3937 inch

**BORT Immobilisation Splint
with Patella Recess**

No. 145 000 KIDS



Med

PRODUCT-BENEFIT

- Knee brace for immobilisation of the knee joint with a 0° and/or 20° flexion
- 3-piece with open patella
- Median and lateral aluminium splints, which can be freely positioned
- Length: 38 cm (14.9 inches)

COLOUR: ■ black

THE KIDS-BENEFIT

- The support adapts to body shape thanks to a 3-piece design
- Lined with terry cloth on the inside for excellent comfort
- The aluminium splints can be removed for washing

INDICATIONS

Pre-operative, post-operative, post-traumatic, joint, muscle or ligament injuries

SPLINT LENGTH

Fits height

Children _____ 80 – 145 cm

Please indicate 0° or 20° flexion

Adult sizes see page 148

right and left symmetrical

1 cm = 0,3937 inch

**BORT Helix S Spiraldynamik®
Lower Ankle Support**

No. 053 200 KIDS



Med

PRODUCT-BENEFIT

- Ankle bandage to support torsion of the back of the foot and to raise the instep
- Thin, elastic straps
- 3-D contoured knit
- Two silicone pads
- Extra soft, snag proof material edge at the calf

COLOUR: ■ black

THE KIDS-BENEFIT

- Stabilisation of the ankle and activation of foot muscles through compression
- Functional supportive restraint in the case of classical foot deformities (i.e. splayfoot)
- Optimal comfort even during intense activity
- Separate slip-on aid included


INDICATIONS

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of lower ankle joint and the sole of the foot, ligament instability, foot problems due to static anomalies and defects in the shape of the foot due to dysfunction of the lower ankle joint e. g. skew and splayfoot

SIZE


Circumference of ankle

Children _____ 17 – 19 cm


 Adult sizes see page 159

please indicate left or right

1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Active ankle support for compression of soft tissues
- **SOFTflex knitted fabric** for more elasticity, stability and very comfortable to wear
- Two silicone pads for intermittent massage
- Compression class II

COLOUR: ■ blue

THE KIDS-BENEFIT

- Velvety soft surface structure
- Pressure relief at the base of the ankle



PRODUCT-BENEFIT

- Combined Air/Gel brace and foam pad to stabilise the ankle, while limiting pronation and supination

COLOURS: ■ black □ white

BORT MalleoStabil® SOFT Brace No. 100 520 KIDS

- With memory foam padding only

BORT MalleoStabil® Air/Gel Brace No. 100 540 KIDS

- With Air/Gel padding only

THE KIDS-BENEFIT

- No impairment when walking
- Easy to put on thanks to circular Velcro fastening

INDICATIONS

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of lower ankle joint and the sole of the foot, sprain, ligament weakness, after removal of plaster cast, ligament injuries

INDICATIONS

Post-operative protection after ligament suture, conservative treatment of ankle distortions of any severity, functional treatment of ligament ruptures, ligament injuries, supination/pronation trauma, fibular ligament rupture

SIZE Circumference of ankle

Children _____ 17 – 19 cm

Adult sizes see page 157

right and left symmetrical

1 cm = 0,3937 inch



SIZE Fits height

Size 1/Children _____ up to 155 cm

Adult sizes see page 165

right and left symmetrical

1 cm = 0,3937 inch

BORT TalofX®

No. 054 800 KIDS



PRODUCT-BENEFIT

- Soft ankle brace for stabilisation of the lower and upper ankle joints
- Bilateral shells padded with soft memory foam
- Cross restraint featuring Dynamic-Force-Control (DFC) technology
- Elastic material that carries away moisture

COLOUR: ■ black

THE KIDS-BENEFIT

- Slim design
- Can be worn without difficulty with sports and walking shoes
- Limitation of talus prolapse with optimal stabilisation
- For positioning the ankle in pronation or supination

INDICATIONS

Post-operative protection after ligament suture, conservative treatment of ankle distortions and ligament ruptures, immobilisation in the case of inflammatory conditions, rheumatism, arthritis, tendinitis, ligament injuries, supination / pronation trauma

BORT Heel Lift made of Silicone

No. 930 200 KIDS



PRODUCT-BENEFIT

- Silicone heel lift to compensate differences in leg-length
- Made of durable silicone
- Heights: 3 mm, 5 mm, 8 mm, 10 mm (0.1, 0.2, 0.3 and 0.4 inches)
- Pack: Single item

THE KIDS-BENEFIT

- Non-slip
- Fits any kind of shoe
- Washable
- Hygienic

INDICATIONS

Compensation of differences in leg-length, achilles tendon disorders (Achillodynia)

SIZE Circumference of ankle



Children _____ 17 – 19 cm

Adult sizes see page 161

please indicate left or right

1 cm = 0,3937 inch

SIZE Shoe size UK

Children _____ 1 – 2


Adult sizes see page 175

right and left symmetrical

1 cm = 0,3937 inch

Keep your posture:
Stabilisation and relief of
cervical spine and shoulder.





PRODUCT-BENEFIT

- Anatomical cervical support for stabilisation of the cervical spine
- **Stiff**, open-pored foam, for daytime wear
- With recess for larynx
- Adjustable Velcro fastening




COLOUR:  skin-tone

THE MED-BENEFIT

- Particularly rounded edges
- Skin-friendly 100% cotton cover

PRODUCT-BENEFIT

- Anatomical cervical support for stabilisation of cervical spine
- **Soft**, open-pored foam
- With recess for larynx
- Adjustable Velcro fastening
- **For children too**, see Page 24

COLOURS:  white (No. 127 260)
 skin-tone (No. 127 360)
 blue/black (No. 127 660)

THE MED-BENEFIT

- Can be worn day and night
- Particularly rounded edges
- Skin-friendly 100% cotton cover

INDICATIONS

Pain in the cervical spine area, minor whiplash injuries, sprain, cervical syndrome, post-operative

SIZES	Circumference of neck	SIZES	Circumference of neck
	small _____ up to 38 cm	Children (only 127 260, blue) – up to 33 cm	
	medium _____ 38 – 43 cm	small _____ up to 38 cm	
	large _____ 43 – 48 cm	medium _____ 38 – 43 cm	
		large _____ 43 – 48 cm	

Chin heights 5.5 cm / 7.5 cm / 9.0 cm / 11.0 cm **Children 5.0 cm / 6.0 cm**

1 cm = 0,3937 inch

BORT StabiloNec®
with Reinforcement Strap
 No. 127 560



PRODUCT-BENEFIT

- Anatomical cervical support for stabilisation of cervical spine
- **Integrated reinforcement splint**
- Stiff, open-pored foam
- With recess for larynx
- Adjustable Velcro fastening

THE MED-BENEFIT

- Optimal support
- Particularly rounded edges
- Skin-friendly 100% cotton cover

COLOUR: skin-tone

INDICATIONS

Severe painful conditions, minor whiplash injuries, cervical syndrome, torticollis (rheumatic), osteochondrosis/spondylarthrosis of the cervical spine, post-operative

SIZES

Circumference of neck



small	_____	up to 38 cm
medium	_____	38 – 43 cm
large	_____	43 – 48 cm

Chin heights

5.5 cm / 7.5 cm / 9.0 cm / 11.0 cm

1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Anatomical cervical support for stabilisation of cervical spine
- With recess for larynx
- **Soft**, open-pored foam
- Adjustable Velcro fastening
- Length: 56 cm (22.05 inches)

COLOURS: ■ skin-tone ■ blue/black

THE ECO-BENEFIT

- It can be worn day and night
- Particularly rounded edges
- Skin-friendly 100% cotton cover



PRODUCT-BENEFIT

- Anatomical cervical support for stabilisation of cervical spine
- **Integrated reinforcement splint**
- With recess for larynx
- Stiff open-pored foam
- Adjustable Velcro fastening
- Length: 56 cm (22.05 inches)

COLOURS: ■ skin-tone ■ blue/black

THE ECO-BENEFIT

- Optimal support
- Particularly rounded edges
- Skin-friendly 100% cotton cover

INDICATIONS

Pain in the cervical spine area, minor whiplash injuries, sprain, post-operative

INDICATIONS

Severe painful conditions, minor whiplash injuries, cervical syndrome, post-operative


SIZES	Chin height
Size 0	5.5 cm
Size 1	7.5 cm
Size 2	9.0 cm
Size 3	11.0 cm
One length only (56 cm/22.1 inches), to fit neck circumference from 32 to 45 cm (12.6 to 17.7 inches).	

1 cm = 0,3937 inch

SIZES	Chin height
Size 0	5.5 cm
Size 1	7.5 cm
Size 2	9.0 cm
Size 3	11.0 cm
One length only (56 cm/22.1 inches), to fit neck circumference from 32 to 45 cm (12.6 to 17.7 inches).	

BORT Clavicle Support

No. 121 700



Med

PRODUCT-BENEFIT

- Backpack support
- Soft padded straps
- Three Velcro fastenings, individually adjustable, ensure an ideal blockage
- For children too, see Page 25

COLOURS: ☐ white ☒ blue

THE MED-BENEFIT

- Conforms ideally to the body's contours, especially in the critical armpit area
- Skin-friendly
- Can be worn day and night

INDICATIONS

Treatment of clavicle fractures

BORT OmoFX

No. 121 100

NEW
DESIGN


Med

PRODUCT-BENEFIT

- Shoulder/arm abduction cushion for immobilisation and relief of the shoulder joint
- Abduction: **ca. 30°**
- Stabilisation straps for the shoulder and elbow
- Skin-friendly Velour cover

COLOUR: ☒ grey

NEW




THE MED-BENEFIT

- Easy to handle
- Comfortable arm rest

INDICATIONS

Pre-operative, post-operative, post-traumatic, for example in the case of periarticular irritations and disorders in the glenohumeral joint area, acromioclavicular joint dislocation (Tosny Classification Grades I + II), supraspinatus syndrome

SIZES Underbust circumference

	Children (only blue)	_____	up to 65 cm
	small	_____	65 – 80 cm
	medium	_____	80 – 95 cm
	large	_____	95 – 110 cm
	x-large	_____	110 – 125 cm


1 cm = 0,3937 inch

SIZES Fits height

Size 1	_____	up to 165 cm
Size 2	_____	over 165 cm

please indicate left or right

1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Shoulder/arm abduction brace for the immobilisation and relief of the shoulder joint
- Freely positionable straps for fixing the arm
- Ball-shaped hand rest for training the forearm muscles
- Skin-friendly Velour cover



THE MED-BENEFIT

- **Reversible cushion: positioning in 30° or 45° abduction**
- High therapeutic safety due to the anti-slip anatomic cut of the cushion
- Comfortable to wear

INDICATIONS


Pre-operative, post-operative, post-traumatic, for example in the case of impingement syndrome, rotator cuff injury, implantation of shoulder prosthesis

SIZES Circumference of torso

small	_____	up to 104 cm
large	_____	over 104 cm



right and left symmetrical
1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Inflatable shoulder/arm abduction cushion for immobilisation and relief of shoulder joint
- Abduction: ca. 50 - 60°
- Hand rest with ball for training forearm muscles
- Skin-friendly Velour cover



THE MED-BENEFIT

- 2 part air chamber system for prevention of undesired anteversion of the arm when lying down
- Arm is blocked with additional arm padding
- Extra lightweight design


INDICATIONS

Pre-operative, post-operative, post-traumatic, for example in the case of impingement syndrome, rotator cuff injury, implantation of shoulder prosthesis

SIZE

One size only

right and left symmetrical

BORT OmoBasic®acc. to Gilchrist
No. 121 350

PRODUCT-BENEFIT

- Open shoulder/arm adduction brace for immobilisation of the shoulder joint
- Stabilising strapping system ensures secure fit and immobilisation of the arm in the adduction position
- Soft, tension-resistant straps
- For children too, see Page 24

COLOURS: ■ blue □ white




THE BASIC-BENEFIT

- Flexible hand and arm adjustment
- Easy access to the surgical area to allow for wound care
- Skin-friendly

INDICATIONS

Pre-operative, post-operative, post-traumatic, sprain, dislocation, sub-capital humerus fractures, scapular fracture, rotator cuff injury

BORT OmoBasic®**Enclosed Style**
No. 121 500

PRODUCT-BENEFIT

- Closed shoulder/arm adduction brace for immobilisation of shoulder joint
- Highly elastic hose material
- Elastic support band at the base of the thumb
- Stabilising strapping system ensures secure fit and immobilisation of the arm in the adduction position

COLOURS: ■ blue □ white



THE BASIC-BENEFIT

- Stabilisation of the shoulder while allowing free movement of the thumb and fingers
- Elastic border in the shoulder area
- Skin-friendly

INDICATIONS


Pre-operative, post-traumatic, post-operative, sprain, dislocation, sub-capital humerus fractures, scapular fracture, rotator cuff injury

SIZES Length Olecranon to metacarpal knuckle point V

	Children (blue)	— 22 – 25 cm
	Size 0 (small)	— 25 – 28 cm
	Size 1 (medium)	— 28 – 31 cm
	Size 2 (large)	— 31 – 34 cm
	Size 3 (x-large)	— 34 – 38 cm

right and left symmetrical
1 cm = 0,3937 inch

SIZES Arm circumference Arm length

	0/S	— 26 – 32 cm	48 – 52 cm
	1/M	— 28 – 34 cm	50 – 56 cm
	2/L	— 32 – 40 cm	54 – 60 cm
	3/XL	— 38 – 48 cm	58 – 68 cm

right and left symmetrical
1 cm = 0,3937 inch



PRODUCT-BENEFIT


- Shoulder/arm adduction brace as a so-called **Desault's Bandage** for the immobilisation of the shoulder joint in a defined position
- Cap-shaped recess in the shoulder area
- Fastens around the torso area
- Perforated composite velour material

COLOURS:  blue  white

THE MED-BENEFIT



- Surface fixation of the shoulder and arm allowing free movement of the hand



PRODUCT-BENEFIT

- Shoulder/arm adduction brace **with pouch** for immobilisation of the shoulder joint
- Hand rest with ball for training forearm muscles

COLOUR:  grey 


THE MED-BENEFIT



- Easy and comfortable to use
- Waist belt for safe immobilisation

INDICATIONS

Pre-operative, post-traumatic, post-operative, contusion, dislocation, subcapital humerus fractures, injury to the rotator cuff

SIZES	Circumference of waist
 Size 0	up to 75 cm
Size 1	75 – 95 cm
Size 2	95 – 125 cm
Size 3	125 – 145 cm

right and left symmetrical
1 cm = 0,3937 inch

INDICATIONS


Pre-operative, post-traumatic, post-operative, sprain, contusion, subcapital humerus fractures, signs of paralysis

SIZES	Fits height
Size 1	up to 160 cm
Size 2	160 – 180 cm
Size 3	over 180 cm

right and left symmetrical
1 cm = 0,3937 inch

BORT OmoXpress

No. 120 900

NEW
DESIGN

Med

PRODUCT-BENEFIT

- Shoulder/arm adduction brace for immobilisation of shoulder joint
- Stable fixation of the shoulder and arm
- Open design for easy wound inspection
- Skin-friendly material

COLOUR: ■ grey

NEW



THE MED-BENEFIT

- Hand strap can be fully opened
- Easy to handle by rotating buckle, which allows an easy early-functional post-treatment
- Comfortable to wear

BORT OmoXpress light

No. 120 920

NEW
DESIGN

Med

PRODUCT-BENEFIT

- Shoulder/arm adduction brace for immobilisation of shoulder joint
- Stable fixation of the shoulder and arm
- Open design for easy wound inspection
- Skin-friendly material

COLOUR: ■ grey

NEW



THE MED-BENEFIT

- Very easy adaptation and safe handling due to the Velcro fastening
- Freely positionable arm slings
- Slight, visually hardly noticeable care, is comfortable to wear

INDICATIONS

Pre-operative, post-traumatic, post-operative, sprain, contusion, dislocation, humerus fractures

INDICATIONS

Pre-operative, post-traumatic, post-operative, sprain, contusion, dislocation, humerus fractures

SIZES Underbust circumference



Size 0 _____ up to 80 cm
Size 1 _____ 80 – 120 cm
Size 2 _____ 120 – 160 cm


right and left symmetrical
1 cm = 0,3937 inch

SIZES Underbust circumference




Size 0 _____ up to 80 cm
Size 1 _____ 80 – 120 cm
Size 2 _____ 120 – 160 cm

right and left symmetrical
1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Arm support cushion for immobilisation and relief of the shoulder joint with outward rotation at 15° or 30° and 15° abduction
- Strap system with Velcro fastening
- Hand rest with ball for training forearm muscles



COLOUR: ■ grey 

ACCESSORIES: No. 121 960
Wedge for positioning of outward rotation at 30°

THE MED-BENEFIT

- Safe stability and positioning through anatomical cushion
- Easy to handle
- Skin-friendly Velour cover

INDICATIONS

Pre-operative, post-traumatic, post-operative, for example in the case of impingement syndrome, rotator cuff rupture, implantation of shoulder prosthesis, subacromial defilee dilation (acc. to Neer), subacromial bursitis

SIZES	Fits height
Size 1	up to 165 cm
Size 2	over 165 cm

right and left symmetrical
1 cm = 0,3937 inch


BORT OmoARS with 3-Point-Strap

No. 121 970

NEW
DESIGN

PRODUCT-BENEFIT

- Arm support cushion for immobilisation and relief of shoulder joint in outward rotation at 15° or 30° and 15° abduction
- 3-point strap system with Velcro fastening
- Hand rest with ball for training the forearm muscles

COLOUR: ■ grey **NEW**

THE MED-BENEFIT

- 3-point strap system ensures optimal fit and stability
- Easy to handle
- Skin-friendly Velour cover

ACCESSORIES: No. 121 960

Wedge for positioning of outward rotation at 30°


BORT OmoARS short

No. 121 030

NEW
DESIGN

PRODUCT-BENEFIT

- Arm support cushion for immobilisation and relief of the shoulder joint in outward rotation at 15° and in abduction at 15°
- Strap system with Velcro fastening
- Hand rest with ball for training the forearm muscles

COLOUR: ■ grey **NEW**

THE MED-BENEFIT

- Safe stability and positioning through anatomical cushion
- Skin-friendly Velour cover

INDICATIONS


Pre-operative, post-traumatic, post-operative, for example in the case of impingement syndrome, rotator cuff rupture, implantation of shoulder prosthesis, subacromial defilee dilation (acc. to Neer), subacromial bursitis

SIZES	Fits height
Size 1 _____	up to 165 cm
Size 2 _____	over 165 cm

right and left symmetrical
1 cm = 0,3937 inch

SIZES	Length Olecranon to metacarpal knuckle point V
Size 1 _____	up to 34 cm
Size 2 _____	35 – 39 cm
Size 3 _____	over 39 cm

right and left symmetrical
1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Shoulder Joint compression support
- Additional Velcro straps to restrict outward rotation
- Silicone pad, which can be freely positioned
- Anatomical, muscle-activating contoured knit

COLOUR: ■ black

THE MED-BENEFIT

- Improved guidance of the glenohumeral joint
- Length of the straps can be adjusted individually
- Comfortable to wear even under clothing



PRODUCT-BENEFIT

- Arm traction brace
- Four weights of 250 g each for subacromial decompression
- Recess in the bending area
- Adjustable Velcro fastening

COLOUR: ■ grey



THE MED-BENEFIT


- Non-slip fit, anatomical shape
- Comfortable to wear
- Easy to handle

INDICATIONS


Chronic, post-traumatic or post-operative irritations of soft tissues in the shoulder area, arthrosis

INDICATIONS

Impingement syndrome, periarthropathia humero-scapularis (PHS)

SIZES	Circumference of upper arm	
	Size 1	24 – 28 cm
	Size 2	28 – 32 cm
	Size 3	32 – 36 cm
	Size 4	36 – 40 cm

right and left symmetrical
1 cm = 0,3937 inch

SIZES	Circumference of upper arm	
	Size 1	18 – 26 cm
	Size 2	27 – 35 cm

right and left symmetrical
1 cm = 0,3937 inch

BORT OmoTex Traction

No. 120 800



PRODUCT-BENEFIT

- **Modular shoulder joint support with separately useable weight band** for subacromial decompression
- **SOFTflex knitted fabric** for more elasticity, stability and very comfortable to wear
- Four 250 g weights for individual traction weighting
- Improved adaptability due to the flexibly adjustable belt guidance plate

THE MED-BENEFIT

- Improved joint guidance and relief of pressure in the humeral head region
- Weight band removable from the support and adjustable with Velcro tape in accordance with the course of therapy
- Comfortable to wear even under clothing

COLOUR: ■ black




INDICATIONS

Impingement syndrome, chronic, post-traumatic or post-operative irritations of soft tissues in the shoulder area, periarthropathia humero-scapularis (PHS)

SIZES

Circumference of upper arm



Size 1	24 – 28 cm
Size 2	28 – 32 cm
Size 3	32 – 36 cm
Size 4	36 – 40 cm

right and left symmetrical

1 cm = 0,3937 inch

Remain mobile: Relief and Stabilisation for the elbow



Overview of indications for selected products	select® EpiPlus® 122 700 P. 48	EpiBasic Sport 122 600 SP P. 48	EpiBasic 122 600 EpiBasic Short 122 600 K P. 49	Kubi Tal® 055 100 P. 50	EpiContur® Sport 022 500 SP P. 53	Stabilo® Epi Strap 122 580 P. 55	Elbow Support 104 200 P. 51	KubiFX 055 300 KubiFX light 055 350 P. 56	KubiFX long 055 400 KubiFX long light 055 450 P. 57	Elbow ROM Brace 055 500 P. 55
Page										
Post-traumatic			●					●	●	●
Post-operative							●	●	●	●
Dislocation										●
Sulcus-ulnaris syndrome				●				●		
Tendopathies						●			●	
Epicondylitis	●	●	●		●	●			●	
Irritations	●	●		●	●					●
Tendomyopathy			●							
Sprain							●			
Arthrosis							●			




PRODUCT-BENEFIT

- High quality elbow brace with additional Stabilo® Strap
- Two structured silicone pads for intermittent massage
- Lengthwise-elastic stretch zone in the elbow area ensures an ideal fit at all bending angles
- Compression class II
- Length: 20 cm (7.9 inches)

COLOURS:  skin-tone  silver
 black

THE SELECT-BENEFIT

- **SOFTflex Plus knit fabric:**
 - BREEZE: Luxury fibre with a cooling effect
 - Soft-velvety surface structure
- Extra-soft, snag proof material edge
- Very comfortable to wear, ergonomic and slip-free



PRODUCT-BENEFIT

- Elbow compression support with additional Stabilo® Strap
- **SOFTflex 3D knitted fabric** with a velvety soft surface structure for more elasticity, stability and very comfortable to wear
- Two structured silicone pads for intermittent massage
- Length: 20 cm (7.9 inches)

COLOUR:  black/green

THE SPORT-BENEFIT

- **COOLMAX®:** wicks moisture away from the body
- Easy fixing and length setting of the Stabilo® Strap using the **easyClick system**
- **Latex-free knitted fabric**
- The pressure on the pads can be individually adjusted through Stabilo® Strap

INDICATIONS

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of the elbow joint, **Epicondylitis** (acute and chronic), tendomyopathy

SIZES

Circumference of forearm 5 cm (2 inches) below elbow joint



xx-small	20 – 22 cm
x-small	22 – 24 cm
small	24 – 26 cm
medium	26 – 28 cm
large	28 – 30 cm
x-large	30 – 32 cm
xx-large	32 – 34 cm

right and left symmetrical
1 cm = 0,3937 inch

BORT EpiBasic

No. 122 600

**BORT EpiBasic Short**

No. 122 600 K



PRODUCT-BENEFIT

- Elbow compression support
- **SOFTflex knitted fabric** for more elasticity, stability and very comfortable to wear
- Two structured silicone pads for intermittent massage

COLOURS: skin-tone silver
 black

THE BASIC-BENEFIT

- Velvety soft surface structure
- Special elastic bending zone, also protects the antecubital fossa when worn for longer periods

PRODUCT-BENEFIT

- Elbow compression support
- **SOFTflex knitted fabric** for more elasticity, stability and very comfortable to wear
- Two structured silicone pads for intermittent massage

COLOURS: skin-tone silver
 black

THE BASIC-BENEFIT

- Velvety soft surface structure
- Special elastic bending zone, also protects the antecubital fossa when worn for longer periods
- **Reduced upper arm length**

*Length information only denotes approximate measures and may vary by size.

INDICATIONS

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of the elbow joint, **Epicondylitis** (acute and chronic), tendomyopathy

SIZES

Circumference of forearm 5 cm (2 inches) below elbow joint

	xx-small	20 – 22 cm
	x-small	22 – 24 cm
	small	24 – 26 cm
	medium	26 – 28 cm
	large	28 – 30 cm
	x-large	30 – 32 cm
	xx-large	32 – 34 cm

right and left symmetrical

1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Extensive padded elbow support
- **SOFTflex knitted fabric** for more elasticity, stability and very comfortable to wear
- Extensive silicone padding with recess at olecranon for dimensional stability during flexion
- Compression class II

THE MED-BENEFIT

- Velvety soft surface structure
- Extensive protection of entire elbow area

COLOURS:

skin-tone

blue

INDICATIONS

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of the elbow joint, for example in the case of compression of the nervus ulnaris at elbow (**sulcus ulnaris syndrome/cubital tunnel syndrome**), bursitis

SIZES		Circumference of forearm 5 cm (2 inches) below elbow joint
	x-small	22 – 24 cm
	small	24 – 26 cm
	medium	26 – 28 cm
	large	28 – 30 cm
	x-large	30 – 32 cm
	xx-large	32 – 34 cm
right and left symmetrical		
1 cm = 0,3937 inch		

BORT activemed Elbow Support No. 220 200



PRODUCT-BENEFIT

- Circular knitted elbow support
- Seamless, microfiber fabric
- Ulnar and radial silicone pads
- Length: 20 cm (7.9 inches)

COLOURS: skin-tone
 mineral grey

THE MED-BENEFIT

- Prevents pressure irritation in the inner elbow area thanks to an elastic comfort zone

BORT Elbow Support No. 104 200



PRODUCT-BENEFIT

- Elbow compression support
- Anatomically contoured knit for optimal fit
- Compression class II

COLOUR: skin-tone

THE MED-BENEFIT

- Comfortable to wear
- Not bulky

INDICATIONS

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of the elbow joint, **Epicondylitis** (acute and chronic), tendomyopathy

INDICATIONS

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of the elbow joint, articular effusion, sprain, arthrosis, contusion

SIZES Circumference of forearm
5 cm (2 inches) below elbow joint



x-small	_____	22 – 24 cm
small	_____	24 – 26 cm
medium	_____	26 – 28 cm
large	_____	28 – 30 cm
x-large	_____	30 – 32 cm

right and left symmetrical
1 cm = 0,3937 inch

SIZES Circumference of elbow



small	_____	24 – 26 cm
medium	_____	26 – 28 cm
large	_____	28 – 30 cm

right and left symmetrical
1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Elbow compression support with **lateral splints**
- Rounded knitted fabric, seamless
- Elastic in longitudinal and lateral directions
- Compression class II

COLOUR:  silver

THE MED-BENEFIT

- Comfortable to wear
- Lateral splints for extra support



PRODUCT-BENEFIT

- Circular knitted elbow compression support
- Elastic in longitudinal and lateral directions
- Compression class II


COLOUR:  skin-tone

THE MED-BENEFIT


- Comfortable to wear thanks to seamless knitted fabric

INDICATIONS

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of the elbow joint, articular effusion, sprain, arthrosis, contusion

SIZES	Circumference of elbow	
	small	24 – 26 cm
	medium	26 – 28 cm
	large	28 – 30 cm
	x-large	30 – 32 cm

right and left symmetrical
1 cm = 0,3937 inch

SIZES	Circumference of elbow	
	x-small	22 – 24 cm
	small	24 – 26 cm
	medium	26 – 28 cm
	large	28 – 30 cm
	x-large	30 – 32 cm

right and left symmetrical
1 cm = 0,3937 inch

BORT EpiContur® Sport

No. 022 500 SP



PRODUCT-BENEFIT

- Epicondylitis soft clasp for relieving the radial or ulnar muscle origin
- **Standard size adjustable to individual circumferences using Velcro fastening**

COLOUR: ■ black/green

THE SPORT-BENEFIT

- Fitted cut brace design enables slip-free fit for sport activities
- Brace pressure individually metered
- Specially shaped friction pad, freely positionable
- Modern, sporty design (with a reflecting label)

INDICATIONS

Epicondylitis, **chronic or recurrent epicondylopathy**, tendopathy, overstraining caused by work and sport, 'tennis/golf elbow'


SIZE

Circumference of forearm
5 cm (2 inches) below elbow joint





One size only _____ 21 – 35 cm

right and left symmetrical
1 cm = 0,3937 inch




PRODUCT-BENEFIT

- Epicondylitis soft clasp for relieving the radial or ulnar muscle origin
- **Standard size adjustable to individual circumferences using Velcro fastening**
- Freely positionable friction pad

COLOURS:  skin-tone  blue


THE MED-BENEFIT

- Brace pressure individually metered
- Fitted cut brace design enables slip-free fit



PRODUCT-BENEFIT

- Soft epicondylitis strap to relieve the muscle origins at the ulnar and/or radial epicondyle
- Ulnar and radial friction pad, which can be freely positioned

COLOURS:  skin-tone  blue


THE MED-BENEFIT

- Strap pressure is individually adjustable

INDICATIONS


Epicondylitis, **chronic or recurrent epicondylopathy**, tendopathies, overstraining, 'tennis/golf elbow'

Circumference of forearm	
SIZES	5 cm (2 inches) below elbow joint
One size only 21 – 35 cm	



right and left symmetrical
1 cm = 0,3937 inch


Circumference of forearm	
SIZES	5 cm (2 inches) below elbow joint
Size 0	21 – 25 cm
Size 1	25 – 30 cm
Size 2	30 – 35 cm



right and left symmetrical
1 cm = 0,3937 inch

BORT Stabilo® Epicondylitis Strap with Ulnar Relief

No. 122 580



PRODUCT-BENEFIT

- Epicondylitis strap to relieve the muscle origins at the ulnar and / or radial epicondyle
- Special anatomical design allows to individually meter strap pressure

COLOUR: ■ silver

THE MED-BENEFIT


- Elastic Velcro fastening, which results in no pressure on the elbow
- Simple design
- Padded inside

INDICATIONS

Epicondylitis, chronic or recurrent epicondylopathy, tendopathies, overstraining, tennis elbow, golf elbow

BORT Elbow ROM Brace

No. 055 500



PRODUCT-BENEFIT

- Elbow brace for mobilisation with simultaneous immobilisation of the proximal radio ulnar joint
- Limitation of extension: from 0° - 90°
- Limitation of flexion from 0° - 120°
- Fixation of the joint at 0°, 15°, 30°, 45°, 60°, 75°, 90°
- Avoidance of supination and pronation in the lower arm with a hand rest

COLOUR: ■ black

THE MED-BENEFIT

- Comfortable to wear due to extensive contact surfaces on the upper arm and lower arm
- Individual length setting
- Extra lightweight brace construction

INDICATIONS

Early functional postoperative or posttraumatic mobilisation with restricted movement of the elbow joint and immobilisation of the proximal radial ulnar joint, i.e conservative/post-operative care, i.e radial head fractures, elbow luxation, distal humerus fractures, after prosthesis implant

Circumference of forearm
SIZES 5 cm (2 inches) below elbow joint



Size 1	22 – 23 cm
Size 2	24 – 25 cm
Size 3	26 – 27 cm
Size 4	28 – 30 cm
Size 5	31 – 34 cm

right and left symmetrical
1 cm = 0,3937 inch

Circumference of upper arm



Size 1	25 – 34 cm
Size 2	34 – 43 cm

please indicate left or right
1 cm = 0,3937 inch




PRODUCT-BENEFIT

- Brace for immobilisation of the elbow joint
- Two removable aluminium splints
- Non-elastic material
- Inner side lined with terry cloth
- Length: ca. 39 cm (15.4 inches)

COLOUR: ■ black

THE MED-BENEFIT

- Secure support
- Individually mouldable from 0° to 90°
- Extra light and open frame construction



PRODUCT-BENEFIT

- Brace for immobilisation of elbow joint
- Two removable aluminium splints
- Non-elastic material
- Inner side lined with terry cloth
- Freely adjustable elbow pad
- Four reversible straps
- Length: ca. 37 cm (14.6 inches)

COLOUR: ■ black

THE MED-BENEFIT

- Reliable support
- Individually shaped from 0° to 90°

INDICATIONS

Pre-operative, post-operative, post-traumatic, for example in the case of proximal sulcus ulnaris syndrome, insertional tendinopathy, post-bursectomy conditions

SIZE

One size only

right and left symmetrical
1 cm = 0,3937 inch

SIZE

One size only

right and left symmetrical
1 cm = 0,3937 inch

BORT KubiFX Long Elbow Brace No. 055 400



PRODUCT-BENEFIT

- Brace for immobilisation of the elbow joint and the proximal radioulnar joint
- **Volar splint extending to palm area**
- Two removable aluminium splints
- Non-elastic material
- Inner side lined with terry cloth
- Freely positionable elbow pad
- Five reversible straps
- Length: ca. 53 cm (20.9 inches)

COLOUR: ■ black

THE MED-BENEFIT

- Reliable support
- Individually shaped from 0° to 90°

BORT KubiFX Long Light Elbow Brace No. 055 450



PRODUCT-BENEFIT

- Brace for immobilisation of the elbow joint and the proximal radioulnar joint
- Inner side lined with terry cloth
- Hand rest which can be used on both sides
- Freely adjustable elbow pad
- Standard size due to stepless setting of the lower arm length (olecranon to metacarpal knuckle point V) from 28 cm (11.0 inches) to 38 cm (15.0 inches)


COLOUR: ■ black

THE MED-BENEFIT

- Reliable support
- Individually shaped from 0° to 90°
- **Extra light and open frame construction**

INDICATIONS

Pre-operative, post-operative, post-traumatic, for example in the case of proximal sulcus ulnar syndrome, insertional tendinopathy, epicondylitis, post-bursectomy conditions

SIZES	Length Olecranon to metacarpal knuckle point V
	Size 1 _____ 28 – 31 cm
	Size 2 _____ 31 – 34 cm
	Size 3 _____ 34 – 38 cm

please indicate left or right
1 cm = 0,3937 inch

SIZE


One size only

right and left symmetrical
1 cm = 0,3937 inch

Fit for everyday commitments:
Pain reduction, more stability and
functionality for the fingers and thumb.



Overview of indications for the selected products	MobiDig® Flexion 112 560	MobiDig® Extension 112 510	DigiSoft® Finger Brace 112 770	SOFT Thumb Splint Long 112 720	SOFT Thumb Splint Plus 112 780	Generation Thumb Splint 215 300	Rhizo Flex® 105 100	Sella- Xpress Short 112 840	Sella- Xpress Plus 112 860	Sella- Xpress Long 112 880	SellaFlex Sport 112 030 SP	SellaTex® Classic 112 050 C SellaTex® 112 050	Generation Thumb-Wrist Brace 215 200
	P. 60	P. 60	P. 62	P. 65	P. 65	P. 68	P. 69	P. 66	P. 67	P. 67	P. 63	P. 74 / 75	P. 73
Immobilisation				●	●	●		●	●	●	●	●	●
				●	●	●	●	●	●	●	●	●	●
									●				
Flexion DIP/PIP joints	●												
Extension PIP joints		●											
Positioning													
Ulnar/radial lesion of ligaments			●			●					●	●	●
Rhizarthrosis				●	●	●	●	●	●	●	●	●	●
Tendopathy					●				●	●		●	●
Arthrosis						●	●					●	●
Conditions after fractures										●			
Additional wrist immobilisation												●	●




PRODUCT-BENEFIT

- Dynamic Quengel splint for mobilisation of interphalangeal joints
- Ergonomic, durable padding

COLOUR: ☐ white

THE MED-BENEFIT

- Linear spring force
- Tension can be increased or reduced by carefully bending the spring



PRODUCT-BENEFIT

- Dynamic Quengel splint for mobilisation of interphalangeal joints
- Proximal padding in the palm area
- Ergonomic, durable padding

COLOUR: ☒ grey 

THE MED-BENEFIT


- Linear spring force
- Tension can be increased or reduced by carefully bending the spring
- Plus sizes

INDICATIONS

Post-operative, post-traumatic flexion of PIP (L₁) or DIP (L₂) joints


INDICATIONS

Post-operative, post-traumatic extension of PIP joints

SIZES	Length
 small	3.5 – 5.0 cm
medium	5.0 – 5.8 cm
large	5.8 – 6.5 cm
x-large	over 6.5 cm

L₁ = Flexion PIP joint or
L₂ = Flexion DIP joint

right and left symmetrical
1 cm = 0,3937 inch


SIZES	Length ¹⁾	Finger width ²⁾
 x-small	5.0– 6.3	1.3–1.5
x-small plus	5.0– 6.3	1.5–1.7
small	6.3– 7.6	1.5–1.7
small plus	6.3– 7.6	1.7–1.9
medium	7.6– 8.9	1.7–1.9
medium plus	7.6– 8.9	1.9–2.1
large	8.9–10.2	1.9–2.1
large plus	8.9–10.2	2.1–2.3

¹⁾ Distal fold of palm to distal finger fold
²⁾ at PIP joint

right and left symmetrical
1 cm = 0,3937 inch

BORT Capener Splint Extension

No. 112 520



PRODUCT-BENEFIT

- Dynamic Quengel splint for mobilisation of interphalangeal joints
- Two lateral spiral springs

COLOUR: ☐ white

THE MED-BENEFIT

- Velcro strip for optimal fixing
- Soft padding
- Progressive spring force

INDICATIONS

Post-operative, post-traumatic extension of PIP joints

SIZES

Length¹⁾


	small	6.3 – 7.0 cm
	medium	7.0 – 7.6 cm
	large	7.6 – 8.0 cm

¹⁾ Length (L) = Distal fold of palm to distal finger fold

right and left symmetrical
1 cm = 0.3937 inch

BORT Stack Finger Splint

No. 470 100



PRODUCT-BENEFIT

- Plastic finger splint for hyperextension of the last damaged phalanx

COLOURS: ☐ skin-tone
☐ transparent

THE MED-BENEFIT


- Skin-tone colour: Can be thermo-plastically remoulded

INDICATIONS

Fracture, dislocation, ligament injuries, extensor tendon rupture

SIZES

Finger width
at DIP joint

	Size 1	1.4 – 1.6 cm
	Size 2	1.5 – 1.8 cm
	Size 3	1.7 – 2.1 cm
	Size 4	1.8 – 2.2 cm
	Size 5	1.9 – 2.3 cm
	Size 5.5	2.3 – 2.6 cm
	Size 6	2.0 – 2.8 cm
	Size 7	2.6 – 2.9 cm

right and left symmetrical
1 cm = 0.3937 inch



PRODUCT-BENEFIT

- Finger brace for immobilisation of interphalangeal joints
- Adjustable aluminium splint with soft padding
- Also for children, see Page 26


COLOUR:  grey 

THE MED-BENEFIT

- Immobilisation of either one or two fingers
- Velcro fasteners, on metacarpals and fingers can be positioned as required

INDICATIONS

Post-operative, post-traumatic, i. e. for redressement, positioning, Immobilisation of the interphalangeal joint

SIZES	Hand length
	Children — 11 – 14 cm
Size 1	13 – 16 cm
Size 2	16 – 19 cm
Size 3	19 – 22 cm

right and left symmetrical
1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Compression wrist support with thumb inlay to support the metacarpophalangeal joint and saddle joint of the thumb and to stabilise the wrist
- New, slim anatomically pre-moulded plastic stabilising element
- Anatomically contoured knit
- Compression class II


COLOURS:  skin-tone  black

THE MED-BENEFIT

- Stabilo® Strap for additional stabilisation of the wrist
- Special lightweight and slim structure for good wearing comfort


INDICATIONS

Chronic, post-traumatic, post-operative, irritations in the area of saddle joint of the thumb, i.e. late stage post-fracture, arthritis, ulnar and radial ligament lesions

SIZES	Circumference of wrist
	small — up to 15 cm
	medium — 15 – 17 cm
	large — 17 – 19 cm
	x-large — 19 – 21 cm

right and left symmetrical
1 cm = 0,3937 inch

BORT SellaFlex Sport Thumb Support No. 112 030 SP



Sport

COOLMAX

PRODUCT-BENEFIT

- Flexible wrist support to stabilise the metacarpophalangeal joint and saddle joint of the thumb using integrated reinforcing elements and additional Stabilo® Strap
- **SOFTflex 3D knitted fabric** with a velvety soft surface structure very comfortable to wear
- Anatomically contoured knit
- Supportive in load situations

COLOUR: ■ black/green

THE SPORT-BENEFIT

- **COOLMAX®**: wicks moisture away from the body
- Soft, snag proof material edge
- **Latex-free knitted fabric**
- Easy fixation using the **easyClick-System**

BORT SellaFlex Thumb Support No. 112 030



Med

PRODUCT-BENEFIT

- Flexible wrist support to stabilise the metacarpophalangeal joint and saddle joint of the thumb
- Anatomically contoured knit
- Additional stabilising elements at the thumb

COLOURS: ■ skin-tone ■ black

THE MED-BENEFIT

- Soft, snag proof material edge
- Flat design, not bulky
- Very comfortable to wear

INDICATIONS

Chronic, post-traumatic or post-operative irritations in the area of saddle joint of the thumb, in the case of ulnar and radial ligament instability, **early stages of rhizarthrosis**

SIZES

Circumference of wrist



x-small	___ up to 13 cm
small	___ 13 – 15 cm
medium	___ 15 – 17 cm
large	___ 17 – 19 cm
x-large	___ 19 – 21 cm

right and left symmetrical
1 cm = 0,3937 inch

SIZES

Circumference of wrist



x-small	___ up to 13 cm
small	___ 13 – 15 cm
medium	___ 15 – 17 cm
large	___ 17 – 19 cm
x-large	___ 19 – 21 cm

right and left symmetrical
1 cm = 0,3937 inch

BORT SOFT Thumb Splint Short
No. 112 710

NEW
DESIGN



PRODUCT-BENEFIT

- Soft brace for immobilisation of the metacarpophalangeal joint of the thumb
- Soft inner padding
- Integrated, anatomically pre-moulded stabilising element
- Length: approx. 11 cm (4.3 inches)

COLOURS:  skin-tone  grey  black

THE MED-BENEFIT

- Free movement of the fingers
- Breathable material provides optimal comfort

INDICATIONS

Post-operative, post-traumatic, degenerative, inflammatory conditions, such as in case of arthritis, ulnar/radial ligament injury to the thumb, late stage after sinew reconstruction, **rhizarthrosis**

SIZES

Circumference of wrist



x-small	_____	13 – 15 cm
small	_____	15 – 17 cm
medium	_____	17 – 19 cm
large	_____	19 – 21 cm

right and left symmetrical
1 cm = 0,3937 inch

BORT SOFT Thumb Splint Long

No. 112 720

NEW
DESIGN

Med

PRODUCT-BENEFIT

- Soft brace for immobilisation of the metacarpophalangeal joint and saddle joint of the thumb
- Soft inner padding
- Integrated, anatomically pre-moulded stabilising element extending to the wrist
- Length: approx. 14 cm (5.5 inches)

COLOURS: skin-tone grey black

NEW

THE MED-BENEFIT

- Free movement of the fingers
- Breathable material provides optimal comfort

INDICATIONS

Post-operative, post-traumatic, degenerative, inflammatory conditions, such as in case of arthritis, ulnar/radial ligament injury to the thumb, late stage after sinew reconstruction, **rhizarthrosis**

SIZES

Circumference of wrist



x-small	13 – 15 cm
small	15 – 17 cm
medium	17 – 19 cm
large	19 – 21 cm

right and left symmetrical
1 cm = 0,3937 inch

BORT SOFT Thumb Splint Plus

No. 112 780

NEW
DESIGN

Med

PRODUCT-BENEFIT

- Soft brace for immobilisation of the metacarpophalangeal joint of the thumb, saddle joint and distal interphalangeal joint of the thumb
- Soft inner padding
- Long integrated, anatomically pre-moulded stabilising element extending to the wrist area

COLOUR: grey

NEW

THE MED-BENEFIT

- Extended thumb element with Velcro fastening
- Free movement of the fingers
- Breathable material provides optimal comfort

INDICATIONS

Post-operative, post-traumatic, i.e. following tendon reconstruction, immobilisation in the case of degenerative processes such as de Quervain's Disease, **rhizarthrosis**

SIZES

Circumference of wrist



x-small	13 – 15 cm
small	15 – 17 cm
medium	17 – 19 cm
large	19 – 21 cm

right and left symmetrical
1 cm = 0,3937 inch

BORT

SellaXpress Short

No. 112 840

NEW DESIGN



Med

PRODUCT-BENEFIT

- Brace for immobilisation of the metacarpophalangeal joint and saddle joint of the thumb
- Integrated aluminium stabilisation element
- Extremely thin lightweight structure
- Latex-free

COLOUR:  grey 

THE MED-BENEFIT

- Customisable, cold-mouldable
- Extra soft padding provides optimal comfort
- Skin-friendly Velour

INDICATIONS

Post-traumatic, post-operative, i. e. following ligament lesions, sprain, immobilisation in the case of rheumatic-inflammatory processes i. e. **rhizarthrosis**

SIZES

Circumference of wrist



x-small	_____	13 – 15 cm
small	_____	15 – 17 cm
medium	_____	17 – 19 cm
large	_____	19 – 21 cm

please indicate left or right
1 cm = 0,3937 inch

BORT SellaXpress Long

No. 112 880

NEW
DESIGN

PRODUCT-BENEFIT

- Brace for partial immobilisation of the thumb saddle joint and metacarpophalangeal joint
- Integrated aluminium stabilisation element
- Extremely lightweight structure
- Soft padding
- Special guide on the thumb tube
- Latex-free

COLOUR: ■ grey 

THE MED-BENEFIT

- Long splint spanning the wrist joint
- Individual, cold-mouldable interlocking system
- Skin-friendly velour, comfortable to wear

INDICATIONS

Post-traumatic, post-operative, after thumb fractures, arthritis, arthrosis of the thumb saddle and/or metacarpophalangeal joints, ulnar/radial ligament injury to the thumb, late stage after sinew reconstruction, tendopathy, for controlled mobilisation in the case of rhizarthrosis, ski thumb

SIZES

Circumference of wrist



x-small	13 – 15 cm
small	15 – 17 cm
medium	17 – 19 cm
large	19 – 21 cm

please indicate left or right

1 cm = 0,3937 inch

BORT SellaXpress Plus

No. 112 860

NEW
DESIGN

PRODUCT-BENEFIT

- Brace for immobilisation of the metacarpophalangeal joint and saddle joint of the thumb plus the distal interphalangeal joint
- Integrated aluminium stabilisation element
- Extremely lightweight structure
- Soft padding
- Latex-free

COLOUR: ■ grey 

THE MED-BENEFIT

- Stabilisation including the distal interphalangeal joint of the thumb
- Customisable, cold-mouldable
- Extra soft padding provides optimal comfort
- Skin-friendly Velour

INDICATIONS

Post-traumatic, post-operative, i. e. following tendon reconstruction, immobilisation in the case of degenerative processes such as de Quervain's Disease, for inflammatory and rheumatic diseases, such as rhizarthrosis, tendopathy

SIZES

Circumference of wrist



x-small	13 – 15 cm
small	15 – 17 cm
medium	17 – 19 cm
large	19 – 21 cm

please indicate left or right

1 cm = 0,3937 inch

BORT Generation Thumb Splint
No. 215 300



PRODUCT-BENEFIT

- Soft brace for immobilisation of the metacarpophalangeal joint and saddle joint of the thumb
- Soft inner padding
- Removable, anatomically pre-moulded stabilising elements extending to the wrist
- Length: approx. 14 cm (5.5 inches)

COLOUR: ■ blue

THE GENERATION-BENEFIT

- Practical lacing enables easy, safe handling
- Variably adjustable thumb strap with a reversal loop for an extra safe hold
- Suitable for use together with a walking aid/wheeled walker
- Attractive, functional design
- Can also be worn at night

BORT Generation:
For better mobility for the aged
Product range with special handling advantages

INDICATIONS

Post-operative, post-traumatic, **rhizarthrosis, degenerative, inflammatory**, i. e. arthritis, ulnar/radial ligament injury of the thumb

SIZES

Circumference of wrist



small	_____	13 – 16 cm
medium	_____	16 – 19 cm
large	_____	19 – 22 cm

please indicate left or right
1 cm = 0,3937 inch

BORT RhizoFlex® Thumb Ring Brace

No. 105 100



Thumb

PRODUCT-BENEFIT

- Thumb ring brace for the stabilisation and relief of the thumb saddle joint
- Positioning the thumb in functional position
- Integrated, cold adjustable stabilisation element
- Optimal wearing comfort
- Easy to put on

COLOUR: ■ grey

THE MED-BENEFIT

- Extremely slim and functional brace design
- Free movement of the fingers and the wrist
- **Washable, highly durable**

INDICATIONS

Post-operative, post-traumatic or degenerative irritations in the area of saddle joint of the thumb, e.g. arthrosis of the CMC-I joint (thumb saddle joint), post-operatively following arthroplasty of the CMC-I joint, instability of the CMC-I joint

SIZES

Metacarpus circumference



Size 0	_____	15 – 17 cm
Size 1	_____	17 – 19 cm
Size 2	_____	19 – 22 cm
Size 3	_____	22 – 26 cm

please indicate left or right

1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Thermoplastic brace for immobilisation of the metacarpophalangeal joint of the thumb
- Padded straps that can be individually positioned
- Adjustable Velcro fastening
- Inside fully lined with textile flocking
- Extensive perforation

COLOUR:  skin-tone

THE MED-BENEFIT

- Free movement of the wrist (short design)
- Thermoplastically mouldable at approx. 60°C
- Washable



PRODUCT-BENEFIT

- Thermoplastic brace for immobilisation of the metacarpophalangeal joint and saddle joint of the thumb
- Padded strap that can be individually positioned
- Circumferential Velcro fastening
- Inside fully lined with textile flocking
- Perforated

COLOUR:  skin-tone

THE MED-BENEFIT


- Individually adjustable
- Thermoplastically mouldable at approx. 70°C
- Washable

INDICATIONS


Post-traumatic, post-operative, sprain, contusion

INDICATIONS

Post-traumatic, post-operative, i. e. sprain, contusion, ligament lesions (skiers thumb), immobilisation in the case of rheumatic-inflammatory processes (rhizarthrosis)

SIZES	Circumference of wrist	
	x-small	up to 15 cm
	small	15 – 17 cm
	medium	17 – 19 cm
	large	19 – 21 cm
	x-large	21 – 23 cm

please indicate left or right
1 cm = 0,3937 inch

SIZES	Circumference of wrist	
	x-small	up to 15 cm
	small	15 – 17 cm
	medium	17 – 19 cm
	large	19 – 21 cm
	x-large	21 – 23 cm

please indicate left or right
1 cm = 0,3937 inch

BORT SellaFix® N Black Edition

No. 112 040 B



Thumb

PRODUCT-BENEFIT

- Thermoplastic brace for immobilisation of the metacarpophalangeal joint and saddle joint of the thumb
- Innovative fastening system for optimum adaptation to the patient's circumstances
- Slim structure

COLOUR: ■ black

THE MED-BENEFIT

- Individually adjustable
- Thermoplastically adjustable at approx. 90°C
- Comfortable to wear due to the use of innovative materials
- Washable

INDICATIONS

Post-traumatic, post-operative, i.e. sprain, contusion, immobilisation in the case of rheumatic-inflammatory processes (i.e. rhizarthrosis), ligament lesions (skiers thumb)

SIZES

Metacarpus circumference



x-small	_____	18 – 20 cm
small	_____	20 – 22 cm
medium	_____	22 – 24 cm
large	_____	24 – 26 cm
x-large	_____	26 – 28 cm

please indicate left or right

1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Brace for immobilisation of the metacarpophalangeal joint and saddle joint of the thumb
- Special raised edge in the critical thumb forefinger area (commissure)
- Extensive perforation

COLOUR: ■ grey

THE MED-BENEFIT

- Soft textile flocking ensures optimal comfort
- Washable



PRODUCT-BENEFIT

- Brace for immobilisation of the metacarpophalangeal joint and saddle joint of the thumb
- Special raised edge in the critical thumb forefinger area (commissure)
- Elastic sleeve with downward traction


COLOUR: ■ anthracite

THE MED-BENEFIT


- Extremely durable brace
- Extensive perforation for comfortable wear
- Washable

INDICATIONS

Post-traumatic, post-operative, sprain, contusion, ligament lesions (skiers thumb), immobilisation in the case of rheumatic-inflammatory processes, rhizarthrosis

SIZES	Circumference of wrist	
	small	15 – 17 cm
	medium	17 – 19 cm
	large	19 – 21 cm

please indicate left or right
1 cm = 0,3937 inch

SIZES	Circumference of wrist	
	small	15 – 17 cm
	medium	17 – 19 cm
	large	19 – 21 cm

please indicate left or right
1 cm = 0,3937 inch

BORT Generation Thumb-Wrist Brace

No. 215 200



Generation



Thumb

PRODUCT-BENEFIT

- Wrist brace with additional immobilisation of the metacarpophalangeal joint and saddle joint of the thumb
- Soft inner padding
- Splints individually adjustable and removable

COLOUR: ■ blue

THE GENERATION-BENEFIT

- Practical lacing enables easy, safe handling
- The specially padded zone in the palm area enables use combined with walking aids/a wheeled walker
- Variably adjustable thumb wrap with a reversal loop for an extra safe hold
- Attractive, functional design




BORT Generation:
For better mobility for the aged
 Product range with special handling advantages

INDICATIONS

Pre-operative, post-traumatic, post-operative, i. e. after thumb ligament reconstruction, immobilisation in the case of rheumatic-inflammatory processes, i. e. **rhizarthrosis**, in the case of degenerative diseases such as **carpal tunnel syndrome**

SIZES

Circumference of wrist

	x-small	_____	up to 15 cm
	small	_____	15 – 17 cm
	medium	_____	17 – 19 cm
	large	_____	19 – 21 cm
	x-large	_____	21 – 23 cm

please indicate left or right

1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Wrist brace with additional immobilisation of the metacarpophalangeal joint and saddle joint of the thumb
- Palm and dorsal splints to stabilise wrist and metacarpus
- Dorsal thumb splint for additional thumb stabilisation
- Length: approx. 20 cm (7.9 inches)
- Also for children, see Page 26

THE MED-BENEFIT

- Stabilo® Strap for additional stabilisation of the wrist
- Extremely flat design
- Breathable material
- Comfortable to wear


COLOUR:

grey

NEW

INDICATIONS

Pre-operative, post-traumatic, post-operative, i.e. after thumb ligament reconstruction, immobilisation in the case of rheumatic-inflammatory processes, rhizarthrosis, in the case of degenerative diseases such as de Quervain's Disease, carpal tunnel syndrome

SIZES		Circumference of wrist
	Children	13 – 15 cm
	x-small	up to 15 cm
	small	15 – 17 cm
	medium	17 – 19 cm
	large	19 – 21 cm
	x-large	21 – 23 cm

please indicate left or right

1 cm = 0,3937 inch

BORT SellaTex® No. 112 050

NEW
DESIGN

PRODUCT-BENEFIT

- Wrist brace with additional immobilisation of the metacarpophalangeal joint and saddle joint of the thumb
- Palm and dorsal splints to stabilise wrist and metacarpus
- Can be fully opened
- Length: approx. 19 cm (7.5 inches)

COLOUR: ■ grey

THE MED-BENEFIT

- Adjustable thumb element, which can be individually positioned, completely removable
- Comfortable to wear
- Stabilo® Strap

INDICATIONS

Pre-operative, post-traumatic, post-operative, i.e. after thumb ligament reconstruction, immobilisation in the case of rheumatic-inflammatory processes, rhizarthrosis, in the case of degenerative diseases such as carpal tunnel syndrome

SIZES

Circumference of wrist



x-small	—	up to 15 cm
small	—	15 – 17 cm
medium	—	17 – 19 cm
large	—	19 – 21 cm
x-large	—	21 – 23 cm

please indicate left or right

1 cm = 0,3937 inch

BORT SellaTex® Plus No. 112 160

NEW
DESIGN

PRODUCT-BENEFIT

- Wrist brace for stabilisation of the metacarpophalangeal joint, saddle joint and distal interphalangeal joint of the thumb
- Palm and dorsal splints to stabilise wrist and metacarpus
- Can be fully opened
- Length approx. 19 cm (7.5 inches)

COLOUR: ■ grey

THE MED-BENEFIT

- Long, adjustable thumb element, which can be individually positioned, completely removable
- Comfortable to wear
- Stabilo® Strap

INDICATIONS

Pre-operative, post-traumatic, post-operative, i.e. after thumb ligament reconstruction, immobilisation in the case of rheumatic-inflammatory processes, rhizarthrosis, in the case of degenerative diseases such as de Quervain's Disease, carpal tunnel syndrome

SIZES

Circumference of wrist



x-small	—	up to 15 cm
small	—	15 – 17 cm
medium	—	17 – 19 cm
large	—	19 – 21 cm
x-large	—	21 – 23 cm

please indicate left or right

1 cm = 0,3937 inch

All in hand:
More functional stability
for the wrist.



Overview of indications for selected products	select® ManuZip® volar 112 220 P. 78	ManuBasic® Sport 112 920 SP P. 79	ManuBasic® 112 920 ManuBasic® Plus 112 940 P. 80/81	Wrist Support with Alu- minium Splint 103 300 P. 83	ManuStabil® Short 112 730 ManuStabil® 102 200 P. 86/87	SOFT Volar CTS Splint 112 750 P. 88	Generation Wrist Brace 215 100 P. 85	ManuCarpal® Combi 102 400 P. 91	Dorsal Intrinsic Plus Splint (D.I.P.S.) 102 500 P. 90	ManuDisc® 102 300 P. 89
Irritations	●	●	●	●						
Inflammatory rheumatic							●			
Degenerative							●			
Sprain	●	●	●	●		●				
Carpal arthrosis	●	●	●	●			●			●
Carpal tunnel syndrome (CTS)				●	●	●	●	●		●
Tenosynovitis	●	●	●	●	●		●	●		
Disc injury										●
(Post-operative) immobilisation					●	●		●		●
Peripheral nerve paralysis				●	●			●		
Metacarpal and metacarpophalangeal joint fractures									●	
Ligament injuries of the long fingers									●	



PRODUCT-BENEFIT

- High quality compression wrist support to stabilise the wrist
- Aluminium palm splint, anatomically premoulded
- Stabilo® Strap, freely positionable
- Viscoelastic silicone pad for intermittent massage
- Anatomically contoured knit

THE SELECT-BENEFIT

- **Zip extends** to the thumb making the product very comfortable for the patient to put on
- Extra-soft, snag proof material edge
- Integrated slip-on aid

COLOURS: ■ silver ■ black

INDICATIONS

Chronic, post-traumatic or post-operative irritations in the wrist area, e. g. sprain, carpal arthrosis, tenosynovitis

SIZES

Circumference of wrist



small	_____	13 – 15 cm
medium	_____	15 – 17 cm
large	_____	17 – 19 cm
x-large	_____	19 – 21 cm

please indicate left or right
1 cm = 0,3937 inch

BORT ManuBasic® Sport

No. 112 920 SP



Sport

COOLMAX®

PRODUCT-BENEFIT

- Compressing wrist support with additional Stabilo® Strap
- **SOFTflex 3D knitted fabric** with a velvety soft surface structure for more elasticity, stability and very comfortable to wear
- Anatomically contoured knit
- Anatomically shaped, internal palm splint in thermoplastic material
- Viscoelastic silicone pad

COLOUR: ■ black/green

THE SPORT-BENEFIT

- **COOLMAX®**: wicks moisture away from the body, for dryness and comfort
- Ideal fit for sports activities
- **Latex-free knitted fabric**
- Extra elastic in the wrist and forearm area, making it very easy to put on
- Easy fixing and length setting of the Stabilo® Strap using the **easyClick-System**
- Modern, sporty design (with a reflecting label)

Hand

INDICATIONS

Chronic, post-traumatic or post-operative irritations in the wrist area, e. g. sprain, carpal arthrosis, tenosynovitis

SIZES

Circumference of wrist



x-small	_____	up to 13 cm
small	_____	13 – 15 cm
medium	_____	15 – 17 cm
large	_____	17 – 19 cm
x-large	_____	19 – 21 cm

please indicate left or right

1 cm = 0,3937 inch

BORT ManuBasic®
No. 112 920



PRODUCT-BENEFIT

- Compression wrist support to stabilise the wrist
- **SOFTflex knitted fabric** for more elasticity, stability and very comfortable to wear
- Anatomically pre-moulded, thermoplastic palm splint
- Stabilo® Strap, freely positionable
- Silicone pad with hollow at the ulnar styloid process

COLOURS:  skin-tone  silver
 black


THE BASIC-BENEFIT

- Extra elastic in the wrist and forearm area, making it easier to put on
- Especially soft material edge in the forearm area
- Plus-size system results in a good fit

INDICATIONS

Chronic, post-traumatic or post-operative irritations in the wrist area, e. g. sprain, carpal arthrosis, tenosynovitis

SIZES

		c1-measure c2-measure		c1-measure c2-measure	
	x-small	up to 13	14 – 16 cm	large	17 – 19 20 – 22 cm
	small	13 – 15	16 – 18 cm	large plus	17 – 19 22 – 24 cm
	small plus	13 – 15	18 – 20 cm	x-large	19 – 21 22 – 24 cm
	medium	15 – 17	18 – 20 cm	x-large plus	19 – 21 24 – 26 cm
	medium plus	15 – 17	20 – 22 cm		

please indicate left or right
1 cm = 0,3937 inch

BORT ManuBasic® Plus Wrist Support with a Thumb Inlay No. 112 940



Med

PRODUCT-BENEFIT

- Compressing wrist support with additional stabilisation of the metacarpophalangeal joint and saddle joint of the thumb
- SOFTflex knitted fabric for more elasticity, stability and very comfortable to wear
- Anatomically pre-moulded, thermoplastic palm splint
- Silicone pad with hollow at the ulnar styloid process

COLOUR: ■ black

THE BASIC-BENEFIT

- Velvety soft surface structure
- Extra-soft, snag proof material edge
- Stabilo® Strap for additional stabilisation of the wrist

INDICATIONS

Chronic post-traumatic or post-operative irritation in the wrist and thumb saddle joint areas, sprain, carpal arthrosis, tenosynovitis, in the case of ulnar and radial ligament instability, **early stages of rhizarthrosis**

SIZES Circumference of wrist



small	_____	13 – 15 cm
medium	_____	15 – 17 cm
large	_____	17 – 19 cm
x-large	_____	19 – 21 cm

please indicate left or right
1 cm = 0,3937 inch

BORT ManuZip® Eco No. 112 620



Med



PRODUCT-BENEFIT

- Compression wrist support to stabilise the wrist
- Anatomically shaped palm splint
- Stabilo® Strap
- Silicone pad
- Anatomically contoured knit

COLOURS: ■ skin-tone ■ black

THE ECO-BENEFIT

- **Zip extends** to the thumb making the product very comfortable for the patient to put on

INDICATIONS

Chronic, post-traumatic or post-operative irritations in the wrist area, e. g. sprain, carpal arthrosis

SIZES Circumference of wrist



small	_____	13 – 15 cm
medium	_____	15 – 17 cm
large	_____	17 – 19 cm
x-large	_____	19 – 21 cm

please indicate left or right
1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Compression wrist support to stabilise the wrist
- Anatomically pre-moulded, thermoplastic palm splint
- Stabilo® Strap
- Silicone pad with hollow at the ulnar styloid process

COLOURS: skin-tone
 mineral grey

THE MED-BENEFIT

- Seamless expansion zone, perfectly surrounding the base of the thumb and optimally adjusting to anatomical form



PRODUCT-BENEFIT

- Compression wrist support to stabilise the wrist
- Anatomically pre-moulded, thermoplastic palm splint
- Stabilo® Strap, for additional stabilisation of the wrist
- Silicone pad with hollow at the ulnar styloid process

COLOUR: skin-tone

THE MED-BENEFIT


- Ganglion pad, freely positionable

INDICATIONS


Chronic, post-traumatic or post-operative irritations in the wrist area, e. g. sprain, carpal arthrosis, tenosynovitis

INDICATIONS

Chronic, post-traumatic or post-operative irritations in the wrist area, e. g. sprain, carpal arthrosis, intracapsular ganglion

SIZES	Circumference of wrist	
	small	13 – 15 cm
	medium	15 – 17 cm
	large	17 – 19 cm
	x-large	19 – 21 cm

please indicate left or right
1 cm = 0,3937 inch

SIZES	Circumference of wrist	
	small	13 – 15 cm
	medium	15 – 17 cm
	large	17 – 19 cm
	x-large	19 – 21 cm

please indicate left or right
1 cm = 0,3937 inch

BORT Wrist Support with Aluminium Splint (Medianus Splint)

No. 103 300



Med

PRODUCT-BENEFIT

- Classical wrist support
- Individually adjustable aluminium palm splint, in a durable pocket
- Length: 16 cm (6.3 inches)

COLOURS:  skin-tone  blue
 black

THE MED-BENEFIT

- Integrated slip-on aid and a Velcro strip on the back of the hand make it easy to apply/take off, even if gripping ability is restricted

Hand

BORT Wrist Support with Aluminium Splint and Additional Strap

No. 103 301

- Version as No. 103 300

- Additional stabilisation thanks to Stabilo® Strap

COLOURS:  skin-tone  blue

INDICATIONS

Chronic, post-traumatic or post-operative irritations in the wrist area, e. g. sprain, carpal arthrosis, carpal tunnel syndrome (CTS), rheumatism, tenosynovitis

SIZES

Circumference of wrist



x-small	_____	up to 15 cm
small	_____	15 – 17 cm
medium	_____	17 – 19 cm
large	_____	19 – 21 cm
x-large	_____	21 – 23 cm


please indicate left or right

1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Classical arm and wrist support to stabilise the wrist
- Individually adjustable aluminium palm splint in durable pocket
- Velcro fastening on back of hand
- Length: 22 cm (8.7 inches)

COLOURS:  skin-tone  blue

THE MED-BENEFIT

- Longer aluminium splint resulting in better stabilisation
- Easy to put on/remove, even if gripping ability is restricted, thanks to integrated slip-on aid and dorsal Velcro fastening



PRODUCT-BENEFIT

- Classical arm and wrist support to stabilise the wrist
- Individually adjustable aluminium palm splint in durable pocket
- Velcro fastening on back of hand
- Length: 22 cm (8.7 inches)


COLOURS:  skin-tone  blue

THE MED-BENEFIT


- **Additional stabilisation thanks to the Stabilo® Strap**
- Longer aluminium splint resulting in better stabilisation
- Easy to put on/remove, even if gripping ability is restricted, thanks to integrated slip-on aid and dorsal Velcro fastening

INDICATIONS

Chronic, post-traumatic or post-operative irritations in the wrist area, e. g. sprain, carpal arthrosis, carpal tunnel syndrome (CTS), tenosynovitis

SIZES	Circumference of wrist	
	x-small	up to 15 cm
	small	15 – 17 cm
	medium	17 – 19 cm
	large	19 – 21 cm
	x-large	21 – 23 cm

please indicate left or right
1 cm = 0,3937 inch

SIZES	Circumference of wrist	
	x-small	up to 15 cm
	small	15 – 17 cm
	medium	17 – 19 cm
	large	19 – 21 cm
	x-large	21 – 23 cm

please indicate left or right
1 cm = 0,3937 inch

BORT Generation Wrist Brace

No. 215 100



Generation



PRODUCT-BENEFIT

- Padded wrist brace with palm and dorsal reinforcement
- Splints individually adjustable and removable
- Breathable material provides optimum comfort

COLOUR: ■ blue

THE GENERATION-BENEFIT

- Can be put on with one hand due to innovative strapping fastening
- The specially padded zone in the palm area enables use combine with a walking aid / wheeled walker
- Attractive, functional design
- Can also be worn at night



BORT Generation:
For better mobility for the aged
Product range with special handling advantages

INDICATIONS

Pre-operative, post-operative, post-traumatic, **carpal arthrosis, inflammatory-rheumatic, degenerative diseases**, e. g. carpal tunnel syndrome (CTS), tenosynovitis, peripheral nerve paralysis

SIZES

Circumference of wrist



x-small	_____	up to 15 cm
small	_____	15 – 17 cm
medium	_____	17 – 19 cm
large	_____	19 – 21 cm
x-large	_____	21 – 23 cm

please indicate left or right

1 cm = 0,3937 inch

BORT ManuStabil® Short
No. 112 730

NEW
DESIGN



PRODUCT-BENEFIT

- Rigid wrist brace with palm and dorsal reinforcement splints, individually shaped
- Short design
- Length: approx. 18 cm (7.1 inches)
- Also for children, see page 27

THE MED-BENEFIT

- Stabilo® Strap for effective stabilisation of the wrist
- Breathable material provides optimal comfort
- Can also be worn at night

COLOURS:  skin-tone  grey  black



INDICATIONS

Pre-operative, post-operative, post-traumatic, e. g. sprain, degenerative, e. g. **carpal tunnel syndrome (CTS)**, tenosynovitis, peripheral nerve paralysis

SIZES

Circumference of wrist



Children (only grey)	_____	13 – 15 cm
x-small	_____	up to 15 cm
small	_____	15 – 17 cm
medium	_____	17 – 19 cm
large	_____	19 – 21 cm
x-large	_____	21 – 23 cm

please indicate left or right
1 cm = 0,3937 inch

BORT ManuStabil® No. 102 200

NEW
DESIGN

Med

PRODUCT-BENEFIT

- Rigid wrist brace with palm and dorsal reinforcement splints
- Breathable material
- Can be fully opened
- Can be used on the left or right
- Long design
- Length: approx 22 cm (8.7 inches)

COLOUR: ■ grey



THE MED-BENEFIT

- Individually mouldable
- Circular strap for effective stabilisation of the wrist
- Can also be worn at night

INDICATIONS

Pre-operative, post-operative, post-traumatic, degenerative, carpal tunnel syndrome (CTS), tenosynovitis, peripheral nerve paralysis

SIZES

Circumference of wrist



x-small	_____ up to 15 cm
small	_____ 15 – 17 cm
medium	_____ 17 – 19 cm
large	_____ 19 – 21 cm
x-large	_____ 21 – 23 cm

right and left symmetrical
1 cm = 0,3937 inch



BORT ManuStabilPro Open Style No. 182 500



Med

PRODUCT-BENEFIT

- Wrist brace with palm and dorsal reinforcing splint
- Breathable, skin-friendly material
- Extra-long design
- Length: approx 25 cm (9.8 inches)

COLOUR: ■ black

THE MED-BENEFIT

- Easy to shape individually
- Excellent stabilisation of the wrist and forearm
- Can be fully opened
- Easy to put on, even following surgery
- Can also be worn at night

INDICATIONS

Post-operative, post-traumatic, degenerative, e. g. carpal tunnel syndrome (CTS), tenosynovitis


SIZES

Circumference of wrist



x-small	_____ up to 16 cm
small	_____ 16 – 18 cm
medium	_____ 18 – 20 cm
large	_____ 20 – 22 cm
x-large	_____ 22 – 24 cm

please indicate left or right
1 cm = 0,3937 inch




PRODUCT-BENEFIT

- Soft palm rest splint for immobilisation of wrist
- Extra soft padded splint and fastening straps
- Latex-free

COLOUR: grey

THE MED-BENEFIT

- Easily mouldable
- Lightweight and comfortable to wear
- Can also be worn at night



PRODUCT-BENEFIT

- Soft dorsal rest splint for immobilising the wrist
- Palm splint
- Extra soft padded splint and fastening straps
- Three point Velcro fastening
- Latex-free

COLOUR: grey

THE MED-BENEFIT

- Maintained gripping ability
- Easily mouldable
- Lightweight and comfortable to wear
- Can also be worn at night

INDICATIONS

Post-operative, post-traumatic, e. g. carpal tunnel syndrome (CTS), sprain

SIZES	Circumference of wrist	
	x-small	13 – 15 cm
	small	15 – 17 cm
	medium	17 – 19 cm
	large	19 – 21 cm

please indicate left or right
1 cm = 0,3937 inch

SIZES	Circumference of wrist	
	small	15 – 17 cm
	medium	17 – 19 cm
	large	19 – 21 cm

please indicate left or right
1 cm = 0,3937 inch

BORT ManuCarpal® Brace No. 103 200



PRODUCT-BENEFIT

- Thermoplastically adjustable planthand brace
- Individual adjustment at 90°C
- Full padding
- Large contact surface with bilateral support of the forearm
- Extra lightweight design
- Washable

COLOUR: ■ black

THE MED-BENEFIT

- Easy for the patient to position and use
- Can be fully opened
- Comfortable to wear

INDICATIONS

Pre-operative, post-operative, post-traumatic, degenerative, e. g. **carpal tunnel syndrome (CTS)**, tenosynovitis, peripheral nerve paralysis

SIZES Circumference of wrist



x-small	—	up to 15 cm
small	—	15 – 17 cm
medium	—	17 – 19 cm
large	—	19 – 21 cm
x-large	—	21 – 23 cm

please indicate left or right

1 cm = 0,3937 inch

BORT ManuDisc® acc. to Dr. Schütz No. 102 300

NEW
DESIGN



PRODUCT-BENEFIT

- Wrist brace for multidirectional stabilisation of wrist
- Ulnar support and bilateral forearm rest
- Individually adjustable
- Soft padding on all inside surfaces

COLOUR: ■ grey

THE MED-BENEFIT

- Lightweight design with T-shaped palm splint
- Can be fully opened
- Comfortable to wear

INDICATIONS

Pre-operative, post-operative, post-traumatic, e. g. **carpal tunnel syndrome (CTS)**, disc injury, carpal arthrosis

SIZES Circumference of wrist



small	—	15 – 17 cm
medium	—	17 – 19 cm
large	—	19 – 21 cm
x-large	—	21 – 23 cm

please indicate left or right

1 cm = 0,3937 inch

BORT Dorsal Intrinsic Plus Splint (D.I.P.S.)
No. 102 500



PRODUCT-BENEFIT

- Wrist splint for immobilisation with finger fixation
- Ideal pre-forming of the splint in Intrinsic Plus position
 - 25° wrist extension
 - 70° metacarpophalangeal flexion
- Two elastic Velcro fasteners in the lower arm region, respectively a non-elastic Velcro fastening in the metacarpophalangeal joint area of the hand as well as the phalanx area
- Continuously softly padded on both sides, ideally comfortable to wear
- Easy post-operative use due to anatomical pre-forming

THE MED-BENEFIT

- Immobilising dorsal attachment
- Post-operative fast and safe attachment
- Avoidance of contractures of the long fingers
- Early mobilisation of the long fingers possible via fist closure
- Stretching the long fingers in the central and interphalangeal joints with Velcro for wearing at night-time

COLOUR: ■ black

INDICATIONS

Pre-operative, post-operative, post-traumatic, e. g. metacarpus and proximal finger phalanx in the area of the long fingers (e.g. a boxer fracture), ligament lesions in the long fingers in the proximity of the metacarpophalangeal joints, injuries and diseases which make the immobilisation of long fingers necessary (for example infections)

SIZES

Hand length



- | | | |
|--------|-------|------------|
| Size 1 | _____ | 14 – 17 cm |
| Size 2 | _____ | 17 – 22 cm |

BORT ManuCarpal® Combi

No. 102 400

NEW
DESIGN**PRODUCT-BENEFIT**

- Hand brace with removable finger rest
- Individually adjustable finger rest
- Internal palm padding
- Can be fully opened
- Can be used on the left or right

COLOUR: ■ grey **THE MED-BENEFIT**

- **Modular orthotic system**
- Additional relief in CTS area when using the finger rest
- Special cushioning of critical carpal tunnel area
- Circular strap for effective stabilisation of the wrist

Hand

INDICATIONS

Pre-operative, post-operative, post-traumatic, degenerative, e. g. **carpal tunnel syndrome (CTS)**, tenosynovitis, peripheral nerve paralysis

SIZES

Circumference of wrist



x-small	_____	up to 15 cm
small	_____	15 – 17 cm
medium	_____	17 – 19 cm
large	_____	19 – 21 cm
x-large	_____	21 – 23 cm




right and left symmetrical 

1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Palm splint for immobilising the wrist while fingers are blocked
- Long, one-piece aluminium splint, individually adjustable
- Integrated slip-on aid
- Breathable material provides optimal comfort


COLOURS:  skin-tone  grey 

THE MED-BENEFIT


- Easy to put on, even if gripping ability is restricted
- Can be worn day and night

INDICATIONS

Pre-operative, post-operative, post-traumatic, degenerative, e. g. carpal tunnel syndrome (CTS), tenosynovitis

SIZES	Circumference of wrist	
	x-small	up to 15 cm
	small	15 – 17 cm
	medium	17 – 19 cm
	large	19 – 21 cm
	x-large	21 – 23 cm

please indicate left or right
1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Soft palm splint for immobilising the wrist while fingers and thumb are blocked
- Wide finger rest and thumb surround, individually adjustable
- Four Velcro fastenings to block the splint securely
- Extra-soft padding material
- Latex-free


COLOUR:  grey 

THE MED-BENEFIT

- Simple immobilisation of the thumb
- Safe and comfortable at day and night

INDICATIONS

Pre-operative, post-operative, post-traumatic, e. g. carpal tunnel syndrome (CTS), paralysis, for inflammatory and rheumatic diseases

SIZES	Circumference of wrist	
	x-small	13 – 15 cm
	small	15 – 17 cm
	medium	17 – 19 cm
	large	19 – 21 cm

please indicate left or right
1 cm = 0,3937 inch

BORT Volar Resting Splint No. 103 100



PRODUCT-BENEFIT

- Palm rest splint for immobilising the wrist while fingers and thumb are blocked
- Adjustable stabilising element
- Soft foam padding

COLOUR: ■ anthracite

THE MED-BENEFIT

- Individually adjustable without using tools
- Extra lightweight structure
- Lower arm hand rest with skin-friendly textile cover

INDICATIONS

Pre-operative, post-operative, post-traumatic, degenerative, e. g. carpal tunnel syndrome (CTS), for paralysis, for inflammatory and rheumatic diseases, sprain

SIZES


Circumference of wrist



small	_____	15 – 17 cm
medium	_____	17 – 19 cm
large	_____	19 – 21 cm

please indicate left or right
1 cm = 0,3937 inch

BORT Paresis Splint No. 103 700



PRODUCT-BENEFIT

- Palm rest splint for immobilising the wrist while fingers and thumb are blocked
- Thermoplastic material
- Inside fully lined with textile flocking
- Washable

COLOUR: ■ skin-tone

THE MED-BENEFIT

- For positioning the hand in the intrinsic-plus position
- Easy, individual adjustment using a hair-dryer at 60°C
- Can be applied by the patient

INDICATIONS

Pre-operative, post-operative, post-traumatic, e. g. carpal tunnel syndrome (CTS), **hand paresis**, for inflammatory and rheumatic diseases

SIZES

Circumference of wrist



x-small	_____	up to 15 cm
small	_____	15 – 17 cm
medium	_____	17 – 19 cm
large	_____	19 – 21 cm
x-large	_____	21 – 23 cm

please indicate left or right
1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Compression wrist support for easy stabilisation of the wrist
- Optimised fit in the palm area
- Anatomically contoured knit
- Hemmed thumb opening

COLOURS:  skin-tone  blue

THE MED-BENEFIT

- Easy to use
- Can be worn on the left or right



PRODUCT-BENEFIT

- Compression wrist support for stabilisation of the wrist and metacarpus
- Accessory: Removable and freely positionable Stabilo® Strap


COLOURS:  skin-tone  blue

THE MED-BENEFIT


- Easy to use
- Can be worn on the left or right

INDICATIONS

Minor sprains, strains, after removal of plaster cast, tenosynovitis

SIZES	Circumference of wrist	
	x-small	up to 13 cm
	small	13 – 15 cm
	medium	15 – 17 cm
	large	17 – 19 cm
	x-large	19 – 21 cm

right and left symmetrical
1 cm = 0,3937 inch

SIZES	Circumference of wrist	
	x-small	up to 13 cm
	small	13 – 15 cm
	medium	15 – 17 cm
	large	17 – 19 cm
	x-large	19 – 21 cm

right and left symmetrical
1 cm = 0,3937 inch

BORT Wrist Support No. 112 110



PRODUCT-BENEFIT

- Elastic wrist support
- Velcro fastening
- Width: 8 cm (3.2 inches)

COLOURS:  skin-tone  black
 multicoloured

THE MED-BENEFIT




- Easy to put on

BORT Stabilo® Wrist Support No. 112 010



PRODUCT-BENEFIT

- Elastic wrist support with Stabilo® Strap
- Velcro fastening
- Width: 8 cm (3.2 inches)

COLOURS:  skin-tone  black
 white with blue stripe

THE MED-BENEFIT

- Additional stabilisation thanks to the Stabilo® Strap
- Easy to put on

USE

Gentle Wrist Support

SIZES Circumference of wrist



Size 1 _____ up to 17 cm
 Size 2 _____ 17 – 19 cm
 Size 3 _____ 19 – 22 cm

right and left symmetrical
 1 cm = 0,3937 inch

SIZES Circumference of wrist




Size 1 _____ up to 17 cm
 Size 2 _____ 17 – 19 cm
 Size 3 _____ 19 – 22 cm

right and left symmetrical
 1 cm = 0,3937 inch

For a strong back:
Correction, stabilisation, relief.



Overview of indications for selected products	Stabilo Fix® 104 600	select® Stabilo® 104 710 104 740	Stabilo- Basic Sport 104 690 SP 104 680 SP	Generation Back Support 215 850	Stabilo® Sacroiliac Brace 104 750	Vario- Plus special Width 112 690	Lumbo Xpress Back Support 112 600	select® Lumbar Spine Brace 104 700	Generation Back Brace 215 800	Stabilo® FA Back Brace 180 500	Stabilo® Lumbar Back Brace 180 100	Stabilo- Plus Back Support 104 730	Stabilo- Plus Bridging Brace 104 760	DorsoFX Back Brace 180 400	Generation Osteoporosis Brace 215 870
Page	P. 98	P. 101	P. 102	P. 108	P. 107	P. 105	P. 106	P. 99	P. 108	P. 109	P. 112	P. 110	P. 111	P. 113	P. 115
Thoracic spine posture correction	●														
ISG (sacroiliac joint) syndrome					●										
Lumbago		●	●	●		●									
Spondylose		●				●	●		●						
Lumbar spine syndrome		●				●	●								
Low back pain			●	●			●		●						
Ischialgia			●						●						
Facet joint arthrosis										●					
Prolaps (conservative/ post-operative)								●	●		●	●	●	●	
Osteoporosis											●				●
Vertebral lumbar/ thoracic spine fractures											●		●	●	
Post-operative														●	
Degenerative, lumbar change				●					●						




PRODUCT-BENEFIT

- Corset for correcting the posture of the thoracic spine
- Soft, non-elastic strap material
- For children too, see Page 25

COLOURS: ☒ blue ☐ white

THE MED-BENEFIT

- Fastening in the chest area makes it possible to put on/remove the item independently
- Moulds itself perfectly to the body contours, especially in the critical underarm area
- Fitted with tucks to enable straps to be shortened easily
- Latex-free



PRODUCT-BENEFIT

- Elastic corset for active straightening of the thoracic spine
- Breathable, skin-friendly material
- Slim design, not bulky

COLOUR: ☐ white

THE MED-BENEFIT

- Individual adjustment thanks to flexible elastic belts
- Active straightening of the spine

INDICATIONS

Straightening and relief of the thoracic spine, e. g. misalignment of thoracic spine


INDICATIONS

Straightening and relief of the thoracic spine, e. g. misalignment of thoracic spine, back pain (pain in the region of the thoracic spine)

SIZES Underbust circumference



SIZES Underbust circumference



BORT select® Lumbar Spine Brace for Mobilisation
No. 104 700



PRODUCT-BENEFIT

- Back brace for functional mobilisation of the lumbar spine
- Two anatomically pre-moulded splints in the dorsal area
- Easy to dose support and relief effect via inelastic belt tensioning
- Height: ventral 17 cm (6.7 inches), dorsal 26 cm (10.2 inches)

THE SELECT-BENEFIT

- **Removable restraint** for adaptation to the course of therapy
- Easy to put on using the hand straps provided


COLOURS:  skin-tone  silver

INDICATIONS

Moderate lumboischialgia with mild deficits with herniated discs/prolapse, moderate radicular, pseudo-radicular lumbar syndrome, grade I spondylolisthesis with low back pain, moderate lumbar deformity (recurrent treatment, complicated history), i. e. after inter-vertebral disc surgery with mild to moderate prolapse

SIZES

Circumference of waist

	Size 1	_____	70 – 80 cm
	Size 2	_____	80 – 90 cm
	Size 3	_____	90 – 100 cm
	Size 4	_____	100 – 110 cm
	Size 5	_____	110 – 120 cm
	Size 6	_____	120 – 130 cm
	Size 7	_____	130 – 140 cm

1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Back support for muscular stabilisation of the lumbar spine
- Active support with anatomically-contoured knit
- Four supporting splints
- High-quality viscoelastic sacral pad, which can be freely positioned
- Height: ventral 17 cm (6.7 inches), dorsal 23 cm (9.1 inches)

PRODUCT-BENEFIT

- Back support for muscular stabilisation of the lumbar spine
- Active support with anatomically-contoured knit
- Fitted at the waist: specific for women
- Four supporting splints
- High-quality viscoelastic sacral pad, which can be freely positioned
- Height: ventral 17 cm (6.7 inches), dorsal 24 cm (9.4 inches)

COLOURS:  silver
 black



COLOURS:  silver
 black



THE SELECT-BENEFIT


- Hand straps allow support to be fitted easily and without effort
- Comfortable to wear

THE SELECT-BENEFIT


- Hand straps allow support to be fitted easily and without effort
- Comfortable to wear

INDICATIONS

Pain in the lumbar spine, recurrent low back pain, osteochondrosis, spondylarthrosis, spondylolysis without spondylolisthesis, mild lumbar deformity (degenerative spondylolisthesis), dorsalgia

SIZES	Circumference of waist
	Size 1 _____ 70 – 80 cm Size 2 _____ 80 – 90 cm Size 3 _____ 90 – 100 cm Size 4 _____ 100 – 110 cm Size 5 _____ 110 – 120 cm Size 6 _____ 120 – 130 cm Size 7 _____ 130 – 140 cm

1 cm = 0,3937 inch

SIZES	Circumference of waist
	Size 0 _____ 60 – 70 cm Size 1 _____ 70 – 80 cm Size 2 _____ 80 – 90 cm Size 3 _____ 90 – 100 cm Size 4 _____ 100 – 110 cm Size 5 _____ 110 – 120 cm Size 6 _____ 120 – 130 cm

1 cm = 0,3937 inch


BORT select® Stabilo® Back Support with Pad No. 104 710



PRODUCT-BENEFIT

- Back brace to stabilise the lumbar spine
- Three anatomically pre-moulded splints in the dorsal area
- Viscoelastic lumbo-sacral pad with cross strap system
- Height: ventral 17 cm (6.7 inches), dorsal 26 cm (10.2 inches)

COLOURS: skin-tone
 silver
 black



THE SELECT-BENEFIT

- **Cross strap system** for additional stabilisation and relief
- Easy to put on using the hand straps provided


BORT select® Stabilo® Lady Back Support with Pad No. 104 740



PRODUCT-BENEFIT

- Back brace to stabilise the lumbar spine
- Fitted at the waist: specific for women
- Three anatomically pre-moulded splints in the dorsal area
- Viscoelastic lumbo-sacral pad with cross strap system
- Height: ventral 17 cm (6.7 inches), dorsal 26 cm (10.2 inches)

COLOURS: skin-tone
 silver
 black



THE SELECT-BENEFIT

- **Cross strap system** for additional stabilisation and relief
- Easy to put on using the hand straps provided

INDICATIONS

Treatment of increased pain in the lumbar spine area, severe low back pain, recurrent low back pain, osteochondrosis, spondylarthrosis, lumbosciatica, mild radicular syndrome, spondylolysis without spondylolisthesis, mild to moderate lumbar deformity (degenerative spondylolisthesis), after intervertebral disc surgery

SIZES

Circumference of waist



Size 1	70 – 80 cm
Size 2	80 – 90 cm
Size 3	90 – 100 cm
Size 4	100 – 110 cm
Size 5	110 – 120 cm
Size 6	120 – 130 cm
Size 7	130 – 140 cm

1 cm = 0,3937 inch

SIZES

Circumference of waist



Size 0	60 – 70 cm
Size 1	70 – 80 cm
Size 2	80 – 90 cm
Size 3	90 – 100 cm
Size 4	100 – 110 cm
Size 5	110 – 120 cm
Size 6	120 – 130 cm

1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Back support for muscular stabilisation of the lumbar spine
- Active support with anatomically-contoured knit
- Four supporting splints
- Viscoelastic sacral pad, which can be freely positioned
- Height: ventral 17 cm (6.7 inches), dorsal 23 cm (9.1 inches)

COLOUR: ■ black/green



THE SPORT-BENEFIT

- **COOLMAX®**: wicks moisture away from the body, for dryness and comfort
- **Latex-free knitted fabric**
- Practical hand straps allow support to be fitted easily and without effort

PRODUCT-BENEFIT

- Back support for muscular stabilisation of the lumbar spine
- Active support with anatomically-contoured knit
- Fitted at the waist: specific for women
- Four supporting splints
- Viscoelastic sacral pad, which can be freely positioned
- Height: ventral 17 cm (6.7 inches), dorsal 24 cm (9.4 inches)

COLOUR: ■ black/green



THE SPORT-BENEFIT

- **COOLMAX®**: wicks moisture away from the body, for dryness and comfort
- **Latex-free knitted fabric**
- Practical hand straps allow support to be fitted easily and without effort

INDICATIONS

Pain in the lumbar spine, recurrent low back pain, osteochondrosis, spondylarthrosis, spondylolysis without spondylolisthesis, mild lumbar deformity (degenerative spondylolisthesis), ischialgia, dorsalgia


SIZES	Circumference of waist	
	Size 1	70 – 80 cm
	Size 2	80 – 90 cm
	Size 3	90 – 100 cm
	Size 4	100 – 110 cm
	Size 5	110 – 120 cm
	Size 6	120 – 130 cm
	Size 7	130 – 140 cm

1 cm = 0,3937 inch

SIZES	Circumference of waist	
	Size 0	60 – 70 cm
	Size 1	70 – 80 cm
	Size 2	80 – 90 cm
	Size 3	90 – 100 cm
	Size 4	100 – 110 cm
	Size 5	110 – 120 cm
	Size 6	120 – 130 cm

1 cm = 0,3937 inch

BORT StabiloBasic Back Support with Pad No. 104 690



PRODUCT-BENEFIT

- Back support for muscular stabilisation of the lumbar spine
- Active support with anatomically-contoured knit
- Four supporting splints
- Viscoelastic sacral pad, which can be freely positioned
- Height: ventral 17 cm (6.7 inches), dorsal 23 cm (9.1 inches)


COLOURS: skin-tone
 silver
 black



THE BASIC-BENEFIT

- Hand straps allow support to be fitted easily and without effort
- Comfortable to wear

BORT StabiloBasic Lady Back Support with Pad No. 104 680



PRODUCT-BENEFIT

- Back support for muscular stabilisation of the lumbar spine
- Active support with anatomically-contoured knit
- Cut fitted at the waist: especially for women
- Four supporting splints
- Viscoelastic sacral pad, which can be freely positioned
- Height: ventral 17 cm (6.7 inches), dorsal 24 cm (9.4 inches)

COLOURS: skin-tone
 silver
 black/lilac




THE BASIC-BENEFIT

- Hand straps allow support to be fitted easily and without effort
- Comfortable to wear

INDICATIONS


Pain in the lumbar spine, recurrent low back pain, osteochondrosis, spondylarthrosis, spondylolysis without spondylolisthesis, mild lumbar deformity (degenerative spondylolisthesis), dorsalgia

SIZES Circumference of waist

	Size 1	70 – 80 cm
	Size 2	80 – 90 cm
	Size 3	90 – 100 cm
	Size 4	100 – 110 cm
	Size 5	110 – 120 cm
	Size 6	120 – 130 cm
	Size 7	130 – 140 cm

1 cm = 0,3937 inch

SIZES Circumference of waist

	Size 0	60 – 70 cm
	Size 1	70 – 80 cm
	Size 2	80 – 90 cm
	Size 3	90 – 100 cm
	Size 4	100 – 110 cm
	Size 5	110 – 120 cm
	Size 6	120 – 130 cm

1 cm = 0,3937 inch



PRODUCT-BENEFIT


- Active support for muscular stabilisation of lumbar spine
- Two splints in the dorsal area
- Viscoelastic sacral pad, which can be freely positioned
- Height: ventral 17 cm (6.7 inches), dorsal 23 cm (9.1 inches)

COLOURS: skin-tone
 mineral grey



THE MED-BENEFIT

- Practical hand straps allow support to be fitted easily and without effort
- Comfortable to wear



PRODUCT-BENEFIT

- Elastic support for stabilisation of lumbar spine
- Sacral pad, which can be freely positioned
- Four anatomical pre-molded splints in the dorsal area
- Elastic cross straps
- Height: ventral 17 cm (6.7 inches), dorsal 26 cm (10.2 inches)

COLOURS: white
 black




THE BASIC-BENEFIT

- Grip-straps allow easy fitting and removal of support
- Comfortable to wear


BORT VarioBasic Back Support without Pad
No. 112 650

INDICATIONS

Treatment of pain in the lumbar spine, low back pain, recurrent low back pain, osteochondrosis, spondylarthrosis, spondylolysis without spondylolisthesis, mild lumbar deformity (degenerative spondylolisthesis)


SIZES	Circumference of waist
	Size 1 70 – 80 cm Size 2 80 – 90 cm Size 3 90 – 100 cm Size 4 100 – 110 cm Size 5 110 – 120 cm Size 6 120 – 130 cm Size 7 130 – 140 cm

1 cm = 0,3937 inch

SIZES	Circumference of waist
	x-small 63 – 74 cm small 75 – 89 cm medium 90 – 105 cm large 106 – 120 cm x-large 121 – 140 cm

1 cm = 0,3937 inch

BORT Vario Back Support with Pad No. 112 610




PRODUCT-BENEFIT

- Non-elastic back support for stabilisation of lumbar spine
- Non-elastic zone in the lumbar area made of breathable textile
- Lumbo-sacral friction pad
- Two supporting splints in the dorsal area
- Height: ventral 17 cm (6.7 inches), dorsal 26 cm (10.2 inches)

COLOURS: skin-tone
 black

THE MED-BENEFIT

- Improved fit due to fitted cut
- Additional elastic pull straps from dorsal to ventral area



BORT VarioPlus Back Support special Width No. 112 690




PRODUCT-BENEFIT

- Non-elastic back support for larger sizes
- Semi-elastic body belt with lateral Velcro fastening
- Fitted cut in the lumbar spinal area
- Lumbo-sacral friction pad
- Three supporting splints in the dorsal area
- Height:
 - ventral 24 cm (9.4 inches)
 - dorsal 30 cm (11.8 inches, sizes 1+2)
 - ventral 28 cm (11.0 inches)
 - dorsal 33 cm (13.0 inches, sizes 3+4)

COLOUR: black

THE MED-BENEFIT

- Very comfortable to wear, even for larger torsos
- Additional elastic strap to lift abdomen for relief of the lumbar spine
- Robust design




INDICATIONS

Treatment of increased pain in the lumbar spine area, severe low back pain, recurrent low back pain, osteochondrosis, spondylarthrosis, mild radicular syndrome, spondylolysis without spondylolisthesis, mild to moderate lumbar deformity (degenerative spondylolisthesis), after intervertebral disc surgery

SIZES Circumference of waist

	small	—	up to 89 cm
	medium	—	89 – 104 cm
	large	—	104 – 119 cm
	x-large	—	119 – 134 cm

SIZES Circumference of torso

	Size 1	—	105 – 120 cm
	Size 2	—	120 – 135 cm
	Size 3	—	135 – 150 cm
	Size 4	—	150 – 165 cm

BORT LumboXpress Back Support with Double Fastening
No. 112 600



PRODUCT-BENEFIT

- Elastic support for stabilisation of lumbar spine
- Four anatomical pre-molded splints in the dorsal area
- Two Velcro belts for safe stabilisation
- Two additional elastic pull straps from dorsal to ventral area enable individual dosage of the pressure
- Height: ventral 16 cm (6.3 inches), dorsal 26 cm (10.2 inches)

COLOUR: ■ silver

THE MED-BENEFIT

- Specially constructed fastening system simplifies putting on
- Very flat construction, therefor hardly visible under clothing

BORT Back Support for Pregnant Woman
No. 104 580



PRODUCT-BENEFIT

- Elastic back support for relief of lumbar spine and to support the abdomen
- Individual adjustment of width using flexible adjustable straps
- Thin spring rods for stabilisation of the lumbar region
- Dorsal height: 21 cm (8.3 inches)

COLOUR: □ white

THE MED-BENEFIT

- No risk to baby and abdomen
- Adaptable to varying abdominal circumferences

INDICATIONS

Treatment of increased pain in the lumbar spine area, severe low back pain, recurrent low back pain, osteochondrosis, spondylarthrosis, mild radicular syndrome, spondylolysis without spondylolisthesis, mild to moderate lumbar deformity (degenerative spondylolisthesis), after intervertebral disc surgery

INDICATIONS

For stabilisation of lumbar spine with metered pressure during pregnancy, for the treatment of pain in the lumbar area

SIZES Circumference of waist



Size 1	70 – 80 cm
Size 2	80 – 90 cm
Size 3	90 – 100 cm
Size 4	100 – 110 cm
Size 5	110 – 120 cm
Size 6	120 – 130 cm
Size 7	130 – 140 cm

1 cm = 0,3937 inch

SIZES Circumference of hip



Size 0	75 – 91 cm
Size 1	91 – 105 cm
Size 2	105 – 119 cm
Size 3	119 – 140 cm

1 cm = 0,3937 inch

BORT Symphysis Belt No. 102 700



PRODUCT-BENEFIT

- Non-elastic support for stabilisation of pelvis
- Pressure on the symphysis can be individually adjusted
- Anatomically designed cut
- Padded on the inside
- Height: ventral 8,5 cm (3.3 inches), dorsal 15 cm (5.9 inches)

COLOUR: ■ silver

THE MED-BENEFIT

- Circular stabilisation and relief of pelvis and symphysis
- Comfortable to wear thanks to narrow closure section (for example when sitting)

INDICATIONS

Symphysis insufficiency, symphysis rupture, structural loosening in the area of the iliosacral joints/pubis symphysis, after birth

BORT Stabilo® Sacroiliac Brace No. 104 750



PRODUCT-BENEFIT

- Brace for the stabilisation of sacroiliac joints
- Special combination of elastic and non-elastic knit fabric
- Two-part sacroiliac pad to stimulate the circulation and encourage the healing process
- Circular tension adjustable with additional tension straps

COLOUR: ■ silver


THE MED-BENEFIT

- Circular stabilisation and relief of pelvis and symphysis
- Very comfortable to wear thanks to narrow closure section (for example when sitting)


INDICATIONS

Sacroiliac joint arthrosis, sacroiliac joint instability, sacroiliac joint blockage, myalgias and tendopathies in the pelvic region, structural abnormalities after lumbar spinal spondylodeses

SIZES Circumference of hip

	small	85 – 95 cm
	medium	95 – 105 cm
	large	105 – 115 cm
	x-large	115 – 125 cm
	xx-large	125 – 135 cm

SIZES Circumference of iliac crest

	Sizes 1	70 – 85 cm
	Sizes 2	85 – 100 cm
	Sizes 3	100 – 115 cm
	Sizes 4	115 – 130 cm
	Sizes 5	130 – 145 cm

1 cm = 0,3937 inch


1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Back brace for stabilisation of lumbar spine
- Integrated stabilisation elements
- Elastic, breathable carrier material
- Viscoelastic pads which can be freely positioned
- Height: ventral 17 cm (6.7 inches), dorsal 21 cm (8.3 inches)

COLOUR: ■ blue



THE GENERATION-BENEFIT

- Wide grip straps on the latch plate
- Innovative light brace construction
- **Energy-saving pulley system**
- Simple and easy to put on

INDICATIONS

Degenerative lumbar change, mild to moderate lumboischialgia, slipped disc with mild deficits, spondylosis, spondylolisthesis, spinal canal stenoses, radicular irritation syndrome, condition after inter-vertebral disc surgery

SIZES	Circumference of waist
Size 1	70 – 80 cm
Size 2	80 – 90 cm
Size 3	90 – 100 cm
Size 4	100 – 110 cm
Size 5	110 – 120 cm
Size 6	120 – 130 cm
Size 7	130 – 140 cm

1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Spinal brace for relief and stabilisation of the lumbar spine
- Integrated plastic pelvic girdle
- Elastic, breathable carrier material
- Back plate heights: 22 cm (8.7 inches), 28 cm (11.0 inches)

COLOUR: ■ blue

THE GENERATION-BENEFIT

- Wide grip straps on the latch plate
- Innovative buckle shape supports own fixation of the brace when putting it on
- **Energy-saving pulley system**
- Simple and easy to put on

INDICATIONS

Facet syndrome, lumbosciatica, intervertebral disc protrusions/prolapse, radicular irritation syndrome, degenerative lumbar change, condition after inter-vertebral disc surgery, stable fractures, spondylosis, spondylolisthesis, spinal canal stenoses

SIZES	Circumference of waist
Size 1	70 – 80 cm
Size 2	80 – 90 cm
Size 3	90 – 100 cm
Size 4	100 – 110 cm
Size 5	110 – 120 cm
Size 6	120 – 130 cm
Size 7	130 – 140 cm

please indicate height 22 cm (8.7 inches) or 28 cm (11.0 inches)
1 cm = 0,3937 inch



BORT Stabilo® FA Back Brace

No. 180 500

NEW



PRODUCT-BENEFIT

- Lumbar brace for partial relief of the lumbar spine and **facet joints**
- Two lateral supporting splints
- Height: ventral 15,5 cm (6.1 inches), dorsal 30 cm (11.8 inches)

COLOUR: ■ grey

THE MED-BENEFIT

- Elastic belt system to reduce lordosis of the lumbar spine
- Hand straps allow easy, effortless fitting
- Adjustable width using ergonomic latch plates


Back

INDICATIONS

Facet joint arthrosis, intervertebral disc protrusions/prolapse, severe low back pain, lumboischialgia, lumbago, lumbago with sacroiliac joint involvement, chronic pain in the lumbosacral area, spondylolisthesis, muscular imbalance, segment instability, lumbar spinal canal stenoses

SIZES

Circumference of waist



Size 1	_____	70 – 80 cm
Size 2	_____	80 – 90 cm
Size 3	_____	90 – 100 cm
Size 4	_____	100 – 110 cm
Size 5	_____	110 – 120 cm
Size 6	_____	120 – 130 cm
Size 7	_____	130 – 140 cm

1 cm = 0,3937 inch


BORT StabiloPlus Back Support
No. 104 730



PRODUCT-BENEFIT

- Lumbar brace for bridging and for relief of the lumbar spine
- Stimulation of abdominal muscles and straightening of pelvis
- Semi-elastic stabilisation strap, individually adjustable
- Height: ventral 17 cm (6.7 inches), dorsal 32 cm (12.6 inches)

COLOUR: ■ silver



THE MED-BENEFIT

- Anatomical cut for precise fit, flexible section in the hip-bending area
- Improved function and fit thanks to rotating strap guide loops
- The degree of stabilisation can be adjusted to suit the particular course of therapy
- Easy to put on using the hand straps provided

INDICATIONS

Moderate lumboschialgia with mild deficits with herniated discs/prolapse, moderate radicular, pseudo-radicular lumbar syndrome, grade I spondylolisthesis with low back pain, moderate lumbar deformity (recurrent treatment, complicated history) i. e. after inter-vertebral disc surgery with mild to moderate prolapse, spinal canal stenoses, lumbar instability syndrome

SIZES

Circumference of waist



Size 1	70 – 80 cm
Size 2	80 – 90 cm
Size 3	90 – 100 cm
Size 4	100 – 110 cm
Size 5	110 – 120 cm
Size 6	120 – 130 cm
Size 7	130 – 140 cm

1 cm = 0,3937 inch

BORT StabiloPlus Bridging Brace

No. 104 760



PRODUCT-BENEFIT

- Lumbar brace for bridging and for relief of the lumbar spine
- Classical brace-shaped stable plastic bridging frame, anatomically pre-moulded, thermoplastically adjustable
- Non-elastic stabilisation strap, individually adjustable
- Height: ventral 17 cm (6.7 inches), dorsal 32 cm (12.6 inches)

THE MED-BENEFIT

- High degree of stabilisation, can be gradually adjusted during the course of therapy (mobility function)
- Additional ventral abdominal reinforcing plate
- Basic support with an anatomic cut for optimum fitting
- Easy to put on using the hand straps provided

COLOUR: ■ black

INDICATIONS

Severe lower back pain with muscular weakness and significant hernia protrusion/prolapse, without absolute indication for surgery, severe lumbar radicular or pseudo-radicular syndrome (refractory to conservative treatment), grade 2 spondylolisthesis with recurring lower back pain, severe lumbar deformity with facet syndrome/arthrosis, post-surgery for herniated disks without resulting paralysis with removal of the hernia, decompression of the spinal canal, single level, without deformity, single level herniated disk surgery, short stenosis of the spinal canal without significant paralysis, vertebral fractures resulting from minor trauma, superior end plate deformities in the lumbar spine

SIZES

Circumference of waist



Size 1	70 – 80 cm
Size 2	80 – 90 cm
Size 3	90 – 100 cm
Size 4	100 – 110 cm
Size 5	110 – 120 cm
Size 6	120 – 130 cm
Size 7	130 – 140 cm

1 cm = 0,3937 inch

BORT Stabilo® Lumbar Back Brace
No. 180 100

NEW
DESIGN



PRODUCT-BENEFIT

- Lumbar brace with integrated plastic pelvic girdle to secure the LV 1 – LV 5 segments
- Continuously adjustable stabilisation of the lumbar spine using Boa-technology (strap system)
- Integrated lumbar support with abdominal closing plate

THE MED-BENEFIT

- Rotation of the lumbar spine is effectively restricted
- Modular design: can be downgraded during the course of therapy
- Easily adjustable: Setting lumbar lordosis possible via cold deformation

COLOUR:  grey 

INDICATIONS

Major lumboschialgia without absolute indication for surgery, when surgery is contraindicated, when surgery is rejected, major radicular and pseudo radicular lumbar syndrome (cautionary, therapy-resistant), spondylolisthesis with severe lumboschialgia, severe lumbar deformity with hypermobility associated with Facet Syndrome / Osteoarthritis, condition after spine surgery, i.e. dorso-ventral/dorso-lateral fusion, spinal column compression (single level with deformity and multi-level), thoracic- surgery with severe prolapse/lumbar instability, spinal canal stenosis with paresis, post-traumatic vertebral fractures with significant spinal damage

SIZES

Circumference of waist



Size 1	_____	70 – 80 cm
Size 2	_____	80 – 90 cm
Size 3	_____	90 – 100 cm
Size 4	_____	100 – 110 cm
Size 5	_____	110 – 120 cm
Size 6	_____	120 – 130 cm
Size 7	_____	130 – 140 cm

BORT DorsoFX Back Brace

No. 180 400

**PRODUCT-BENEFIT**

- Immobilisation brace for safe stabilisation of the thoracic spine up to TH 8
- Modular back support with reclination brace for immobilisation
- Easy height adjustment using sliding mechanism

THE MED-BENEFIT

- Individual brace design (Modular-brace framework, back support)
- Divided pelvis buckle prevents unwanted pressure on iliac spine
- Very comfortable to wear due to full-surface padding

COLOUR: ■ grey/black

INDICATIONS

Conservative treatment of fresh, traumatic, stable fractures of the thoracic spine (Th 8-12) and lumbar spine (L 1-5), conservative treatment of osteoporotic compression fractures of the thoracic spine (T 8-12) and lumbar spine (L 1-5), immobilization in the case of acute and/or chronic spinal diseases (in the above mentioned section), if immobilization is indicated according to the physician's instructions

SIZES (Thoracic strap)	Circumference of chest	Circumf. of Hip (H) (Pelvic corset)	Circumf. of Waist (T) (Support)
Size B0 _____	70 – 78 cm	H0 _ 86 – 94 cm	T0 _ 60 – 70 cm
Size B1 _____	78 – 86 cm	H1 _ 94 – 102 cm	T1 _ 70 – 80 cm
Size B2 _____	86 – 94 cm	H2 _ 102 – 110 cm	T2 _ 80 – 90 cm
Size B3 _____	94 – 102 cm	H3 _ 110 – 118 cm	T3 _ 90 – 100 cm
Size B4 _____	102 – 110 cm	H4 _ 118 – 126 cm	T4 _ 100 – 110 cm
Size B5 _____	110 – 118 cm	H5 _ 126 – 134 cm	T5 _ 110 – 120 cm
SIZES (Back plate)	Back Length (L)		
small _____	35 – 42 cm		
medium _____	42 – 50 cm		
large _____	50 – 58 cm		

1 cm = 0,3937 inch

BORT DynamicFX Back Brace
No. 180 300



PRODUCT-BENEFIT

- Back brace for osteoporosis characterised by dynamic reclination
- Therapeutic spring tension can be individually adjusted for permanent stabilisation of the achieved straightening
- With flexible, lightweight lumbar support

COLOUR: ■ grey

THE MED-BENEFIT

- Activation of the torso muscles
- Cold-formable back plate made of ultra light honeycomb aluminium
- Increase in mobility and decrease in the risk of falling
- The intensity of the dynamics can be selected by means of five different spring elements depending upon the patient's requirements

INDICATIONS

Active straightening for relief and correction of the lumbar/thoracic spine associated with restricted movement in the sagittal plane, stable, osteoporotic vertebral fractures, Osteoporosis, cautionary treatment of secondary kyphosis (stable fractures, tumors), muscular insufficiency, pre-operative and post-operative stabilisation

SIZES Support	Circumf. of waist	SIZES	Length back plate
Size 0	60 – 70 cm	small	32 cm
Size 1	70 – 80 cm	medium	40 cm
Size 2	80 – 90 cm	large	48 cm
Size 3	90 – 100 cm	x-large	56 cm
Size 4	100 – 110 cm	The length of the back plate is determined by the required stabilisation level	
Size 5	110 – 120 cm		
Size 6	120 – 130 cm		

Please ensure that the support and back plate are always ordered together (e.g. size 3 / medium)!

1 cm = 0,3937 inch

BORT Generation Osteoporosis Brace

No. 215 870



Generation



PRODUCT-BENEFIT

- Activating and spine-straightening back brace
- Stable aluminium back splint, anatomically pre-shaped and individual, without tools, adjustable
- Relief and stabilisation of the lumbar spine
- Active straightening by activating the back muscles
- Increasing mobility

COLOUR: ■ blue

THE GENERATION-BENEFIT

- Strap system according to the rucksack principle
- Guide loops in the chest area for easy, stepless straightening of the upper body
- Lightweight and slim structure for good wearing comfort
- Easy to handle
- Fastener with a large finger hold and hand straps for effortless fitting




BORT Generation:
For better mobility for the aged
Product range with special handling advantages

INDICATIONS

Active straightening for relief and correction of the lumbar/thoracic spine associated with restricted movement in the sagittal plane, stable, osteoporotic vertebral fractures, Osteoporosis, cautionary treatment of secondary kyphosis (stable fractures, tumors), muscular insufficiency, pre-operative and post-operative stabilisation

SIZES

Back length




x-small	_____	up to 43 cm ^{NEW}
small	_____	44 – 49 cm
medium	_____	50 – 55 cm
large	_____	56 – 61 cm
x-large	_____	62 – 67 cm

1 cm = 0,3937 inch

Active and secure for everyday life:
Straighten, compress, stabilise.



Overview of Indications for selected products	Rib Belt for Men 102 800 Rib Belt for Women 102 900 S. 119 / 118	Scar Rupture Support 104 050 S. 120	Umbilical Hernia Support 104 070 S. 120	Stoma Support 104 090 S. 121	Vario Stoma Support 104 010 S. 121	PostOban® 104 150 PostOban® Special Width 104 180 S. 122 / 123	PostOban® SOFT Thorax-Abdominal- Support 104 100 S. 122	Abdominal Support for Pregnant Women 104 620 S. 123
Page								
Rib fracture/contusion	●							
Incisional hernia		●						
Umbilical hernia			●					
Weakness of the connective tissues						●		
Paralysis of abdominal wall				●	●			
Condition after ileostomy and colostomy				●	●			
Relief of abdomen								●
Pendulous abdomen							●	
Post-operative		●				●		




PRODUCT-BENEFIT

- Non-elastic circular support for stabilisation of the thorax
- **Non-elastic** side and front section with Velcro fastening, elastic back section
- Breathable fabric with Velcro fastening
- Feminine cut fitted at the waist
- Height: dorsal 16 cm (6.3 inches), ventral 12 cm (4.7 inches)

COLOUR: ☐ white

THE MED-BENEFIT

- Hardly noticeable, invisible under clothing
- Very comfortable to wear due to the use of cotton satin and fine terry cloth material



PRODUCT-BENEFIT

- Elastic circular support for stabilisation of the thorax
- **Elastic**, breathable fabric with Velcro fastening
- Feminine cut fitted at the waist
- Height: dorsal 16 cm (6.3 inches), ventral 12 cm (4.7 inches)


COLOUR: ☐ white

THE MED-BENEFIT


- Hardly noticeable, not visible under clothing
- Good wearing comfort

INDICATIONS

Rib fracture, trauma consequences in the rib area (e. g. contusion, sprain)


SIZES	Underbust circumference	
	small	up to 66 cm
	medium	66 – 76 cm
	large	76 – 89 cm
	x-large	89 – 102 cm
	xx-large	102 – 125 cm

1 cm = 0,3937 inch

SIZES	Underbust circumference	
	small	65 – 78 cm
	medium	78 – 91 cm
	large	91 – 106 cm
	x-large	106 – 125 cm

1 cm = 0,3937 inch

**BORT Rib Belt
for Men**
No. 102 800



PRODUCT-BENEFIT

- Non-elastic circular support for stabilisation of the thorax
- **Non-elastic** side and front section, elastic back section
- Breathable fabric with Velcro fastening
- Straight cut
- Height: 16 cm (6.3 inches)

COLOUR: ☐ white

THE MED-BENEFIT

- Hardly noticeable, not visible under clothing
- Very comfortable to wear due to the use of cotton satin and fine terry cloth material

**BORT Elastic Rib Belt
for Men**
No. 102 850



PRODUCT-BENEFIT

- Elastic circular support for stabilisation of the thorax
- **Elastic**, breathable fabric with Velcro fastening
- Straight cut
- Height: 16 cm (6.3 inches)

COLOUR: ☐ white


THE MED-BENEFIT

- Hardly noticeable, not visible under clothing
- Good wearing comfort

INDICATIONS


Rib fracture, trauma consequences in the rib area (e. g. contusion, sprain)

SIZES Circumference of chest

	small	—	up to 82 cm
	medium	—	82 – 96 cm
	large	—	96 – 112 cm
	x-large	—	112 – 128 cm
	xx-large	—	128 – 144 cm

1 cm = 0,3937 inch

SIZES Circumference of chest

	small	—	65 – 78 cm
	medium	—	78 – 91 cm
	large	—	91 – 106 cm
	x-large	—	106 – 125 cm
	xx-large	—	125 – 145 cm

1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Elastic compression support for incisional hernia
- Round silicone pad reinforced with plastic, for individual positioning
- Two supporting splints
- Elastic, breathable fabric with Velcro fastening
- Heights: 13 cm (5.1 inches), 21 cm (8.3 inches), 26 cm (10.2 inches)

COLOUR: ☐ white



Schematic representation of pad

THE MED-BENEFIT

- Removable pad cover made of soft terrycloth
- Height 13 cm (5.1 inches): round pad with 10 cm (3.94 inches) diameter
- Height 21 cm (8.3 inches): pad 18 (7.09) × 10 cm (3.94 inches), oval
- Height 26 cm (10.2 inches): pad 22 (8.7) × 10 cm (3.94 inches), oval

PRODUCT-BENEFIT

- Elastic compression support for umbilical hernia
- Anatomically concaved umbilical hernia pad, with adhesive layer for individual positioning
- Two supporting splints
- Elastic, breathable fabric with Velcro fastening
- Height: 13 cm (5.1 inches)

COLOUR: ☐ white



Schematic representation of pad

THE MED-BENEFIT

- Removable pad cover made of soft terry cloth


INDICATIONS

Scar rupture, scar insufficiency, abdominal muscle infirmity/paralysis

INDICATIONS

Umbilical hernia, reducible umbilical hernias

SIZES Circumference of waist

	Size 1	_____	60 – 95 cm
	Size 2	_____	95 – 125 cm
	Size 3	_____	125 – 150 cm

please indicate height 13 cm (5.1 inches),
21 cm (8.3 inches) or 26 cm (10.2 inches)
1 cm = 0,3937 inch

SIZES Circumference of waist

Size 1	_____	60 – 95 cm
Size 2	_____	95 – 125 cm
Size 3	_____	125 – 150 cm

1 cm = 0,3937 inch

BORT Stoma Support

No. 104 090



PRODUCT-BENEFIT

- Circular abdominal support with a variable stoma recess
- Non-elastic composite velour material, elastic insert in the back section
- Two supporting splints
- Plastic reinforcing plate to be inserted
- Heights: 15 cm (5.9 inches), 22 cm (8.7 inches)

COLOUR: ☐ white

THE MED-BENEFIT

- The variable stoma recess width can be adjusted using a Velcro fastening
- Large, skin-friendly lining in stoma area
- Easy to put on and position the stoma recess using the Velcro fastener on the right and left in the abdominal region
- Washable at 90°C

ACCESSORIES:

No. 104 099 | BORT Stoma Protector

INDICATIONS

Paralysis of abdominal wall, hernia of abdominal wall, condition after ileostomy and colostomy with accompanying weakness of connective tissue

SIZES Circumference of waist

small	70 – 90 cm
medium	90 – 115 cm
large	115 – 150 cm

Please indicate height 15 cm (5.9 inches) or 22 cm (8.7 inches)
1 cm = 0,3937 inch

BORT Vario Stoma Support

No. 104 010

NEW



PRODUCT-BENEFIT

- Circular abdominal support with a variable stoma recess
- Non-elastic composite velour material, elastic insert in the back section
- Two supporting splints
- The size, shape and positioning of the stoma recess can be individually determined
- Heights: 15 cm (5.9 inches), 22 cm (8.7 inches)

COLOUR: ☐ white

THE MED-BENEFIT

- Stencil provided for easy, individual size determination of the stoma recess
- Large, skin-friendly lining in stoma area
- Washable at 90°C

SIZES Circumference of waist

small	70 – 90 cm
medium	90 – 115 cm
large	115 – 150 cm

Please indicate height 15 cm (5.9 inches) or 22 cm (8.7 inches)
1 cm = 0,3937 inch




PRODUCT-BENEFIT

- Elastic circular support for stabilisation of thoracic and abdominal regions
- Elastic, soft fabric
- Wide, two-piece Velcro fastening
- Heights: 21 cm (8.3 inches), 26 cm (10.2 inches)

COLOUR: ☒ skin-tone

THE MED-BENEFIT

- Individually adjustable
- Comfortable to wear



PRODUCT-BENEFIT

- Elastic circular support for stabilisation of thoracic and abdominal regions
- Elastic, breathable fabric
- Wide, two-part Velcro fastening
- Heights: 21 cm (8.3 inches), 26 cm (10.2 inches)

COLOUR: ☐ white

THE MED-BENEFIT

- Individually adjustable
- Two-part Velcro fastening makes the support easy to put on and ensures a better fit

INDICATIONS

After abdominal surgery, weakness of the connective tissue in the abdominal area, pendulous abdomen

INDICATIONS


After abdominal/thoracic surgery, weakness of the connective tissue in the abdominal area, pendulous abdomen

SIZES Circumference of waist

	small	—	up to 87 cm
	medium	—	87 – 102 cm
	large	—	102 – 117 cm
	x-large	—	117 – 132 cm

please indicate height 21 cm (8.3 inches) or 26 cm (10.2 inches)
1 cm = 0,3937 inch

SIZES Circumference of waist

	Size 1	—	75 – 85 cm
	Size 2	—	85 – 95 cm
	Size 3	—	95 – 105 cm
	Size 4	—	105 – 115 cm
	Size 5	—	115 – 125 cm
	Size 6	—	125 – 135 cm
	Size 7	—	135 – 145 cm

please indicate height 21 cm (8.3 inches) or 26 cm (10.2 inches)
1 cm = 0,3937 inch

BORT PostOban® Special Width Thorax-Abdominal-Support No. 104 180



PRODUCT-BENEFIT

- Elastic circular support for stabilisation of thoracic and abdominal regions for larger sizes
- Two-piece elastic fastener
- Elastic, breathable fabric
- Heights: 21 cm (8.3 inches), 26 cm (10.2 inches)

COLOUR: ☐ white

THE MED-BENEFIT

- Very comfortable to wear, even for larger torsos
- Individually adjustable
- Integrated dressing aids
- Easy to put on


INDICATIONS

After abdominal surgery, weakness of the connective tissue in the abdominal area, pendulous abdomen

SIZES

Circumference of torso

Size 1	105 – 120 cm
Size 2	120 – 135 cm
Size 3	135 – 150 cm
Size 4	150 – 165 cm



BORT Abdominal Support for Pregnant Women No. 104 620



PRODUCT-BENEFIT

- Elastic bandage to support the abdomen
- Elastic, breathable fabric
- Four supporting rods which can be removed, if required
- Dorsal height: 24 cm (9.5 inches)

COLOUR: ☐ white

THE MED-BENEFIT

- Back relief due to additional elastic strap
- Adjustable to fit varying abdominal circumferences

INDICATIONS

Lumbar back pain in pregnancy, obesity, multiple pregnancy

SIZES

Circumference of iliac crest

Size 0	75 – 90 cm
Size 1	90 – 105 cm
Size 2	105 – 120 cm
Size 3	120 – 135 cm



please indicate height 21 cm
(8.3 inches) or 26 cm (10.2 inches)
1 cm = 0,3937 inch

1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Single-sided, rounded pressure pad
- Non-elastic hip belt with approx. 50 cm (19.7 inches) long spring made of spring steel


COLOUR:  skin-tone



PRODUCT-BENEFIT

- Single-sided, anatomically-shaped pad
- Non-elastic hip belt with approx. 50 cm (19.7 inches) long spring made of spring steel
- Thigh strap prevents pad from slipping

COLOUR:  skin-tone



PRODUCT-BENEFIT

- Single-sided, rounded pressure pad
- Elastic hip belt with short spring 20 cm (7.9 inches) made of spring steel

COLOUR:  skin-tone




PRODUCT-BENEFIT

- Single-sided, anatomically-shaped pad
- Elastic hip belt with short spring 20 cm (7.9 inches) made of spring steel
- Thigh strap prevents pad from slipping


COLOUR:  skin-tone

INDICATIONS

Reducible inguinal rupture

SIZES	Circumference of iliac crest	
	85 _____	83 – 88 cm
	90 _____	88 – 93 cm
	95 _____	93 – 98 cm
	100 _____	98 – 103 cm
	105 _____	103 – 108 cm

please indicate left or right
1 cm = 0,3937 inch

SIZES	Circumference of iliac crest	
	85 _____	83 – 88 cm
	90 _____	88 – 93 cm
	95 _____	93 – 98 cm
	100 _____	98 – 103 cm
	105 _____	103 – 108 cm

please indicate left or right
1 cm = 0,3937 inch

BORT HerniaFix® Elastic Spring Truss, Double-Sided
No. 109 500




PRODUCT-BENEFIT

- Double-sided, rounded pressure pad
- Elastic hip belt with approx. 20 cm (7.9 inches) long spring made of spring steel

COLOUR: skin-tone

BORT HerniaFix® Special Spring Truss with Adjustable Anatomically-Shaped Pad No. 110 300



PRODUCT-BENEFIT

- Special spring truss, no pressure or spring, single sided
- Anatomically adjustable pad
- Elastic hip belt
- Thigh strap

COLOUR: skin-tone

BORT HerniaFix® Special Spring Truss with Adjustable Rounded Pad No. 110 100

PRODUCT-BENEFIT


- Special spring truss, no pressure or spring, single sided
- Rounded pad, which can be adjusted
- Elastic hip belt

COLOUR: skin-tone

INDICATIONS

Reducible inguinal rupture

SIZES Circumference of iliac crest

	85	_____	83 – 88 cm
	90	_____	88 – 93 cm
	95	_____	93 – 98 cm
	100	_____	98 – 103 cm
	105	_____	103 – 108 cm

right and left symmetrical
1 cm = 0,3937 inch

SIZES Circumference of iliac crest

	_____	85 – 94 cm
	_____	95 – 105 cm
	_____	106 – 115 cm

please indicate left or right
1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Inguinal hernia belt
- Double sided, with skin patch
- Two detachable pads
- Double sided thigh strap

COLOUR: ☐ white

THE MED-BENEFIT

- Easy, safe to handle
- Can be used with just one pad on one side



PRODUCT-BENEFIT

- Special ring cushion with recess in the coccyx area
- **Inflatable**, via valve
- Supplied with manual air pump
- Size: approx. 45 cm × 45 cm, (17.7 × 17.7 inches)
- Inner diameter approx. 16 cm (6.3 inches)

COLOUR: ☒ blue

THE MED-BENEFIT


- Suitable for all types of seating
- Inflatable, therefore easy to carry and can be used anywhere

INDICATIONS

Reducible, single and double sided inguinal rupture

INDICATIONS

Coccygodynia, post coccyx contusion and fracture conditions, perianal relief of haemorrhoids, post low rectal surgery, periproctitis abscess, post episiotomy condition

SIZES	Circumference of hip	
	small	— up to 89 cm
	medium	— 89 – 104 cm
	large	— 104 – 119 cm

SIZE
One length only

BORT Generation Hip Protector No. 215 900



PRODUCT-BENEFIT

- Hip protection belt with a shock absorbing effect
- Protectors made of soft visco-elastic foam can be removed
- Strap with an anatomically designed cut can be fully opened

COLOUR: ■ blue

THE GENERATION-BENEFIT

- Easy to put on/remove thanks to Velcro fastening
- Wide frontal grip straps
- Safe handling and comfortable to wear
- Ideal for temporary use during phases of activity
- Attractive, functional design

ACCESSORIES:

- Hip Protector set (No. 201 120)

INDICATIONS

Tendency to fall for various reasons, danger of falling

BORT StabiloHip Hip Protector No. 201 100



PRODUCT-BENEFIT

- Hip protection pants with shock absorbing effect
- Protectors made of soft visco-elastic foam can be removed

COLOUR: □ white

THE MED-BENEFIT

- Highly elastic waistband and leg edge
- Protectors can be easily removed
- High-quality functional pants made of cotton and Lycra


ACCESSORIES:

- Hip Protector set (No. 201 120)
- Pants (No. 201 200)

INDICATIONS

Tendency to fall for various reasons

SIZES	Circumference of iliac crest	
	small	— 80 – 90 cm
	medium	— 90 – 100 cm
	large	— 100 – 110 cm
	x-large	— 110 – 120 cm
	xx-large	— 120 – 130 cm
	xxx-large	— 130 – 140 cm



1 cm = 0,3937 inch

SIZES	Dress Size	Women USA	Men USA
small	_____	38	46
medium	_____	40	48
large	_____	42	50
x-large	_____	44	52
xx-large	_____	46/48	54/56
xxx-large	_____	50/52	58/60

Walking into the future:
functional stabilisation with
maximum mobility.



Overview of indications for selected products	select® Stabilo- Gen® 114 420	Stabilo- Gen® Sport 114 520 SP	Patella Tendon Strap Sport 114 010 SP	Stabilo Gen®Eco 114 520	Generation Knee Support 215 420	Asymme- tric® 114 700	Osgood- Schlatter Knee Support 114 510	Stabilo® Knee Support with Hinge 150 140	Knee Support with Patella Recess 114 150	Immob Splint 145 000 Adjustable Immobi- lisation Splint 145 080	Generation Knee Brace 215 400	StabiloPro® 182 150 StabiloPro® Open Style 182 300	Genu- Xpress 100 700	OA- Xpress 101 600
Page	P. 130/131	P. 132	P. 139	P. 134	P. 133	P. 141/142	P. 137	P. 142	P. 136	P. 148	P. 145	P. 144	P. 147	P. 146
Irritations	●	●		●	●		●							
Patellar lateralisation						●								
Patellar (sub) luxation						●								
Patellar chondropathy			●						●					
Patellofemoral arthrosis			●											
Meniscopathy	●	●		●	●	●			●					
Ligament instability	●	●		●				●			●	●	●	
Gonarthrosis	●	●		●	●						●	●		
Varus/Valgus Gonarthrosis														●
Complete immobilisation										●				
Osgood-Schlatter disease							●							
ACL/PCL ruptur													●	●



PRODUCT-BENEFIT

- High-quality knee support for compression of soft tissues with silicone pad for stabilisation and relief of the knee joint
- Lengthwise-elastic stretch zone above the silicone pad for ideal fit at all bending angles
- Pad for the popliteal fossa ensures the best possible comfort for wearing
- One spiral stay on each side for lateral support
- Compression class II

COLOURS:  skin-tone  silver
 black


THE SELECT-BENEFIT

- **SOFTflex Plus knit fabric:**
 - BREEZE: Luxury fibre with a cooling effect
 - Soft-velvety surface structure for top wearing comfort
- Extra-soft, snag proof material edge
- Plus-size system results in a good fit

BORT INNOVATION:
- Woven silicone bonding zone in the thigh area

INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, gonarthrose, arthritis, e. g. ligament instability, meniscopathy

SIZES		d-measure f-measure		d-measure f-measure	
	Size 1 (small)	26 – 29	32 – 36 cm	Size 4+ (x-large plus)	37 – 41 48 – 52 cm
	Size 2 (medium)	29 – 33	36 – 40 cm	Size 5 (xx-large)	41 – 45 52 – 56 cm
	Size 2+ (medium plus)	29 – 33	40 – 44 cm	Size 5+ (xx-large plus)	41 – 45 56 – 60 cm
	Size 3 (large)	33 – 37	40 – 44 cm	Size 6 (xxx-large)	45 – 49 60 – 64 cm
	Size 3+ (large plus)	33 – 37	44 – 48 cm	Size 6+ (xxx-large plus)	45 – 49 64 – 68 cm
	Size 4 (x-large)	37 – 41	44 – 48 cm		

right and left symmetrical
1 cm = 0,3937 inch

BORT select® GenuZip®

No. 114 390



select®

BREEZE



PRODUCT-BENEFIT

- High-quality knee support for compression of soft tissues with silicone pad for stabilisation and relief of the knee joint
- **Side zip fastener**
- Lengthwise-elastic stretch zone above the silicone pad ensures an ideal fit at all bending angles
- **With an highlyly elastic pad for the popliteal fossa, providing the best possible comfort for wearing**
- One spiral stay on each side for lateral support
- Compression class II

COLOURS:  skin-tone  silver

BORT INNOVATION:

- Woven silicone bonding zone in the thigh area

THE SELECT-BENEFIT

- **SOFTflex Plus knit fabric:**
 - BREEZE: Luxury fibre with a cooling effect
 - Soft-velvety surface structure for top wearing comfort
- Side zip fastener makes it easy to put on/remove
- Extra-soft, snag proof material edge
- Silicone adhesive border for a secure fit
- Very comfortable to wear
- Plus-size system results in a very good fit

INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, gonarthrose, arthritis, e. g. ligament instability, meniscopathy

SIZES	d-measure	f-measure	d-measure	f-measure
Size 1 (small)	26 – 29	32 – 36 cm	Size 4+ (x-large plus)	37 – 41 48 – 52 cm
Size 2 (medium)	29 – 33	36 – 40 cm	Size 5 (xx-large)	41 – 45 52 – 56 cm
Size 2+ (medium plus)	29 – 33	40 – 44 cm	Size 5+ (xx-large plus)	41 – 45 56 – 60 cm
Size 3 (large)	33 – 37	40 – 44 cm	Size 6 (xxx-large)	45 – 49 60 – 64 cm
Size 3+ (large plus)	33 – 37	44 – 48 cm	Size 6+ (xxx-large plus)	45 – 49 64 – 68 cm
Size 4 (x-large)	37 – 41	44 – 48 cm		

please indicate left or right

1 cm = 0,3937 inch

BORT StabiloGen® Sport
No. 114 520 SP





Anatomically and functionally optimised pad form

Sport

COOLMAX®

PRODUCT-BENEFIT

- Knee support for compression of soft tissues for stabilisation and relief of the knee joint
- **SOFTflex 3D knitted fabric** with a velvety soft surface structure for more elasticity, stability and very comfortable to wear
- One lateral spiral stay on each side
- Plus-size system results in a good fit
- Highly elastic, soft flexion region in the hollow of the knee
- Compression class II


THE SPORT-BENEFIT

- Woven silicone bonding zone in the thigh area
- **COOLMAX®**: wicks moisture away from the body, for dryness and comfort
- Ideal fit for sports activities
- **Latex-free knitted fabric**
- **Anatomically and functionally optimised pad form** for safe guidance of the patella during intense activity
- Modern, sporty design (with a reflecting label)

COLOUR: ■ black/green

INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, e. g. gonarthrose, arthritis, ligament instability, meniscopathy

SIZES		d-measure f-measure		d-measure f-measure	
	Size 1 (small)	26 – 29	32 – 36 cm	Size 4+ (x-large plus)	37 – 41 48 – 52 cm
	Size 2 (medium)	29 – 33	36 – 40 cm	Size 5 (xx-large)	41 – 45 52 – 56 cm
	Size 2+ (medium plus)	29 – 33	40 – 44 cm	Size 5+ (xx-large plus)	41 – 45 56 – 60 cm
	Size 3 (large)	33 – 37	40 – 44 cm	Size 6 (xxx-large)	45 – 49 60 – 64 cm
	Size 3+ (large plus)	33 – 37	44 – 48 cm	Size 6+ (xxx-large plus)	45 – 49 64 – 68 cm
	Size 4 (x-large)	37 – 41	44 – 48 cm		

right and left symmetrical
1 cm = 0,3937 inch

BORT Generation Knee Support No. 215 420

NEW



Generation



BORT Generation Knee Support Plus No. 114 490 G



Generation



PRODUCT-BENEFIT

- Knee support with adjustable thigh width for stabilisation of the knee joint
- One lateral spiral stay each side
- Functional knitted fabric, very comfortable to wear
- Anatomic silicone pad
- **Integrated plus size system**

COLOUR: ■ blue

THE GENERATION-BENEFIT

- **Easy to put on/remove** thanks to closed design on the thigh and open design on the shank
- Grip straps allow safe and energy-saving handling, also in case of restricted gripping ability
- Simple and easy to put on
- Attractive, functional design

PRODUCT-BENEFIT

- Knee support with adjustable thigh width for stabilisation of the knee joint
- Circular closure strap for improved stabilisation on the thigh
- One lateral spiral stay each side
- Functional knitted fabric, very comfortable to wear
- Anatomic silicone pad
- **Integrated plus size system**

COLOUR: ■ blue

THE GENERATION-BENEFIT

- **Easy to put on/remove** thanks to closed design on the thigh and open design on the shank
- Grip straps allow safe and energy-saving handling, also in case of restricted gripping ability
- Simple and easy to put on
- Attractive, functional design

INDIKATIONEN

Gonarthrosis, recurrent articular effusion, arthritis, slight ligament instability

INDICATIONS

Gonarthrosis, recurrent articular effusion, arthritis, slight ligament instability

SIZES



		d-measure	f-measure
Size 1 (small)	_____	26 – 29 cm	32 – 36 cm
Size 2 (medium)	_____	29 – 33 cm	36 – 40 cm
Size 3 (large)	_____	33 – 37 cm	40 – 48 cm
Size 4 (x-large)	_____	37 – 41 cm	44 – 52 cm
Size 5 (xx-large)	_____	41 – 45 cm	52 – 60 cm
Size 6 (xxx-large)	_____	45 – 49 cm	60 – 68 cm

right and left symmetrical

1 cm = 0,3937 inch

BORT StabuloGen® Eco
No. 114 520



Med

PRODUCT-BENEFIT

- Classical knee support for compression of soft tissues with silicone pad for stabilisation and relief of the knee joint
- **SOFTflex knitted fabric** for more elasticity, stability and very comfortable to wear
- One lateral spiral stay each side
- Plus-size system resulting in very good fitting
- Compression class II
- For children too, see Page 28


COLOURS:  skin-tone  silver
 black  denim

THE ECO-BENEFIT

- Velvety soft surface structure
- Highly elastic bending zone, also protects the popliteal fossa when worn for longer periods
- Silicone adhesive border from 5 size

INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, e. g. gonarthrose, arthritis, ligament instability, meniscopathy

SIZES	d-measure	f-measure		d-measure	f-measure
 children (blue)	23 – 26	28 – 32 cm	Size 4 (x-large)	37 – 41	44 – 48 cm
Size 1 (small)	26 – 29	32 – 36 cm	Size 4+ (x-large plus)	37 – 41	48 – 52 cm
Size 2 (medium)	29 – 33	36 – 40 cm	Size 5 (xx-large)	41 – 45	52 – 56 cm
Size 2+ (medium plus)	29 – 33	40 – 44 cm	Size 5+ (xx-large plus)	41 – 45	56 – 60 cm
Size 3 (large)	33 – 37	40 – 44 cm	Size 6 (xxx-large)	45 – 49	60 – 64 cm
Size 3+ (large plus)	33 – 37	44 – 48 cm	Size 6+ (xxx-large plus)	45 – 49	64 – 68 cm

right and left symmetrical
1 cm = 0,3937 inch

BORT StabiloGen® Eco with a Grip Top No. 114 520 H



Med

PRODUCT-BENEFIT

- Classical knee support for compression of soft tissues with silicone pad
- Stabilises and relieves the knee joint
- **With additional silicone grip top in all sizes**
- **SOFTflex knitted fabric** for more elasticity, stability and very comfortable to wear
- One lateral spiral stay each side
- Plus-size system resulting in very good fitting

COLOURS: ■ skin-tone ■ silver

THE ECO-BENEFIT

- Optimal adhesion and pressure distribution on the edge of the thigh
- Velvety soft surface structure
- Highly elastic bending zone, also protects the popliteal fossa when worn for longer periods

INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, gonarthrosis, arthritis, ligament instability, meniscopathy

BORT StabiloGen® Latex-Free No. 114 250



Med

PRODUCT-BENEFIT

- Latex-free knee support for compression of soft tissues
- Stabilisation and relief of the knee joint
- **SOFTflex knitted fabric** for more elasticity, stability and very comfortable to wear
- Highly elastic, soft flexion region in the hollow of the knee
- Silicone pad
- One spiral stay on each side
- Compression class II

COLOURS: ■ black ■ silver

THE MED-BENEFIT

- **Latex-free knitted fabric**
- Velvety soft surface structure
- Plus-size system resulting in very good fitting
- Silicone adhesive border from 5 size

INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, e. g. arthritis, gonarthrose, ligament instability

SIZES

d-measure f-measure

d-measure f-measure



Size 1 (small)	26 – 29	32 – 36 cm
Size 2 (medium)	29 – 33	36 – 40 cm
Size 2+ (medium plus)	29 – 33	40 – 44 cm
Size 3 (large)	33 – 37	40 – 44 cm
Size 3+ (large plus)	33 – 37	44 – 48 cm
Size 4 (x-large)	37 – 41	44 – 48 cm

Size 4+ (x-large plus)	37 – 41	48 – 52 cm
Size 5 (xx-large)	41 – 45	52 – 56 cm
Size 5+ (xx-large plus)	41 – 45	56 – 60 cm
Size 6 (xxx-large)	45 – 49	60 – 64 cm
Size 6+ (xxx-large plus)	45 – 49	64 – 68 cm

right and left symmetrical

1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Circular knitted knee support for compression of soft tissues, stabilisation and relief of the knee joint
- Anatomically contoured knit
- Silicone pad
- One lateral spiral stay each side
- Compression class II

COLOURS: skin-tone
 mineral grey

THE MED-BENEFIT

- Seamless round knitted fabric
- Plus-size system resulting in very good fitting



PRODUCT-BENEFIT

- Knee support with patella recess for stabilisation and relief of the knee joint
- Highly elastic, soft flexion region in the hollow of the knee
- Silicone pad
- One lateral spiral stay each side
- Silicone adhesive border from 5 size

COLOURS: skin-tone silver

THE MED-BENEFIT

- **SOFTflex knitted fabric** for more elasticity, stability and very comfortable to wear
- Complete relief of the pressure on patella
- Plus-size system resulting in very good fitting

INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, e. g. arthritis, gonarthrose, ligament instability

INDICATIONS

Patellar chondropathy, meniscopathy, arthritis, bursitis

SIZES

d-measure f-measure

d-measure f-measure



Size 1 (small)	26 – 29	32 – 36 cm
Size 2 (medium)	29 – 33	36 – 40 cm
Size 2+ (medium plus)	29 – 33	40 – 44 cm
Size 3 (large)	33 – 37	40 – 44 cm
Size 3+ (large plus)	33 – 37	44 – 48 cm
Size 4 (x-large)	37 – 41	44 – 48 cm

Size 4+ (x-large plus)	37 – 41	48 – 52 cm
Size 5 (xx-large)	41 – 45	52 – 56 cm
Size 5+ (xx-large plus)	41 – 45	56 – 60 cm
Size 6 (xxx-large)	45 – 49	60 – 64 cm
Size 6+ (xxx-large plus)	45 – 49	64 – 68 cm

right and left symmetrical
1 cm = 0,3937 inch

BORT Osgood-Schlatter Knee Support

No. 114 510

**PRODUCT-BENEFIT**

- Knee support for local decompression of tibial tuberosity
- **SOFTflex knitted fabric** for more elasticity, stability and very comfortable to wear
- Anatomically contoured knit
- Multi-zone silicone pad
- One lateral spiral stay on each side
- **For children too**, see Page 27

THE MED-BENEFIT

- Velvety soft surface structure
- Plus-size system resulting in very good fitting
- Adjustable setting using V-shaped strap correction

COLOUR: ■ silver

INDICATIONS

Facilitation of the movement of the patella, **Osgood-Schlatter disease**, tendinosis of the tibial tuberosity ligament, gonarthrosis, arthritis, articular effusion, irritations

SIZES

d-measure f-measure

d-measure f-measure

**children (blue)**

Size 1 (small)

Size 2 (medium)

Size 2+ (medium plus)

Size 3 (large)

Size 3+ (large plus)

23 – 26 28 – 32 cm

26 – 29 32 – 36 cm

29 – 33 36 – 40 cm

29 – 33 40 – 44 cm

33 – 37 40 – 44 cm

33 – 37 44 – 48 cm

Size 4 (x-large)

Size 4+ (x-large plus)

Size 5 (xx-large)

Size 5+ (xx-large plus)

37 – 41 44 – 48 cm

37 – 41 48 – 52 cm

41 – 45 52 – 56 cm

41 – 45 56 – 60 cm

right and left symmetrical

1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Knee support for compression of soft tissues with V-shaped stabilisation bands for stabilisation and relief of the knee joint
- Silicone pad
- Silicone non-slip border in the proximal and distal area
- Long femoral part (22 cm, 8.7 inches)
- Length: 37 cm (14.6 inches)

COLOUR: ■ silver

THE MED-BENEFIT

- No lateral splints, there is no pressure in this area
- Relax zone in the knee bend
- Plus-size system resulting in very good fitting

INDICATIONS

Chronic post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, arthritis, arthrosis



PRODUCT-BENEFIT

- Active thigh support with circular straps and cross restraints
- Anterior/posterior application is possible
- Anatomically contoured knit
- Relief zones in the injury region
- Length: 23 cm (9.1 inches)


COLOUR: ■ black/green

THE SPORT-BENEFIT


- COOLMAX®: wicks moisture away from the body
- Its operating principle consists of compression, conventional tape and Kinesio tape
- Proprioceptive effect (muscular detonisation, pain reduction)

INDICATIONS

Preventive: For highly strained muscle groups, **Acute:** Swollen contusions with haematoma, muscle injury, **Chronic:** following muscle injury, protection against new muscle lesions, reduction of the burden on damaged muscle structures (e.g. scarring after old injuries)

SIZES		d-measure	f-measure
	Size 1 (small)	26-29	32-36 cm
	Size 2 (medium)	29-33	36-40 cm
	Size 2+ (medium plus)	29-33	40-44 cm
	Size 3 (large)	33-37	40-44 cm
	Size 3+ (large plus)	33-37	44-48 cm
	Size 4 (x-large)	37-41	44-48 cm
	Size 4+ (x-large plus)	37-41	48-52 cm
	Size 5 (xx-large)	41-45	52-56 cm
	Size 5+ (xx-large plus)	41-45	56-60 cm
	Size 6 (xxx-large)	45-49	60-64 cm
	Size 6+ (xxx-large plus)	45-49	64-68 cm


right and left symmetrical
1 cm = 0,3937 inch

SIZES		Circumference of thigh
	small	_____ up to 43 cm
	medium	_____ 43 - 50 cm
	large	_____ 50 - 56 cm
	x-large	_____ 56 - 61 cm
	xx-large	_____ 61 - 66 cm

right and left symmetrical
1 cm = 0,3937 inch

BORT Patella Tendon Strap Sport

No. 114 010 SP



PRODUCT-BENEFIT

- So-called Kassel Patellar Tendon Strap for target relief of the patella tendon
- Sports version **without a guide loop**
- Optimised construction for high activity levels
- Soft, elastic strap material

COLOUR: ■ black/green

THE SPORT-BENEFIT


- Also suitable for contact sport types
- Comfortable to wear during physical activity
- Individually controllable pad pressure
- Anti-slip positioning of the silicone pad via direct skin contact
- Modern, sporty design

INDICATIONS

Patellar chondropathy, femoropatellar arthrosis, patellar tendonitis

BORT Patella Tendon Strap

No. 114 010



PRODUCT-BENEFIT

- So-called Kassel Patellar Tendon Strap for target relief of the patella tendon
- Short stretch elastic material
- Anatomical silicone proprioceptive pad
- Soft bordering
- Reversible strap with Velcro fastening


COLOUR: ■ black


THE MED-BENEFIT

- Pad pressure adjustable without graduation
- Improved joint coordination
- Pain relief
- No pressure in the hollow of the knee

INDICATIONS


Patellar chondropathy, patellar tendonitis, femoropatellar arthrosis

SIZES		Circumference of knee joint below patella	
	Size 1	_____	up to 34 cm
	Size 2	_____	34 – 41 cm
	Size 3	_____	41 – 48 cm

SIZES		Circumference of knee joint below patella	
	Size 1	_____	up to 34 cm
	Size 2	_____	34 – 41 cm
	Size 3	_____	41 – 48 cm

right and left symmetrical
1 cm = 0,3937 inch

right and left symmetrical
1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Knee support with side joint splints, which can be fully opened
- Four Velcro fastenings
- Ideal in case of restrictions of movement and injuries that do not allow it to be slipped on
- Non-slip fit

COLOUR:  silver

THE MED-BENEFIT

- The support can be opened completely thanks to its Velcro fastenings
- **Support can be put on/removed easily (ideal for use in the elderly)**



PRODUCT-BENEFIT

- Knee support with adjustable thigh-width for guidance and stabilisation of the knee joint
- With two lateral spiral stays on each side and lateral silicone pads
- Patella recess
- Two Velcro fastening straps for optional adjustment of upper thigh width

COLOUR:  skin-tone

THE MED-BENEFIT


- Particularly suitable for patients with broader thighs
- Support can be put on/removed easily
- Non-slip fit

INDICATIONS


Mild to moderate lateral ligament instability of the knee joint, arthrosis, arthritis

INDICATIONS

Recurrent articular effusion, arthrosis, arthritis, ligament instability

SIZES	Circumference of knee	
	Size 1 (small)	up to 34 cm
	Size 2 (medium)	34 – 37 cm
	Size 3 (large)	37 – 40 cm
	Size 4 (x-large)	40 – 43 cm
	Size 5 (xx-large)	43 – 46 cm
	Size 6 (xxx-large)	46 – 49 cm

right and left symmetrical
1 cm = 0,3937 inch

SIZES	Circumference of knee	
	Size 1	34 – 37 cm
	Size 2	37 – 40 cm
	Size 3	40 – 43 cm
	Size 4	43 – 46 cm
	Size 5	46 – 49 cm
	Size 6	49 – 52 cm
	Size 7	52 – 55 cm

right and left symmetrical
1 cm = 0,3937 inch

BORT Asymmetric®

No. 114 700

**PRODUCT-BENEFIT**

- Knee support to facilitate movement of the patella
- One lateral spiral stay on each side
- DynamicFlex functional knitted fabric:
 - Relax zone in the knee bend
 - Lengthwise-elastic stretch zone in the lateral area of the silicone pad improves function and fit at all bending angles
- Compression class II

THE MED-BENEFIT

- Individual guidance of the patella in the medial and distal direction by means of elastic straps, supported by lateral silicone half-ring
- Plus size system, resulting in very good fitting
- Silicone adhesive border for a secure fit

COLOUR: ■ silver

INDICATIONS

Patellofemoral pain syndrome, high riding patella, patellar lateralisation, condition after lateral release surgery, recurrent or chronic patellar instability, meniscopathy

SIZES

d-measure f-measure

d-measure f-measure



Size 1 (small)	26 – 29	32 – 36 cm
Size 2 (medium)	29 – 33	36 – 40 cm
Size 2+ (medium plus)	29 – 33	40 – 44 cm
Size 3 (large)	33 – 37	40 – 44 cm
Size 3+ (large plus)	33 – 37	44 – 48 cm
Size 4 (x-large)	37 – 41	44 – 48 cm
Size 4+ (x-large plus)	37 – 41	48 – 52 cm

Size 5 (xx-large)	41 – 45	52 – 56 cm
Size 5+ (xx-large plus)	41 – 45	56 – 60 cm
Size 6 (xxx-large)	45 – 49	60 – 64 cm
Size 6+ (xxx-large plus)	45 – 49	64 – 68 cm

please indicate left or right

1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Soft brace made of textile with joints to facilitate movement of the patella
- Lateral silicone half ring
- Side joint splints
- Circular fastening straps on the thigh and lower leg
- Compression class II

COLOUR: ■ silver

THE MED-BENEFIT

- Vector correction filaments: Innovative strap system which enables individual gradual guidance of the patella
- Individually adjustable joint splints
- Plus size system, resulting in very good fitting



PRODUCT-BENEFIT

- Knee support with side joint splints for guidance of the knee joint
- Extension limit at 0°
- Circular fastening straps on the thigh and lower leg
- Pressure-reduced zones in the edge and popliteal fossa areas
- Plus-size system resulting in very good fitting
- For children too, see Page 28

COLOUR: ■ silver

THE MED-BENEFIT

- Easily anatomically adjustable joint splints with effective side stabilisation
- Slim joint splint design

INDICATIONS

High riding patella, patellar lateralisation, patellofemoral pain syndrome, patellar luxation or subluxation combined with lateral ligament instability, post-operative

INDICATIONS

Mild to moderate lateral ligament instability of the knee joint, gonarthrosis, arthritis

SIZES

d-measure f-measure

d-measure f-measure

	*Children	23 – 26	28 – 32 cm	Size 4 (x-large)	37 – 41	44 – 48 cm
	Size 1 (small)	26 – 29	32 – 36 cm	Size 4+ (x-large plus)	37 – 41	48 – 52 cm
	Size 2 (medium)	29 – 33	36 – 40 cm	Size 5 (xx-large)	41 – 45	52 – 56 cm
	Size 2+ (medium plus)	29 – 33	40 – 44 cm	Size 5+ (xx-large plus)	41 – 45	56 – 60 cm
	Size 3 (large)	33 – 37	40 – 44 cm	Size 6 (xxx-large)	45 – 49	60 – 64 cm
	Size 3+ (large plus)	33 – 37	44 – 48 cm	Size 6+ (xxx-large plus)	45 – 49	64 – 68 cm

please indicate left or right
1 cm = 0,3937 inch

*Children size only available for No. 150 140

BORT Knee Support with Adjustable Hinge

No. 114 460



PRODUCT-BENEFIT

- Knee support with side joint splints for guidance of the knee joint
- **Extension limit at 0° and 10°**
- **Flexion limitation at 30°, 60° and 90°**
- Circular fastening straps on the thigh and lower leg
- Plus-size system resulting in very good fitting

THE MED-BENEFIT

- Adjustable joint for limiting extension and flexion
- Easily anatomically adjustable joint splints with effective side stabilisation
- Slim joint splint design

COLOUR: ■ silver

INDICATIONS

Meniscus injuries, mild to moderate lateral ligament instability of the knee joint, gonarthrosis, arthritis, bursitis

SIZES



	d-measure	f-measure		d-measure	f-measure
Size 1 (small)	26 – 29	32 – 36 cm	Size 4+ (x-large plus)	37 – 41	48 – 52 cm
Size 2 (medium)	29 – 33	36 – 40 cm	Size 5 (xx-large)	41 – 45	52 – 56 cm
Size 2+ (medium plus)	29 – 33	40 – 44 cm	Size 5+ (xx-large plus)	41 – 45	56 – 60 cm
Size 3 (large)	33 – 37	40 – 44 cm	Size 6 (xxx-large)	45 – 49	60 – 64 cm
Size 3+ (large plus)	33 – 37	44 – 48 cm	Size 6+ (xxx-large plus)	45 – 49	64 – 68 cm
Size 4 (x-large)	37 – 41	44 – 48 cm			

please indicate left or right

1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Soft brace with limitation of mobility for guidance and relief of the knee joint
- With open patella
- Polycentric aluminium sidebars, medial and lateral stabilisation
- Adjustable extension: 0°, 10°, 15°, 20°, 30°, 40°, 45°
- Adjustable flexion: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°
- Breathable, elastic material for optimal comfort

COLOUR: ■ grey

THE MED-BENEFIT

- With two circular Velcro straps for additional support
- Plus-size system resulting in very good fitting

INDICATIONS

Mild to moderate ligament instability of the knee joint, meniscus injuries, gonarthrosis



PRODUCT-BENEFIT

- Soft brace with movement limitation for guidance and relief of the knee joint
- Polycentric aluminium sidebars, medial and lateral stabilisation
- Adjustable extension: 0°, 10°, 15°, 20°, 30°, 40°, 45°
- Adjustable flexion: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°
- Breathable, elastic material for optimal comfort

COLOUR: ■ grey

THE MED-BENEFIT

- With four circular Velcro straps for particularly secure support
- **It can be fully opened** so is very easy to put on
- Plus-size system resulting in very good fitting


INDICATIONS

Moderate ligament instability of the knee joint, meniscus injuries, after-care in arthroscopic procedures involving the knee, chronic, postoperative and post-traumatic inflammations

SIZES

d-measure f-measure

d-measure f-measure



Size 1 (small)	26 – 29	32 – 36 cm
Size 2 (medium)	29 – 33	36 – 40 cm
Size 2+ (medium plus)	29 – 33	40 – 44 cm
Size 3 (large)	33 – 37	40 – 44 cm
Size 3+ (large plus)	33 – 37	44 – 48 cm
Size 4 (x-large)	37 – 41	44 – 48 cm
Size 4+ (x-large plus)	37 – 41	48 – 52 cm

Size 5 (xx-large)	41 – 45	52 – 56 cm
Size 5+ (xx-large plus)	41 – 45	56 – 60 cm
Size 6 (xxx-large)	45 – 49	60 – 64 cm
Size 6+ (xxx-large plus)	45 – 49	64 – 68 cm

right and left symmetrical

1 cm = 0,3937 inch

BORT Generation Knee Brace

No. 215 400

NEW
DESIGN

Generation

**PRODUCT-BENEFIT**

- Soft brace with movement limitation for guidance and relief of the knee joint
- Polycentric aluminium sidebars, medial and lateral stabilisation
- Adjustable extension: 0°, 10°, 15°, 20°, 30°, 40°, 45°
- Adjustable flexion: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°
- Soft, breathable material, making it comfortable to wear

COLOUR: ■ blue

THE GENERATION-BENEFIT

- Can be fully opened
- Extremely easy handling and putting on the brace
- Wide grip straps allow safe and energy-sparing handling also if gripping ability is restricted
- Popliteal opening facilitates flexion when sitting
- Extra-wide guide loops
- Physiologically optimised slim joint splints
- Attractive, functional design



BORT Generation:
For better mobility for the aged
Product range with special handling advantages

INDICATIONS

Degenerative diseases, e. g. mild to moderate ligament instability of the knee joint, gonarthrosis, meniscus injuries

SIZES

d-measure f-measure

d-measure f-measure



Size 1 (small)	26 – 29	32 – 36 cm
Size 2 (medium)	29 – 33	36 – 40 cm
Size 2+ (medium plus)	29 – 33	40 – 44 cm
Size 3 (large)	33 – 37	40 – 44 cm
Size 3+ (large plus)	33 – 37	44 – 48 cm
Size 4 (x-large)	37 – 41	44 – 48 cm
Size 4+ (x-large plus)	37 – 41	48 – 52 cm

Size 5 (xx-large)	41 – 45	52 – 56 cm
Size 5+ (xx-large plus)	41 – 45	56 – 60 cm
Size 6 (xxx-large)	45 – 49	60 – 64 cm
Size 6+ (xxx-large plus)	45 – 49	64 – 68 cm
Size 7 (xxxx-large)	49 – 53	68 – 72 cm
Size 7+ (xxxx-large plus)	49 – 53	72 – 76 cm

right and left symmetrical

1 cm = 0,3937 inch

BORT OA-Xpress
No. 101 600



PRODUCT-BENEFIT

- Knee relief brace for isolated medial / lateral gonarthrosis up to 5° axial deviation
- Varus / valgus correction freely adjustable
- Adjustable, polycentric joint
- Limitation of extension: 0°, 5°, 10°, 15°, 20° and 25°
- Limitation of flexion: 0° (180°), 30°, 45°, 60°, 75° and 90°
- Brace length: 40 cm

COLOUR: ■ grey

THE MED-BENEFIT

- Partially elastic shells for good fitting to different leg circumferences
- Targeted relief of the medial / lateral compartment
- Lightweight and slim structure for ideal wearing comfort

Application of model A

- Right leg lateral:
Varus gonarthrosis (for valgisation in case of bowleg position)
- Left leg medial:
Valgus gonarthrosis (for varisation in case of knock knee position)

Application of model B

- Right leg medial:
Valgus gonarthrosis (for varisation in case of knock knee position)
- Left leg lateral:
Varus gonarthrosis (for valgisation in case of bowleg position)

INDICATIONS

Partial relief of the medial or lateral compartment in case of level II-III gonarthrosis, Bone bruise (post-traumatic bone marrow edema), post-operative care following tibial plateau fracture

IMPLEMENTATION

Model A _____ varus right/valgus left
Model B _____ valgus right/varus left

One size only

BORT GenuXpress

No. 100 700



Med

PRODUCT-BENEFIT

- Frame knee orthosis by the 4-point principle for guidance and stabilisation of the knee joint
- Adjustable, polycentric joint
- Limitation of extension: 0°, 5°, 10°, 15°, 20° and 25°
- Limitation of flexion: 0° (180°), 30°, 45°, 60°, 75° and 90°

THE MED-BENEFIT

- High frame strength combined with low self-weight thanks to ultra-modern material
- Very good adaptability of brace frame
- Lightweight and slim structure for ideal wearing comfort
- Easy to handle

COLOUR: ■ grey

Knee

INDICATIONS

Serious and/or complex instability of the knee joint, meniscus injuries, genu recurvatum, functional pre-and/or post-operative treatment of ruptured ligaments (ACL, PCL, collateral ligaments)

SIZES

f-measure

Width of knee




small	_____	39 – 47 cm	_____	8.5 – 9.5 cm
medium	_____	47 – 54 cm	_____	9.5 – 10.5 cm
large	_____	54 – 61 cm	_____	10.5 – 11.5 cm
x-large	_____	61 – 68 cm	_____	11.5 – 12.5 cm

After being processed by the manufactures it can be used again.

please indicate left or right

1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Knee brace for immobilisation of the knee joint
- **Adjustable**, 0° or 15° flexion
- With open patella, 1-piece
- Median and lateral aluminium splints
- Length: 55 cm (21.7 inches)

COLOUR: ■ blue/black

THE MED-BENEFIT

- Easy to handle
- Lined with terry cloth on the inside for excellent comfort
- The aluminium splints can be removed for washing

BORT Immob Splint Single-Piece
No. 145 070

- Knee brace for immobilisation of the knee joint 0° or 15° flexion



PRODUCT-BENEFIT

- Knee brace for immobilisation of the knee joint with a 0° - 20° flexion
- 3-piece with open patella
- Median and lateral aluminium splints, which can be freely positioned
- **For children too**, see Page 29

COLOUR: ■ black

THE MED-BENEFIT

- The support adapts to body shape thanks to the 3-piece design
- Lined with terry cloth on the inside for excellent comfort
- The aluminium splints can be removed for washing

INDICATIONS

Pre-operative, post-operative, post-traumatic, joint, muscle or ligament injuries

SIZE

One length only, 55 cm

For No. 145 070
please indicate 0° or 15° flexion

right and left symmetrical
1 cm = 0,3937 inch

SPLINT LENGTH

Fits height

Children (38 cm)	80 – 145 cm
45 cm	145 – 155 cm
50 cm	155 – 165 cm
55 cm	165 – 180 cm
60 cm	over 180 cm
Please indicate 0° or 20° flexion	

right and left symmetrical
1 cm = 0,3937 inch

BORT Knee Support No. 104 400



Med

PRODUCT-BENEFIT

- Simple knee support for easy stabilisation
- Anatomically contoured knit
- Compression class II

COLOURS:  skin-tone  black

THE MED-BENEFIT

- Easy to use
- Not bulky
- Short design

BORT Dual-Tension Knee Support No. 054 400



Med

PRODUCT-BENEFIT

- Simple round knitted knee support for easy stabilisation
- Compression class II

COLOURS:  skin-tone  blue

THE MED-BENEFIT


- Easy to handle
- Not bulky

INDICATIONS

Sports injuries, articular effusion, sprains, contusions

SIZES

Circumference of knee

	small	_____	up to 32 cm
	medium	_____	32 – 37 cm
	large	_____	37 – 42 cm
	x-large	_____	42 – 47 cm


right and left symmetrical
1 cm = 0,3937 inch

INDICATIONS

Sports injuries, articular effusion, sprains, contusions

SIZES

Circumference of knee

	small	_____	up to 31 cm
	medium	_____	31 – 37 cm
	large	_____	37 – 42 cm
	x-large	_____	42 – 48 cm
	xx-large	_____	48 – 54 cm

right and left symmetrical
1 cm = 0,3937 inch

Sure footed: with all-round protection and stability at every step.



Application area	Upper ankle joint							Lower ankle joint	Foot				
Overview of indications for selected products	select® TaloStabil® 054 700	select® Achillo Stabil® Plus 053 100	TaloStabil® Sport 054 600 SP	Achillo Stabil® Plus Sport 053 100 SP	Generation Ankle Support 215 600	TaloStabil® 054 600	Achillo Stabil® Eco 054 900	Stabilo® Ankle Support 114 100	TaloFX® 054 800 Talo Xpress 053 900	Helix S <i>Spiraldynamik</i> 053 200	Malleo Stabil® Brace 100 500	Malleo Xpress Active 100 580 Malleo Xpress 100 560	Foot Levator Orthotic Device 054 300
	P. 152/153	P. 152	P. 154/155	P. 154	P. 156	P. 157	P. 157	P. 160	P. 161/164	P. 159	P. 165	P. 163/162	P. 166
	●		●		●	●		●	●		●	●	
	●		●		●	●		●	●		●	●	
		●		●			●						
	●		●		●	●					●	●	
								●					
													●
													●
	Foot problems due to static anomalies and defects in the shape of the foot due to dysfunction of the lower ankle joint, e. g. skew and splayfoot										●		



PRODUCT-BENEFIT

- High quality active ankle support
- Highly elastic stretch zone in the sensitive instep area simplifies putting on and removal
- Two structured silicone pads for intermittent massage

COLOURS: skin-tone silver
 black

THE SELECT-BENEFIT

- **SOFTflex Plus knit fabric**
 - Soft-velvety surface structure for top wearing comfort
- Extra soft, snag proof material edge at the calf
- Separate slip-on aid included

INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, arthrosis, ligament weakness, after removal of plaster cast, ligament injuries, distortion of upper ankle joint, tendomyopathies

SIZES Circumference of ankle

	small	_____	up to 21 cm
	medium	_____	21 – 23 cm
	large	_____	23 – 25 cm
	x-large	_____	25 – 27 cm

please indicate left or right
1 cm = 0,3937 inch



PRODUCT-BENEFIT

- High-quality **active Achilles tendon support** for compression of soft tissue
- Highly elastic stretch zone in the instep area
- Extended, structured silicone pad for intermittent massage of the Achilles tendon and the musculotendon transitional area towards the calf muscle

COLOURS: silver black

THE SELECT-BENEFIT

- Its half-open structure makes it easier to put on and take off
- Easy and simple fixing of the support using three Velcro fasteners
- Heel wedge to relieve the Achilles tendon during the acute phase

INDICATIONS

Chronic or post-traumatic irritation or tenderness in the Achilles tendon area, Haglund syndrome heel irritation, achilodynia, irritation and pain of the Achilles tendon and surrounding areas, after surgery, e. g. Achilles tendon rupture, bursitis subachillea

SIZES Circumference of ankle

	small	_____	up to 21 cm
	medium	_____	21 – 23 cm
	large	_____	23 – 25 cm
	x-large	_____	25 – 27 cm

right and left symmetrical
1 cm = 0,3937 inch

BORT select®
TaloStabil® Plus
 No. 054 100



BREEZE

PRODUCT-BENEFIT

- High quality active ankle support with Stabilo® Strap, partially elastic, it can be freely positioned
- Highly elastic stretch zone in the sensitive instep area simplifies putting on and removal
- Treatment-supporting positioning of the foot in the pronation / supination / neutral position
- Two structured silicone pads for intermittent massage

COLOURS: skin-tone silver
 black

THE SELECT-BENEFIT

- **SOFTflex Plus knit fabric:**
 - BREEZE: Luxury fibre with a cooling effect
 - Soft-velvety surface structure for top wearing comfort
- Extra soft, snag proof material edge at the cal
- Separate slip-on aid included

INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, arthrosis, ligament weakness, after removal of plaster cast, ligament injuries, distortion of upper ankle joint, supination prophylaxis

Ankle

SIZES

Circumference of ankle

	small	_____	up to 21 cm
	medium	_____	21 – 23 cm
	large	_____	23 – 25 cm
	x-large	_____	25 – 27 cm

please indicate left or right

1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Ankle support for compression of soft tissues
- **SOFTflex 3D knitted fabric** with a velvety soft surface structure for more elasticity, stability and very comfortable to wear
- Two silicone pads for intermittent massage
- Compression class II


COLOUR: ■ black/green

THE SPORT-BENEFIT

- **COOLMAX®**: wicks moisture away from the body, for dryness and comfort
- Ideal fit for sports activities
- **Latex-free knitted fabric**
- Modern, sporty design (with a reflecting label)

INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, sprain, ligament weakness, after removal of plaster cast, ligament injuries, distortion of upper ankle joint, arthrosis

SIZES	Circumference of ankle	
	small	_____ up to 21 cm
	medium	_____ 21 – 23 cm
	large	_____ 23 – 25 cm
	x-large	_____ 25 – 27 cm

right and left symmetrical
1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Active Achilles tendon support
- **SOFTflex 3D knitted fabric** with a velvety soft surface structure very comfortable to wear
- Extended structured silicone pads for intermittent massage of the Achilles tendon
- Extra soft, snag proof material edge at the calf


COLOUR: ■ black/green

THE SPORT-BENEFIT

- **COOLMAX®**: wicks moisture away from the body
- Heel wedge to relieve the Achilles tendon during the acute phase
- Integrated comfort zone thanks to Dynamic-Flex-fabric
- Separate slip-on aid included
- **Latex-free knitted fabric**

INDICATIONS

Chronic or post-traumatic irritation or tenderness in the Achilles tendon area, in the case of Haglund syndrome heel irritation, achillodynia, after surgery, e. g. Achilles tendon rupture, bursitis subachillea

SIZES	Circumference of ankle	
	small	_____ up to 21 cm
	medium	_____ 21 – 23 cm
	large	_____ 23 – 25 cm
	x-large	_____ 25 – 27 cm

right and left symmetrical
1 cm = 0,3937 inch

BORT TaloStabil® Plus Sport

No. 054 100 SP



COOLMAX®

PRODUCT-BENEFIT

- Active ankle support for compression of soft tissues, with Stabilo® Strap, partially elastic can be freely positioned
- **SOFTflex 3D knitted fabric** with a velvety soft surface structure for more elasticity, stability and very comfortable to wear
- Two silicone pads for intermittent massage
- Compression class II

COLOUR: ■ black/green

THE SPORT-BENEFIT

- **COOLMAX®**: wicks moisture away from the body, for dryness and comfort
- Treatment-supporting positioning of the foot in the pronation/supination/neutral position
- Ideal fit for sports activities
- **Latex-free knitted fabric**
- Modern, sporty design (with a reflecting label)

INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, arthrosis, ligament weakness, after removal of plaster cast, ligament injuries, distortion of upper ankle joint, supination prophylaxis

Ankle

SIZES

Circumference of ankle



small	_____	up to 21 cm
medium	_____	21 – 23 cm
large	_____	23 – 25 cm
x-large	_____	25 – 27 cm

please indicate left or right

1 cm = 0,3937 inch

BORT Generation Ankle Support
No. 215 600



Generation

PRODUCT-BENEFIT

- High-quality active ankle support for compression of soft tissue
- Highly elastic stretch zone in the sensitive instep area simplifies putting on and removal
- Functional knitted fabric, very comfortable to wear
- Two structured silicone pads for intermittent massage

COLOUR: ■ blue

THE GENERATION-BENEFIT

- Easy to put on/remove thanks to open design, also with restricted movement
- Easy and simple fixing of the support using two Velcro fasteners
- Attractive, functional design



BORT Generation:
For better mobility for the aged
Product range with special handling advantages

INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, sprain, ligament weakness, after removal of plaster cast, ligament injuries, arthrosis

SIZES

Circumference of ankle



small	_____	up to 21 cm
medium	_____	21 – 23 cm
large	_____	23 – 25 cm
x-large	_____	25 – 27 cm

right and left symmetrical
1 cm = 0,3937 inch

BORT TaloSabil® Eco

No. 054 600



PRODUCT-BENEFIT

- Active ankle support for compression of soft tissues
- **SOFTflex knitted fabric** for more elasticity, stability and very comfortable to wear
- Two silicone pads for intermittent massage
- Compression class II
- **Also for children**, see Page 30

COLOURS: skin-tone silver
 black

THE ECO-BENEFIT

- Velvety soft surface structure
- Pressure relief at the base of the ankle

INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, sprain, ligament weakness, after removal of plaster cast, ligament injuries, arthrosis

SIZES Circumference of ankle

	Children (blue)	17 – 19 cm
	small	up to 21 cm
	medium	21 – 23 cm
	large	23 – 25 cm
	x-large	25 – 27 cm
	xx-large	27 – 29 cm

right and left symmetrical
 1 cm = 0,3937 inch

BORT AchilloSabil® Eco

No. 054 900



PRODUCT-BENEFIT

- Achilles tendon support
- **SOFTflex knitted fabric** for more elasticity, stability and very comfortable to wear
- Viscoelastic pad alongside the Achilles tendon for intermittent massage
- Heel wedge to relieve the Achilles tendon during the acute phase
- Compression class II

COLOURS: skin-tone silver
 black

THE ECO-BENEFIT

- Velvety soft surface structure
- Reduced pressure in the ankle area

INDICATIONS

Chronic or post-traumatic irritation or tenderness in the Achilles tendon area, Haglund syndrome heel irritation, achillodynia, irritation and pain of the Achilles tendon and surrounding areas, after surgery, e. g. Achilles tendon rupture, bursitis subachillea

SIZES Circumference of ankle

	small	up to 21 cm
	medium	21 – 23 cm
	large	23 – 25 cm
	x-large	25 – 27 cm

right and left symmetrical
 1 cm = 0,3937 inch




PRODUCT-BENEFIT

- Anatomically knitted ankle joint support for soft tissue compression
- Two silicone pads for intermittent massage
- Microfiber knitted fabric, soft to the touch
- Compression class II

COLOURS: skin-tone
 mineral grey

THE MED-BENEFIT

- Reduced pressure in the ankle area
- Anatomical knitted heel for optimal fit



PRODUCT-BENEFIT

- Achilles tendon support
- Viscoelastic pad alongside the Achilles tendon for intermittent massage
- Heel wedge to relieve the Achilles tendon during the acute phase
- Compression class II

COLOURS: skin-tone
 mineral grey

THE MED-BENEFIT

- Reduced pressure in the ankle area
- Anatomical knitted heel for optimal fit

INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, sprain, ligament weakness, after removal of plaster cast, ligament injuries, arthrosis

INDICATIONS

Chronic or post-traumatic irritation or tenderness in the Achilles tendon area, Haglund syndrome heel irritation, achillodynia, irritation and pain of the Achilles tendon and surrounding areas, after surgery, e. g. Achilles tendon rupture, bursitis subachillea

SIZES



small	_____	up to 21 cm
medium	_____	21 – 23 cm
large	_____	23 – 25 cm
x-large	_____	25 – 27 cm

right and left symmetrical
1 cm = 0,3937 inch

Circumference of ankle

BORT Helix S Spiraldynamik® Lower Ankle Support

No. 053 200

**PRODUCT-BENEFIT**

- Ankle bandage to support torsion of the back of the foot and to raise the instep
- Ligament separation thanks to longitudinal cuts for an anatomical fit and optimal function
- Thin, elastic straps
- 3-D contoured knit
- Two silicone pads
- Extra soft, snag proof material edge at the calf
- Also for children, see page 29

COLOUR: ■ black

THE MED-BENEFIT

- Stabilisation of the ankle joint and activation of foot muscles through compression
- Functional supportive restraint in the case of classical foot deformities (i.e. skew and splayfoot)
- Optimal comfort even during intense activity
- Intermittent massage
- Simplified handling
- Separate slip-on aid included
- Pack: single item

INDICATIONS

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of lower ankle joint and the sole of the foot, e. g. arthrosis, ligament instability, foot problems due to static anomalies and defects in the shape of the foot due to dysfunction of the lower ankle joint e. g. skew and splayfoot

SIZES

Circumference of ankle

	x-small/children	17 – 19 cm
	small	19 – 21 cm
	medium	21 – 23 cm
	large	23 – 25 cm
	x-large	25 – 27 cm

please indicate left or right

1 cm = 0,3937 inch

BORT Stabilo® Ankle Support
No. 114 100

NEW
DESIGN



PRODUCT-BENEFIT

- Ankle brace **with laces** for stabilisation of the lower and upper ankle joints
- Safe handling due to variable hook and eyelet elements
- Individually adjustable stabilisation level
- Plastic U-shaped stabilisation element
- Non-elastic breathable material

THE MED-BENEFIT

- Stabilo® Strap for additional stabilisation and limitation of talus prolapse
- Can be worn with sports and walking shoes


COLOURS: ■ grey 

INDICATIONS

Post-operative protection after ligament suture/reconstruction, conservative treatment of ankle distortions and ligament ruptures, immobilisation in the case of inflammatory conditions, rheumatism, arthritis, tendinitis, ligament injuries

SIZES

Circumference of ankle

	x-small	_____	up to 19 cm
	small	_____	19 – 21 cm
	medium	_____	21 – 23 cm
	large	_____	23 – 25 cm
	x-large	_____	25 – 27 cm

right and left symmetrical
1 cm = 0,3937 inch

BORT TaloFX®

No. 054 800



Med

PRODUCT-BENEFIT

- Soft ankle brace for stabilisation of the lower and upper ankle joints
- Bilateral shells padded with soft memory foam
- Cross restraint featuring Dynamic-Force-Control (DFC) technology
- Elastic material that carries away moisture
- Also for children, see Page 31

THE MED-BENEFIT

- **Slim design**
- Can be worn without difficulty with sports and walking shoes
- Limitation of talus prolapse with optimal stabilisation
- For positioning the ankle in pronation or supination

COLOUR: ■ black

INDICATIONS

Post-operative protection after ligament suture, conservative treatment of ankle distortions and ligament ruptures, immobilisation in the case of inflammatory conditions, rheumatism, arthritis, tendinitis, ligament injuries, supination/pronation trauma

Ankle

SIZES

Circumference of ankle



Children

small

medium

large

x-large

17 – 19 cm

19 – 21 cm

21 – 23 cm

23 – 25 cm

25 – 27 cm

please indicate left or right

1 cm = 0,3937 inch



PRODUCT-BENEFIT

- Soft ankle joint brace for stabilisation of the ankle joint
- Modular structure, can be downgraded during the course of therapy
- Reduction of the talus prolapse
- Bilateral shells padded softly with memory foam
- **Additional stabilisation by retrievable soleplate in the acute phase**

THE MED-BENEFIT

- **Can be fully opened, therefore easy to put on**
- Can be worn without difficulty with sports and walking shoes
- Easy to handle and stabilise without graduation using the Boa fastening system

COLOUR: ■ black/bordeaux

INDICATIONS

Post-operative, post-traumatic, e. g. conservative treatment of ankle distortions (level II-III) and ligament ruptures, immobilisation in the case of inflammatory conditions, post-operative protection after ligament suture/reconstruction, chronic instabilities of the upper and lower ankle joints, acute and chronic capsular ligament instabilities

SIZES

Circumference of ankle



small	_____	19 – 21 cm
medium	_____	21 – 23 cm
large	_____	23 – 25 cm
x-large	_____	25 – 27 cm

please indicate left or right

1 cm = 0,3937 inch

BORT MalleoXpress Active

No. 100 580



PRODUCT-BENEFIT

- Soft ankle joint brace for stabilisation of the ankle joint
- Reduction of the talus prolapse by intelligent track control
- Bilateral shells padded softly with memory foam

COLOUR: black

THE MED-BENEFIT

- **Can be fully opened, therefore easy to put on**
- Easy, safe handling
- Easy to handle and stabilise without graduation using the Boa fastening system
- Can be worn without difficulty with sports and walking shoes

INDICATIONS

Post-operative protection after ligament suture, conservative treatment of ankle distortions and ligament ruptures, immobilisation in the case of inflammatory conditions, rheumatism, arthritis, tendinitis, ligament injuries, supination/pronation trauma

SIZES

Circumference of ankle



small	_____	19 – 21 cm
medium	_____	21 – 23 cm
large	_____	23 – 25 cm
x-large	_____	25 – 27 cm

please indicate left or right

1 cm = 0,3937 inch

BORT TaloXpress
No. 053 900



PRODUCT-BENEFIT

- Ankle joint brace for stabilisation of the upper and lower ankle joint
- Stabilisation degree adjustable using restraints
- Very thin and slim construction
- Can be worn in shoes without any problems
- Effective stabilisation of the ankle joint and very comfortable to wear

COLOUR: ■ grey

THE MED-BENEFIT

- Can be fully opened, therefore easy to put on and simplest handling
- High stabilisation combined with low self-weight thanks to ultra-modern material
- Lightweight and slim structure for good wearing comfort

INDICATIONS

Conservative treatment of ankle distortions and ligament ruptures, post-operative protection after ligament suture/reconstruction, immobilisation in the case of inflammatory conditions, rheumatism, arthritis, tendinitis, supination/pronation trauma

SIZES

h measure



small

medium

large

27 – 31 cm

31 – 34 cm

34 – 40 cm

Please indicate left or right

1 cm = 0,3937 inch

BORT MalleoStabil® Brace

No. 100 500



Med

PRODUCT-BENEFIT

- Combined Air/Gel brace and foam pad to stabilise the ankle, while limiting pronation and supination
- Anatomically shaped plastic shells
- Two sets of padded variants:
 - Pre-coolable Air/Gel pad to be used during the acute phase
 - Memory foam padding for comfortable wear during the rehabilitation phase

- Also for children, see Page 30

COLOURS: ☒ black ☐ white

THE MED-BENEFIT

- No impairment when walking
- Easy to put on thanks to circular Velcro fastening

Other product versions:

BORT MalleoStabil®
SOFT Brace No. 100 520

- With Memory foam padding only

BORT MalleoStabil®
Air/Gel Brace No. 100 540

- With Air/Gel padding only

INDICATIONS

Post-operative protection after ligament suture, conservative treatment of ankle distortions of any severity, functional treatment of ligament ruptures, ankle fractures (Weber A), chronic ligament instability, shortening of time spent in plaster, ligament injuries, distortion of upper ankle joint, fibular ligament rupture, pronation and/or supination traumas

SIZES

Fits height

Size 1/Children _____ up to 155 cm
Size 2 _____ over 155 cm

right and left symmetrical

1 cm = 0,3937 inch

BORT Foot Levator Orthotic Device

No. 054 300



PRODUCT-BENEFIT

- Dynamic textile Tibialis Anterior Muscle brace
- 8-strap system
- Non-elastic lateral strapping
- Supports flexion of the foot in the upper ankle joint
- Stabilisation of the lower ankle joint
- Pronation effect

COLOUR: ■ black

THE MED-BENEFIT

- Easy to put on thanks to labelling of the bracing elements
- Can be applied using one hand
- Two support lengths towards the forefoot for optimal flexion

INDICATIONS

Peroneal paralysis of any severity, especially for mild to medium paralysis of the tibialis anterior muscle and the flexor hallucis longus (degree of severity 2-4)

SIZES

Circumference of ankle




small	_____	19 – 21 cm
medium	_____	21 – 23 cm
large	_____	23 – 25 cm
x-large	_____	25 – 27 cm

short	_____	up to shoe size 6
long	_____	over shoe size 6

Indicate right/left, short/long

1 cm = 0,3937 inch

BORT Peroneal Splint No. 470 620



PRODUCT-BENEFIT

- Polypropylene tibialis anterior muscle brace to flex the foot during the swing phase
- $\frac{3}{4}$ sole up to the metacarpophalangeal joints of the big toe

COLOUR: ☐ white

THE MED-BENEFIT

- Thermoplastic reprocessing
- Can be used in closed shoes
- Highly flexible semi-shell in the calf area
- Velcro fastening for individual adjustment to suit different calf circumferences
- Light weight
- Easy to clean

INDICATIONS

Flaccid paralysis of the tibialis anterior muscle, reduced mobility of the patient, in the case of peroneal paralysis, lumbar nerve root damage, polio, stroke

SIZES	Shoe size UK
small	3 – 5.5
medium	6 – 8.5
large	9 – 13

please indicate left or right

BORT Peroneal Splint, long sole No. 470 630



PRODUCT-BENEFIT

- Polypropylene tibialis anterior muscle brace to flex the foot during the swing phase
- Long sole design

COLOUR: ☐ white

THE MED-BENEFIT

- Thermoplastic reprocessing
- Can be used in closed shoes
- Velcro fastening for individual adjustment to suit different calf circumferences
- Light weight
- Easy to clean

INDICATIONS

Flaccid paralysis of the tibialis anterior muscle, reduced mobility of the patient, in the case of peroneal paralysis, lumbar nerve root damage, polio, stroke

SIZES	Shoe size UK
small	2 – 4
medium	4 – 5.5
large	5.5 – 7
x-large	7 – 9

please indicate left or right



PRODUCT-BENEFIT

- Orthopaedic shoe, closed design
 - Anti-slip sole
 - Breathable materials
- Pack: Single item

COLOUR: ■ black

THE MED-BENEFIT

- The foot bedding and EVA rubber sole can be easily adjusted and re-adjusted to suit any individual
- Can be fully opened
- Easy to put on
- Comfortable to wear



PRODUCT-BENEFIT

- Orthopaedic shoe, open
 - Anti-slip sole
 - Breathable materials
- Pack: Single item

COLOUR: ■ black

THE MED-BENEFIT

- The foot bedding and EVA rubber sole can be easily adjusted and re-adjusted to suit any individual
- Can be fully opened
- Easy to put on
- Comfortable to wear

INDICATIONS

Post-operative and other wounds, even when treated with extensive and possibly padded bandaging, where the healing process takes place over an extended period of time and weight can be placed partially on the foot whilst walking, even outdoors. An additional insole is not required

	Shoe size UK
	1.5 – 2
	3 – 4
	4.5 – 5.5
	6 – 7
	7.5 – 8.5
	9 – 10
	11 – 11.5
right and left symmetrical	

	Shoe size UK
	1.5 – 2
	3 – 4
	4.5 – 5.5
	6 – 7
	7.5 – 8.5
	9 – 10
	11 – 11.5
right and left symmetrical	

BORT Therapeutic Shoe Comfort No. 930 320



PRODUCT-BENEFIT

- Orthopaedic shoe
- Extra-wide fit at the forefoot and toe area
- Zipper at the heel
- Flexible, anti-slip sole

- Pack: Pair or single item (No. 930 320 R or 930 320 L)

COLOUR: ■ black

THE MED-BENEFIT

- Secure hold and comfort for sensitive feet
- Width can be individually adjusted thanks to Velcro fastening
- Can be fully opened, comfortable walking

INDICATIONS

Post-operative and other wounds, also being treated with large and possibly padded bandages whose healing process take place over an extended period of time and weight can be placed partially on the foot whilst walking, even outdoors. An additional insole is not required

BORT Diabetic Orthopaedic Shoe No. 930 400 R/930 400 L



PRODUCT-BENEFIT

- Orthopaedic shoe, open with an optionally removable forefoot/toe cap
- Removable soft pad cushioning with a honeycomb structure
- Reinforced anti-slip sole with foot rolling motion aid

- Pack: Single item

COLOUR: ■ black

THE MED-BENEFIT

- Soft pad cushioning with individually removable honeycomb elements for local pressure relief
- Additional protection by the removable forefoot/toe cap
- No inner seams

INDICATIONS

Wounds in the foot and ankle region due to diabetes mellitus, peripheral artery occlusive disease (PAOD), further wounds in the foot and ankle region, post-traumatic, post-operative

Shoe size UK

	3	7.5
	4	8.5
	4.5	9
	5.5	10
	6	11
	7	

Supplied as pairs or single item

please indicate left or right

Shoe size UK

	2 – 3
	3.5 – 4.5
	5 – 6
	6.5 – 7.5
	8 – 9
	9.5 – 10.5
	11 – 12



PRODUCT-BENEFIT

- Forefoot relief shoe
- Anti-slip sole
- With additional protective toecap
- Breathable materials

- Pack: Single item

COLOUR: ■ black

THE MED-BENEFIT

- Redistribution of pressure towards the middle or back of the foot (heel)
- The foot bedding and EVA rubber sole can be easily adjusted and re-adjusted to suit any individual



PRODUCT-BENEFIT

- Heel relief shoe
- Anti-slip sole
- Concave heel edge
- Breathable materials

- Pack: Single item

COLOUR: ■ black

THE MED-BENEFIT

- Redistribution of pressure towards the middle or forefoot
- Improved stability during walking
- The foot bedding and EVA rubber sole can be easily adjusted and re-adjusted to suit any individual

INDICATIONS

Post-operative, wounds/injuries to the rear foot area at the sole of the foot, post-operative following corrective osteotomy in the forefoot area

INDICATIONS

Post-operative, wounds/injuries to the rear foot area and to the heel

	Shoe size UK
	1 – 2
	3 – 4.5
	5.5 – 7
	7.5 – 9
	10 – 11.5
right and left symmetrical	

	Shoe size UK
	2 – 4
	4.5 – 6
	7 – 9
	9 – 11
right and left symmetrical	

BORT Metatarsal Support with Pad

No. 112 070



Med

PRODUCT-BENEFIT

- Elastic bandage support with pad for splayfoot
- Pad inserted in a slit pocket
- Support and lifting of the transverse arch of the metatarsal bones

- Pack: Pair

COLOUR:  skin-tone

THE MED-BENEFIT

- Very easy to use
- Removable pad

BORT Dual-Tension Ankle Support

No. 054 500



Med

PRODUCT-BENEFIT

- Simple ankle support
- Anatomical design
- Closed heel
- Compression class II

COLOURS:  skin-tone  blue

THE MED-BENEFIT

- Very thin, comfortable to wear
- Seamless circular knitted fabric

INDICATIONS

Chronic fallen arch disorders, fallen arches and splayfoot, metatarsalgia

INDICATIONS

Ankle joint injuries, after removal of plaster cast, tendon irritations

Circumference in cm around forefoot



_____ 16, 17, 18, 19, 20,
_____ 21, 22, 23, 24, 25, 26

SIZES

Circumference of ankle



small _____ up to 21 cm
medium _____ 21 – 23 cm
large _____ 23 – 25 cm
x-large _____ 25 – 27 cm

Supplied as pairs

right and left symmetrical




PRODUCT-BENEFIT

- Textile Hallux Valgus Brace for corrective or postoperative stabilisation of the large toe and the metatarso-phalangeal joint
- Side reinforcement splint and plantar splay foot pad
- Can be fully opened using the Vecro fastener
- Additional toe straps for correction of the 2nd - 5th toes

COLOUR: ■ black

THE MED-BENEFIT

- Easy to adapt, also on bandages
- Very comfortable to wear due to soft padded material
- Low material thickness and individual adjustment options
- Can also be worn during the daytime and at night



PRODUCT-BENEFIT

- Hallux Valgus Brace for physiological alignment of the metatarsophalangeal joint of the big toe
- Anatomical fit
- Three adjustable Velcro straps
- For use during periods of rest

- Pack: Single item


COLOUR: □ white/grey 

THE MED-BENEFIT

- Simple and precise adjustment of the leverage force
- Inner surface of the splint entirely lined with foam
- Ventilation hole in the big toe area

INDICATIONS

Hallux valgus, after hallux surgery to ensure the success of treatment, postoperative capsule relief

SIZES	Circumference in cm around forefoot	
 small	_____	16 – 19 cm
medium	_____	19 – 22 cm
large	_____	22 – 25 cm

please indicate left or right
1 cm = 0,3937 inch

SIZES	Shoe size UK	
small	_____	up to 4.5
medium	_____	4.5 – 7
large	_____	over 7

please indicate left or right

BORT Generation Hallux Valgus Splint No. 215 930



PRODUCT-BENEFIT

- Soft brace for corrective or postoperative stabilisation of the large toe and the metatarsophalangeal joint
- Stepless setting of the corrective pressure using the laterally positioned corrective tape
- Relief pad on the same level as the metatarsophalangeal joint

COLOUR: ■ blue

THE GENERATION-BENEFIT

- Low material thickness and individual adjustment options
- Can also be worn during the daytime and at night
- Comfortable to wear due to thin, soft material
- Simple and easy to put on
- Attractive, functional design



BORT Generation:
For better mobility for the aged
Product range with special handling advantages

BORT Silicone Heel Cushion No. 930 050



PRODUCT-BENEFIT

- Silicone heel cushion for shock absorption
- Reduction of impact load on the ankle, knee and hip joints and the spine
- Shell-shaped heel
- Height: 6 mm (0.2 inches)
- Pack: Pair

THE MED-BENEFIT

- Non-slip
- Fits any kind of shoe
- Washable
- Hygienic

INDICATIONS

Hallux valgus, after hallux surgery to ensure the success of treatment, post-operative capsule relief

INDICATIONS

Achillodynia, Haglund's deformity, heel disorders, shock absorption in the case of arthrosis of the ankle, knee or hip joint

SIZES	Shoe size UK
Size 0	4 – 6
Size 1	7 – 9
Size 2	10 – 11

please indicate left or right

SIZES	Shoe size UK
small	up to 4.5
medium	5.5 – 7
large	7.5 – 9
x-large	10 – 11.5

Supplied as pairs



PRODUCT-BENEFIT

- Viscoelastic spur cushion with removable padding module
- Reduction of impact load on the ankle, knee and hip joints and at the spine
- Upper side is lined, under side is self-adhesive
- Height: 6 mm (0.2 inches)

- Pack: Pair

THE MED-BENEFIT

- Non-slip
- Fits any kind of shoe
- Washable
- Hygienic



PRODUCT-BENEFIT

- Silicone heel spur padding with SoftSpot insole
- Reduction of impact load on the ankle, knee and hip joints and the spine
- Shell-shaped heel
- Soft anatomical insole for the heel spur
- Height: 6 mm (0.2 inches)

- Pack: Pair

THE MED-BENEFIT

- Non-slip
- Fits any kind of shoe
- Washable
- Hygienic

INDICATIONS

Heel spur, achillodynia, Haglund's deformity, axial shock absorption in the case of leg joint arthrosis

INDICATIONS

Heel spur, achillodynia, Haglund's deformity, shock absorption in the case of arthrosis of the ankle, knee or hip joint

SIZES	Shoe size UK
A _____	0 – 2
B _____	3 – 6
C _____	7 – 9
D _____	10 – 13

Supplied as pairs

SIZES	Shoe size UK
small _____	up to 4.5
medium _____	5.5 – 7
large _____	7.5 – 9
x-large _____	10 – 11.5

Supplied as pairs

BORT Heel Lift

No. 950 240



Med

BORT Heel Lift made of Silicone

No. 930 200



Med

PRODUCT-BENEFIT

- Viscoelastic silicone heel lift to compensate differences in leg-length
- Upper side is lined, under side is self-adhesive
- Heights: 5 mm, 10 mm (0.2 and 0.4 inches)

- Pack: Single item

THE MED-BENEFIT

- Non-slip
- Fits any kind of shoe
- Washable
- Hygienic

PRODUCT-BENEFIT

- Silicone heel lift to compensate differences in leg-length
- Made of durable silicone
- Heights: 3 mm, 5 mm, 8 mm, 10 mm (0.1, 0.2, 0.3 and 0.4 inches)
- Also for children, see Page 31

- Pack: Single item

THE MED-BENEFIT

- Non-slip
- Fits any kind of shoe
- Washable
- Hygienic

INDICATIONS

Compensation of differences in leg-length, achilles tendon disorders (Achillodynia)

SIZES Shoe size UK

small _____ up to 6
large _____ over 6

Heights _____ 5 mm / 10 mm

right and left symmetrical

SIZES Shoe size UK

Children _____ 1 – 2
small _____ 3 – 5.5
medium _____ 6 – 8.5
large _____ 9 – 11.5

Heights _____ 3 mm / 5 mm / 8 mm / 10 mm

right and left symmetrical

Quality brands for individual
all-round care.

Products for rehabilitation



BORT Therapeutic Putty Standard

No. 951 500 (80 g) | No. 951 600 (440 g)



Med

PRODUCT-BENEFIT

- 100 % silicone therapy putty
- 6 different types of consistency
- Latex-free
- Box sizes: 80 g, 440 g

THE MED-BENEFIT

- Strengthening of muscles in the fingers, hand and forearm
- Comprehensive, therapeutic finger and hand exercises
- Exercise instructions with many examples available in each box

INDICATIONS


For regaining mobility of the hand

HARDNESS

Colour

extra-soft	_____	light blue
soft	_____	pink
medium-soft	_____	light green
medium	_____	dark blue
medium-hard	_____	violet
hard	_____	dark green

BORT EasyLife® Reachers



PRODUCT-BENEFIT

- High quality gripping aid for picking up large and small objects (bottles, needles, etc.)
- Special rubber grips allow calibrated and targeted grasping of objects
- Locking: Grasped objects are held even after the grip lever is released
- Also available in a collapsible version: Art. 118 100, length 75 cm (29.5 inches), when closed 45 cm (17.7 inches)

THE MED-BENEFIT

- Two possible settings:
 - Short leverage: normal effort with unrestricted gripping ability of the hand
 - Long leverage: less effort with restricted gripping ability of the hand
- Collapsible version, ideal for travel or wheelchair use

INDICATIONS

E.g. paraplegia, restricted mobility

Lengths

	50 cm (No. 118 150)
	75 cm (No. 118 170)
	90 cm (No. 118 190)
	Collapsible 75/45 cm (No. 118 100)

1 cm = 0,3937 inch

BORT Aluminium Grippers

No. 117 170



PRODUCT-BENEFIT

- Standard aluminium gripper
- Special rubber pads with anatomical handle
- Length: 75 cm (29.5 inches)

THE MED-BENEFIT

- For gripping even the smallest objects (needles etc.)

INDICATIONS

E.g paraplegia, restricted mobility

SIZE

One length only

BORT Codyn® Seat Cushion

No. 190 530



PRODUCT-BENEFIT

- Special ring cushion with recess in the coccyx area
- **Inflatable**, via valve
- Supplied with manual air pump
- Size: approx. 45 cm × 45 cm, (17.7 × 17.7 inches)
- Inner diameter approx. 16 cm (6.3 inches)

COLOUR: ■ blue

THE MED-BENEFIT

- Suitable for all types of seating
- Inflatable, therefore easy to carry and can be used anywhere

INDICATIONS

Coccygodynia, post coccyx contusion and fracture conditions, perianal relief of haemorrhoids, post low rectal surgery, periproctitis abscess, post episiotomy condition

SIZE

One length only

BORT. Benefit at your side.

**BORT Pads
for Forearm Crutches**
No. 190 800



PRODUCT-BENEFIT

- Self-adhesive padding for forearm crutches
- Fabric surface, can be cut easily
- Size: 11 x 7.5 cm (4.3 x 3.0 inches)
- Pack: Pair

COLOURS: blue black

THE MED-BENEFIT

- Can be trimmed to size, therefore personal use

**BORT Grip Pads
for Forearm Crutches**
No. 190 820



PRODUCT-BENEFIT

- Covered inside with durable neoprene
- Outside covered with skin-friendly fabric
- Easy to apply thanks to Velcro fastening
- Size: 10 x 10 cm (3.9 x 3.9 inches)
- Pack: Pair

COLOUR: black

THE MED-BENEFIT

- High dimensional stability

APPLICATION

Gripping deficit caused by abnormal exertion in the hand or elbow area

SIZE

One length only

SIZE

One length only

Quality brands for individual all-round care.

OTC Brands



Practical aids for daily support.

BORT Special Detergent

Neat and clean.



Stabilising joint supports for free work and time.



Stabilising joint supports for free time and sport.



The thermal programme for feeling at ease.



Soft foot protection with maximum comfort.



Pressure relieving protection for sensitive feet.

Practical aids for daily support.



EasyLife® 7-Day Tablet Box
No. 108 090



PRODUCT-BENEFIT

- Clearly organised 7 day planner
- Braille lettering
- Size: 16 × 3 × 2 cm
(6.3 × 1.2 × 0.8 inches)

COLOURS: ☐ transparent ☒ blue

EasyLife® Tablet Splitter
No. 108 040



PRODUCT-BENEFIT

- Sharp stainless steel blade

COLOURS: ☐ transparent ☒ blue

THE OTC-BENEFIT

- Safe, controlled taking of medicines
- Easy to open

THE OTC-BENEFIT

- Makes cutting tablets and medicines easy and safe
- Easy to clean

BORT Special Detergent

Neat and Clean.

BORT Special Detergent No. 599 070



OTC

PRODUCT-BENEFIT

- Specifically for knitted medical compression fabrics
- Top purity when hand-washing
- For easy everyday use

Contents: 250 ml

THE OTC-BENEFIT

- Already dissolves skin fat at low temperatures
- Biodegradable
- The transparent bottle shows the filling level

THE CARE-BENEFIT

- Long-term satisfaction when wearing knitted compression fabric
- Environmentally friendly to wash
- High productivity

Stabilising joint supports for work and free time.







PRODUCT-BENEFIT

- Wrist support with adjustable Velcro fastening
- Compression can be individually adjusted
- Very good supporting effect

PRODUCT-BENEFIT

- Compression wrist support for easy stabilisation of the wrist
- Optimized fit in the palm area
- Anatomically contoured knit
- Hemmed thumb opening

COLOURS:  white  red
 skin-tone  blue
 black


COLOURS:  skin-tone  blue
 black


THE ACTIVECOLOR®-BENEFIT

- Relief in the case of heavy overstraining of specific muscle groups
- Preventive application in the case of expected overstraining

INDICATIONS

Prevention of overstraining, strains, sprains, irritations due to fatigue

SIZES	Circumference of wrist
	medium — up to 19 cm large — 19 – 21 cm
right and left symmetrical 1 cm = 0,3937 inch	


SIZES	Circumference of wrist
	small — up to 15 cm medium — 15 – 17 cm large — 17 – 19 cm x-large — 19 – 21 cm
right and left symmetrical 1 cm = 0,3937 inch	

active
COLOR®

ActiveColor®
Wrist Support
No. 1480



ActiveColor®
Elbow Support
No. 1420



PRODUCT-BENEFIT

- Compression wrist support to stabilise the wrist
- Anatomically pre-moulded, thermo-plastic palm splint
- Stabilo® Strap
- Seamless, elastic thumb surround
- Optimal wearing comfort

COLOURS: skin-tone black

PRODUCT-BENEFIT

- Round knitted compression elbow support
- Elastic in longitudinal and lateral directions
- Compression class II

COLOURS: red skin-tone
 blue black

- Very effective stabilisation thanks to compressive knitted fabric
- Can be used for work or free time

INDICATIONS

Prevention of overstraining, strains, sprains, irritations due to fatigue

SIZES Circumference of wrist



small	13 – 15 cm
medium	15 – 17 cm
large	17 – 19 cm
x-large	19 – 21 cm

Please indicate left or right
1 cm = 0,3937 inch

SIZES Circumference of elbow



small	up to 26 cm
medium	26 – 28 cm
large	28 – 30 cm
x-large	30 – 32 cm

right and left symmetrical
1 cm = 0,3937 inch

Stabilising joint supports for work and free time.



PRODUCT-BENEFIT

- Mild active support for muscular stabilisation of the lumbar spine
- 2 splints in the dorsal area
- Height: ventral 17 cm (6.7 inches), dorsal 22 cm (8.7 inches)

COLOURS: skin-tone black



PRODUCT-BENEFIT

- Active thigh support
- Anatomically flat knit
- Compression class II
- Length: 21 cm (8.3 inches)

COLOURS: skin-tone blue

THE ACTIVECOLOR®-BENEFIT

- Relief in the case of heavy overstraining of specific muscle groups
- Preventive application in the case of expected overstraining

INDICATIONS

Back pain, lumbago, strains

SIZES Circumference of waist

	small	_____	65 – 80 cm
	medium	_____	80 – 95 cm
	large	_____	95 – 110 cm
	x-large	_____	110 – 125 cm

1 cm = 0,3937 inch

INDICATIONS

Prevention of overstraining, strains, sprains, irritations due to fatigue

SIZES Circumference of thigh

	small	_____	up to 43 cm
	medium	_____	43 – 50 cm
	large	_____	50 – 56 cm
	x-large	_____	56 – 61 cm

right and left symmetrical
1 cm = 0,3937 inch

active
COLOR®

ActiveColor®
Knee Support
No. 1440



ActiveColor®
Ankle Support
No. 1450



PRODUCT-BENEFIT

- Round knitted knee support for prevention of overstraining
- Supporting effect
- Non-slip fit
- Compression class II

COLOURS: ■ red ■ skin-tone
■ blue ■ black

PRODUCT-BENEFIT

- Circular knitted ankle brace to support the ankle joint
- Anatomical form, reinforced everted heel
- Comfortable edging
- Can be worn with sports and walking shoes
- Compression class II


COLOURS: ■ red ■ skin-tone
■ blue ■ black

- Very effective stabilisation thanks to compressive knitted fabric
- Can be used for work or free time

INDICATIONS


Prevention of overstraining, strains, sprains, irritations due to fatigue

SIZES Circumference of knee

	small	_____ up to 32 cm
	medium	_____ 32 – 37 cm
	large	_____ 37 – 42 cm
	x-large (except red)	_____ 42 – 47 cm
	xx-large (except red)	_____ 47 – 52 cm

right and left symmetrical
1 cm = 0,3937 inch

SIZES Circumference of ankle

	small	_____ up to 21 cm
	medium	_____ 21 – 23 cm
	large	_____ 23 – 25 cm
	x-large (except red)	_____ 25 – 27 cm
	xx-large (except red)	_____ 27 – 29 cm

right and left symmetrical
1 cm = 0,3937 inch

Stabilising joint supports for free time and sport.



PRODUCT-BENEFIT

- Light neoprene support to stabilise the metacarpophalangeal joint and saddle joint of the thumb
- Supportive in load situations
- Additional stabilising elements at the thumb
- Variably adjustable using a Velcro fastener

COLOUR: ■ black/green




PRODUCT-BENEFIT

- Light neoprene support for wrist stabilisation
- Supportive in load situations
- Variably adjustable using a Velcro fastener


COLOUR: ■ black/green

THE ACTIVECOLOUR® SPORT-BENEFIT

- Doesn't influence the gripping function
- Full support for sports activity
- Can be worn in water 

INDICATIONS

Prevention of overstraining, strains, sprains, irritations due to fatigue

SIZES	Circumference of wrist
	small/ medium _____ 13 – 17 cm large/ x-large _____ 17 – 21 cm

please indicate left or right
1 cm = 0,3937 inch

SIZES
One size only

please indicate left or right
1 cm = 0,3937 inch

active
COLOR®SPORT

ActiveColor® Sport
Thumb/Hand Support
No. 1470 SP



AQUARIUS



PRODUCT-BENEFIT

- Compression wrist support for easy stabilisation of the wrist
- Anatomically contoured knit
- Hemmed thumb opening
- Regulates moisture with comfort fibre: **AQUARIUS**

COLOUR: ■ black/green

ActiveColor® Sport
Elbow Support
No. 1420 SP



AQUARIUS



PRODUCT-BENEFIT

- Round knitted compression elbow support
- Elastic in longitudinal and lateral directions
- Compression class II
- Regulates moisture with comfort fibre: **AQUARIUS**

COLOUR: ■ black/green

THE ACTIVECOLOR® SPORT-BENEFIT

- Relief in the case of heavy overstraining of specific muscle groups
- Preventive application in the case of expected overstraining
- Very effective stabilisation thanks to compressive knitted fabric

INDICATIONS

Prevention of overstraining, strains, sprains, irritations due to fatigue

SIZES Circumference of wrist



small	_____	up to 15 cm
medium	_____	15 – 17 cm
large	_____	17 – 19 cm
x-large	_____	19 – 21 cm

right and left symmetrical
1 cm = 0,3937 inch

SIZES Circumference of elbow



small	_____	up to 26 cm
medium	_____	26 – 28 cm
large	_____	28 – 30 cm
x-large	_____	30 – 32 cm

right and left symmetrical
1 cm = 0,3937 inch

Stabilising joint supports for free time and sport.



PRODUCT-BENEFIT

- Active thigh support
- Anatomically flat knit
- Compression class II
- Length: 21 cm (8.3 inches)
- Regulates moisture with comfort fibre: **AQUARIUS**

COLOUR: ■ black/green



PRODUCT-BENEFIT

- Compression calf support
- Anatomically contoured knit
- Compression class II
- Length: 23 cm (9.06 inches)
- Regulates moisture with comfort fibre: **AQUARIUS**


COLOUR: ■ black/green

THE ACTIVECOLOR® SPORT-BENEFIT

- Relief in the case of heavy overstraining of specific muscle groups
- Preventive application in the case of expected overstraining

INDICATIONS


Prevention of overstraining, strains, sprains, irritations due to fatigue

SIZES	Circumference of thigh
 small	_____ up to 43 cm
medium	_____ 43 – 50 cm
large	_____ 50 – 56 cm
x-large	_____ 56 – 61 cm

right and left symmetrical
1 cm = 0,3937 inch

INDICATIONS

After injuries to the calf, for treating sports injuries, e. g. ruptures of the muscle fibre, muscle injury

SIZES	Circumference of calf
 small	_____ up to 33 cm
medium	_____ 33 – 38 cm
large	_____ 38 – 43 cm
x-large	_____ 43 – 47 cm

right and left symmetrical
1 cm = 0,3937 inch

active
COLOR® SPORT

ActiveColor® Sport Knee Support

No. 1440 SP



AQUARIUS



PRODUCT-BENEFIT

- Round knitted knee support for prevention of overstraining
- Supporting effect
- Non-slip fit
- Compression class II
- Regulates moisture with comfort

fibre: **AQUARIUS**

COLOUR: ■ black/green

- Very effective stabilisation thanks to compressive knitted fabric
- Can be used for free time or when practicing sports

INDICATIONS

Prevention of overstraining, strains, sprains, irritations due to fatigue

SIZES

Circumference of knee

	small	_____	up to 32 cm
	medium	_____	32 – 37 cm
	large	_____	37 – 42 cm
	x-large	_____	42 – 47 cm
	xx-large	_____	47 – 52 cm

right and left symmetrical

1 cm = 0,3937 inch

Stabilising joint supports for free time and sport.



PRODUCT-BENEFIT

- Circular knitted ankle support to relieve the ankle joint
- Anatomical form, reinforced everted heel
- Can be worn with sports and walking shoes
- Compression class II
- Regulates moisture with comfort fibre: **AQUARIUS**

COLOUR: ■ black/green

THE ACTIVECOLOUR® SPORT-BENEFIT

- Relief in the case of heavy overstraining of specific muscle groups
- Preventive application in the case of expected overstraining

INDICATIONS

Prevention of overstraining, strains, sprains, irritations due to fatigue

SIZES

Circumference of ankle

	small	_____	up to 21 cm
	medium	_____	21 – 23 cm
	large	_____	23 – 25 cm
	x-large	_____	25 – 27 cm
	xx-large	_____	27 – 29 cm

right and left symmetrical

1 cm = 0,3937 inch

active
COLOR®SPORT

ActiveColor® Sport Compression Stocking
No. 1310 SP



AQUARIUS

PRODUCT-BENEFIT

- Sports compression stocking consisting of open-pored knitted mesh fabric
- Good climate management and high breathability
- Even compression pressure increases blood circulation and venous performance under load
- Reduces the period for muscle regeneration
- Pack: Pair

COLOUR: ■ black/green

- Very effective stabilisation thanks to compressive knitted fabric
- Can be used for free time or when practicing sports

INDICATIONS

Prevention of overstraining, strains, sprains, irritations due to fatigue

SIZES

Circumference of calf



Size II	_____	25 – 31 cm
Size III	_____	31 – 38 cm
Size IV	_____	38 – 44 cm
Size V	_____	44 – 50 cm


The thermal programme for feeling at ease.



PRODUCT-BENEFIT

- Protection from cold and draughts in the shoulder-neck area
- Material composition:
80 % cotton, 20 % angora
- Pack: Single item

COLOUR: ☐ white



PRODUCT-BENEFIT

- Excellent protection from cold in the lumbar spinal area
- High angora content
- Material composition:
44 % angora, 26 % pure new wool,
26 % polyamide, 4 % spandex
- Pack: Single item

COLOURS: ☒ skin-tone ☐ white
☒ silver

THE CLIMACARE®-BENEFIT

- Beneficial heat
- High angora content

SIZES	Dress Size	
	Women	Men
small	8/10	38
medium	12/14	40
large	16/18	42
x-large	20/22	44

SIZES	Circumference of hip	
	small	up to 98 cm
	medium	98 – 108 cm
	large	108 – 118 cm
	x-large (skin-tone, silver)	118 – 128 cm

1 cm = 0,3937 inch



CLIMA Care®

CLIMACare® Joint Warmer No. 106 700



PRODUCT-BENEFIT

- Joint warmer for knee and elbow (only size small)
- Elastic knit, no restriction of movement
- Material composition:
40 % angora, 30 % polyamide,
25 % pure new wool, 5 % spandex
- Pain relief in the elbow and knee areas during sport or other activities
- Pack: Pair

COLOURS: skin-tone white

CLIMACare® Foot Warmer No. 106 800



PRODUCT-BENEFIT

- Ideal for people who have difficulty sleeping due to cold feet
- Non-slip pads provide greater safety when worn around the house
- Material composition:
67 % pure new wool, 25 % angora,
8 % polyamide
- Durable
- Pack: Pair

COLOURS: white blue

- Can be worn unnoticed under clothes
- Very comfortable to wear thanks to the material composition

SIZES	Circumference of joint
small	_____ up to 32 cm
medium	_____ 32 – 38 cm
large	_____ 38 – 44 cm
x-large (skin-tone)	_____ 44 – 50 cm
right and left symmetrical	
1 cm = 0,3937 inch	

SIZES	Shoe size UK
x-small (blue)	_____ 0 – 2
small	_____ 3 – 4.5
medium	_____ 5.5 – 6
large	_____ over 7
1 cm = 0,3937 inch	

Soft foot protection with the highest wearing comfort.

**PediSoft®
Toe Separator**
No. 107 020



Silicone

small, large, x-large

PRODUCT-BENEFIT

- Slim silicone toe separator
- Separates dislocated and tight toes
- Contents: **2 pieces**, includes talcum powder

**PediSoft®
Toe Separator with Ring**
No. 107 650



Gel

small or large

PRODUCT-BENEFIT

- Slim gel toe separator with a ring for safe fixing
- Separates dislocated and tight toes
- Contents: **2 pieces**

**PediSoft®
Silicone Ring**
No. 107 060



Silicone

small or large

PRODUCT-BENEFIT

- Silicone ring for relieving pressure on the fingers and toes
- Contents: **2 pieces**, includes talcum powder

**PediSoft®
Toe/Finger Cap**
No. 107 640



Gel

small or large

PRODUCT-BENEFIT

- Cover made of gel with a ventilation hole
- For relieving pressure on the fingers and toes
- Contents: **2 pieces**

THE PEDISOFT®-BENEFIT

- Very comfortable to wear thanks to soft and smooth Silicone & Gel
- Optimal fit to any foot contour thanks to the "memory effect"
- The material returns to its original shape immediately after wearing



**PediSoft®
Toe Cushion**
No. 107 100



Silicone

left or right

PRODUCT-BENEFIT

- Silicone toe pad
- For toe relief
- Fixation with a toe sleeve
- Contents: **1 piece**, includes talcum powder

**PediSoft®
Forefoot Insole**
No. 107 120



Silicone

one size only

PRODUCT-BENEFIT

- Silicone forefoot insole
- For pressure relief in the forefoot area
- It is inserted directly into the shoe
- Long-lasting, durable and non-slip
- Contents: **1 pair**, includes talcum powder

**PediSoft®
Forefoot Protection**
No. 107 680



Gel

one size only

PRODUCT-BENEFIT

- Forefoot gel protection
- For relief of the ball or arch of the foot
- Fixed to the toe with a gel ring
- Contents: **1 pair**

**PediSoft®
Forefoot Cushion**
No. 107 080



Silicone

one size only

PRODUCT-BENEFIT

- Silicone forefoot pad
- Pressure relief in the forefoot area
- As with a toe spreader, the forefoot pad is fixed to the toe
- Contents: **2 pieces**, includes talcum powder

- Designed for frequent use as it is durable
- Hygienic and long-lasting
- Easy to clean with soap or to disinfect, if required

Soft foot protection with the highest wearing comfort.



PRODUCT-BENEFIT

- Wide, ergonomic silicone toe spreader
- Separates dislocated and tight toes
- Contents: **2 pieces**, includes talcum powder



PRODUCT-BENEFIT

- Wide, ergonomic gel toe spreader with a ring for safe fixing
- Separates dislocated and tight toes
- Contents: **2 pieces**



PRODUCT-BENEFIT

- Double gel toe spreader for safe fixing
- Separates dislocated and tight toes
- Contents: **2 pieces**



PRODUCT-BENEFIT

- Ergonomic gel toe spreader with toe ball protection
- Separates dislocated and tight toes
- Contents: **1 piece**

THE PEDISOFT®-BENEFIT

- Very comfortable to wear thanks to soft and smooth Silicone & Gel
- Optimal fit to any foot contour thanks to the "memory effect"
- The material returns to its original shape immediately after wearing



**PediSoft®
Bunion Shield**
No. 107 030



Silicone

one size only

PRODUCT-BENEFIT

- Silicone protection of the ball of the foot for protection against pressure sores, e.g. in shoes
- Contents: **1 piece**, includes talcum powder

**PediSoft®
Little Toe Shield**
No. 107 130



Silicone

one size only

PRODUCT-BENEFIT

- Little toe silicone shield for protection on the little toe
- Contents: **1 piece**, includes talcum powder

**PediSoft®
Heel Cushion with SoftSpot**
No. 107 070



Silicone

small or large

PRODUCT-BENEFIT

- For shock absorption and relief in case of heel spur
- For complaints in the ankle, knee and hip joints, in the heel area and in the arch of the foot
- Height: 6 mm (0.2 inches)
- Contents: **1 pair**

- Designed for frequent use as it is durable
- Hygienic and long-lasting
- Easy to clean with soap or to disinfect, if required

Protection for sensitive feet – with fabric cover.

PediSoft® TexLine Forefoot Pad No. 137 080



small or large

PRODUCT-BENEFIT

- For pressure relief in the forefoot area
- Thanks to its memory effect, the material reverts to its original shape
- Contents: **1 pair**

PediSoft® TexLine Hallux Pad No. 137 030



small or large

PRODUCT-BENEFIT

- For pressure relief on the metacarpophalangeal joint in shoes
- Very hygienic, skin-friendly in shoes
- Non-slip in shoes
- Contents: **1 piece**

PediSoft® TexLine Forefoot/Hallux Pad No. 137 040



small or large

PRODUCT-BENEFIT

- For the simultaneous relief of pressure in the forefoot and ball area
- Surrounds the ball of the foot and the metacarpophalangeal joint softly
- Contents: **1 pair**

PediSoft® TexLine Metatarsal Bunion Pad No. 137 090



small or large

PRODUCT-BENEFIT

- For pressure relief on the ball of the foot, little toe and in the forefoot area
- Contents: **1 pair**

THE PEDISOFT®-TEXLINE-BENEFIT

- Pads provide optimal pressure distribution
- Protection from additional pressure from outside
- Comfortable to wear thanks to the fabric cover
- Optimal fit to any foot contour thanks to the “memory effect”



PediSoft® TexLine Bunion Cap No. 137 140



one size only

PediSoft® TexLine Forefoot Cap No. 137 130



small or large

PRODUCT-BENEFIT

- For pressure relief in the ball area of the foot
- Fits into all shoes
- Slightly adhesive
- Contents: **2 pieces**

PRODUCT-BENEFIT


- For relieving pressure on the fingers and toes
- Skin-friendly, hygienic material
- Protection against pain and pressure-sensitive toes
- Length: 15 cm (5.9 inches)
- Contents: **1 piece**

PediSoft® TexLine Hammer Toe Correction Pad No. 137 120



small or large

PediSoft® TexLine Toe Pad No. 137 100



left or right

PRODUCT-BENEFIT

- Toe correction pad with adjustable, elastic toe straps
- Corrects dislocation of hammer toes
- Up to three toes can be corrected at the same time
- Contents: **1 piece**

PRODUCT-BENEFIT

- Silicone toe pad
- For toe relief
- Fixation with a toe sleeve
- Contents: **1 piece**

- The material returns to its original shape immediately after wearing
- Very hygienic
- Machine washable at 30 °C (with the exception of the Hammer Toe Correction Pad and the Anti-pressure Protection)

Protection for sensitive feet – with fabric cover.

**PediSoft® TexLine
Pressure Protection**
No. 137 110



small, medium, large

PRODUCT-BENEFIT

- Foam tube for padding pressure sores
- Can be individually cut to the fitting length with scissors
- Length: 12 cm (4.7 inches)
- Packaging unit: **3 or 2 (large size)**

**PediSoft® TexLine
Toe/Finger Protection Pad**
No. 137 010



small, medium, large

PRODUCT-BENEFIT

- Tube bandage with round gel pad inside
- Partial pressure relief on the fingers and toes
- Length: approx. 3.5 cm (1.38 inches)
- Contents: **2 pieces**

**PediSoft® TexLine
Tube Bandage**
No. 137 060



small or large

PRODUCT-BENEFIT

- Tube bandage coated with gel, can be trimmed to size, individually adaptable
- For relieving pressure on the fingers and toes
- Length: 15 cm (5.9 inches)
- Contents: **1 pair**

**PediSoft® TexLine
Toe/Finger Cap**
No. 137 050



small, medium, large

PRODUCT-BENEFIT

- Closed toe cap, fully coated with gel
- For relieving pressure on the fingers and toes
- Contents: **1 piece**, includes talcum powder

THE PEDISOFT®-TEXLINE-BENEFIT

- Pads provide optimal pressure distribution
- Protection from additional pressure from outside
- Comfortable to wear thanks to the fabric cover
- Optimal fit to any foot contour thanks to the “memory effect”



(except 137 150)



PediSoft® TexLine
Shin Protection
No. 137 150



one size only

PRODUCT-BENEFIT

- Prophylactic protection against pressure and friction
- External textile fabric, easy to position
- Internal gel pad, even load distribution
- Contents: **1 piece**

- The material returns to its original shape immediately after wearing
- Very hygienic
- Machine washable at 30 °C (with the exception of the Hammer Toe Correction Pad and the Anti-pressure Protection)



BORT GmbH

P.O. Box 1330 | D-71367 Weinstadt, Germany

Am Schweizerbach 1 | D-71384 Weinstadt

Phone +49 (7151) 99200-0

Fax +49 (7151) 99200-50

www.bort.com | medical@bort.com



Art. M0021 V002



BORT. Benefit at your side.