



WESTMINSTER

INTERNATIONAL LTD

Telephone : +44 (0)1295 756300

Fax : +44 (0)1295 756302

E-Mail : info@wi-ltd.com

Website : www.wi-ltd.com

WG KN95 Particulate Respirator Cone



The WG KN95 Particulate Respirator Cone is ideal for persons working in Agriculture, Clean-rooms, Crime Scene Investigation, General Paint Spraying, Pharmaceutical Industries, Veterinary Services, etc.

Feature

- Description Particulate Respirator Cone
- Non-Valved Single Use
- Standards KN95 (GB2626-2006)
- Filter Performance $\geq 95\%$
- Flow rate 85 L/min
- Latex Free
- Shelf Life 6 Years
- Storage Conditions 5 °C to 45 °C & humidity ≤ 70
- Case Qty 10

Note: Depending on stock availability KN95 equivalents may be substituted.

Please note: The wearing of a mask is an aid to help prevent the spread of airborne infectious allergens, bacteria, dust or contaminants it is not a guarantee and masks should always be used in accordance with official advice.

How to use a Mask - World Health Organisation Guidelines

- Before putting on a mask, clean hands with alcohol-based hand rub or soap and water.
- Cover mouth and nose with mask and make sure there are no gaps between your face and the mask.
- Avoid touching the mask while using it; if you do, clean your hands with alcohol-based hand rub or soap and water.
- Replace the mask with a new one as soon as it is damp and do not re-use single-use masks.
- To remove the mask: remove it from behind (do not touch the front of mask); discard immediately in a closed bin; clean hands with alcohol-based hand rub or soap and water.