

STORIES FROM THE WARD

EDWARD

PICU

Tracking trends in vital signs – COVID-19

As told by Dr Faith Ndebele

Edward¹ was admitted onto the PICU on 27th April². Four days later, Edward had a fever and reported feeling unwell. Dr Faith Ndebele and her team used Oxevision to take measurements of Edward's pulse and breathing rates whilst he was in his bedroom.

“

We were concerned for Edward because he was **experiencing COVID-19 symptoms**. When his temperature spiked, we ran blood tests and a COVID-19 test and continued to regularly measure Edward's pulse and breathing rate using the system. By lunch time, **his MEWS (Modified Early Warning Score) was high enough that we had to transfer him to A&E** to be treated, before his test results came back.


Dr Faith Ndebele, Consultant Psychiatrist

Edward returned to the PICU later that evening. Dr Faith Ndebele and her team continued to use Oxevision to frequently measure Edward's pulse and breathing rates without entering his bedroom and tracked changes over time using the Vital Signs Trends functionality.

We used the system over the night to measure Edward's pulse and breathing rates, whilst he was resting in his bedroom and **to track his vital signs trends**. His pulse and breathing rates were gradually increasing and early the next morning, **Edward's resting pulse rate spiked to 125 bpm** and he was transferred back to A&E.

Using the system to track the trends in his vital signs measurements was useful because we knew as **soon as his MEWS (Modified Early Warning Score) changed and we could act immediately** to get him the care he needed.

Dr Faith Ndebele, Consultant Psychiatrist




¹A pseudonym has been used to maintain patient anonymity

²The date has been changed to maintain patient anonymity

Edward's Vital Signs

1 May - 2 May, 2020

- Pulse Rate
- Breathing Rate



ABOUT OXEHEALTH

Oxehealth's Oxevision supports Mental Health organisations to improve inpatient safety and quality whilst giving clinicians more time to care.

Oxevision uses an optical sensor to pay attention to patients in bedrooms when staff can't be there. The Vital Signs module enables staff to take contact-free vital sign spot checks without entering the room. The system provides real time alerts to highrisk activity and displays activity and behaviour reports.

Oxehealth has won the HSJ Patient Safety Award for "Best HealthTech Solution" and Regional Parliamentary "The Future NHS" Award. Oxevision has been highlighted as an example of outstanding practice by the CQC¹.

FOOTNOTES

The Oxehealth Vitals Signs device is a certified medical device in Europe. It is intended for non-invasive spot measurements of pulse rate and respiratory rate. It is a fixed-installed device for use within single occupancy rooms within hospitals, general care, domestic and secured environments where a framework exists which mandates period checks by a trained professional to ensure subject safety. See Instructions for Use for intended use, contraindications, warnings, cautions, usage directions and maintenance.

The activity detection, fall risk, activity tracker and sleep products are activity trackers and do not have a medical purpose or functionality. They are not patient monitors or vital signs monitors. They cannot be used to make a medical or clinical decision.

¹See <https://oxe.health/cqc-cwpt>

oxehealth

+44 (0) 1865 900 599
info@oxehealth.com
Magdalen Centre North,
The Oxford Science Park,
Oxford OX4 4GA