

Borrelia Burgdorferi

IgG/IgM EIA for Lyme Screening



GSD's Lyme IgG/IgM EIA screen is FDA cleared as a qualitative presumptive (first tier) test for Lyme disease, and is optimized for both superior sensitivity and specificity. The simple, easy to use kit includes ready to use controls and a fast protocol which can easily be automated.



IgG/IgM EIA for Lyme Screening

- **SUPERIOR SENSITIVITY**
Includes antigens from *B. burgdorferi* strains B31 and 2591.
- **OPTIMAL SPECIFICITY**
Contains VlsE immunogenic lipoprotein.
- **POLYVALENT SCREEN**
Both IgG and IgM antibodies are included.
- **CONVENIENT**
Kit includes ready to use controls for ease of use.
- **FAST PROTOCOL**
Protocol includes 1:100 serum dilution with an incubation time of 15/15/15.
- **AUTOMATABLE**
Simple procedure is easy to automate.

LYME IgG/IgM EIA SCREEN IS FDA CLEARED AND CAN BE EASILY
AUTOMATED