



# Radiology

ADAWA CLINICAL RESOURCE

*Clinical Professor Bernard Koong, Dr May Lam, Dr Tom Huang*

- 1 Intraoral Radiography
- 2 Panoramic Radiographs
- 3 Multi-detector & Cone Beam CT
- 4 MRI, Ultrasound & Nuclear Medicine
- 5 Prescription & Radiologic Interpretation

Continued advances in radiology have substantially contributed to diagnosis and management. Increasingly, clinicians need to be up-to-date with all available modalities in order to prescribe the optimal radiological test for their patients.

This publication authored by **Clinical Professor Bernard Koong\***, **Drs May Lam†** and **Tom Huang‡** aims to keep the busy clinician apprised on the latest in dental, oral and maxillofacial diagnostic imaging.

The first four chapters discuss the various modalities, including intraoral radiography, cone beam computed tomography, multislice computed tomography and MRI. The final chapter aims to assist the clinician in the prescription of the appropriate test for the common conditions encountered in practice.

---

*\*Clinical Professor Bernard Koong, BDSc(WAust), MSc(OMR) (Toronto), FICD, FADI, FPFA, Oral and Maxillofacial Radiologist.*

*† Dr May Lam, BDSc (Hons)(UWA), FRACDS (GDP), Registrar, Oral and Maxillofacial Radiology.*

*‡Dr Tom Huang, BDSc (Hons1)(UQ), DClinDent(DMFR)(UQ), Oral and Maxillofacial Radiologist.*

---

**References supplied are available on request from  
[adawa@adawa.com.au](mailto:adawa@adawa.com.au)**

**ADAWA Western Articulator**  
October 2019

Originally published in  
ADA News Bulletin 2017-18

# 1. Intraoral radiography

## – a review

*This chapter focuses on intraoral 2D radiography, namely periapical and bitewing radiographs.*

The table below compares the various intraoral radiography technologies based upon the type of detector.

### Technology

	Film radiography	Direct capture radiography	Indirect capture radiography
		CCD (charged coupled device) CMOS (complementary metal oxide semi-conductor).	PSP (photostimulable phosphor) plates.
<b>Method of image capture</b>	Exposure to an X-ray beam results in a latent image recorded on the emulsion-coated film. Chemical processing generates the radiographic image.	The X-ray beam causes silicon crystals in the sensor to convert photons to electrons. The signal is transmitted to the computer which depicts a real time image.	The surface of the PSP is ionised upon X-ray exposure. When the plate is fed into a scanner, the latent image is transmitted to the computer which depicts the image.
<b>Ease of use</b>		Positioning is often more difficult and uncomfortable due to sensor bulk and rigidity of the cable.	Resembles film radiography and better tolerated by patients.
<b>Radiation dose levels</b>	Effective dose commonly ranges from 5–9 $\mu\text{Sv}$ , although 0.1–9 $\mu\text{Sv}$ have been reported.	Reduced radiation dose due to increased sensitivity of sensor compared to film. With optimal settings, the potential for dose reductions of up to 40–60% have been reported.	Depending on exposure settings, PSP radiation dose is less than film, but more than direct capture radiography.
<b>Risk of overexposure</b>	Dark film will alert the operator to overexposure.	Blooming artefacts (i.e. the pixels are 'burnt out' and appear black) will alert the operator to overexposure.	A useable image will still be produced at high radiation doses. As a result, overexposure and increased patient radiation dose may occur without operator awareness.
<b>Frequency of retakes and errors</b>	Chemical processing of film is technique sensitive, potentially contributing to retakes.	Higher number of retakes (up to 28%) have been quoted due to a smaller active capture area and difficulty in positioning of a thick and rigid sensor.	The larger dynamic range of PSP may decrease retakes, but the potential of unrecognised overexposure must be considered. Latent image quality may be lost due to the spontaneous release of electrons, or exposure to ambient lighting. PSP plates should be scanned no later than 10 minutes after exposure. Erased plates must be kept in light-tight containers.
<b>Image resolution (measured in linepairs per mm lp/mm)</b>	$\geq 20$ lp/mm	7–15 lp/mm Theoretically may achieve $\geq 20$ lp/mm, but this is not often possible in practice.	5–13 lp/mm has been quoted. Lower than film and direct digital radiography. Visibility of small endodontic file tips are potentially suboptimal.
<b>Cost</b>	Comparatively low cost	High initial system cost. Sensor is more durable than PSP, but physical damage can occur.	Regular replacement of PSP is necessary. The lifespans of a PSP have been quoted between 50–200 uses. Damage to the phosphor layer (e.g. scratches, folds) may render the PSP non-useable.
<b>Other advantages</b>		<ul style="list-style-type: none"> <li>Ease of digital storage and transmission of image.</li> <li>No need for darkroom equipment.</li> </ul>	
<b>Other disadvantages</b>	Associated hazardous waste and lead foil from film processing.	<ul style="list-style-type: none"> <li>Ease of use may potentially increase the number of exposures.</li> <li>Infection control may be a problem as the detectors cannot be sterilised.</li> <li>Wiping with an alcohol-impregnated tissue has been suggested, but it is uncertain to what degree the sensors tolerate wiping.</li> </ul>	

# 1. INTRAORAL RADIOGRAPHY

## Reducing Radiation Exposure


*The following summarises some of the ways in which dose minimisation can be achieved:*

- Adopt the principle of ALARA (as low as reasonably achievable).
- Review relevant previous studies prior to further imaging. Request from previous clinician if necessary.
- Use the fastest image receptor compatible with the diagnostic task (F-speed film or digital).
- Benefits associated with the use of lead aprons have been shown to be minimal compared to other methods of dose reduction. As such, this could be considered optional unless required by relevant governing bodies.
- Thyroid shields have been shown to reduce thyroid radiation dose.
- Rectangular collimation has been suggested, substantially reducing dose. However, it is technically more demanding and its use must be balanced against the possibility of repeats, potentially increasing the overall dose delivered.
- Use a film positioning device to aid in the parallel alignment of the film with the collimated beam to minimise geometric distortion and for reproducibility of subsequent radiographs.
- Use optimal settings to minimise exposure times.
- Minimise the number of images necessary to obtain essential diagnostic information.
- Optimal film exposure and processing technique.
- Implement quality assurance protocols to regularly evaluate the film processor, image receptor devices and processing chemistry.



## Application

The justification and selection criteria for intraoral radiographic examination have been thoroughly outlined by various bodies, including the American Dental Association and the European Commission. A radiologic examination should be undertaken only when it provides additional information about a condition or planned procedure and/or influences management or surgical approach. Other factors, including the patient's history and susceptibility to dental diseases should also be considered. The benefit must outweigh the accumulated biological risk to the patient.





The clinician needs to prescribe the optimal technique based on clinical indications. The practice of initially prescribing conventional examinations (e.g. plain film 2D radiographs) before determining whether 3D volumetric information is needed has been suggested. However, the limitations of 2D imaging (geometric distortion and projection errors, overlapping of anatomical structures and image magnification) have been well documented and may lead to difficulty with diagnosis and misinterpretation of results. In situations where 2D radiograph is not sufficiently sensitive nor specific, or multiple exposures are required to visualise an area, more advanced techniques should be considered. In the following subsections, brief references to cone-beam computed tomography (CBCT) and multi detector computed tomography (MDCT) are made where relevant.

A forthcoming article will be dedicated to these volumetric techniques.

## Bitewing radiographs

The bitewing (BW) radiograph remains the most optimal technique for detecting interproximal caries, but may not be necessary in patients who show no evidence of the disease, have open proximal contacts or a low caries risk. Diagnosis should be made in conjunction with a clinical examination as low sensitivity (0.24–0.42) and high specificity (0.70–0.97) have been reported. The frequency of their prescription should be determined based on caries risk assessment and altered as the individual circumstances of the patient changes with time. Several studies have reported that CBCT may be useful for caries detection, but is not recommended as the primary technique for caries diagnosis on account of the potential for higher radiation dose, low specificity and beam hardening artefacts when metallic or radiodense restorations are present.

Furcation defects and the more superficial periodontal bone can be demonstrated in a BW, which is taken at right angles to the X-ray beam.

## Periapical radiographs

For periodontal disease, the periapical (PA) radiograph demonstrates a 2D view of periodontal bone and root morphology, with associated limitations. Foreshortening or elongation of the teeth occur, resulting in distortion and inaccurate assessment of the true horizontal bone height. Image distortion is worse when using a bisecting angle technique. Limited sensitivity for measuring the degree of periodontal bone loss, mapping infrabony and vertical defects and the assessment of furcation involvements have been reported. In addition, the existence and severity of bone loss demonstrated in one imaging episode does not indicate active disease. Therefore, intraoral radiography should be correlated with clinical findings, and compared with previous imaging where relevant. Increasingly, for more complex cases, the use of 3D volumetric imaging has been suggested to allow for more accurate detection, diagnosis, prognosis prediction and treatment planning.

For periapical disease, a PA radiograph correlated with clinical findings is generally sufficient for most cases. However, the limitations of these 2D views must be recognised. It has been shown that MDCT and CBCT are more sensitive for the detection of periapical lesions.

# 1. INTRAORAL RADIOGRAPHY

For the initial diagnostic stages of implant planning, PA radiographs may be helpful as a preliminary guide on bone availability and proximity to anatomical structures. However, PA radiographs are insufficiently accurate for the final planning for implant placement. The 3D morphology and relationship to relevant structures is not demonstrated. In this regard, volumetric imaging techniques are superior.

PA radiographs have been recommended for the measurement of baseline bone levels after implant placement. It is obvious that 2D radiographs only demonstrate the proximal peri-implant bone and optimal paralleling technique is essential. CBCT and MDCT may be useful in some cases, but are not without their limitations, namely the adjacent artefact.

For most cases of dental trauma, at least two PA radiographs with different horizontal and vertical angulations have been recommended. They may be sufficiently accurate for minor dental injuries, but low sensitivity in detecting fractures have been reported. CBCT and MDCT may not be essential for acute management of minor dento-alveolar trauma. These techniques should be considered where 2D radiography is inconclusive. Volumetric imaging could also be considered in the long-term management of trauma cases. Other situations requiring volumetric imaging includes cases with unresolved symptoms and more serious trauma, especially where injuries extend beyond the dento-alveolar complex.

The information provided by a PA radiograph may be insufficient for surgical and/or orthodontic treatment planning for impacted canines and their relationship to adjacent structures. Ultra low dose CBCT is emerging as a modality of choice, as it is able to demonstrate the precise location and morphology of the tooth, the surrounding structures and detection of root resorption of adjacent teeth.

Overall, there are increasingly fewer indications for occlusal views as well as the full mouth series with the advent of CBCT. This will be discussed in a later issue.

## Viewing Conditions

To ensure the highest diagnostic accuracy, it is important to review radiographs under optimal conditions. Low ambient lighting (no more than 50 lux), a bright backlight source and magnification is recommended for film radiographs.

With more steps involved in digital image acquisition and display, any component in this process can influence the quality of the final image.

## Display Monitor

With advances in commercial display monitor technology, these monitors (with correct specifications – see below) are comparable with medical grade monitors for most diagnostic purposes.

<b>Maximum brightness</b>	High luminance profile of $\geq 300$ cd/m <sup>2</sup> (room brightness reduces the contrast ratio and making small contrast differences more difficult to observe)
<b>Calibration</b>	Compliance to DICOM Part 14 GSDF standards
<b>Screen Properties</b>	Size: 20 inches recommended; ranges of 17-24" have been quoted Finish: matt finish to minimise the reflection of ambient lighting
<b>Resolution</b>	Ideally, 2048 x 1536 pixels (3 megapixels), but at least 1920 x 1080 (2MP) is indicated. If incorrect display resolution is chosen, geometric distortion or excessive magnification can result.
<b>Contrast ratio</b>	$\geq 500:1$
<b>Response time</b>	Around 8 milliseconds.
<b>Greyscale bit-depth</b>	$\geq 8$ -bit greyscale (24-bit or 32-bit colour)
<b>Video display interface</b>	Digital (e.g. DVI, displayport)



### Graphics card

- Should be appropriately matched and of high quality to avoid information loss and inferior quality images.

### Image enhancements and manipulation

- Zoom control may aid in radiologic diagnosis of proximal carious lesions. It has been suggested that the optimum magnification should be no more than six times the original.
- Enhancement of density and contrast may improve diagnostic accuracy.
- Inverting greyscale values has not been found to improve diagnostic accuracy. Instead, it may hinder dentinal lesion detectability.
- Pseudo-colour enhancement has not been documented to aid diagnosis.

### Storage and transmission

- The image should be stored in its original, uncompressed format with an automatic backup function. Minimising compression of the radiographic image during transmission is necessary to avoid loss of relevant data.

Printed copies of radiographs compromise diagnostic accuracy and are not recommended. They are highly dependent on the quality of printer and paper, and will not demonstrate the same optical range as film or high quality monitors.

### Interpretation

The importance of applying a methodical approach to radiological interpretation is critical.

*Radiological interpretation is discussed in Chapter 5.*

## 2. Panoramic radiographs

### – a review

#### Technology

Film, direct and indirect digital image receptors are used in panoramic radiography, with digital receptors becoming increasingly popular. There are a few differences with regards to the film and detector technology in relation to those employed in intraoral radiography (Chapter 1):

**Film systems:** an intensifying screen with rare earth elements is used to minimise radiation exposure to the patient. A dose reduction of 50-55% has been quoted in the literature. This is achieved because film is also sensitive to the fluorescent light emitted from the intensifying screen. This light will either be of a blue or green wavelength and the appropriate film must be matched to the screen. Screen film combinations with a speed of 400 or greater are recommended.

**Direct digital systems:** due to the expense of large detectors, three to four CCD sensors are arranged vertically to cover the height of the panoramic X-ray beam, and the image is constructed in increments.

**Indirect digital systems:** storage phosphor plates of appropriate size replace screen films in traditional cassettes.

The radiation dose levels between the three systems are comparable. Although, direct digital radiography can produce lower doses than film and storage phosphor systems if the lowest possible setting appropriate to the patient is chosen. The effective dose has been quoted at a range of 3.85–30  $\mu\text{Sv}$ . However, much higher doses have been reported with older film systems. If the patient is a child or of a smaller stature, protocols should be appropriately adjusted to reduce the radiation dose delivered.

#### Technical Aspects & Limitations

The panoramic radiograph is produced by using the principles of conventional tomography. It is a simple curvi-linear form of tomography where the X-ray source and the image receptor simultaneously rotate around the patient's head, capturing the structures within the focal trough. Structures outside of this zone are significantly distorted, blurred or magnified to the point where they are not recognisable, ideally leaving only the dentition and adjacent structures in clearest view possible.

In addition to the usual disadvantages associated with 2D imaging (geometric distortion and projection errors, overlapping of anatomical structures, image magnification and lack of 3D information), there are additional limitations unique to panoramic imaging, detailed in most radiology texts. Some are listed below:

- Lower resolution compared with intraoral radiography and computed tomography (CT).
- Clinically relevant features may be missed if outside the focal trough.
- Real, double and ghost images are always present and familiarity with their appearances is important in the evaluation of anatomical structures and pathoses.
- Overlapping of the proximal surfaces of teeth is common, usually the premolars.
- Superimposition of the cervical spine over the incisor region is usually present.



- Unequal magnification and distortion is found throughout the image, making linear measurements unreliable. Horizontal magnification is much more unpredictable than vertical magnification.
- Objects located more lingually will be projected superiorly due to the slight craniocaudal orientation of the beam.

The technical aspects are well covered in many texts. Incorrect positioning on the bite-block, and/or rotation of the patient's head are two of the most common positioning errors which can lead to significant geometric distortion and horizontal magnification/minification of the image.

- Magnification of the teeth occur when the jaw has been lingually positioned in relation to the focal trough, and is therefore closer to the X-ray source,

causing the beam to pass through it more slowly. Conversely, objects more buccally placed will appear narrower.

- One method to assist in the identification of horizontal distortion is comparing the width of the mandibular first molars. The smaller side was positioned too close to the receptor, while the larger side would have been too close to the X-ray source.
- If the chin is tilted excessively high (chin up), the mandible will be distorted with a flat or inverted occlusal plane, and the hard palate will be superimposed on the roots of the maxillary teeth. Conversely, if the chin is too low (chin down), both the condyles and symphyseal region of the mandible will not be captured and there will be excessive overlapping of the dentition.

## 2. PANORAMIC RADIOGRAPHS

### Reducing Radiation Exposure

Similar principles to intraoral radiography apply, discussed in Chapter 1.

### Application

For a relatively low radiation dose, the panoramic radiograph remains a useful overview of the dentoalveolar and surrounding structures. However, a lack of understanding of the substantial limitations can contribute to misdiagnosis. On the other hand, even with these limitations in mind, it can be sufficient for many situations and procedures. Where clinically appropriate, it may be a useful initial test, where relevant further imaging could then be considered. It may also be useful for patients who cannot tolerate intraoral radiography, but it is not a substitute for the information that can be obtained from an intraoral radiograph.

Panoramic radiographic imaging is reported to be inadequate for the diagnosis of dental caries. The intraoral bitewing radiograph remains the imaging of choice. Some panoramic machines offer an 'extraoral bitewing radiograph', suggested to help with patients who are unable to tolerate intraoral radiography. Reduced specificity due to ghost artefacts, superimposition of air spaces and overlapping of premolars has been reported. Extraoral bitewing radiographs have been shown to be inferior to intraoral radiographs.

Also, the panoramic radiograph is not the optimal technique for most other common dentoalveolar diseases. Compared to intraoral radiographs, it is less sensitive in its ability to detect periapical lesions, particularly those exhibiting early periapical changes. As well, it often underestimates the extent of periodontal osseous defects. Findings on a panoramic radiograph may need to be supplemented with intraoral radiographs or 3D imaging as appropriate.

The panoramic radiograph has been considered to be appropriate for most cases of tooth extraction, including removal of third molars.

Volumetric imaging should be considered for more complicated or difficult cases, and where plain film appearances suggest that the tooth is in close proximity to critical structures, such as the mandibular canal.

It has been suggested that panoramic imaging could be used for the initial radiologic assessment in evaluation for the dental implant, supplemented by 3D imaging techniques. Given its low radiation dose, low cost and easy accessibility, some authors suggest that panoramic imaging can be used for pre-implant vertical linear measurements in the posterior mandibular region by calculation with the appropriate magnification factor and allowing a 2 mm safety margin superior to the inferior alveolar canal or other significant anatomical structures. However, the degree of magnification and minification is unpredictable, and other important 3D variables such as the prominence of the submandibular fossa or the precise morphology and proximity of the maxillary sinus floor cannot be appreciated in these 2D views. It has been suggested that without 3D imaging, a safety margin of 6 mm from the mental foramen would be required, which may contribute to suboptimal treatment planning. There are many factors which influence the accuracy of this method, such as incorrect patient positioning, distortion, discrepancies between the shape of the dental arch and focal trough, and beam angulation. It has been reported that only 17% of measurements from the alveolar crest to the inferior alveolar canal have errors within one millimetre. Panoramic imaging is inferior to volumetric imaging for visualisation of many critical anatomical structures. In general, pre-implant volumetric imaging is considered essential.

The panoramic radiograph is considered to be particularly useful in the evaluation of the developing dentition and any anomalies. However, it is important to remember that supernumerary teeth or pathologies are likely to be missed if not located within the focal trough.

For the assessment of impacted and/or ectopic canines, panoramic radiography alone is considered to be inadequate due to its limitations including the inability to provide 3D information on the buccal or palatal position of the tooth. Cone beam computed tomography better demonstrates the relationship of the impacted canines and adjacent structures, root position and morphology, and possible associated root resorption. The use of volumetric imaging in treatment planning is considered to improve outcomes.



Even though the temporomandibular joints (TMJ) are seen in most panoramic radiographs, it is considered to be an insufficient test where there are specific indications for a radiologic examination of these joints. These views are unable to depict the precise morphology of these joints. The typically oblique projection of these joints together with variations in condylar angulation contribute to a distorted image. Without the ability to view the articular surface, pathologies such as erosions and osteophytes in the TMJ are detected with low sensitivity. CT and magnetic resonance imaging are the optimal techniques for the bony and soft tissue structures of the TMJ respectively. These techniques will be discussed in later issues. The TMJs are almost always depicted in the panoramic radiograph and, despite the limitations, their appearances should be thoroughly evaluated.

Panoramic imaging is inadequate and unreliable for visualisation of sinus anatomy. There is often misinterpretation of the proximity of the maxillary tooth roots to the sinus due to the 2D nature of the image and the presence of distortion. In situations where 3D information is necessary for accurate and precise planning (such as in sinus grafting procedures, or evaluating the likelihood of an oroantral communication), volumetric imaging is recommended. Sinus disease cannot be fully excluded with the panoramic radiograph. However, the sinuses are demonstrated in these views and due attention should be made as significant disease is sometimes demonstrated.

Although routine panoramic radiographic screening for extra-gnathic occult diseases cannot be recommended, the operator should be aware of these additional findings, including the clinical relevance. For instance, panoramic imaging has a low sensitivity for detecting calcified atheromas in the carotid artery or mandibular cortical erosion in the case of osteoporosis, but the presence of these findings contribute to the diagnosis and management of such conditions. Routine panoramic imaging of the patient at set arbitrary intervals is not considered to be justified.

## Viewing Conditions

The viewing conditions in relation to panoramic imaging are similar to that of intraoral radiographs (Chapter 1).


## Interpretation

The relatively broad coverage and tomographic nature of the panoramic radiograph contributes to challenges in interpretation. It has been reported that dental students and dentists often have difficulty identifying the radiologic anatomy, positioning errors and pathologies/anomalies relating to panoramic images. This highlights the need for advanced and continued education in the interpretation of these radiographs. The prerequisites include a

thorough understanding of the limitations, the radiologic anatomy, orofacial pathology and radiologic features of pathoses. It is obvious that a systematic approach must be taken. All structures included in these radiographs must be appropriately evaluated. The practitioner presiding over the study is responsible for the thorough interpretation of the entire image, not only the region of interest.

Radiologic interpretation is discussed in Chapter 5.

### 3. Multi-detector and Cone Beam CT



The following table compares multidetector computed tomography (MDCT) and cone beam computed tomography (CBCT).

#### Technology

	MDCT	CBCT
<b>Method of capture</b>	Thinly collimated, fan shaped X-ray beams rotate in a helical fashion around the patient to capture multiple image slices. A second collimator reduces the scattered photons before it contacts the multiple rows of detectors to improve image quality. The signal is then relayed to the computer for analysis and image construction.	A divergent cone or pyramidal shaped X-ray beam is directed through the region of interest to an area detector on the opposite side of the rotating gantry. Multiple sequential projections are performed as the platform moves through a single arc of 180 or 360 degrees. The multiple images obtained are computed to generate the final 3D volumetric data.
<b>Detector type</b>	Solid state detectors	Flat panel detector or image intensifier/CCD sensor
<b>Patient positioning</b>	Supine	Standing, sitting or supine
<b>Image reconstruction</b>	A back projection algorithm (most commonly the Feldkamp algorithm) is applied to reconstruct the volumetric images with a computer.	
<b>Radiation dose levels</b>	In imaging of the jaws, effective dose ranges of 280-1,410 $\mu\text{Sv}$ have been quoted. While there is variation between different types of scanners, imaging protocols substantially alter the doses delivered.	Reported range of 5-1,073 $\mu\text{Sv}$ . Not all CBCT scans deliver low radiation doses. It is highly dependent on the type of unit and the protocols employed.
	The difficulty in making comparisons between CBCT units and between the different imaging modalities has been documented. A smaller field of view (FOV) in one unit may not necessarily deliver less radiation than a larger FOV in another unit. While ultra-low dose CBCT units, with appropriate low dose protocols, can deliver doses comparable to some panoramic radiographs, other types of CBCT units may deliver doses greater than a well optimised low-dose MDCT.	
<b>Scan times</b>	Substantially faster than CBCT in the imaging of orofacial structures.	Scan times are generally substantially longer than MDCT, with a typical range of 10-40 seconds. Scan times can be shorter for low resolution scans.
<b>Spatial resolution</b>	Sub-millimeter imaging (as small as 0.5 mm) is possible to generate high-spatial resolution images.	Usually higher spatial resolution than MDCT due to smaller voxel sizes, ranging from 0.076-0.4 mm. However, image resolution and quality range widely (refer to 'Disadvantages' below), potentially lower than MDCT



## Technical aspects

MDCT machine operation should only be undertaken by an appropriately trained radiographer/radiologist.

The technical aspects of CBCT has been covered in various texts. Adequate training is essential, including a thorough understanding on the influence of protocols on image quality and radiation safety. The basic principles of reducing radiation exposure was discussed in Chapter 1. The following highlights some of the key issues in relation to CBCT:

### Field of view:


- Selecting the smallest FOV for the region of interest (ROI) will improve the image quality and significantly reduce radiation dose.
- Large volume CBCT generally produce moderate resolution scans which may not be sufficient for fine-detailed tasks such as examining nondisplaced tooth fractures or small root canals. High resolution small FOV CBCT may be more appropriate in such cases.
- The FOV is limited by the size of the detector. Some units are capable of scanning a ROI which is larger than the FOV of the machine by stitching multiple scans together. However, this will result in an increased radiation dose and acquisition time.

**Voxel size:** Smaller voxel sizes capture less X-ray photons, which results in images with more noise. Most CBCT units compensate for this by increasing the radiation dose.

**Frame rate:** The more basis projection images captured, the more data available for reconstruction. This improves spatial and contrast resolution, decreases noise and reduces metallic artefacts. However, it increases patient dose, and also increases scan times with associated potential for motion artefact.

**Trajectory arc:** Reduction in the rotation arc from 360 to 180 degrees will reduce the radiation dose by 50%. However, this will have a corresponding decrease in image quality.

**X-ray generation:** Using a pulsed X-ray beam to coincide with detector sampling will reduce patient exposure by up to 50%.



# 3. MULTI-DETECTOR AND CONE BEAM CT

## Advantages

3D volumetric imaging overcomes the limitations of 2D imaging (superimposition, geometric distortion and magnification). The data from a single scan can be reformatted to be viewed along the axial, coronal, sagittal or any other plane (multiplanar reformatted images). Surface rendering facilitates viewing of the data as 3D virtual “models”.

Compared with 2D plain films, CBCT and MDCT have both been found to be more accurate in measurements in any dimension. Several studies have found CBCT and MDCT linear measurements (over lengths comparable to a tooth) to be accurate within 1mm.

### CBCT potential advantages:

- Easier accessibility, usually lower cost and smaller physical footprint than MDCT.
- Potentially lower radiation dose than MDCT. It must be noted CBCT scans can deliver relative large radiation dose levels, potentially higher than a low dose MDCT scan (refer to ‘The technology’ above).
- High osseous spatial resolution due to smaller voxel size, which can better depict bony trabeculae, small bony/root structures, soft tissue calcifications, etc.

### MDCT potential advantages:

- A more powerful and flexible imaging modality.
- Visualisation of soft tissues and scans with intravenous contrast may be critical in some cases.
- High contrast resolution, such that tissues with even 1% difference in physical density can be identified.
- Improved image quality due to better signal-to-noise ratio.

## Disadvantages

CBCT is subject to more detrimental effects on image quality than MDCT:

- **Increased noise compared to MDCT:** a result of Compton scattering. The amount of scattered radiation is proportional to the volume of tissues in the X-ray beam, and can be reduced by decreasing the FOV.

- **Poor soft tissue contrast:** the contrast resolution is limited to bony or calcified structures. CBCT has less overall contrast resolution than MDCT.
- **Beam hardening:** preferential absorption of lower energy photons through dense objects results in distortion of objects, appearance of bands or streaks (extinction or missing value artefacts) which can contribute to misinterpretation of the scan. The size of the patient’s head as well as the density of the object also contributes to the degree of beam hardening.
- **Motion artefacts:** occurs if the patient is unable to keep still for the entire scan time. This is a common problem as CBCT scan times are generally substantially longer than with MDCT. Selection of a shorter scan time may reduce this, but a decrease in scan time reduces image quality.
- **Metal streaking artefacts:** the absorption of nearly the entire X-ray beam energy by metallic structures results in an opaque streak. This may obscure key structures.
- **Cone-beam effect:** the peripheries of the image are subject to more image distortion, streaking artefacts and greater noise.
- **Moire artefact:** too few basis projection images occur resulting in undersampling of the object.
- **Scanner related artefacts:** appears as a circular streak, often due to poor calibration or imperfections in scanner detection.
- **Inaccurate bone density estimation:** the grey values in CBCT cannot be quantified as Hounsfield units, which measure the relative density of body tissues based on a calibrated grey-level scale. This is possible with MDCT.

The main limitations of MDCT are the potential for higher radiation doses, and reduced accessibility. In some cases, the MDCT image quality of teeth and osseous structures is comparable or better than CBCT scans. However, high resolution small FOV CBCT scan for persons who are not particularly large (with no motion artefact) produces higher image quality images of the jaws and teeth than MDCT scans.



## Application

In general, the prescription of 3D imaging should only be carried out in cases where lower dose imaging techniques are not able to provide the information required for diagnosis and treatment planning. Without specific indications, 3D imaging for routine screening is not recommended. In recent years, there has been an increase in the use of CBCT in dentistry, primarily related to improved understanding of the application, accessibility and dose considerations. While a few authors consider that it has largely replaced MDCT for dentistry, the limitations of CBCT require consideration and it is not the optimal modality for many cases. Understanding of the strengths and limitations of both CBCT and MDCT is necessary for selection of the most appropriate radiologic test. In some instances, other modalities such as MRI, ultrasound or nuclear medicine may be optimal.

The indications for a CBCT/MDCT image has been outlined by various bodies. These position statements and guidelines will alter as the research in this field continues. The prescription of volumetric imaging, in relation to intraoral radiography (IOR) and OPG has been discussed in preceding articles of this series. The following discusses the more common application of CBCT and MDCT in dentistry:

- Despite its increased sensitivity for the detection of dental caries, CBCT has low specificity and is associated with a higher radiation dose

than bitewing radiography. As such, it is not recommended as the primary imaging technique for dental caries. Also, metallic artefacts and beam hardening from adjacent restorations render CBCT inadequate for caries detection.


- MDCT and CBCT are more accurate than IORs and OPGs at identifying the severity of periodontal bone loss and infrabony defects, detecting furcation involvements and providing information on root morphology. With the associated contribution to diagnosis and prognosis, it has been suggested that treatment time and cost may be reduced with the better ability to make appropriate treatment decisions on whether to extract or maintain periodontally compromised teeth. However, it remains important that 3D imaging for periodontal disease should only be applied when clinically indicated.
- CBCT and MDCT have been shown to be more sensitive in the detection of periapical lesions. However, IOR may suffice in many cases. Several authors consider that limited FOV, high resolution CBCT should only be considered in cases where there are contradictory clinical and radiographic signs and symptoms, non-specific or persistent pain, complex root morphology or extra canals, surgical planning and suspected vertical root fracture not detectable with 2D radiography. The application of MDCT and MRI must also be considered.

### 3. MULTI-DETECTOR AND CONE BEAM CT

- Potentially more complicated or difficult extraction cases may require volumetric imaging, particularly when plain film suggests the tooth is in close proximity to critical structures, such as the mandibular canal or maxillary sinus. While the presence of any of Rood and Shehab's criteria may warrant further radiologic examination, it has been shown that the absence of these radiologic signs does not preclude a close relationship with the inferior alveolar nerve. Diagnostic imaging, including 3D imaging, should only be performed where it has the potential to contribute to diagnosis, prognosis, treatment planning or influences a planned procedure.
- Most authors consider that 3D imaging should be used in implant planning. Where required, virtual planning and fabrication of computer-generated surgical guides can be performed, based on the scan data sets.
- When indicated, 3D imaging has been recommended for the assessment of sinus anatomy or pathology, and the preoperative assessment for sinus augmentation or grafting procedures. Plain film imaging for paranasal sinus disease is insufficient. Presently, MDCT remains the most commonly used modality for paranasal sinus disease, considered to be the optimal technique by most clinicians and surgeons.
- MDCT and MRI are appropriate for the evaluation of the bony and soft tissue structures of the TMJ respectively. Some consider that CBCT may be more suitable if soft tissue imaging is not required due to its high bony spatial resolution and reported lower radiation, but factors affecting the image quality must be considered (refer to 'Disadvantages' above), especially beam hardening and the relatively low signal-to-noise ratio associated with CBCT scans. Inappropriate selection of the imaging modality can lead to misdiagnosis and a delay in appropriate management.
- MDCT has been recommended over CBCT for more complex, serious or significant cases, especially where there is potential soft tissue involvement. The list includes benign or malignant tumours, cysts, undiagnosed pain, facial fractures, cellulitis, osteomyelitis or osteonecrosis of the jaw. Other tests, including MRI, ultrasound and nuclear medicine must also be considered.
- CBCT use in orthodontics for diagnosis and treatment planning is becoming increasingly popular, but due consideration should be given to the increased lifetime risk of radiation exposure in children/adolescents. Not all CBCT scans are low dose (refer to 'The Technology' above). The appropriate low dose CBCT unit together with appropriate protocols must be employed. Common applications include assessment of impacted canines, root resorption, tooth position and morphology, and craniofacial anomalies.
- It has been reported that volumetric imaging is more accurate than 2D plain radiographs in the evaluation for obstruction in the upper airways and craniofacial structures. Radiological tests, including 3D imaging, should be prescribed only where it has the potential to contribute to diagnosis, prognosis, treatment planning or influence a planned procedure.

#### Interpretation

The interpretation of CBCT or MDCT scans require a higher level of knowledge and skill sets than plain film. Prerequisite knowledge includes radiologic anatomy, orofacial pathology, radiologic features of pathoses and a thorough understanding of the limitations of the technique. Appropriate image interrogation is critical. Some authors recommend that these scans should be interpreted by persons with sufficient advanced training. Interpretation is discussed in Chapter 5.





## 4. MRI, Ultrasound and Nuclear Medicine

*B Koong, M Lam, T Huang.*

*With contributing authors: Dr Jerry Moschilla, Radiologist and Nuclear Medicine Specialist and Clin A/Prof Michael Bynevelt, Neuroradiologist.*

### The Technology

#### Magnetic resonance imaging

The basis of MRI imaging modality is a magnetic field generated around a patient which aligns spinning nuclei, the majority of which are protons, electromagnetically with the field flux. A radiofrequency pulse is then applied, causing some of the nuclei to absorb energy (resonate), where protons spin in a higher energy state. When the pulse is turned off, relaxation occurs, protons return to the lower energy state and the stored energy is released, which is detected by a receiver coil in the MRI system. These signals are then used to generate the image.

#### Ultrasonography

This application uses ultrasonic energy (in the range of 7.5–30 MHz) to obtain a topographical map of tissue interfaces at differing depths in the body. A transducer

converts electrical energy into sonic energy using piezoelectric crystals which is held against the body part of interest. This ultrasonic beam interacts with the various tissues which all have different acoustic impedance. Some of the sonic waves will reflect (echo) back to the transducer, generating an electric signal which is used to produce the diagnostic image. Each tissue has a characteristic echo pattern, allowing detection of tissue boundaries or pathological changes within the tissue. Real-time imaging is possible, as the processing of these echoes occur at a rapid enough rate to allow perception of motion.

#### Nuclear medicine

A functional imaging technique, which detects abnormal metabolic processes in the body, rather than anatomical/morphological changes, which may not be discernible in the early stages of some diseases.

## 4. MRI, ULTRASOUND AND NUCLEAR MEDICINE

- **Radionuclide imaging:** evaluates tissue function by utilising radioactive atoms or molecules (radionuclides) which emit gamma rays (e.g. Technetium 99m). These radionuclides are combined with a pharmaceutical to form a radiotracer which is distributed to various parts of the body based on their chemical properties. A gamma camera captures the emitted photons, converting them to light, then into a voltage signal for image reconstruction.
- **Single photon emission computed tomography (SPECT):** Is a method for acquiring tomographic slices through a patient, where a single or dual headed gamma camera rotates around the patient detected emitted gamma rays. This data is processed via filtered back-projection or iterative reconstruction algorithms to form axial slices, similar to a CT image.
- **Positron emission tomography (PET):** Is an imaging technique using positron emitting radionuclides (e.g.  $^{11}\text{C}$ ,  $^{13}\text{N}$ ,  $^{15}\text{O}$ ,  $^{18}\text{F}$ ) which are usually combined with pharmaceuticals such as glucose or amino acids, to assess metabolic processes in the body. After a set period of time, positron emission decay occurs, and two photons are produced which travel in opposite directions. A PET camera has a ring of detectors which can map the photons that arrive at the same time and this information is used to produce a functional image of organs and tissues.
- There are now hybrid imaging systems where the nuclear medicine images are coregistered with a CT or MRI images (i.e. SPECT/CT, PET/CT and PET/MRI) allowing for combined morphological and functional imaging.

### Advantages

#### MRI

- No ionising radiation.
- Excellent soft tissue contrast compared to X-ray based techniques due to the higher water content in soft tissues. Certain anatomical and pathological

structures with greater vascularity and permeability can be enhanced by intravenous paramagnetic contrast agents such as gadolinium.

#### Ultrasonography

- No ionising radiation.
- Good soft tissue discrimination and sensitivity for superficial mass lesions.
- Colour Doppler sonography for evaluation of blood flow is possible.
- May be a useful alternative for patients who are contraindicated for MRI.




#### Nuclear medicine

- Evaluates physiologic alterations of tissues.
- Identify early changes of some diseases not demonstrated in other techniques.
- PET has very high spatial resolution and is able to detect very small lesions.

### Disadvantages

#### MRI

- Ferromagnetic objects may move, overheat and therefore injure the patient when in the vicinity of the magnetic field. Therefore, this modality is contraindicated for some patients with some implanted metallic objects or medical devices.
- Metals used in dentistry will not move but may distort the image in its vicinity. Titanium implants only cause minor degradation of the image.
- Longer scan times.
- May not be suitable for claustrophobic patients.
- The use of gadolinium-based contrast media must be used with caution in those with renal impairment as this has been associated with nephrogenic systemic fibrosis. Gadolinium deposition within regions of the brain has recently been discovered and is currently being investigated.



## Ultrasonography

- Difficulty in imaging deeper structures and structures obscured by bone.

## Nuclear medicine

- Associated with ionising radiation.

## Application in the Orofacial Region

### MRI

- Evaluating of soft tissue anatomy and pathology, characterisation and extent of lesions e.g. evaluating for perineural spread of tumours.
- Additional characterisation of soft tissue components of bone lesions.
- Considered the gold standard in the assessment of the soft tissues of the temporomandibular joint, particularly the articular disc position. Also demonstrates joint effusions, synovitis, marrow oedema, and changes in the adjacent masticatory muscles.
- **Implant dentistry:** Identifying the location of inferior alveolar neurovascular bundle where multislice CT or cone beam CT is not able to demonstrate the location of the mandibular canal.

## Ultrasonography

- Commonly used in the orofacial region for the evaluation of salivary gland, cervical lymph nodes and neck lumps.
- Ultrasound guided fine needle aspiration and core biopsies.
- Other applications in this region include evaluation of the thyroid glands and carotid vessels.


## Nuclear medicine

- Osteoblastic metastatic neoplasms involving bone.
- SPECT: has been used to assess mandibular growth in patients with asymmetry. The limited specificity of these studies must be considered. Correlation with CT and/or MRI is useful.
- The extent of medication-induced osteonecrosis of the jaw (MRONJ).
- Gallium and radiolabelled white cell radionuclide and SPECT/CT imaging is useful in diagnosing base of skull osteomyelitis.
- PET/CT imaging: is useful for skeletal imaging for assessment of primary bone tumours, locating metastases in bone and detecting osteomyelitis. It is often correlated with post contrast CT or MRI scans.
- PET/CT is particularly useful for staging squamous cell carcinoma and other head and neck malignancies.

## Interpretation

Radiologists, neuro-radiologists, maxillofacial radiologists and nuclear medicine specialists perform interpretation of these studies.

## 5. Prescription and radiologic interpretation



This chapter summarises the prescription of the various imaging techniques in relation to the more common conditions/clinical situations encountered in dentistry. An introduction to the key responsibilities and principles involved in radiological interpretation is also included.

**Abbreviations:** IOR – intraoral radiography; OPG – panoramic radiographs; CBCT – cone beam computed tomography; MDCT – multidetector (multislice) computed tomography; MRI – magnetic resonance imaging.

### Prescription of the optimal imaging technique - a summary:

The various imaging techniques, including the associated radiation dose levels delivered, have been discussed in the preceding chapters. The following summarises the indications in relation to the conditions and clinical scenarios which are more commonly encountered in dentistry:

#### **Dental caries:**

- Bitewing radiographs remain the optimal technique – the limitations, especially in relation to sensitivity, are noted and clinical correlation is essential.
- OPGs are considered inadequate – caries cannot be fully excluded.
- CBCT is associated with increased sensitivity, but decreased specificity. Artefact related to restorations is another limitation.

#### **Periodontal disease:**

- OPG provides a good overview. The associated limitations compared to the IOR are recognised.
- IOR demonstrates periodontal bone levels and root morphology, but is less accurate than CBCT and MDCT. Limited sensitivity in relation to vertical defects and furcations has also been demonstrated.

- CBCT and MDCT provide 3D information, which allows for more accurate diagnosis and treatment planning. Could be considered for cases with moderate to severe periodontal disease.

#### **Periapical inflammatory disease:**

- A periapical radiograph and clinical findings are sufficient for most cases.
- OPG is less sensitive for detecting periapical lesions.
- CBCT and MDCT are more sensitive, should be considered in cases with contradictory findings or non-specific/unresolved pain.
- In rare instances, MRI can be considered where the clinical suspicion for periapical inflammatory lesion is high, but this is not demonstrated with MDCT/CBCT.

#### **Other conditions affecting the jaws:**

- For example, cysts, tumours, osteomyelitis, etc.
- MDCT is usually the technique of choice. CBCT may suffice for some lesions.

#### **Orofacial pain:**

- Dentoalveolar inflammatory disease is a common cause – clinical findings combined with IOR and/or OPG is sufficient for diagnosis in most cases.
- Dentoalveolar inflammatory lesions contributing to orofacial pain cannot be fully excluded with IOR and OPG.
- MDCT or CBCT should be considered if the cause of pain is not identified clinically and with 2D imaging.
- Volumetric imaging should be considered if other causes for the orofacial pain is suspected, e.g. sinus disease.
- MRI should be considered if diagnosis is not made following volumetric imaging.

### **Dental implants:**

- IOR and OPG may be useful for initial assessment.
- Pre-implant MDCT or CBCT must be considered.
- MRI may be used to identify the location of the inferior alveolar nerve where the mandibular canal borders are not demonstrated with MDCT or CBCT.

### **Dental extractions:**

- IOR and OPG are appropriate for most cases.
- Volumetric imaging could be considered for more complicated cases or where the tooth is potentially in close relationship to significant anatomic structures, e.g. the mandibular canal.
- Third molar and related morphology including the relationships with the mandibular canal are better demonstrated with volumetric imaging. The indication for CBCT or MDCT is based upon clinical parameters, potential intra-surgical implications and the OPG appearances.
- Ultra-low dose CBCT should be considered for impacted/ectopic canines.

### **Dental trauma:**

- PA radiographs with different horizontal and vertical angulations has been recommended.
- CBCT and MDCT should be considered where 2D imaging is inconclusive, for complex/severe cases of dentoalveolar trauma and where jaw fractures are suspected.

### **Paranasal sinuses:**

- OPGs are inadequate – sinus disease cannot be fully excluded with this technique.
- Volumetric imaging (particularly MDCT) is recommended for the assessment of sinus disease.

### **TMJ:**

- OPGs are inadequate – arthropathy cannot be fully excluded with this technique.
- MDCT is optimal for the evaluation of the bony structures. MDCT also demonstrates the soft tissues and the articular disc may be visualised.
- CBCT demonstrates the bony structures – image degradation related to beam hardening must be considered.
- MRI is also optimal in the evaluation for effusion, synovitis and marrow oedema. MRI is also optimal in the evaluation for effusion, synovitis and marrow oedema. Bony structures are also demonstrated although the lower spatial resolution is recognised.

### **Soft tissue lesions:**

- CBCT does not demonstrate the soft tissues sufficiently well.
- MDCT or MRI could be considered. Ultrasound may be useful, especially for more superficial lesions



## 5. PRESCRIPTION & RADIOLOGIC INTERPRETATION



### Radiological Interpretation

All structures included in any imaging must be appropriately evaluated. The practitioner presiding over the radiologic study is responsible for the thorough interpretation of the entire image, not only the region of interest. Should the entirety of the lesion not be included in the initial scan, preliminary interpretation should still be carried out, which can be useful in deciding the optimal imaging technique for further evaluation. Studies have demonstrated the challenges that dentists encounter in relation to the interpretation of OPGs and CBCT scans.

Volumetric data requires a different level of knowledge and skill for interpretation. Several authors recommend that CBCT scans should be interpreted by persons with sufficient advanced training.

The following publications introduce the key principles of radiologic interpretation:

Koong B. The basic principles of radiological interpretation. *Aust Dent J*. 2012;57 Suppl 1:33–9.

Koong B. Diagnostic imaging of the periodontal and implant patient. In Lang PL & Lindhe J (eds.), *Clinical periodontology and implant dentistry* (6th edition., 574–608). UK: Wiley Blackwell.

Some of the key points, based upon these two publications, are summarised:

#### **Prerequisites:**

**Radiologic anatomy:** An in-depth knowledge of anatomy and their normal variants is critical. Knowledge of the appearances of all normal structures will aid in identifying the presence of pathology.

**Pathology:** Knowledge of the pathology which may occur in all regions included in the field of view or scan, as well as the radiologic appearances of these lesions, is obviously of importance.

**Imaging modality:** Understanding the strengths and limitations of the various imaging techniques impacts on the interpretation of the studies. Examples include the tomographic nature of OPGs and beam hardening seen in CBCT scans.

**Viewing conditions:** The lighting conditions and display quality of the monitor can have a significant influence on the accuracy of the radiologic interpretation (refer to Chapter 1).

#### **Identifying the presence of disease:**

- A methodical approach during evaluation of the entire image or dataset is critical.
- The interrogation of volumetric data requires a different skill set to that for plain 2D films.

#### **Radiologic evaluation of a lesion(s)**

An algorithm should be followed to identify the relevant radiologic features. Below is an example:

1. Location
2. Shape and contour
3. Border
4. Internal appearances
5. Adjacent anatomical structures

Lesions will not always present classically, nor will they necessarily demonstrate all the typical features. As a result, it is necessary to weigh the identified features. For example, a lucency at the apex of a tooth is often inflammatory in nature. However, a malignant lesion may also present as a lucent lesion apically and weight must be given to the marginal appearances.

## Conclusion

As advances in diagnostic imaging continue, it is increasingly essential that clinicians remain up-to-date in order to prescribe the optimal radiological test for their patients. The practitioner presiding over the radiologic study is responsible for the examination in its entirety and should ensure that it is interpreted by appropriately skilled person.




# Mediref

The national standard in  
secure medical communication

**NOW FREE  
for ADAWA  
Members!**

## Get started in 3 easy steps:

1. Register an account at [mediref.com.au](http://mediref.com.au)
2. Download the Mediref Printer
3. To send any document or photo - simply click on the Print button and select your recipient!



## See the benefits below:



Send large files  
and images



Keep a log of  
transferred files



Reduce postage  
costs



Legal and secure  
messaging for  
practitioners and  
patients



Works alongside  
CRM systems



Built-in directory  
of state  
practitioners

## See what existing clients say:

It is intuitive and even  
works on my Macs!

Support is very good [with] after  
hours response and constant  
development, they are always willing  
to listen to suggestions. Mediref  
is a great Australian product,

**I recommend it unreservedly.**

Dr. Kevin Morris  
Specialist Periodontist  
Northern Periodontics & Implants,  
Victoria

I compile a high volume of medical  
reports for a four surgeon practice.

Posting large quantities of letters  
was a very time-consuming  
and costly process.

Mediref has provided a quick  
and secure way of emailing  
correspondence that is easy too.

Erin  
Administration Officer/Medical Typist  
West Perth Oral & Maxillofacial  
Surgery, Western Australia

Mediref is so quick and easy to  
use and it is reassuring to know  
patient referrals are sent securely

**It is great to have such  
comprehensive backup and support.**

Dr. Andrew Welsch  
Bannockburn Dental, Victoria



**Mediref**  
[mediref.com.au](http://mediref.com.au)

For further information contact Andrea Andrys [mediref@adawa.com.au](mailto:mediref@adawa.com.au)