

# POWERFUL + PORTABLE

Mobile Blood Clot Prevention





The VenaPro™ is intended to be an easy to use portable system, prescribed by a physician, for use in the home or clinical setting to help prevent the onset of DVT in patients by stimulating blood flow in the extremities (stimulating muscle contractions).

### CIRCUL8 PRO Can Be Used To:

- Aid in the prevention of DVT
- Enhance blood circulation
- Diminish post-operative pain and swelling
- Reducing wound healing time



TUBELESS &  
CORDLESS  
OPERATION

EASY TO USE

UP TO 20 HOUR  
RECHARGEABLE  
BATTERY

PORTABLE &  
LIGHTWEIGHT

USB COMPLIANCE  
MONITORING

## HOW IT WORKS




### Place on Calf

Open wrap and place around your calf with inflatable part of wrap on calf.



### Turn On

Press and hold power button on both units for 3 seconds and release. Device will start.



### Pressure & Charging

The device will gently fill with air creating pressure on your calf. Recharge device when needed. Scan QR to see product tutorial.



Ortho8 Inc.  
2217 Plaza Drive  
Rocklin, CA 95765  
800-604-2487



2140 Hutton Drive  
Czrrollton, TX 75006

888-963-6265  
info@pmpmed.com

in Precision Medical Products  
@pmpmed



Enhancing Life  
Through Innovation