

ENFit[®] ENTERAL FEEDING SYSTEM



Ameritus ENFit® ENTERAL FEEDING SYSTEM

- Ameritus ENFit® feeding tubes with ENFit® connector hubs are specifically designed to help prevent misconnections
- The ENFit® connector hub/tip are ONLY compatible with other ENFit® enteral devices
- Orange color coding provides easy visual recognition and continuity of feeding tubes, extension sets, and syringes. ONE COLOR SYSTEM FROM PATIENT TO PUMP
- ENFit® enteral system offers the clinician a technically advanced, safe enteral system
- Sterile, Single packed, BPA/Latex/DEHP free

Ameritus ENFit® FEEDING TUBES

- A full line of polyurethane nasogastric feeding tubes from the micro-preemie to the pediatric patient
- Feeding tube identification including French size, length and lot number printed on tube
- Unique, open-ended distal tip promotes greater fluid flow, effective flushing and minimizes pockets of fluid build-up
- Second large distal lateral hole placed only 5.5 mm from the distal open-ended tip reducing the risk of high feeds in the esophageal tract
- ENFit® feeding tubes come in Dual and Single proximal connector hubs to fit clinical needs
- ENFit® feeding tubes also offered with Lateral Access Port. Two proximal ports: an administration port and medication port. Designed to eliminate disconnection when administering medications/fluids

Ameritus ENFit® EXTENSION SETS

- ENFit® Extension set comes in two lengths: 36" and 60" in straight and bifurcated configurations (bifurcated set offers additional side medication port eliminating the need to disconnect the feeding line)
- Medication port is spliced directly into the distal connector. Dead space of the distal hub is reduced to eliminate fluid waste

Ameritus ENFit® SYRINGES

- One piece, sterile, ENFit® syringe with orange plunger, black metric gradients, and cap
- ENFit® Low Dose Tip: 1 mL, 3 mL, 5 mL
- Standard ENFit® 10 mL, 20 mL, 30 mL, 60 mL (off-set tip to maximize caloric intake)

ENFit® ENTERAL FEEDING SYSTEM



Ameritus ENFit Enteral Feeding Tubes		
Catalog Number	Description	Unit of Measure
ENF-40P-40	4 Fr, 40 cm ENFit Enteral Feeding Tube	10/bx
ENF-50P-40	5 Fr, 40 cm ENFit Enteral Feeding Tube	10/bx
ENF-Y50P-40	5 Fr, 40 cm Dual ENFit Enteral Feeding Tube	10/bx
ENF-60P-40	6 Fr, 40 cm ENFit Enteral Feeding Tube	10/bx
ENF-Y60P-40	6 Fr, 40 cm Dual ENFit Enteral Feeding Tube	10/bx
ENF-80P-40	8 Fr, 40 cm ENFit Enteral Feeding Tube	10/bx
ENF-Y80P-40	8 Fr, 40 cm Dual ENFit Enteral Feeding Tube	10/bx
ENF-40P-50	4 Fr, 50 cm ENFit Enteral Feeding Tube	10/bx
ENF-50P-50	5 Fr, 50 cm ENFit Enteral Feeding Tube	10/bx
ENF-Y50P-50	5 Fr, 50 cm Dual ENFit Enteral Feeding Tube	10/bx
ENF-60P-50	6 Fr, 50 cm ENFit Enteral Feeding Tube	10/bx
ENF-Y60P-50	6 Fr, 50 cm Dual ENFit Enteral Feeding Tube	10/bx
ENF-80P-50	8 Fr, 50 cm ENFit Enteral Feeding Tube	10/bx
ENF-Y80P-50	8 Fr, 50 cm Dual ENFit Enteral Feeding Tube	10/bx
ENF-50P-80	5 Fr, 80 cm ENFit Enteral Feeding Tube	10/bx
ENF-Y50P-80	5 Fr, 80 cm Dual ENFit Enteral Feeding Tube	10/bx
ENF-60P-80	6 Fr, 80 cm ENFit Enteral Feeding Tube	10/bx
ENF-Y60P-80	6 Fr, 80 cm Dual ENFit Enteral Feeding Tube	10/bx
ENF-80P-80	8 Fr, 80 cm ENFit Enteral Feeding Tube	10/bx
ENF-Y80P-80	8 Fr, 80 cm Dual ENFit Enteral Feeding Tube	10/bx
ENF-50P-90	5 Fr, 90 cm ENFit Enteral Feeding Tube	10/bx
ENF-Y60P-90	6 Fr, 90 cm Dual ENFit Enteral Feeding Tube	10/bx
ENF-Y80P-90	8 Fr, 90 cm Dual ENFit Enteral Feeding Tube	10/bx
ENF-Y100P-90	10 Fr, 90 cm Dual ENFit Enteral Feeding Tube	10/bx
ENF-100P-110	10 Fr, 110 cm ENFit Enteral Feeding Tube	10/bx

Ameritus ENFit Enteral Extension Sets		
Catalog Number	Description	Unit of Measure
ENF-EXT-036	36" ENFit Extension Set	100/ca
ENF-EXT-060	60" ENFit Extension Set	100/ca
ENF-EXT-Y060	60" Bifurcated ENFit Extension Set	100/ca

Ameritus ENFit Enteral Syringes		
Catalog Number	Description	Unit of Measure
ENF1-KMS-LD	1mL ENFit Enteral Syringe (Low Dose Tip)	200/ca
ENF3-KMS-LD	3mL ENFit Enteral Syringe (Low Dose Tip)	100/ca
ENF5-KMS-LD	5mL ENFit Enteral Syringe (Low Dose Tip)	100/ca
ENF10-KMS	10mL ENFit Enteral Syringe	100/ca
ENF20-KMS	20mL ENFit Enteral Syringe	100/ca
ENF30-KMS	30mL ENFit Enteral Syringe	50/ca
ENF60-KMS	60mL ENFit Enteral Syringe	20/ca

Ameritus ENFit Enteral Accessories		
Catalog Number	Description	Unit of Measure
ENF-STRAW	5" ENFit Enteral Milk Straw	100/bx

ENFit® is a federally registered trademark of GEDSA, used with permission. Ameritus ENFit® feeding tubes and extension sets manufactured for Kentec Medical, Inc. Ameritus ENFit® syringes and milk straw distributed by Kentec Medical, Inc. Ameritus is a registered trademark of Kentec Medical, Inc.

Manufactured for and Distributed by
Kentec Medical, Inc.

17871 Fitch Irvine, CA (949) 863-0810 FAX (949) 833-9730
TOLL FREE - USA Only (800) 825-5996
sales@kentecmedical.com | www.kentecmedical.com



Visit us on Facebook, YouTube and LinkedIn!