

Safety Oral ENTERAL FEEDING SYSTEM



AMERITUS
Product Brand

AMERITUS SAFETY ORAL FEEDING SYSTEM

- Enteral safety feeding tube, extension set and syringe prevent connection to Luer devices and IV lines
- Orange color coding provides easy visual recognition and continuity of feeding tubes, extension sets and syringes. ONE COLOR SYSTEM FROM PUMP TO PATIENT
- Sterile, single packed
- Does not contain natural latex rubber
- Does not contain di(2-ethyl hexyl)phthalate (DEHP)
- Does not contain bisphenol A (BPA)

AMERITUS SAFETY ORAL FEEDING TUBES

- A full line of polyurethane safety enteral feeding tubes serving micro-preemie to pediatric patient needs
- Multiple sizes and lengths
- Feeding tubes available in dual and single proximal connector hubs to fit clinical needs
- The dual hub feeding connector is designed with two proximal ports: an administration port and medication port. This design eliminates disconnection when administering medications/ fluids
- Feeding tube identification including french size, length and lot number pad printed on tube
- Unique, open-ended distal tip promotes greater fluid flow, effective flushing and minimizes pockets of fluid build-up
- Second large distal lateral hole placed only 5mm from the distal open-ended tip reducing the risk of high feeds in the esophageal tract



AMERITUS SAFETY ORAL EXTENSION SETS

- Male and female connectors are designed to prevent misconnections
- Extension sets come in two lengths: 36" and 60" in straight and bifurcated configurations (bifurcated set offers additional side medication port eliminating the need to disconnect the feeding line)



AMERITUS SAFETY ORAL SYRINGES

- Sterile, one piece, oral syringe with orange plunger. Black metric gradients and labeled: ORAL/ENTERAL ONLY
- Self-righting orange cap
- Multiple syringe sizes including: 1mL, 3mL, 5mL, 10mL, 20mL, 30mL, 60mL
- Compatible with most syringe pumps



Ameritus Safety Oral Enteral Feeding Tubes		
Catalog Number	Description	Unit of Measure
AEG-40P-40	4 Fr, 40 cm Single Hub Enteral Feeding Tube	10/bx
AEG-40P-50	4 Fr, 50 cm Single Hub Enteral Feeding Tube	10/bx
AEG-50P-40	5 Fr, 40 cm Single Hub Enteral Feeding Tube	10/bx
AEG-50P-50	5 Fr, 50 cm Single Hub Enteral Feeding Tube	10/bx
AEG-50P-80	5 Fr, 80 cm Single Hub Enteral Feeding Tube	10/bx
AEG-60P-40	6 Fr, 40 cm Single Hub Enteral Feeding Tube	10/bx
AEG-60P-50	6 Fr, 50 cm Single Hub Enteral Feeding Tube	10/bx
AEG-60P-80	6 Fr, 80 cm Single Hub Enteral Feeding Tube	10/bx
AEG-80P-40	8 Fr, 40 cm Single Hub Enteral Feeding Tube	10/bx
AEG-80P-50	8 Fr, 50 cm Single Hub Enteral Feeding Tube	10/bx
AEG-80P-80	8 Fr, 80 cm Single Hub Enteral Feeding Tube	10/bx
AEG-Y50P-40	5 Fr, 80 cm Dual Hub Enteral Feeding Tube	10/bx
AEG-Y50P-50	5 Fr, 80 cm Dual Hub Enteral Feeding Tube	10/bx
AEG-Y50P-80	5 Fr, 90 cm Dual Hub Enteral Feeding Tube	10/bx
AEG-Y60P-40	6 Fr, 40 cm Dual Hub Enteral Feeding Tube	10/bx
AEG-Y60P-50	6 Fr, 50 cm Dual Hub Enteral Feeding Tube	10/bx
AEG-Y60P-80	6 Fr, 80 cm Dual Hub Enteral Feeding Tube	10/bx
AEG-Y80P-40	8 Fr, 40 cm Dual Hub Enteral Feeding Tube	10/bx
AEG-Y80P-50	8 Fr, 50 cm Dual Hub Enteral Feeding Tube	10/bx
AEG-Y80P-80	8 Fr, 80 cm Dual Hub Enteral Feeding Tube	10/bx

Ameritus Safety Oral Enteral Extension Sets		
Catalog Number	Description	Unit of Measure
OC-ENT-036	36" Enteral Extension Set	100/ca
OC-ENT-060	60" Enteral Extension Set	100/ca
YOC-ENT-60	60" Bifurcated Enteral Extension Set	100/ca

Ameritus Safety Oral Enteral Syringes		
Catalog Number	Description	Unit of Measure
AEG1-KMS	1mL Enteral Syringe	200/ca
AEG3-KMS	3mL Enteral Syringe	100/ca
AEG5-KMS	5mL Enteral Syringe	100/ca
AEG10-KMS	10mL Enteral Syringe	100/ca
AEG20-KMS	20mL Enteral Syringe	100/ca
AEG30-KMS	30mL Enteral Syringe	50/ca
AEG60-KMS	60mL Enteral Syringe	20/ca

Ameritus Safety Oral Enteral Accessories		
Catalog Number	Description	Unit of Measure
AEG-STRAW	5" Enteral Milk Straw	100/bx



Ameritus feeding tubes and extension sets manufactured for Kentec Medical, Inc.
 Ameritus syringes and milk straw distributed by Kentec Medical, Inc.
 Ameritus is a registered trademark of Kentec Medical, Inc.

Manufactured for and Distributed by

Kentec Medical, Inc.

Irvine, CA (949) 863-0810

sales@kentecmedical.com | www.kentecmedical.com



Visit us on Facebook, YouTube and LinkedIn!