


### Set Includes

Before your first use, ensure that you have all components listed below.



- a. Foot model with swivel stand ···· 1 piece
- b. Toe nail A ······ 20 pcs  
(thickened ingrown nail)
- c. Toe nail B ······ 20 pcs  
(thickened nail with ringworm)
- d. Callosity ······ 10 pcs
- e. Corn ······ 10 pcs

### Consumables and replacement parts

Product Code	Part Name
11273-030	Toe nail A (thickened ingrown nail) 20 pcs
11273-040	Toe nail B (thickened nail with ringworm) 20 pcs
11273-050	Callosity 10 pcs
11273-060	Corn 10 pcs



Toe nail A  
(thickened ingrown nail)



Toe nail B  
(thickened nail with ringworm)



Callosity



Corn

### ⚠ DOs and DON'Ts

#### DOs

- Handle with care. The materials for the models are a special composition of soft resin which could be damaged in case of fall. Please use this model carefully and refrain from subjecting it to unnecessary shock.
- Store the model at room temperature, away from heat, moisture and direct sunlight. Bad storage condition may cause damage.
- Clean the surface of the foot model with a wet, soft cloth and mildly soapy water or diluted detergent. Put talcum powder before storage to prolong the model life.

#### DON'Ts

- Do not write on the foot model with pens, or let printed matters such as newspapers come in contact with the surface of the model. The resin of the model absorbs ink, and those stains cannot be cleaned off.
- Never use organic solvent such as paint thinner to clean the model surface.
- The foot model can remain usable even though the color of the model may fade over time.