


Set includes

Before your first use, ensure that you have all components listed below.



① Lower torso manikin	1
● Case holder	
② Inner frame	1
③ Outer frame	1
④ Skin	1
● Interchangeable rectum units	
⑤ Normal	1
⑥ Rectal cancers	1
⑦ Rectal cancers	1
⑧ Rectum with polyps	1
⑨ Case module	3
⑩ Knob nut	1
⑪ Spring	1
⑫ Washer	1
⑬ Pillow	1
⑭ Lubricant gel	1
⑮ Gloves	1
Storage case (no image)	



*Product is packaged as shown on the left.

⚠ DOs and DON'Ts

DO s

Handle with care

The materials for the models are special compositions of soft resin. Please handle them with the utmost care at all times.

Use vaseline

Please coat your finger(s) with lubricant gel when performing the training for digital examination.

Storage

Store the training set at room temperature, away from heat, moisture and direct sunlight.

DON'Ts

- Never wipe the simulator skin or pads with thinner or organic solvent.
- Do not mark on the phantom and models with pens or leave any printed materials in contact with the surface. Ink marks on the models are not removable.