








A complete range of DSI Sutures

THE BRAND YOU CAN TRUST

	 <p>Nylon Sutures Non-absorbable, Non-coated Monofilament of Polyamide SURGICAL SUTURES</p>	 <p>Silk Sutures Non-Absorbable Braided & Coated SURGICAL SUTURES</p>	 <p>Chromic Gut Absorbable Chromic Gut Natural Monofilament SURGICAL SUTURES</p>	 <p>Plain Catgut Absorbable Plain Catgut Natural Monofilament SURGICAL SUTURES</p>	 <p>PGA Sutures Absorbable, Braided & Coated of Polyglycolic Acid SURGICAL SUTURES</p>	 <p>PDO Sutures Synthetic Absorbable, Monofilament Polydioxanone SURGICAL SUTURES</p>
	Monofilament Polyamide	Braided & Coated Silk	Natural Monofilament	Natural Monofilament	Polyglycolic acid	Monofilament Polydioxanone
Colour	Blue	Black	Brown	Ivory	Violet	Violet
Absorbable	Non-absorbable	Non-absorbable	Absorbable	Absorbable	Absorbable	Absorbable
Wound support			18-21 days	12-15 days	21-28 days	50 days
Complete Absorption			50-70 days	40-50 days	60-75 days	180 days
Tensile strength	High	Very High	High	Medium	Very High	Very High
Tissue reaction	Minimal	Minimal	Low	Medium	Minimal	Very minimal
Wicking behaviour / Capillarity	No	Limited (coated)	Limited (coated)	Yes	Limited (coated)	No
Shape memory	Low	No	Yes	Yes	No	Low
Knot safety	Good	Excellent	Good	Good	Excellent	Very good
Smooth	Good	Very good	Low	Low	Very good	Very good
Sawing effect	No	Limited (coated)	Limited (coated)	Yes	Limited (coated)	No