


DSI Degranulation Bur Kit

Ref: DG-BUR-K



DSI Degranulation Kit are is used specifically to tear out all the remnants of soft tissue which are attached to the bone without harming the bone but quickly and effectively. The diamond burs are more precise and do not grab the bone. DSI burs permit optimal debridement of the inflammatory tissue from the bone which cannot be removed manually by a surgical spoon curette because it is penetrating into the bone and it is firmly attached to it.



The kit is supplied with four specially designed burs (1.0, 2.5, 3.0, 3.5 mm) a cleaning brush and a suction tube.

ADVANTAGES:

- Optimal debridement of the inflammatory tissue is achieved, without the use of surgical spoon curette.
- Multi-task Burs can be used in implantology, periodontal surgery, endodontic surgery, and root removal.
- Maximum safety when working near anatomical structures.
- Multiple use burs.