



| Practitioner: | | | | | | |
|---------------------|----------------------------|--|--|--|--|--|
| Facility: | | | | | | |
| Address: | | | | | | |
| City/ State/ Zip: _ | | | | | | |
| Phone: | | | | | | |
| Email: | | | | | | |
| Patient Name: | | | | | | |
| Left | Right Bilateral | | | | | |
| Clear Black | Inner Boot 1/8" 3/16" | | | | | |
| Valgus | F/H Metatarsal Heel Height | | | | | |

Malleo-Loc

Ankle Gauntlet

Code Options (When Clinically Justified)

L1950, L2250, L2280, L2820, L2755

The orthosis is comprised of an Ortho-Flex molded inner boot and a

Pre-Preg Carbon Fiber external structure.

The inner boot is adjustable and washable.

The external structure has two uprights that spiral across the ankle and a foot plate that extends to the

metatarsal heads.

The uprights have a unique "Oblique" orientation to the ankle. This orientation combined with Pre-Preg Carbon fiber provides a small amount of natural motion with a quickly increasing resistance.

Functional Benefit:

With this small amount of controlled motion the patient is supported while allowing just enough motion to Balance and walk with a reasonable gait.

The ramp-up of the resistance is the key.

The orthosis is:

Washable, Adjustable and Easy to Donn

Special Instructions