TENSION BASEPLATE KITS & ACCESSORY PARTS

Add 3D tissue culture, strain in high through-put plates, or equibiaxial or uniaxial tension to your Flexcell® Tension System.

- Increase your strain potential with an additional baseplate kit:
 - o BioFlex[®] kit for applying equibiaxial tension to cells in monolayer culture.
 - o HT BioFlex[®] kit for applying equibiaxial tension to cells in monolayer culture and in 24-well plates.
 - o Tissue Train[®] kit for creating 3D tissues.
 - o UniFlex® kit for applying uniaxial tension to cells in monolaver culture.
- All tension accessory parts such as baseplates, gaskets, and tubing can be replaced as needed. Parts are sold separately or as a part of the baseplate kits.
- Single plate baseplate kit includes one baseplate, 1 gasket, 1 loading station, 1 cell seeder, and two sample plates.



- 4 Gaskets
- **Corresponding Loading Stations**
- 4 Cell Seeders
- 4 Sample Plates
- Acrylic Window
- Software update, if needed

Figure 27. Single plate baseplate.

ORDERING INFORMATION (Please contact distributor for pricing)		
Catalog Number	Product/Item	Price
BFBK-4000	BioFlex® Baseplate Kit	
HTBK-4000	24-well HT BioFlex® Baseplate Kit	
TTBK-4000	Tissue Train® Baseplate Kit	
UFBK-4000	UniFlex® Baseplate Kit	
SPBK-1000	Single Plate Baseplate Kit	
Tension Baseplate Accessory Parts		
BPG-3000T	BioFlex® Baseplate and Gasket Set	
B-3000T	BioFlex® Baseplate Only	
G-3000T	BioFlex® Tension Gaskets Only (Set of 4)	
HTBPG-3000T	24-well HT Baseplate and Gasket Set	
HTB-3000T	24-well HT Baseplate Only	
HTG-3000T	24-well High Throughput Gaskets (Set of 4)	
TTB-3000T	Tissue Train® Baseplate Only	
TAK-4000	Tubing and Adaptor Kit	
LS-LUB	Loading Station™ Grease	
BP-3000W	Baseplate Acrylic Window	