

# Products

---

	Page
White Label .....	02
THC .....	03
CBD .....	06
THC:CBD .....	09
Primary packaging .....	12

---

## Valcon Medical A/S

Industrivej 12  
3540 Lyngø - Denmark

Phone: +45 22 95 80 19  
Mail: [sales@valcon-medical.com](mailto:sales@valcon-medical.com)

CVR: 40086714



# White Label

## What we offer

Valcon Medical offers supply of white label products to partners, who wish to sell and distribute medical cannabis oil under their own brand.

As your supplier, we will support you in the product approval process. This includes compilation of all regulatory documentation required for a smooth introduction of your white label products to the market.

Our promise is to deliver the highest quality of products while assuring consistent and scalable supply.

## Starting material

The starting material – Cannabis Sativa L – used for extraction and manufacturing is cultivated in accordance with Good Agricultural and Collection Practice (GACP) for Starting Materials on Herbal Origin and required provisions of the EU GMP guideline. Starting material is cultivated in Denmark by Schroll Medical to meet the requirements of the European Pharmacopoeia and the Danish Cannabis Flower Monograph.

## Compliance

White label products are manufactured under EU-GMP and tested in accordance with the German Monograph for cannabis extracts, DAB Eingestellter Cannabisextract, and the European Pharmacopoeia, Ph. Eur. We only apply well

established accredited laboratories performing the required range of analysis, including chemical, physical, and microbiological analysis of starting materials, intermediates and finished products. All according to validated analytical methods.

White label products are subject to stability according to ICH guidelines Q1A, Q1B, and EU-GMP guidelines.

## Packaging

Valcon Medical offers a selection of primary- and secondary packaging solutions. Whether you require a lid with- or without child resistance, or encapsulation for a dosing syringe, we can manufacture accordingly. The product can be delivered labeled or unlabeled and packaged in outer carton with dosing syringe and leaflet.

## GDP

Valcon Medical's supply chain complies with Good Distribution Practice (GDP). Qualified shipping solutions to customers include temperature monitoring according to individual product specifications.



# THC

Available now

---



## Full-spectrum THC cannabis oil

- Full-spectrum CO<sub>2</sub> – extracted
- 30 mL standardized oil solution
- THC: 25 mg/mL; CBD: <1 mg/mL  
(2.5% THC (w/w), <0.1% CBD (w/w))
- MCT carrier oil
- EU GACP/GMP cannabis floss input, origin Denmark

### Select preferred potency

We offer THC cannabis oil in one formulation of 25 mg/mL according to the table below. Special formulations can be ordered under manufacturing service request.

THC	X% w/w	X% w/w	X% w/w	2.5% w/w	X% w/w
	10 mg/mL	XX mg/mL	XX mg/mL	25 mg/mL	XX mg/mL
CBD	<0.1% w/w	<0.1% w/w	<0.1% w/w	<0.1% w/w	<0.1% w/w
	<1 mg/mL	<1 mg/mL	<1 mg/mL	<1 mg/mL	<1 mg/mL

Products are manufactured under EU-GMP and tested in accordance with the German Monograph for cannabis extracts, DAB Eingestellter Cannabisextract, and the European Pharmacopoeia, Ph. Eur. The product is composed of a soft extract from Cannabis Sativa L. 'White Widow' flower (floss) standardized with MCT (medium-chain triglycerides) oil to 25 mg/mL THC and <1 mg/mL CBD. Input: EU GACP/GMP cannabis floss, origin Denmark. The extraction solvent is carbon dioxide.



# THC

## Product Specification Table

<b>Product name</b>	To be determined by customer
<b>Product description</b>	Standardized full-spectrum cannabis oil
<b>Product form</b>	Oral fluid
<b>Product composition</b>	THC: 25 mg/mL; CBD <1 mg/mL (THC 2,5% w/w; CBD 0,1% w/w)
<b>Input material</b>	Cannabis Sativa L. (White Widow), Schroll Medical, Denmark
<b>Storage conditions</b>	2-8 °C or 15-25 °C
<b>Primary Packaging</b>	Amber glass bottle 30 mL, type III; Screw cap PP/HD-PE/LD-PE
<b>Product ID</b>	060518
<b>Product monograph</b>	DAB Eingestellter Cannabisextract, and the European Pharmacopoeia, Ph. Eur.
<b>Application</b>	For oral use
<b>Carrier</b>	MCT, CAS: 73398-61-5
<b>Extraction</b>	Carbon dioxide
<b>Shelf life / Hold life</b>	Stability according to ICH guidelines Q1A, Q1B, and EU-GMP
<b>Dosing</b>	To be determined by customer
<b>Secondary packaging Optional:</b>	<ol style="list-style-type: none"><li>1) Label with branding, product information, expiry and batch number placed on the primary packaging.</li><li>2) Dosing syringe/pipette separate or in outer carton.</li><li>3) Outer carton with printed product information, expiry and batch number.</li><li>4) Leaflet / preparation method information sheet placed in outer carton.</li></ol>



# CBD

Available Q1 2023

---



## Full-spectrum mixed CBD cannabis oil

- Full-spectrum CO<sub>2</sub> – extracted
  - 30 or 50 mL standardized oil solution
  - CBD: Formulation between of 2,5-10% w/w (25 to 100mg/mL); THC: less than 0,2% w/w (THC < 2 mg/mL)
  - MCT carrier oil
  - EU GACP/GMP cannabis floss input, origin Denmark
- 

### Select preferred potency

CBD cannabis oil is under development. Products can be formulated by request according to the range in the specification table below in a 30 or 50 mL bottle.

CBD	2.5% w/w	X% w/w	X% w/w	X% w/w	X% w/w	10% w/w
	25 mg/mL	XX mg/mL	XX mg/mL	XX mg/mL	XX mg/mL	100 mg/mL
THC	< 0.2% w/w	< 0.2% w/w	< 0.2% w/w	< 0.2% w/w	< 0.2% w/w	< 0.2% w/w
	< 2 mg/mL	< 2 mg/mL	< 2 mg/mL	< 2 mg/mL	< 2 mg/mL	< 2 mg/mL

---

Products are manufactured under EU-GMP and tested in accordance with the German Monograph for cannabis extracts, DAB Eingestellter Cannabisextract, and the European Pharmacopoeia, Ph. Eur. The product is composed of a soft extract from Cannabis Sativa L. “Nightingale” flower (floss) standardized with MCT oil (medium-chain triglycerides) to the required formulation. Input: EU GACP/GMP cannabis floss/origin Denmark. The extraction solvent is carbon dioxide.



# CBD

## Product Specification Table

<b>Product name</b>	To be determined by customer
<b>Product description</b>	Standardized full-spectrum cannabis oil
<b>Product form</b>	Oral fluid
<b>Product composition</b>	CBD: Formulation between of 2,5-10% w/w (25 to 100mg/mL); THC: less than 0,2% w/w (THC< 2 mg/mL)
<b>Input material</b>	Cannabis Sativa L. (Nightingale), Schroll Medical, Denmark
<b>Storage conditions</b>	2-8 °C or 15-25 °C
<b>Primary Packaging</b>	Amber glass bottle 30 mL, type III; Screw cap PP/HD-PE/LD-PE
<b>Product ID</b>	360118
<b>Product monograph</b>	DAB Eingestellter Cannabisextract, and the European Pharmacopoeia, Ph. Eur.
<b>Application</b>	For oral use
<b>Carrier</b>	MCT, CAS: 73398-61-5
<b>Extraction</b>	Carbon dioxide
<b>Shelf life / Hold life</b>	Stability according to ICH guidelines Q1A, Q1B, and EU-GMP
<b>Dosing</b>	To be determined by customer
<b>Secondary packaging Optional:</b>	<ol style="list-style-type: none"><li>1) Label with branding, product information, expiry and batch number placed on the primary packaging.</li><li>2) Dosing syringe/pipette separate or in outer carton.</li><li>3) Outer carton with printed product information, expiry and batch number.</li><li>4) Leaflet / preparation method information sheet placed in outer carton.</li></ol>





# THC:CBD

Available Q4 2022

---



## Full-spectrum mixed THC:CBD cannabis oil

- Full-spectrum CO<sub>2</sub> – extracted
- 30 mL standardized oil solution
- THC: Formulation between of 1-3% w/w (10 to 30mg/mL);  
CBD: Formulated between 1-3% w/w (10 to 30mg/mL).
- MCT carrier oil
- EU GACP/GMP cannabis floss input, origin Denmark

### Select preferred potency

THC:CBD cannabis oil can be formulated in 25 combinations as THC/CBD dominant or balanced in the preferred ratios. Potency options are listed in the specification table below.

THC	1% w/w	1.5% w/w	2% w/w	2.5% w/w	3% w/w
	10 mg/mL	15 mg/mL	20 mg/mL	25 mg/mL	30 mg/mL
CBD	1% w/w	1.5% w/w	2% w/w	2.5% w/w	3% w/w
	10 mg/mL	15 mg/mL	20 mg/mL	25 mg/mL	30 mg/mL

Products are manufactured under EU-GMP and tested in accordance with the German Monograph for cannabis extracts, DAB Eingestellter Cannabisextract, and the European Pharmacopoeia, Ph. Eur. The product is composed of a soft extract from Cannabis Sativa L. “White Widow” & “Nightingale” flower (floss) standardized with MCT oil (medium-chain triglycerides) to the required formulation. Input: EU GACP/GMP cannabis floss/origin Denmark. The extraction solvent is carbon dioxide.



# THC:CBD

## Product Specification Table

<b>Product name</b>	To be determined by customer
<b>Product description</b>	Standardized full-spectrum cannabis oil
<b>Product form</b>	Oral fluid
<b>Product composition</b>	THC: Formulation between of 1-3% w/w (10 to 30mg/mL); CBD: Formulated between 1-3% w/w (10 to 30mg/mL)
<b>Input material</b>	Cannabis Sativa L (White Widow & Nightingale), Schroll Medical, Denmark
<b>Storage conditions</b>	2-8 °C or 15-25 °C
<b>Primary Packaging</b>	Amber glass bottle 30 ml, type III; Screw cap PP/HD-PE/LD-PE
<b>Product ID</b>	060518/360118
<b>Product monograph</b>	DAB Eingestellter Cannabisextract, and the European Pharmacopoeia, Ph. Eur.
<b>Application</b>	For oral use
<b>Carrier</b>	MCT, CAS: 73398-61-5
<b>Extraction</b>	Carbon dioxide
<b>Shelf life / Hold life</b>	Stability according to ICH guidelines Q1A, Q1B, and EU-GMP
<b>Dosing</b>	To be determined by customer
<b>Secondary packaging Optional:</b>	<ol style="list-style-type: none"><li>1) Label with branding, product information, expiry and batch number placed on the primary packaging.</li><li>2) Dosing syringe/pipette separate or in outer carton.</li><li>3) Outer carton with printed product information, expiry, and batch number.</li><li>4) Leaflet / preparation method information sheet placed in outer carton.</li></ol>



# Primary packaging

Amber glass bottle 30 ml,  
type III



Valcon Medical perform technical release of secondary packaging material compliant with all applicable laws and regulations and approved by customer.



**Lid 1** ID: PM-0089

- Tamper evident
  - Child resistant
- Screw cap  
PP/HD-PE/LD-PE



**Lid 2** ID: PM-0070

- Tamper evident
  - Child resistant
  - Insert for dosing pipette
- Screw cap  
PP/HD-PE/LD-PE

Insert opening for dosing pipette Ø4,8mm





## Lid 3: ID: PM-0113

---

- Tamper evident
- Insert for dosing pipette

Screw cap

PP/HD-PE/LD-PE

Insert opening for dosing pipette Ø4,8mm



## Dosing pipette:

---

- 1 mL oral dosing syringe
- Fits opening Ø4,9mm.
- Individually packed

Usage: Dosage and administration of liquid medicine, intended for human use.

Type: Medical device, class I with measuring function according to EU Directive 93/42/CEE.

Body: PP/ Polypropylene, natural.

Plunger: Low-density Polyethylene, natural.

